

STELLAR

45° of tilt manual



- Ideal for hand propulsion with accessible rear wheels
- Knees remain low in tilt—wheel under desks, into vans and access other environments with ease
- Short wheelbase provides excellent maneuverability
- Features adjustable seat height, seat depth and back angle
- Available in Heavy Duty (HD) extra-wide sizes
- Transit-ready option
- ISO tested



Stellar HD



WHEELCHAIRS FOR FUNCTION AND INDEPENDENCE

WWW.PDGMOBILITY.COM



Nine standard frame colors to choose from and custom colors also available. Contact us for details.



Multi-adjustable



Adjustable push bar option



Easily adjustable axle position



Actuator adjusts for optimum size



HD option

Seat width (standard) 14" – 22"
Seat width (HD) 22" – 32"
Seat depth (adjustable) 16" – 22"
Seat height (adjustable) 14" – 20"
Back angle (adjustable) 90° – 120°
Back height 18" – 26"
Dynamic tilt 45°

Overall width
 12" & 16" rear wheels seat width + 9.5"
 20" & 24" rear wheels seat width + 10.5"

Overall length 30" – 37" (without front rigging)
Weight capacity 250 lbs, HD to 450 lbs
Chair weight 40 lbs (standard); 50 lbs (HD)
Transit option available ANSI/RESNA WC-19 (Standard frame only)
Warranty PDG's exclusive TGIF warranty (Thank Goodness It Fits) Lifetime warranty on the frame.

PDG recommends all wheelchairs be fitted to the consumer under the guidance of a health care professional knowledgeable of the particular needs of the consumer and the limitations of the product. US Patent # 6 447064 B1

USA
 PDG Product Design Group
 PMB # 173 - 1124 Fir Avenue
 Blaine, WA
 98230, USA

CANADA
 PDG Product Design Group
 Unit 103 - 318 East Kent Avenue South
 Vancouver, BC
 V5X 4N6

T: +1 604 323 9220
 F: +1 604 323 9097
 Toll Free (in N. America): 888-858-4422
 Toll Free Fax (in N. America): 877-323-9093
www.pdgmobility.com