

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877-323-9093.

### Order Notes


### Frame Type

- |   |         |                                       |          |
|---|---------|---------------------------------------|----------|
| <input checked="" type="checkbox"/>     | 2101000 | Stellar HD 45° Manual Tilt Wheelchair | \$ 3,995 |
| Extended delivery lead times may apply. |         |                                       |          |

### Frame Options

#### Heavy Duty Base and Seat Frame

- |                          |         |                                      |        |
|--------------------------|---------|--------------------------------------|--------|
| <input type="checkbox"/> | 2102045 | 450 lbs Weight Capacity. 45° of Tilt | \$ 600 |
| <input type="checkbox"/> | 2102060 | 600 lbs Weight Capacity. 25° of Tilt | \$ 900 |

#### User Weight Range

Required for selection of Gas Struts.

- |                          |         |  |  |
|--------------------------|---------|--|--|
| <input type="checkbox"/> | 2102212 | Up to 250 lbs (2 x 1200 Newton Gas Struts)     |  |
| <input type="checkbox"/> | 2102308 | 250 lbs - 350 lbs (3 x 800 Newton Gas Struts)  |  |
| <input type="checkbox"/> | 2102310 | 350 lbs - 425 lbs (3 x 1000 Newton Gas Struts) |  |
| <input type="checkbox"/> | 2102312 | 425 lbs - 600 lbs (3 x 1200 Newton Gas Struts) |  |

### Seat Width

#### Standard Seat Widths

- |                          |         |                |        |
|--------------------------|---------|----------------|--------|
| <input type="checkbox"/> | 2103022 | 22" Seat Width | \$ 325 |
| <input type="checkbox"/> | 2103024 | 24" Seat Width | \$ 325 |
| <input type="checkbox"/> | 2103026 | 26" Seat Width | \$ 325 |
| <input type="checkbox"/> | 2103028 | 28" Seat Width | \$ 625 |
| <input type="checkbox"/> | 2103030 | 30" Seat Width | \$ 625 |
| <input type="checkbox"/> | 2103032 | 32" Seat Width | \$ 625 |

### Seat Depth

#### Standard Seat Depths (Adjustable)

- |                          |         |                |        |
|--------------------------|---------|----------------|--------|
|                          |         |                | \$ 325 |
| <input type="checkbox"/> | 2104018 | 18" Seat Depth |        |
| <input type="checkbox"/> | 2104019 | 19" Seat Depth |        |
| <input type="checkbox"/> | 2104020 | 20" Seat Depth |        |
| <input type="checkbox"/> | 2104021 | 21" Seat Depth |        |
| <input type="checkbox"/> | 2104022 | 22" Seat Depth |        |

### Frame Colour

#### PDG Primary Colours

- |                          |         |                    |
|--------------------------|---------|--------------------|
| <input type="checkbox"/> | 2105000 | Tracer Red         |
| <input type="checkbox"/> | 2105001 | Gloss Black        |
| <input type="checkbox"/> | 2105002 | Silver             |
| <input type="checkbox"/> | 2105003 | Black Shimmer      |
| <input type="checkbox"/> | 2105004 | Orange Fusion      |
| <input type="checkbox"/> | 2105005 | Karl's Imagination |
| <input type="checkbox"/> | 2105006 | Blue Shimmer       |
| <input type="checkbox"/> | 2105007 | Gloss White        |
| <input type="checkbox"/> | 2105008 | Dark Titanium      |

#### PDG Premium Colours

\$ 150

- |                          |         |               |
|--------------------------|---------|---------------|
| <input type="checkbox"/> | 2105009 | Geysir Green  |
| <input type="checkbox"/> | 2105010 | Lollipop Blue |
| <input type="checkbox"/> | 2105011 | Lollipop Red  |

#### Custom Colour

\$ 300

Contact PDG for specifications.

### Arm Type

- |                          |         |                                       |        |
|--------------------------|---------|---------------------------------------|--------|
| <input type="checkbox"/> | 2106000 | No Armrest                            |        |
| <input type="checkbox"/> | 2106214 | Dual Post Adjustable Height (10"-14") | \$ 255 |

### Armrests

#### Regular Armpads

- |                          |         |                               |  |
|--------------------------|---------|-------------------------------|--|
| <input type="checkbox"/> | 2107100 | Desk Length - 2" x 10" (pair) |  |
| <input type="checkbox"/> | 2107101 | Full Length - 2" x 14" (pair) |  |

#### Wide Armpads

\$ 115

May interfere with wheels in some tilt configurations.

- |                          |         |                               |  |
|--------------------------|---------|-------------------------------|--|
| <input type="checkbox"/> | 2107110 | Desk Length - 3" x 10" (pair) |  |
| <input type="checkbox"/> | 2107111 | Full Length - 3" x 12" (pair) |  |

#### Gel Armpads

\$ 115

May interfere with wheels in some tilt configurations.

- |                          |         |                                 |  |
|--------------------------|---------|---------------------------------|--|
| <input type="checkbox"/> | 2107120 | Full Length - 2" x 12" (pair)   |  |
| <input type="checkbox"/> | 2107121 | Full Length - 3.5" x 12" (pair) |  |

#### Contour Arm Trough

\$ 125

May interfere with wheels in some tilt configurations.

- |                          |         |                                  |  |
|--------------------------|---------|----------------------------------|--|
| <input type="checkbox"/> | 2107130 | Right (4.5" x 18") with Hand Pad |  |
| <input type="checkbox"/> | 2107131 | Left (4.5" x 18") with Hand Pad  |  |

#### Clothing Guards

\$ 165

### Seat Pan

- |                                     |         |                |        |
|-------------------------------------|---------|----------------|--------|
| <input checked="" type="checkbox"/> | 1208001 | Solid Seat Pan | \$ 120 |
|-------------------------------------|---------|----------------|--------|

Cushion		
<input type="checkbox"/>	2109000	<b>No Cushion</b>
<b>3" Basic Foam (Nylon Cover)</b>		
<input type="checkbox"/>	2109200	Up to 28" Width \$ 484
<input type="checkbox"/>	2109201	Above 28" Width \$ 605
<b>3" Contour Foam (Nylon Cover)</b>		
<input type="checkbox"/>	2109210	Up to 28" Width \$ 836
<input type="checkbox"/>	2109211	Above 28" Width \$ 1,089
<b>3" Contour Pressure Reduction (Nylon Cover)</b>		
<input type="checkbox"/>	2109220	Up to 28" Width \$ 1,485
<input type="checkbox"/>	2109221	Above 28" Width \$ 1,903
<b>Spare Cushion Cover</b>		
<input type="checkbox"/>	1109230	Basic \$ 391
<input type="checkbox"/>	1109231	Contour \$ 490
<b>Cushion Size</b>		
	2109400	Specify Size (Width" x Depth") _____

Seat Height	
<b>Seat to Floor Height</b>	
<input type="checkbox"/>	2110016 Set at 16" STF
<input type="checkbox"/>	2110017 Set at 17" STF
<input type="checkbox"/>	2110018 Set at 18" STF

Casters	
<input checked="" type="checkbox"/>	2111050 <b>5" Caster With Heavy Duty Steel Fork</b>

Rear Wheels	
<input checked="" type="checkbox"/>	2112224 <b>24" Heavy Duty Spoke</b> \$ 350

Axle Type	
<input checked="" type="checkbox"/>	2113010 <b>Fixed Rear Axles</b>

Rear Tires	
<input checked="" type="checkbox"/>	2114000 <b>Urethane</b>

Hand Rims	
<input checked="" type="checkbox"/>	2115010 <b>Standard</b>

Wheel Locks	
<input type="checkbox"/>	2116010 <b>Push to Lock</b>
<input type="checkbox"/>	2116020 <b>Pull to Lock</b>
<input type="checkbox"/>	2116100 <b>Wheel Lock Extensions</b> \$ 60

Anti-Tippers	
<input checked="" type="checkbox"/>	2117010 <b>Rear Anti-Tippers</b> \$ 60 Required Safety Option.



Scan this code or visit [www.pdgmobility.com](http://www.pdgmobility.com) to learn more about the Stellar HD Manual Tilt Wheelchair.

Back Post Style	
<input checked="" type="checkbox"/>	2118010 <b>Push Bar Back</b>
<b>Back Post Height</b>	
Seat Pan to Handle	
<input type="checkbox"/>	2118425 25" Height
<input type="checkbox"/>	2118431 31" Height

Back Options	
<input checked="" type="checkbox"/>	2119090 <b>Seat to Back Angle set at 90°</b>
<input type="checkbox"/>	2119240 <b>Dynamic Rocker Back Assembly</b> \$ 2,995 Not compatible with 30° Recline. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller.

Back Upholstery	
<input type="checkbox"/>	2120000 <b>No Back Upholstery</b>
<b>Adjustable Tension</b> \$ 250	
Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	2120321 Set at 21"
<input type="checkbox"/>	2120322 Set at 22"
<input type="checkbox"/>	2120323 Set at 23"

Backrest	
<input checked="" type="checkbox"/>	2121000 <b>No Backrest</b>

Tilt/Recline Options	
<b>Attendant Tilt only</b>	
<input type="checkbox"/>	2123011 Installed on Right Back Post
<input type="checkbox"/>	2123012 Installed on Left Back Post
<b>Tilt Limit</b>	
<input type="checkbox"/>	2123220 Set at 20°
<input type="checkbox"/>	2123225 Set at 25° (Maximum for 600 lbs chair)
<input type="checkbox"/>	2123230 Set at 30°
<input type="checkbox"/>	2123235 Set at 35°
<input type="checkbox"/>	2123240 Set at 40°
<input type="checkbox"/>	2123245 Set at 45°
<input type="checkbox"/>	2123410 <b>Reclining Backrest (Manual)</b> \$ 830 Maximum recline = 30°. Not compatible with Tension Adjustable Upholstery and 20"-22" Seat Depth.

Headrest		
<input type="checkbox"/>	2124000	<b>No Headrest</b>
<input type="checkbox"/>	2124010	<b>Headrest Mounting Plate</b> \$ 75 <small>Not required if mounting to a solid Back.</small>
<input type="checkbox"/>	2124100	<b>nxt MultiFit Head Support Hardware</b> \$ 155
		<b>nxt MultiFit Head Support Pad</b> \$ 195
<input type="checkbox"/>	2124105	Regular (9" x 5")
<input type="checkbox"/>	2124106	Large (10" x 5")
<input type="checkbox"/>	2124107	X-Large (14" x 5")

Front Rigging		
<input checked="" type="checkbox"/>	2125001	<b>Sector Blocks only</b>
<input type="checkbox"/>	2125130	<b>Swing Away 60°</b>
<input type="checkbox"/>	2125200	<b>Swing Away Elevating Legrest</b> \$ 355

Foot Plates		
<input checked="" type="checkbox"/>	2126020	<b>Angle Adjustable</b> \$ 140

Footrest/Legrest Accessories		
<input type="checkbox"/>	2127030	<b>Swing Away ELR Gel Knee Pad</b> \$ 100
<input type="checkbox"/>	2127040	<b>Adjustable Length Heel Loops</b> \$ 50

Extension Tubes		
<input checked="" type="checkbox"/>	2128010	<b>Regular Extension Tubes (13"-17")</b>

Positioning Straps		
<b>Bodypoint - Hip Belts</b>		
<input type="checkbox"/>	2129010	Center-Pull Auto Style \$ 105
<input type="checkbox"/>	2129015	Center-Pull Auto Style, Padded \$ 155

Miscellaneous Options		
<input type="checkbox"/>	2130120	<b>O2 Holder</b> \$ 295
<input type="checkbox"/>	2130170	<b>Adhesive Back Hook Strip, 2" x 30"</b> \$ 10 <small>For Cushion mounting.</small>