

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877-323-9093.

### Order Notes


### Frame Type

- |   |  |              |          |
|---|--|--------------|----------|
| <input checked="" type="checkbox"/> 2101000 | <b>Stellar HD 45° Manual Tilt Wheelchair</b> | <i>E1161</i> | \$ 2,995 |
| Extended delivery lead times may apply.     |  |              |          |

### Frame Options

- |   |  |  |        |
|---|--|--|--------|
| <b>Heavy Duty Base and Seat Frame</b> <i>K0108</i>                |  |  |        |
| <input type="checkbox"/> 2102045                                  | 450 lbs Weight Capacity. 45° of Tilt           |  | \$ 600 |
| <input type="checkbox"/> 2102060                                  | 600 lbs Weight Capacity. 25° of Tilt           |  | \$ 900 |
| <b>User Weight Range</b><br>Required for selection of Gas Struts. |  |  |        |
| <input type="checkbox"/> 2102212                                  | Up to 250 lbs (2 x 1200 Newton Gas Struts)     |  |        |
| <input type="checkbox"/> 2102308                                  | 250 lbs - 350 lbs (3 x 800 Newton Gas Struts)  |  |        |
| <input type="checkbox"/> 2102310                                  | 350 lbs - 425 lbs (3 x 1000 Newton Gas Struts) |  |        |
| <input type="checkbox"/> 2102312                                  | 425 lbs - 600 lbs (3 x 1200 Newton Gas Struts) |  |        |

### Seat Width

- |                                  |                |              |        |
|----------------------------------|----------------|--------------|--------|
| <b>Standard Seat Widths</b>      |                |              |        |
| <input type="checkbox"/> 2103022 | 22" Seat Width | <i>E2201</i> | \$ 325 |
| <input type="checkbox"/> 2103024 | 24" Seat Width | <i>E2202</i> | \$ 525 |
| <input type="checkbox"/> 2103026 | 26" Seat Width | <i>E2202</i> | \$ 525 |
| <input type="checkbox"/> 2103028 | 28" Seat Width | <i>K0108</i> | \$ 925 |
| <input type="checkbox"/> 2103030 | 30" Seat Width | <i>K0108</i> | \$ 925 |
| <input type="checkbox"/> 2103032 | 32" Seat Width | <i>K0108</i> | \$ 925 |

### Seat Depth

- |  |                |              |        |
|--|----------------|--------------|--------|
| <b>Standard Seat Depths (Adjustable)</b> |                |              |        |
| <input type="checkbox"/> 2104018         | 18" Seat Depth |              | \$ 325 |
| <input type="checkbox"/> 2104019         | 19" Seat Depth |              | \$ 325 |
| <input type="checkbox"/> 2104020         | 20" Seat Depth | <i>E2203</i> | \$ 475 |
| <input type="checkbox"/> 2104021         | 21" Seat Depth | <i>E2203</i> | \$ 475 |
| <input type="checkbox"/> 2104022         | 22" Seat Depth | <i>E2204</i> | \$ 750 |

### Frame Color

#### PDG Primary Colors

- 2105000 Tracer Red
- 2105001 Gloss Black
- 2105002 Silver
- 2105003 Black Shimmer
- 2105004 Orange Fusion
- 2105005 Karl's Imagination
- 2105006 Blue Shimmer
- 2105007 Gloss White
- 2105008 Dark Titanium

#### PDG Premium Colors

- 2105009 Geyser Green
- 2105010 Lollipop Blue
- 2105011 Lollipop Red

#### 2105200 Custom Color

Contact PDG for specifications.

\$ 150

\$ 300

### Arm Type

#### 2106000 No Armrest

- 2106214 **Dual Post Adjustable Height (10"-14")** *E0973*

\$ 255

### Armrests

#### Regular Armpads

- 2107100 Desk Length - 2" x 10" (pair)
- 2107101 Full Length - 2" x 14" (pair)

#### Wide Armpads

May interfere with wheels in some tilt configurations.

- 2107110 Desk Length - 3" x 10" (pair)
- 2107111 Full Length - 3" x 12" (pair)

\$ 115

#### Gel Armpads

May interfere with wheels in some tilt configurations.

- 2107120 Full Length - 2" x 12" (pair)
- 2107121 Full Length - 3.5" x 12" (pair)

\$ 115

#### Contour Arm Trough *E2209*

May interfere with wheels in some tilt configurations.

- 2107130 Right (4.5" x 18") with Hand Pad
- 2107131 Left (4.5" x 18") with Hand Pad

\$ 125

- 2107320 **Clothing Guards** *K0108*

\$ 165

### Seat Pan

- 1208001 **Solid Seat Pan** *E2231*

\$ 120

Cushion			
<input type="checkbox"/>	2109000	<b>No Cushion</b>	
<b>Matrx VI Heavy Duty</b>			
<input type="checkbox"/>	2109115	18"-30" width, 18"-22" depth <i>E2608</i>	\$ 1,150
<b>3" Basic Foam (Nylon Cover) <i>K0108</i></b>			
<input type="checkbox"/>	2109200	Up to 28" Width	\$ 485
<input type="checkbox"/>	2109201	Above 28" Width	\$ 605
<b>3" Contour Foam (Nylon Cover) <i>K0108</i></b>			
<input type="checkbox"/>	2109210	Up to 28" Width	\$ 835
<input type="checkbox"/>	2109211	Above 28" Width	\$ 1,090
<b>3" Contour Pressure Reduction (Nylon Cover) <i>K0108</i></b>			
<input type="checkbox"/>	2109220	Up to 28" Width	\$ 1,485
<input type="checkbox"/>	2109221	Above 28" Width	\$ 1,900
<b>Spare Cushion Cover (Nylon)</b>			
<input type="checkbox"/>	2109230	Basic	\$ 390
<input type="checkbox"/>	2109231	Contour	\$ 490
<b>Cushion Size</b>			
	2109400	Specify Size (Width" x Depth") _____	

Seat Height	
<b>Seat to Floor Height</b>	
<input type="checkbox"/>	2110016 Set at 16" STF
<input type="checkbox"/>	2110017 Set at 17" STF
<input type="checkbox"/>	2110018 Set at 18" STF

Casters	
<input checked="" type="checkbox"/>	2111050 <b>5" Caster With Heavy Duty Steel Fork</b>

Rear Wheels	
<input checked="" type="checkbox"/>	2112224 <b>24" Heavy Duty Spoke</b> \$ 350

Axle Type	
<input checked="" type="checkbox"/>	2113010 <b>Fixed Rear Axles</b>

Rear Tires	
<input checked="" type="checkbox"/>	2114000 <b>Urethane</b>

Hand Rims	
<input checked="" type="checkbox"/>	2115010 <b>Standard</b>

Wheel Locks	
<input type="checkbox"/>	2116010 <b>Push to Lock</b>
<input type="checkbox"/>	2116020 <b>Pull to Lock</b>
<input type="checkbox"/>	2116100 <b>Wheel Lock Extensions</b> <i>E0961</i> \$ 60

Anti-Tippers	
<input checked="" type="checkbox"/>	2117010 <b>Rear Anti-Tippers</b> <i>E0971</i> \$ 60 Required Safety Option.



Scan this code or visit [www.pdgmobility.com](http://www.pdgmobility.com) to learn more about the Stellar HD Manual Tilt Wheelchair.

Back Post Style	
<input checked="" type="checkbox"/>	2118010 <b>Push Bar Back</b>
<b>Back Post Height</b> Seat Pan to Handle	
<input type="checkbox"/>	2118425 25" Height
<input type="checkbox"/>	2118431 31" Height

Back Options	
<input checked="" type="checkbox"/>	2119090 <b>Seat to Back Angle set at 90°</b>
<input type="checkbox"/>	2119240 <b>Dynamic Rocker Back Assembly</b> <i>K0108</i> \$ 2,995 Not compatible with 30° Recline. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller.

Back Upholstery	
<input type="checkbox"/>	2120000 <b>No Back Upholstery</b>
<b>Adjustable Tension</b> \$ 350 Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	2120321 Set at 21"
<input type="checkbox"/>	2120322 Set at 22"
<input type="checkbox"/>	2120323 Set at 23"

Backrest	
<input type="checkbox"/>	2121000 <b>No Backrest</b>
<b>Matrx Elite</b>	
<input type="checkbox"/>	2121112 HD - 22"-26" width, 16"-20" height <i>E2621</i> \$ 1,095

Tilt/Recline Options	
<b>Attendant Tilt only</b>	
<input type="checkbox"/>	2123011 Installed on Right Back Post
<input type="checkbox"/>	2123012 Installed on Left Back Post
<b>Tilt Limit</b>	
<input type="checkbox"/>	2123220 Set at 20°
<input type="checkbox"/>	2123225 Set at 25° (Maximum Tilt for 600 lbs chair)
<input type="checkbox"/>	2123230 Set at 30°
<input type="checkbox"/>	2123235 Set at 35°
<input type="checkbox"/>	2123240 Set at 40°
<input type="checkbox"/>	2123245 Set at 45°
<input type="checkbox"/>	2123410 <b>Reclining Backrest (Manual)</b> <i>E1225</i> \$ 830 Maximum recline = 30°. Not compatible with Tension Adjustable Upholstery and 20"-22" Seat Depth.

Headrest			
<input type="checkbox"/>	2124000	<b>No Headrest</b>	
<input type="checkbox"/>	2124010	<b>Headrest Mounting Plate</b>	\$ 75
Not required if mounting to a solid Back.			
<input type="checkbox"/>	2124100	<b>nxt MultiFit Head Support Hardware</b> <i>E1028</i>	\$ 155
		<b>nxt MultiFit Head Support Pad</b> <i>E0955</i>	\$ 195
<input type="checkbox"/>	2124105	Regular (9" x 5")	
<input type="checkbox"/>	2124106	Large (10" x 5")	
<input type="checkbox"/>	2124107	X-Large (14" x 5")	
		<b>Stealth Headrest Hardware</b> <i>E1028</i>	
<input type="checkbox"/>	2124110	Multi-Axis HMO475	\$ 150
<input type="checkbox"/>	2124111	Link TWB480	\$ 325
		<b>Stealth Headrest Pad</b> <i>E0955</i>	\$ 200
<input type="checkbox"/>	2124115	10" width	
<input type="checkbox"/>	2124116	14" width	

Front Rigging			
<input type="checkbox"/>	2125001	<b>Sector Blocks only</b>	
<input type="checkbox"/>	2125130	<b>Swing Away 60°</b>	
<input type="checkbox"/>	2125200	<b>Swing Away Elevating Legrest</b> <i>E0990</i>	\$ 355

Foot Plates			
<input checked="" type="checkbox"/>	2126020	<b>Angle Adjustable</b> <i>K0040</i>	\$ 140

Footrest/Legrest Accessories			
<input type="checkbox"/>	2127030	<b>Swing Away ELR Gel Knee Pad</b>	\$ 100
<input type="checkbox"/>	2127040	<b>Adjustable Length Heel Loops</b> <i>E0951</i>	\$ 50

Extension Tubes			
<input checked="" type="checkbox"/>	2128010	<b>Regular Extension Tubes (13"-17")</b>	

Positioning Straps			
		<b>Bodypoint - Hip Belt</b> <i>E0978</i>	
<input type="checkbox"/>	2129010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	2129015	Center-Pull Auto Style, Padded	\$ 155

Miscellaneous Options			
<input type="checkbox"/>	2130120	<b>O2 Holder</b> <i>E2208</i>	\$ 295
<input type="checkbox"/>	2130170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	\$ 10
For Cushion mounting.			