



PDG Product Design Group Inc. | 318 East Kent Ave South | Vancouver BC V5X 4N6, Canada | Phone 604-323-9220 | 888-858-4422 | Email info@pdgmobility.com

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

### Frame Type

Stellar LEAP Manual Posterior & Anterior Tilt Wheelchair		
<input checked="" type="checkbox"/>	2201000 LEAP Manual Tilt Wheelchair Base Wheelchair with up to 20° of Posterior Tilt.	\$ 2,995
<input type="checkbox"/>	2201101 LEAP Anterior Tilt with Linked Recline Up to 30° of Anterior Tilt.	\$ 1,250
<input type="checkbox"/>	2201102 Dynamic Back Recline Up to 30° of Recline.	\$ 500

### Frame Options

<input checked="" type="checkbox"/>	2202000 <b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.
<input type="checkbox"/>	2202300 <b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

### Seat Width

Standard Seat Widths		
<input type="checkbox"/>	2203014 14" Seat Width	
<input type="checkbox"/>	2203015 15" Seat Width	
<input type="checkbox"/>	2203016 16" Seat Width	
<input type="checkbox"/>	2203017 17" Seat Width	
<input type="checkbox"/>	2203018 18" Seat Width	
<input type="checkbox"/>	2203019 19" Seat Width	
<input type="checkbox"/>	2203020 20" Seat Width <i>E2201</i>	\$ 325
<input type="checkbox"/>	2203021 21" Seat Width <i>E2201</i>	\$ 325
<input type="checkbox"/>	2203022 22" Seat Width <i>E2201</i>	\$ 325

### Seat Depth

Standard Seat Depths (Adjustable)		
<input type="checkbox"/>	2204016 16" Seat Depth	
<input type="checkbox"/>	2204017 17" Seat Depth	
<input type="checkbox"/>	2204018 18" Seat Depth	
<input type="checkbox"/>	2204019 19" Seat Depth	
<input type="checkbox"/>	2204020 20" Seat Depth <i>E2203</i>	\$ 475
Custom Seat Depths		
<input type="checkbox"/>	2204021 21" Seat Depth <i>E2203</i>	\$ 475
<input type="checkbox"/>	2204022 22" Seat Depth <i>E2204</i>	\$ 750

### Frame Color

PDG Primary Colors		
<input type="checkbox"/>	2205000 Tracer Red	
<input type="checkbox"/>	2205001 Gloss Black	
<input type="checkbox"/>	2205002 Silver	
<input type="checkbox"/>	2205003 Black Shimmer	
<input type="checkbox"/>	2205004 Orange Fusion	
<input type="checkbox"/>	2205005 Karl's Imagination	
<input type="checkbox"/>	2205006 Blue Shimmer	
<input type="checkbox"/>	2205007 Gloss White	
<input type="checkbox"/>	2205008 Dark Titanium	
PDG Premium Colors		\$ 150
<input type="checkbox"/>	2205009 Geyser Green	
<input type="checkbox"/>	2205010 Lollipop Blue	
<input type="checkbox"/>	2205011 Lollipop Red	
<input type="checkbox"/>	2205200 <b>Custom Color</b>	\$ 300
Contact PDG for specifications.		

### Arm Type

<input type="checkbox"/>	2206000 <b>No Armrest</b>	
<input type="checkbox"/>	2206100 <b>Single Post Fixed Height (12")</b> Full Length Armpad only.	
	<b>Single Post Adjustable Height</b> <i>E0973</i>	\$ 255
<input type="checkbox"/>	2206111 8"-11"	
<input type="checkbox"/>	2206114 10"-14"	
<input type="checkbox"/>	2206200 <b>Dual Post Fixed Height (12")</b> Full Length Armpad only.	
	<b>Dual Post Adjustable Height</b> <i>E0973</i>	\$ 255
<input type="checkbox"/>	2206211 8"-11"	
<input type="checkbox"/>	2206214 10"-14"	
<input type="checkbox"/>	2206216 12"-16"	
	<b>Locking Flip-up</b> <i>E0973</i>	\$ 255
Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.		
<input type="checkbox"/>	2206310 Set at 10"	
<input type="checkbox"/>	2206311 Set at 11"	
<input type="checkbox"/>	2206312 Set at 12"	
<input type="checkbox"/>	2206313 Set at 13"	
<input type="checkbox"/>	2206314 Set at 14"	

Armrests	
<b>Regular Armrests</b>	
<input type="checkbox"/> 2207100	Desk Length - 2" x 10" (pair)
<input type="checkbox"/> 2207101	Full Length - 2" x 14" (pair)
<b>Wide Armrests</b> \$ 115	
May interfere with wheels in some tilt configurations.	
<input type="checkbox"/> 2207110	Desk Length - 3" x 10" (pair)
<input type="checkbox"/> 2207111	Full Length - 3" x 12" (pair)
<b>Gel Armrests</b> \$ 115	
May interfere with wheels in some tilt configurations.	
<input type="checkbox"/> 2207120	Full Length - 2" x 12" (pair)
<input type="checkbox"/> 2207121	Full Length - 3.5" x 12" (pair)
<b>Contour Arm Trough</b> E2209 \$ 125	
May interfere with wheels in some tilt configurations.	
<input type="checkbox"/> 2207130	Right (4.5" x 18") with Hand Pad
<input type="checkbox"/> 2207131	Left (4.5" x 18") with Hand Pad
<b>Armrest Extensions</b> \$ 60	
Not compatible with Flip-up or Single Post Armrests.	
<input type="checkbox"/> 2207400	3" Extension - Forward
<input type="checkbox"/> 2207401	3" Extension - Rearward

Seat Pan	
<input type="checkbox"/> 2208000	<b>No Seat Pan</b>
<input type="checkbox"/> 2208001	<b>Solid Seat Pan</b> E2231 \$ 120

Cushion	
<input type="checkbox"/> 2209000	<b>No Cushion</b>
<b>nxt BioFit</b>	
<input type="checkbox"/> 2209101	14"-20" width, 16"-20" depth E2607 \$ 550
<input type="checkbox"/> 2209102	22" width, 16"-20" depth E2608 \$ 600
<b>Matrx VI</b>	
<input type="checkbox"/> 2209111	14"-20" width, 16"-22" depth E2607 \$ 550
<input type="checkbox"/> 2209112	22" width, 18"-22" depth E2608 \$ 600
<input type="checkbox"/> 2209115	HD - 18"-22" width, 18"-22" depth E2608 \$ 1,150
<b>ROHO Quadtro Select HP</b>	
<input type="checkbox"/> 2209121	14"-20" width, 16"-22" depth E2624 \$ 650
<input type="checkbox"/> 2209122	22" width, 16"-22" depth E2625 \$ 750
<b>ROHO High Profile</b>	
<input type="checkbox"/> 2209125	14"-20" width, 16"-22" depth E2622 \$ 650
<input type="checkbox"/> 2209126	22" width, 16"-22" depth E2623 \$ 750
<b>Varilite Reflex</b>	
<input type="checkbox"/> 2209131	14"-20" width, 16"-20" depth E2603 \$ 375
<b>Cushion Size</b>	
2209400	Specify Size (Width" x Depth") _____



Scan this code or visit [www.pdgmobility.com](http://www.pdgmobility.com) to learn more about the **STELLAR LEAP** Manual Anterior and Posterior Tilt Wheelchair

Seat Height																																																												
Circle desired Seat-to-Floor Height or enter information below:																																																												
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>20"</td> <td>13"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td></td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td></td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> </tbody> </table>	RW	6" Casters						12"	13"	14"	15"	16"			16"			15"	16"	17"	18"	20"	13"	14"	15"	16"			22"		14"	15"	16"	17"		24"			15"	16"	17"	18"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="4">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> <tr> <td>24"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> </tr> </tbody> </table>	RW	8" Casters				16"		15"	16"	17"	18"	24"		15"	16"	17"	18"
RW	6" Casters																																																											
12"	13"	14"	15"	16"																																																								
16"			15"	16"	17"	18"																																																						
20"	13"	14"	15"	16"																																																								
22"		14"	15"	16"	17"																																																							
24"			15"	16"	17"	18"																																																						
RW	8" Casters																																																											
16"		15"	16"	17"	18"																																																							
24"		15"	16"	17"	18"																																																							
Wheel size _____ Caster size _____ STF _____																																																												
STF Height is adjustable within the range of the selected wheel size.																																																												

Casters	
<b>Urethane</b>	
8" Casters are compatible only with 16" and 24" Rear Wheels.	
<input type="checkbox"/> 2211061	6" x 1.25" Caster size
<input type="checkbox"/> 2211063	6" x 1.50" Caster size \$ 85
<input type="checkbox"/> 2211081	8" x 1.25" Caster size
<input type="checkbox"/> 2211084	8" x 2.00" Caster size \$ 135
<b>Pneumatic</b>	
8" Casters are compatible only with 16" and 24" Rear Wheels.	
<input type="checkbox"/> 2211161	6" x 1.25" Caster size \$ 110
<input type="checkbox"/> 2211181	8" x 1.25" Caster size \$ 120
<input type="checkbox"/> 2211184	8" x 2.00" Caster size \$ 150
<input type="checkbox"/> 2211260	<b>6" Soft Roll</b> \$ 95
<input type="checkbox"/> 2211400	<b>Ultra Sport Fork</b> E1015 \$ 395
Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	

Rear Wheels	
<b>Composite</b>	
<input type="checkbox"/> 2212012	12" Wheel size
<input type="checkbox"/> 2212016	16" Wheel size
<input type="checkbox"/> 2212020	20" Wheel size
<input type="checkbox"/> 2212022	22" Wheel size
<input type="checkbox"/> 2212024	24" Wheel size
<b>Spoke</b> \$ 75	
<input type="checkbox"/> 2212120	20" Wheel size
<input type="checkbox"/> 2212122	22" Wheel size
<input type="checkbox"/> 2212124	24" Wheel size
<input type="checkbox"/> 2212224	<b>24" Heavy Duty Spoke</b> \$ 350

Axle Type	
<input type="checkbox"/> 2213000	<b>Quick Release Rear Axles</b>
<input type="checkbox"/> 2213010	<b>Fixed Rear Axles</b>

Rear Tires	
<input type="checkbox"/> 2214000	<b>Urethane</b>
<input type="checkbox"/> 2214010	<b>Pneumatic</b> \$ 60 Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.
<input type="checkbox"/> 2214020	<b>Pneumatic with Flat-Free Inserts</b> E2213 \$ 95 Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.

Hand Rims	
<input type="checkbox"/> 2215000	<b>No Hand Rims</b>
<input type="checkbox"/> 2215010	<b>Standard</b> 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/> 2215020	<b>Plastic Coated</b> 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks	
<input type="checkbox"/> 2216010	<b>Push to Lock</b>
<input type="checkbox"/> 2216020	<b>Pull to Lock</b>
<input type="checkbox"/> 2216030	<b>Attendant Wheel Lock (Push to Lock)</b> \$ 85 Not compatible with Wheel Lock Extensions.
<input type="checkbox"/> 2216040	<b>User &amp; Attendant Wheel Lock</b> \$ 125 Available only on 20", 22", or 24" Rear Wheels.
<input type="checkbox"/> 2216100	<b>Wheel Lock Extensions</b> E0961 \$ 60 Not compatible with Attendant Wheel Lock.
<input type="checkbox"/> 2216200	<b>Hill Holders</b> \$ 350 Requires Pneumatic Rear Tires with Flat-Free Inserts.

Anti-Tippers	
<input checked="" type="checkbox"/> 2217010	<b>Rear Anti-Tippers</b> E0971 \$ 60 Required Safety Option.
<input type="checkbox"/> 2217020	<b>Front Anti-Tippers</b> \$ 300
<input type="checkbox"/> 2217030	<b>Heavy Duty Rear Anti-Tipper</b> \$ 600 Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.

Back Post Style	
<input type="checkbox"/> 2218010	<b>Push Bar Back</b>
<input type="checkbox"/> 2218020	<b>Cane Handle Back</b> Not compatible with Headrest Mounting Plate.
<b>Back Post Height</b> Seat Pan to Handle	
<input type="checkbox"/> 2218420	20" Height
<input type="checkbox"/> 2218425	25" Height
<input type="checkbox"/> 2218431	31" Height

Back Options	
<b>Angle/Height Adjustable Push Bar Handle</b> \$ 150 Not compatible with Headrest Mounting Plate or Cane Handle Back.	
<input type="checkbox"/> 2219010	Button on Outside
<input type="checkbox"/> 2219011	Button on Inside
<b>Seat to Back Angle</b>	
<input type="checkbox"/> 2219090	Set at 90°
<input type="checkbox"/> 2219098	Set at 98°
<input type="checkbox"/> 2219106	Set at 106°
<input type="checkbox"/> 2219114	Set at 114°
<input type="checkbox"/> 2219122	Set at 122°
<input type="checkbox"/> 2219210	<b>Backrest Extension 10" high</b> \$ 120 Not compatible with Angle/Height Adjustable Push Bar Handle.
<input type="checkbox"/> 2219220	<b>Fold Down Back</b>

Back Upholstery	
<input type="checkbox"/> 2220000	<b>No Back Upholstery</b>
<b>Standard</b> \$ 125 Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/> 2220018	Set at 18"
<input type="checkbox"/> 2220019	Set at 19"
<input type="checkbox"/> 2220020	Set at 20"
<input type="checkbox"/> 2220021	Set at 21"
<input type="checkbox"/> 2220022	Set at 22"
<input type="checkbox"/> 2220023	Set at 23"
<b>Adjustable Tension</b> \$ 250 Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/> 2220321	Set at 21"
<input type="checkbox"/> 2220322	Set at 22"
<input type="checkbox"/> 2220323	Set at 23"

Backrest	
<input type="checkbox"/> 2221000	<b>No Backrest</b>
<b>next Optima HA</b>	
<input type="checkbox"/> 2221100	AC - 16"-20" width, 16"-20" height E2620 \$ 690
<input type="checkbox"/> 2221101	PC - 16"-20" width, 16"-20" height E2620 \$ 690
<input type="checkbox"/> 2221105	AC - 22" width, 16"-20" height E2621 \$ 720
<input type="checkbox"/> 2221106	PC - 22" width, 16"-20" height E2621 \$ 720
<b>Matrx Elite</b>	
<input type="checkbox"/> 2221110	14"-20" width, 10"-20" height E2620 \$ 825
<input type="checkbox"/> 2221112	HD - 22" width, 16"-20" height E2621 \$ 1,095
<input type="checkbox"/> 2221116	Deep - 14"-20" width, 12"-20" height E2620 \$ 950
<b>ROHO Agility</b>	
<input type="checkbox"/> 2221120	Min - 14"-20" width, 10"-18" height E2620 \$ 825
<input type="checkbox"/> 2221122	Mid - 14"-20" width, 10"-18" height E2620 \$ 850
<input type="checkbox"/> 2221124	Max - 14"-20" width, 10"-18" height E2620 \$ 950

See next page for more Backrest Options →

Varilite Icon			
<input type="checkbox"/>	2221130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/>	2221132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
Backrest Size			
	2221400	Specify Size (Width" x Height") _____	

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	2223011	Installed on Right Backpost	
<input type="checkbox"/>	2223012	Installed on Left Backpost	
User Tilt only			\$ 100
Full Length Arm Pads only.			
<input type="checkbox"/>	2223021	Installed on Right Armrest	
<input type="checkbox"/>	2223022	Installed on Left Armrest	
Attendant & User Tilt			\$ 195
Full Length Arm Pads only.			
<input type="checkbox"/>	2223031	Backpost mount - Right side	
<input type="checkbox"/>	2223032	Backpost mount - Left side	
<input type="checkbox"/>	2223041	Armrest mount - Right side	
<input type="checkbox"/>	2223042	Armrest mount - Left side	
Anterior Tilt Limit			
<input type="checkbox"/>	2223220	Set at 20°	
<input type="checkbox"/>	2223225	Set at 25°	
<input type="checkbox"/>	2223230	Set at 30°	

Headrest			
<input type="checkbox"/>	2224000	<b>No Headrest</b>	
<input type="checkbox"/>	2224010	<b>Headrest Mounting Plate</b>	\$ 75
Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/>	2224100	<b>nxt MultiFit Head Support Hardware</b>	E1028 \$ 155
		<b>nxt MultiFit Head Support Pad</b>	E0955 \$ 195
<input type="checkbox"/>	2224105	Regular (9" x 5")	
<input type="checkbox"/>	2224106	Large (10" x 5")	
<input type="checkbox"/>	2224107	X-Large (14" x 5")	
Stealth Headrest Hardware			E1028
<input type="checkbox"/>	2224110	Multi-Axis HMO475	\$ 150
<input type="checkbox"/>	2224111	Link TWB480	\$ 325
Stealth Headrest Pad			E0955 \$ 200
<input type="checkbox"/>	2224115	10" width	
<input type="checkbox"/>	2224116	14" width	

Front Rigging			
<input type="checkbox"/>	2225000	<b>No Front Rigging</b>	
<input type="checkbox"/>	2225100	<b>Easy Out Plug In 70°</b>	
<input type="checkbox"/>	2225120	<b>Plug In Adjustable Knee Angle</b>	\$ 975

Foot Plates			
<input type="checkbox"/>	2226010	<b>Composite</b>	
<input type="checkbox"/>	2226020	<b>Angle Adjustable</b>	K0040 \$ 140
<input type="checkbox"/>	2226030	<b>One Piece Fixed</b>	\$ 350
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	2226040	<b>One Piece Flip-up</b>	\$ 395
<input type="checkbox"/>	2226050	<b>Interlocking Angle Adjustable</b>	\$ 195
Seat Widths of 14"-21" only.			

Footrest/Legrest Accessories			
<input type="checkbox"/>	2227040	<b>Adjustable Length Heel Loops</b>	E0951 \$ 50
<input type="checkbox"/>	2227050	<b>Bodypoint - Padded Calf Strap</b>	K0038 \$ 105
<input type="checkbox"/>	2227060	<b>Miller's Dynamic Footrest Mounts</b>	\$ 190

Extension Tubes			
Regular Extension Tubes			
<input type="checkbox"/>	2228010	Standard (13"-17")	
<input type="checkbox"/>	2228020	Long (15"-19")	\$ 75

Positioning Straps			
Bodypoint - Hip Belt			E0978
<input type="checkbox"/>	2229010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	2229015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	2229020	4-Point, Padded	\$ 175
<input type="checkbox"/>	2229030	<b>Bodypoint - Chest Strap</b>	E0960 \$ 85

Miscellaneous Options			
<input type="checkbox"/>	2230110	<b>IV Pole</b>	E1029 \$ 295
Not compatible with Attendant Wheel Locks or TTDS.			
<input type="checkbox"/>	2230120	<b>O2 Holder</b>	E2208 \$ 295
Not compatible with Attendant Wheel Locks or TTDS.			
<input type="checkbox"/>	2230170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	\$ 10
For Cushion mounting.			

Order Notes			