

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

Frame Type
Stellar LEAP Manual Posterior & Anterior Tilt Wheelchair
<input checked="" type="checkbox"/> 2201000 20° Manual Tilt Wheelchair \$ 3,995 <input checked="" type="checkbox"/> 2201101 30° Anterior Tilt Feature \$ 1,250 <input type="checkbox"/> 2201102 20° Dynamic Recline Backrest Tilt Feature \$ 500

Frame Options
<input checked="" type="checkbox"/> 2202000 Standard Base and Seat Frame 250 lbs Weight Capacity. <input type="checkbox"/> 2202300 Transit Tie Down System (TTDS) TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.

Seat Width
Standard Seat Widths
<input type="checkbox"/> 2203014 14" Seat Width <input type="checkbox"/> 2203015 15" Seat Width <input type="checkbox"/> 2203016 16" Seat Width <input type="checkbox"/> 2203017 17" Seat Width <input type="checkbox"/> 2203018 18" Seat Width <input type="checkbox"/> 2203019 19" Seat Width <input type="checkbox"/> 2203020 20" Seat Width <input type="checkbox"/> 2203021 21" Seat Width \$ 215 <input type="checkbox"/> 2203022 22" Seat Width \$ 215

Seat Depth
Standard Seat Depths (Adjustable)
<input type="checkbox"/> 2204016 16" Seat Depth <input type="checkbox"/> 2204017 17" Seat Depth <input type="checkbox"/> 2204018 18" Seat Depth <input type="checkbox"/> 2204019 19" Seat Depth <input type="checkbox"/> 2204020 20" Seat Depth

Custom Seat Depths	\$ 325
<input type="checkbox"/> 2204021 21" Seat Depth <input type="checkbox"/> 2204022 22" Seat Depth	

Frame Colour	
PDG Primary Colours	
<input type="checkbox"/> 2205000 Tracer Red <input type="checkbox"/> 2205001 Gloss Black <input type="checkbox"/> 2205002 Silver <input type="checkbox"/> 2205003 Black Shimmer <input type="checkbox"/> 2205004 Orange Fusion <input type="checkbox"/> 2205005 Karl's Imagination <input type="checkbox"/> 2205006 Blue Shimmer <input type="checkbox"/> 2205007 Gloss White <input type="checkbox"/> 2205008 Dark Titanium	
PDG Premium Colours	\$ 150
<input type="checkbox"/> 2205009 Geyser Green <input type="checkbox"/> 2205010 Lollipop Blue <input type="checkbox"/> 2205011 Lollipop Red	
Custom Colour	\$ 300
Contact PDG for specifications.	

Arm Type	
<input type="checkbox"/> 2206000 No Armrest <input type="checkbox"/> 2206100 Single Post Fixed Height (12") Full Length Armpad only.	
Single Post Adjustable Height	\$ 255
<input type="checkbox"/> 2206111 8"-11" <input type="checkbox"/> 2206114 10"-14" <input type="checkbox"/> 2206200 Dual Post Fixed Height (12") Full Length Armpad only.	
Dual Post Adjustable Height	\$ 255
<input type="checkbox"/> 2206211 8"-11" <input type="checkbox"/> 2206214 10"-14" <input type="checkbox"/> 2206216 12"-16"	
Locking Flip-up	\$ 255
Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.	
<input type="checkbox"/> 2206310 Set at 10" <input type="checkbox"/> 2206311 Set at 11" <input type="checkbox"/> 2206312 Set at 12" <input type="checkbox"/> 2206313 Set at 13" <input type="checkbox"/> 2206314 Set at 14"	

Armrests

Regular Armrests

- 2207100 Desk Length - 2" x 10" (pair)
- 2207101 Full Length - 2" x 14" (pair)

Wide Armrests

May interfere with wheels in some tilt configurations.

- 2207110 Desk Length - 3" x 10" (pair)
- 2207111 Full Length - 3" x 12" (pair)

Gel Armrests

May interfere with wheels in some tilt configurations.

- 2207120 Full Length - 2" x 12" (pair)
- 2207121 Full Length - 3.5" x 12" (pair)

Contour Arm Trough

May interfere with wheels in some tilt configurations.

- 2207130 Right (4.5" x 18") with Hand Pad
- 2207131 Left (4.5" x 18") with Hand Pad

Armrest Extensions

Not compatible with Flip-up or Single Post Armrests.

- 2207400 3" Extension - Forward
- 2207401 3" Extension - Rearward

Seat Pan

- 2208000 **No Seat Pan**
- 2208001 **Solid Seat Pan** \$ 120

Cushion

- 2209000 **No Cushion**
- next BioFit**
- 2209101 14"-20" width, 16"-20" depth \$ 595
- 2209102 22"-24" width, 16"-20" depth \$ 625

Cushion Size

Available in even number sizes only.

2209400 Specify Size (Width x Depth) _____

Seat Height

Circle desired Seat-to-Floor Height or enter information below:

RW	6" Casters					
12"	13"	14"	15"	16"	17"	18"
16"			15"	16"	17"	18"
20"	13"	14"	15"	16"	17"	
22"		14"	15"	16"	17"	
24"			15"	16"	17"	18"

RW	8" Casters					
16"			15"	16"	17"	18"
24"			15"	16"	17"	18"

Wheel size _____ Caster size _____ STF _____

STF Height is adjustable within the range of the selected wheel size.

Casters

Urethane

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2211061 6" x 1.25" Caster size
- 2211063 6" x 1.50" Caster size \$ 85
- 2211081 8" x 1.25" Caster size
- 2211084 8" x 2.00" Caster size \$ 135



Scan this code or visit www.pdgmobility.com to learn more about the Stellar LEAP Manual Tilt Wheelchair.

Pneumatic

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2211161 6" x 1.25" Caster size \$ 110
- 2211181 8" x 1.25" Caster size \$ 120
- 2211184 8" x 2.00" Caster size \$ 150

6" Soft Roll

- 2211260 \$ 95

- 2211400 **Ultra Sport Fork** \$ 395
Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.

Rear Wheels

Composite

- 2212012 12" Wheel size
- 2212016 16" Wheel size
- 2212020 20" Wheel size
- 2212022 22" Wheel size
- 2212024 24" Wheel size

Spoke

- 2212120 20" Wheel size
- 2212122 22" Wheel size
- 2212124 24" Wheel size
- 2212224 **24" Heavy Duty Spoke** \$ 350

Axle Type

- 2213000 **Quick Release Rear Axles**
- 2213010 **Fixed Rear Axles**

Rear Tires

- 2214000 **Urethane**
- 2214010 **Pneumatic** \$ 60
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.
- 2214020 **Pneumatic with Flat-Free Inserts** \$ 95
Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.

Hand Rims	
<input type="checkbox"/> 2215000	No Hand Rims
<input type="checkbox"/> 2215010	Standard 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/> 2215020	Plastic Coated 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks	
<input type="checkbox"/> 2216010	Push to Lock
<input type="checkbox"/> 2216020	Pull to Lock
<input type="checkbox"/> 2216030	Attendant Wheel Lock (Push to Lock) \$ 85 Not compatible with Wheel Lock Extensions.
<input type="checkbox"/> 2216040	User & Attendant Wheel Lock \$ 125 Available only on 20", 22", or 24" Rear Wheels.
<input type="checkbox"/> 2216100	Wheel Lock Extensions \$ 60 Not compatible with Attendant Wheel Lock.
<input type="checkbox"/> 2216200	Hill Holders \$ 350 Requires Pneumatic Rear Tires with Flat-Free Inserts.

Anti-Tippers	
<input checked="" type="checkbox"/> 2217010	Rear Anti-Tippers \$ 60 Required Safety Option.
<input type="checkbox"/> 2217020	Front Anti-Tippers \$ 300
<input type="checkbox"/> 2217030	Heavy Duty Rear Anti-Tipper \$ 600 Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.

Back Post Style	
<input type="checkbox"/> 2218010	Push Bar Back
<input type="checkbox"/> 2218020	Cane Handle Back Not compatible with Headrest Mounting Plate.
Back Post Height Seat Pan to Handle	
<input type="checkbox"/> 2218420	20" Height
<input type="checkbox"/> 2218425	25" Height
<input type="checkbox"/> 2218431	31" Height

Back Options	
<input type="checkbox"/> 2219000	Angle/Height Adjustable Push Bar Handle \$ 150 Not compatible with Headrest Mounting Plate or Cane Handle Back.
<input type="checkbox"/> 2219010	Button on Outside
<input type="checkbox"/> 2219011	Button on Inside
Seat to Back Angle	
<input type="checkbox"/> 2219090	Set at 90°
<input type="checkbox"/> 2219098	Set at 98°
<input type="checkbox"/> 2219106	Set at 106°
<input type="checkbox"/> 2219114	Set at 114°
<input type="checkbox"/> 2219122	Set at 122°
<input type="checkbox"/> 2219210	Backrest Extension 10" high \$ 120 Not compatible with Angle/Height Adjustable Push Bar Handle.
<input type="checkbox"/> 2219220	Fold Down Back Not compatible with Dynamic Recline Backrest.

Back Upholstery	
<input type="checkbox"/> 2220000	No Back Upholstery
Standard \$ 125 Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/> 2220018	Set at 18"
<input type="checkbox"/> 2220019	Set at 19"
<input type="checkbox"/> 2220020	Set at 20"
<input type="checkbox"/> 2220021	Set at 21"
<input type="checkbox"/> 2220022	Set at 22"
<input type="checkbox"/> 2220023	Set at 23"
Adjustable Tension \$ 250 Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/> 2220321	Set at 21"
<input type="checkbox"/> 2220322	Set at 22"
<input type="checkbox"/> 2220323	Set at 23"

Backrest	
<input type="checkbox"/> 2221000	No Backrest
next Optima HA Active Contour	
<input type="checkbox"/> 2221100	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/> 2221105	22"-24" width, 16"-20" height \$ 720
next Optima HA Plus Contour	
<input type="checkbox"/> 2221101	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/> 2221106	22"-24" width, 16"-20" height \$ 720
Backrest Size Available in even number sizes only.	
2221400	Specify Size (Width" x Height") _____

Tilt/Recline Options	
Attendant Tilt only	
<input type="checkbox"/> 2223011	Installed on Right Backpost
<input type="checkbox"/> 2223012	Installed on Left Backpost
User Tilt only \$ 100 Full Length Arm Pads only.	
<input type="checkbox"/> 2223021	Installed on Right Armrest
<input type="checkbox"/> 2223022	Installed on Left Armrest
Attendant & User Tilt \$ 195 Full Length Arm Pads only.	
<input type="checkbox"/> 2223031	Backpost mount - Right side
<input type="checkbox"/> 2223032	Backpost mount - Left side
<input type="checkbox"/> 2223041	Armrest mount - Right side
<input type="checkbox"/> 2223042	Armrest mount - Left side
Anterior Tilt Limit	
<input type="checkbox"/> 2223220	Set at -20°
<input type="checkbox"/> 2223225	Set at -25°
<input type="checkbox"/> 2223230	Set at -30°

Headrest	
<input type="checkbox"/> 2224000	No Headrest
<input type="checkbox"/> 2224010	Headrest Mounting Plate \$ 75 Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.
<input type="checkbox"/> 2224100	nxt MultiFit Head Support Hardware \$ 155
	nxt MultiFit Head Support Pad \$ 195
<input type="checkbox"/> 2224105	Regular (9" x 5")
<input type="checkbox"/> 2224106	Large (10" x 5")
<input type="checkbox"/> 2224107	X-Large (14" x 5")

Front Rigging	
<input type="checkbox"/> 2225000	No Front Rigging
<input type="checkbox"/> 2225100	Easy Out Plug In 70°
<input type="checkbox"/> 2225120	Plug In Adjustable Knee Angle \$ 975

Foot Plates	
<input type="checkbox"/> 2226010	Composite
<input type="checkbox"/> 2226020	Angle Adjustable \$ 140
<input type="checkbox"/> 2226030	One Piece Fixed \$ 350 Available on Plug In Front Rigging only.
<input type="checkbox"/> 2226040	One Piece Flip-up \$ 395
<input type="checkbox"/> 2226050	Interlocking Angle Adjustable \$ 195 Seat Widths of 14"-21" only.

Footrest/Legrest Accessories	
<input type="checkbox"/> 2227040	Adjustable Length Heel Loops \$ 50
<input type="checkbox"/> 2227050	Bodypoint - Padded Calf Strap \$ 105
<input type="checkbox"/> 2227060	Miller's Dynamic Footrest Mounts \$ 190

Extension Tubes	
	Regular Extension Tubes
<input type="checkbox"/> 2228010	Standard (13"-17")
<input type="checkbox"/> 2228020	Long (15"-19") \$ 75

Positioning Straps	
	Bodypoint - Hip Belt
<input type="checkbox"/> 2229010	Center-Pull Auto Style \$ 105
<input type="checkbox"/> 2229015	Center-Pull Auto Style, Padded \$ 155
<input type="checkbox"/> 2229020	4-Point, Padded \$ 175
<input type="checkbox"/> 2229030	Bodypoint - Chest Strap \$ 85

Miscellaneous Options	
<input type="checkbox"/> 2230110	IV Pole \$ 295 Not compatible with Attendant Wheel Locks or TTDS.
<input type="checkbox"/> 2230120	O2 Holder \$ 295 Not compatible with Attendant Wheel Locks or TTDS.
<input type="checkbox"/> 2230170	Adhesive Back Hook Strip, 2" x 30" \$ 10 For Cushion mounting.