

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877-323-9093.

**Order Notes**

---



---



---

**Frame Type**

**Stellar LEAP Manual Posterior & Anterior Tilt Wheelchair**

<input checked="" type="checkbox"/>	2201000	20° Manual Tilt Wheelchair	\$ 3,995
<input checked="" type="checkbox"/>	2201101	30° Anterior Tilt Feature	\$ 1,250
<input type="checkbox"/>	2201102	20° Dynamic Recline Backrest Tilt Feature	\$ 500

**Frame Options**

<input checked="" type="checkbox"/>	2202000	<b>Standard Base and Seat Frame</b> 250 lbs Weight Capacity.	
<input type="checkbox"/>	2202300	<b>Transit Tie Down System (TTDS)</b> TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.	

**Seat Width**

**Standard Seat Widths**

<input type="checkbox"/>	2203014	14" Seat Width	
<input type="checkbox"/>	2203015	15" Seat Width	
<input type="checkbox"/>	2203016	16" Seat Width	
<input type="checkbox"/>	2203017	17" Seat Width	
<input type="checkbox"/>	2203018	18" Seat Width	
<input type="checkbox"/>	2203019	19" Seat Width	
<input type="checkbox"/>	2203020	20" Seat Width	
<input type="checkbox"/>	2203021	21" Seat Width	\$ 215
<input type="checkbox"/>	2203022	22" Seat Width	\$ 215

**Seat Depth**

**Standard Seat Depths (Adjustable)**

<input type="checkbox"/>	2204016	16" Seat Depth	
<input type="checkbox"/>	2204017	17" Seat Depth	
<input type="checkbox"/>	2204018	18" Seat Depth	
<input type="checkbox"/>	2204019	19" Seat Depth	
<input type="checkbox"/>	2204020	20" Seat Depth	

**Custom Seat Depths** \$ 325

<input type="checkbox"/>	2204021	21" Seat Depth	
<input type="checkbox"/>	2204022	22" Seat Depth	

**Frame Colour**

**PDG Primary Colours**

<input type="checkbox"/>	2205000	Tracer Red	
<input type="checkbox"/>	2205001	Gloss Black	
<input type="checkbox"/>	2205002	Silver	
<input type="checkbox"/>	2205003	Black Shimmer	
<input type="checkbox"/>	2205004	Orange Fusion	
<input type="checkbox"/>	2205005	Karl's Imagination	
<input type="checkbox"/>	2205006	Blue Shimmer	
<input type="checkbox"/>	2205007	Gloss White	
<input type="checkbox"/>	2205008	Dark Titanium	

**PDG Premium Colours** \$ 150

<input type="checkbox"/>	2205009	Geyser Green	
<input type="checkbox"/>	2205010	Lollipop Blue	
<input type="checkbox"/>	2205011	Lollipop Red	

2205200 **Custom Colour** \$ 300  
Contact PDG for specifications.

**Arm Type**

<input type="checkbox"/>	2206000	<b>No Armrest</b> Full Length Armpad only.	
<input type="checkbox"/>	2206100	<b>Single Post Fixed Height (12")</b> Full Length Armpad only.	
		<b>Single Post Adjustable Height</b> <span style="float: right;">\$ 255</span>	
<input type="checkbox"/>	2206111	8"-11"	
<input type="checkbox"/>	2206114	10"-14"	
<input type="checkbox"/>	2206200	<b>Dual Post Fixed Height (12")</b> Full Length Armpad only.	
		<b>Dual Post Adjustable Height</b> <span style="float: right;">\$ 255</span>	
<input type="checkbox"/>	2206211	8"-11"	
<input type="checkbox"/>	2206214	10"-14"	
<input type="checkbox"/>	2206216	12"-16"	
		<b>Locking Flip-up</b> <span style="float: right;">\$ 255</span> Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.	
<input type="checkbox"/>	2206310	Set at 10"	
<input type="checkbox"/>	2206311	Set at 11"	
<input type="checkbox"/>	2206312	Set at 12"	
<input type="checkbox"/>	2206313	Set at 13"	
<input type="checkbox"/>	2206314	Set at 14"	

## Armrests

### Regular Armrests

- 2207100 Desk Length - 2" x 10" (pair)
- 2207101 Full Length - 2" x 14" (pair)

### Wide Armrests

May interfere with wheels in some tilt configurations.

- 2207110 Desk Length - 3" x 10" (pair)
- 2207111 Full Length - 3" x 12" (pair)

### Gel Armrests

May interfere with wheels in some tilt configurations.

- 2207120 Full Length - 2" x 12" (pair)
- 2207121 Full Length - 3.5" x 12" (pair)

### Contour Arm Trough

May interfere with wheels in some tilt configurations.

- 2207130 Right (4.5" x 18") with Hand Pad
- 2207131 Left (4.5" x 18") with Hand Pad

### Armrest Extensions

Not compatible with Flip-up or Single Post Armrests.

- 2207400 3" Extension - Forward
- 2207401 3" Extension - Rearward

## Seat Pan

- 2208000 **No Seat Pan**
- 2208001 **Solid Seat Pan** \$ 120

## Cushion

- 2209000 **No Cushion**
- next BioFit**
- 2209101 14"-20" width, 16"-20" depth \$ 595
- 2209102 22"-24" width, 16"-20" depth \$ 625

### Cushion Size

Available in even number sizes only.

2209400 Specify Size (Width" x Depth") \_\_\_\_\_

## Seat Height

Circle desired Seat-to-Floor Height or enter information below:

RW	6" Casters					
12"	13"	14"	15"	16"	17"	18"
16"			15"	16"	17"	18"
20"	13"	14"	15"	16"	17"	
22"		14"	15"	16"	17"	
24"			15"	16"	17"	18"

RW	8" Casters					
16"			15"	16"	17"	18"
24"			15"	16"	17"	18"

Wheel size \_\_\_\_\_ Caster size \_\_\_\_\_ STF \_\_\_\_\_

STF Height is adjustable within the range of the selected wheel size.

## Casters

### Urethane

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2211061 6" x 1.25" Caster size
- 2211063 6" x 1.50" Caster size \$ 85
- 2211081 8" x 1.25" Caster size
- 2211084 8" x 2.00" Caster size \$ 135



Scan this code or visit [www.pdgmobility.com](http://www.pdgmobility.com) to learn more about the Stellar LEAP Manual Tilt Wheelchair.

### Pneumatic

8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2211161 6" x 1.25" Caster size \$ 110
- 2211181 8" x 1.25" Caster size \$ 120
- 2211184 8" x 2.00" Caster size \$ 150

### 6" Soft Roll

- 2211260 **6" Soft Roll** \$ 95

### Ultra Sport Fork

Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.

- 2211400 **Ultra Sport Fork** \$ 395

## Rear Wheels

### Composite

- 2212012 12" Wheel size
- 2212016 16" Wheel size
- 2212020 20" Wheel size
- 2212022 22" Wheel size
- 2212024 24" Wheel size

### Spoke

- 2212120 20" Wheel size
- 2212122 22" Wheel size
- 2212124 24" Wheel size

### 24" Heavy Duty Spoke

- 2212224 **24" Heavy Duty Spoke** \$ 350

## Axle Type

- 2213000 **Quick Release Rear Axles**
- 2213010 **Fixed Rear Axles**

## Rear Tires

- 2214000 **Urethane**
- 2214010 **Pneumatic** \$ 60  
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.
- 2214020 **Pneumatic with Flat-Free Inserts** \$ 95  
Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.

Hand Rims	
<input type="checkbox"/> 2215000	<b>No Hand Rims</b>
<input type="checkbox"/> 2215010	<b>Standard</b> 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/> 2215020	<b>Plastic Coated</b> 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks	
<input type="checkbox"/> 2216010	<b>Push to Lock</b>
<input type="checkbox"/> 2216020	<b>Pull to Lock</b>
<input type="checkbox"/> 2216030	<b>Attendant Wheel Lock (Push to Lock)</b> \$ 85 Not compatible with Wheel Lock Extensions.
<input type="checkbox"/> 2216040	<b>User &amp; Attendant Wheel Lock</b> \$ 125 Available only on 20", 22", or 24" Rear Wheels.
<input type="checkbox"/> 2216100	<b>Wheel Lock Extensions</b> \$ 60 Not compatible with Attendant Wheel Lock.
<input type="checkbox"/> 2216200	<b>Hill Holders</b> \$ 350 Requires Pneumatic Rear Tires with Flat-Free Inserts.

Anti-Tippers	
<input checked="" type="checkbox"/> 2217010	<b>Rear Anti-Tippers</b> \$ 60 Required Safety Option.
<input type="checkbox"/> 2217020	<b>Front Anti-Tippers</b> \$ 300
<input type="checkbox"/> 2217030	<b>Heavy Duty Rear Anti-Tipper</b> \$ 600 Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.

Back Post Style	
<input type="checkbox"/> 2218010	<b>Push Bar Back</b>
<input type="checkbox"/> 2218020	<b>Cane Handle Back</b> Not compatible with Headrest Mounting Plate.
<b>Back Post Height</b> Seat Pan to Handle	
<input type="checkbox"/> 2218420	20" Height
<input type="checkbox"/> 2218425	25" Height
<input type="checkbox"/> 2218431	31" Height

Back Options	
<input type="checkbox"/>	<b>Angle/Height Adjustable Push Bar Handle</b> \$ 150 Not compatible with Headrest Mounting Plate or Cane Handle Back.
<input type="checkbox"/> 2219010	Button on Outside
<input type="checkbox"/> 2219011	Button on Inside
<b>Seat to Back Angle</b>	
<input type="checkbox"/> 2219090	Set at 90°
<input type="checkbox"/> 2219098	Set at 98°
<input type="checkbox"/> 2219106	Set at 106°
<input type="checkbox"/> 2219114	Set at 114°
<input type="checkbox"/> 2219122	Set at 122°
<input type="checkbox"/> 2219210	<b>Backrest Extension 10" high</b> \$ 120 Not compatible with Angle/Height Adjustable Push Bar Handle.
<input type="checkbox"/> 2219220	<b>Fold Down Back</b>

Back Upholstery	
<input type="checkbox"/> 2220000	<b>No Back Upholstery</b>
<b>Standard</b> \$ 125 Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/> 2220018	Set at 18"
<input type="checkbox"/> 2220019	Set at 19"
<input type="checkbox"/> 2220020	Set at 20"
<input type="checkbox"/> 2220021	Set at 21"
<input type="checkbox"/> 2220022	Set at 22"
<input type="checkbox"/> 2220023	Set at 23"
<b>Adjustable Tension</b> \$ 250 Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/> 2220321	Set at 21"
<input type="checkbox"/> 2220322	Set at 22"
<input type="checkbox"/> 2220323	Set at 23"

Backrest	
<input type="checkbox"/> 2221000	<b>No Backrest</b>
<b>next Optima HA Active Contour</b>	
<input type="checkbox"/> 2221100	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/> 2221105	22"-24" width, 16"-20" height \$ 720
<b>next Optima HA Plus Contour</b>	
<input type="checkbox"/> 2221101	16"-20" width, 16"-20" height \$ 690
<input type="checkbox"/> 2221106	22"-24" width, 16"-20" height \$ 720
<b>Backrest Size</b> Available in even number sizes only.	
2221400	Specify Size (Width" x Height") _____

Tilt/Recline Options	
<b>Attendant Tilt only</b>	
<input type="checkbox"/> 2223011	Installed on Right Backpost
<input type="checkbox"/> 2223012	Installed on Left Backpost
<b>User Tilt only</b> \$ 100 Full Length Arm Pads only.	
<input type="checkbox"/> 2223021	Installed on Right Armrest
<input type="checkbox"/> 2223022	Installed on Left Armrest
<b>Attendant &amp; User Tilt</b> \$ 195 Full Length Arm Pads only.	
<input type="checkbox"/> 2223031	Backpost mount - Right side
<input type="checkbox"/> 2223032	Backpost mount - Left side
<input type="checkbox"/> 2223041	Armrest mount - Right side
<input type="checkbox"/> 2223042	Armrest mount - Left side
<b>Anterior Tilt Limit</b>	
<input type="checkbox"/> 2223220	Set at -20°
<input type="checkbox"/> 2223225	Set at -25°
<input type="checkbox"/> 2223230	Set at -30°

Headrest	
<input type="checkbox"/> 2224000	<b>No Headrest</b>
<input type="checkbox"/> 2224010	<b>Headrest Mounting Plate</b> \$ 75 Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.
<input type="checkbox"/> 2224100	<b>nxt MultiFit Head Support Hardware</b> \$ 155
	<b>nxt MultiFit Head Support Pad</b> \$ 195
<input type="checkbox"/> 2224105	Regular (9" x 5")
<input type="checkbox"/> 2224106	Large (10" x 5")
<input type="checkbox"/> 2224107	X-Large (14" x 5")

Front Rigging	
<input type="checkbox"/> 2225000	<b>No Front Rigging</b>
<input type="checkbox"/> 2225100	<b>Easy Out Plug In 70°</b>
<input type="checkbox"/> 2225120	<b>Plug In Adjustable Knee Angle</b> \$ 975

Foot Plates	
<input type="checkbox"/> 2226010	<b>Composite</b>
<input type="checkbox"/> 2226020	<b>Angle Adjustable</b> \$ 140
<input type="checkbox"/> 2226030	<b>One Piece Fixed</b> \$ 350 Available on Plug In Front Rigging only.
<input type="checkbox"/> 2226040	<b>One Piece Flip-up</b> \$ 395
<input type="checkbox"/> 2226050	<b>Interlocking Angle Adjustable</b> \$ 195 Seat Widths of 14"-21" only.

Footrest/Legrest Accessories	
<input type="checkbox"/> 2227040	<b>Adjustable Length Heel Loops</b> \$ 50
<input type="checkbox"/> 2227050	<b>Bodypoint - Padded Calf Strap</b> \$ 105
<input type="checkbox"/> 2227060	<b>Miller's Dynamic Footrest Mounts</b> \$ 190

Extension Tubes	
	<b>Regular Extension Tubes</b>
<input type="checkbox"/> 2228010	Standard (13"-17")
<input type="checkbox"/> 2228020	Long (15"-19") \$ 75

Positioning Straps	
	<b>Bodypoint - Hip Belt</b>
<input type="checkbox"/> 2229010	Center-Pull Auto Style \$ 105
<input type="checkbox"/> 2229015	Center-Pull Auto Style, Padded \$ 155
<input type="checkbox"/> 2229020	4-Point, Padded \$ 175
<input type="checkbox"/> 2229030	<b>Bodypoint - Chest Strap</b> \$ 85

Miscellaneous Options	
<input type="checkbox"/> 2230110	<b>IV Pole</b> \$ 295 Not compatible with Attendant Wheel Locks or TTDS.
<input type="checkbox"/> 2230120	<b>O2 Holder</b> \$ 295 Not compatible with Attendant Wheel Locks or TTDS.
<input type="checkbox"/> 2230170	<b>Adhesive Back Hook Strip, 2" x 30"</b> \$ 10 For Cushion mounting.