

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

Frame Type
<input checked="" type="checkbox"/> 1201000 Stellar 45° Manual Tilt Wheelchair \$ 3,995

Frame Options
<input type="checkbox"/> 1202000 Standard Base and Seat Frame 250 lbs Weight Capacity.
<input type="checkbox"/> 1202050 Weight Capacity Upgrade \$ 200 Increases Weight Capacity to 300 lbs through addition of a secondary Gas Strut. Available for Seat Width of up to 22".
<input type="checkbox"/> 1202100 Heavy Duty Base and Seat Frame \$ 300 450 lbs Weight Capacity. Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Required for 23" & 24" widths.
<input type="checkbox"/> 1202200 Extended Base with Rear Floor Jacks \$ 600 Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Incompatible with Attendant wheel lock, HD Anti-Tipper, IV Pole, O2 Holder.
<input type="checkbox"/> 1202201 Shortened Strut Bracket \$ 125 Increases heel clearance for foot propelling. Available only for 23" & 24" widths.
<input type="checkbox"/> 1202300 Transit Tie Down System (TTDS) TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and 23" & 24" widths.

Seat Width
Standard Seat Widths
<input type="checkbox"/> 1203014 14" Seat Width
<input type="checkbox"/> 1203015 15" Seat Width
<input type="checkbox"/> 1203016 16" Seat Width
<input type="checkbox"/> 1203017 17" Seat Width
<input type="checkbox"/> 1203018 18" Seat Width
<input type="checkbox"/> 1203019 19" Seat Width
<input type="checkbox"/> 1203020 20" Seat Width
<input type="checkbox"/> 1203021 21" Seat Width \$ 215
<input type="checkbox"/> 1203022 22" Seat Width \$ 215

Custom Seat Widths	\$ 325
Requires 1202100 Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.	
<input type="checkbox"/> 1203023 23" Seat Width	
<input type="checkbox"/> 1203024 24" Seat Width	
Use Stellar HD Order Form for Seat Widths beyond 24"	

Seat Depth
Standard Seat Depths (Adjustable)
<input type="checkbox"/> 1204016 16" Seat Depth
<input type="checkbox"/> 1204017 17" Seat Depth
<input type="checkbox"/> 1204018 18" Seat Depth
<input type="checkbox"/> 1204019 19" Seat Depth
<input type="checkbox"/> 1204020 20" Seat Depth

Custom Seat Depths	\$ 325
<input type="checkbox"/> 1204021 21" Seat Depth	
<input type="checkbox"/> 1204022 22" Seat Depth	

Frame Colour	
PDG Primary Colours	
<input type="checkbox"/> 1205000 Tracer Red	
<input type="checkbox"/> 1205001 Gloss Black	
<input type="checkbox"/> 1205002 Silver	
<input type="checkbox"/> 1205003 Black Shimmer	
<input type="checkbox"/> 1205004 Orange Fusion	
<input type="checkbox"/> 1205005 Karl's Imagination	
<input type="checkbox"/> 1205006 Blue Shimmer	
<input type="checkbox"/> 1205007 Gloss White	
<input type="checkbox"/> 1205008 Dark Titanium	
PDG Premium Colours	\$ 150
<input type="checkbox"/> 1205009 Geyser Green	
<input type="checkbox"/> 1205010 Lollipop Blue	
<input type="checkbox"/> 1205011 Lollipop Red	
Custom Colour	\$ 300
Contact PDG for specifications.	

PDG Product Design Group Inc.
318 East Kent Avenue South
Vancouver BC V5X 4N6 Canada
Phone 604-323-9220 | 888-858-4422



Scan this code or visit www.pdgmobility.com to learn more about the Stellar Manual Tilt Wheelchair.

Arm Type	
<input type="checkbox"/>	1206000 No Armrest
<input type="checkbox"/>	1206100 Single Post Fixed Height (12") Full Length Armpad only.
	Single Post Adjustable Height \$ 255 Not compatible with the Heavy Duty Package and Seat Widths greater than 22".
<input type="checkbox"/>	1206111 8"-11"
<input type="checkbox"/>	1206114 10"-14"
<input type="checkbox"/>	1206200 Dual Post Fixed Height (12") Full Length Armpad only.
	Dual Post Adjustable Height \$ 255
<input type="checkbox"/>	1206211 8"-11"
<input type="checkbox"/>	1206214 10"-14"
<input type="checkbox"/>	1206216 12"-16"
	Locking Flip-up \$ 255 Not recommended for clients with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery.
<input type="checkbox"/>	1206310 Set at 10"
<input type="checkbox"/>	1206311 Set at 11"
<input type="checkbox"/>	1206312 Set at 12"
<input type="checkbox"/>	1206313 Set at 13"
<input type="checkbox"/>	1206314 Set at 14"

Armrests	
	Regular Armpads
<input type="checkbox"/>	1207100 Desk Length - 2" x 10" (pair)
<input type="checkbox"/>	1207101 Full Length - 2" x 14" (pair)
	Wide Armpads \$ 115 May interfere with wheels in some tilt configurations.
<input type="checkbox"/>	1207110 Desk Length - 3" x 10" (pair)
<input type="checkbox"/>	1207111 Full Length - 3" x 12" (pair)
	Gel Armpads \$ 115 May interfere with wheels in some tilt configurations.
<input type="checkbox"/>	1207120 Full Length - 2" x 12" (pair)
<input type="checkbox"/>	1207121 Full Length - 3.5" x 12" (pair)
	Contour Arm Trough \$ 125 May interfere with wheels in some tilt configurations.
<input type="checkbox"/>	1207130 Right (4.5" x 18") with Hand Pad
<input type="checkbox"/>	1207131 Left (4.5" x 18") with Hand Pad
<input type="checkbox"/>	1207220 Armrest Wideners for Dual Post Armrests \$ 375 Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels or Clothing Guards.
	Clothing Guards Not compatible with Flip-up Armrests.
<input type="checkbox"/>	1207310 For Single Post Armrests \$ 110 10"-14" Height range only
<input type="checkbox"/>	1207320 For Dual Post Armrests \$ 165
	Armrest Extensions \$ 60 Not compatible with Flip-up or Single Post Armrests.
<input type="checkbox"/>	1207400 3" Extension - Forward
<input type="checkbox"/>	1207401 3" Extension - Rearward

Seat Pan	
<input type="checkbox"/>	1208000 No Seat Pan
<input type="checkbox"/>	1208001 Solid Seat Pan \$ 120

Cushion	
<input type="checkbox"/>	1209000 No Cushion
	next BioFit
<input type="checkbox"/>	1209101 14"-20" width, 16"-20" depth \$ 595
<input type="checkbox"/>	1209102 22"-24" width, 16"-20" depth \$ 625
	Cushion Size Available in even number sizes only.
1209400	Specify Size (Width" x Depth") _____

Seat Height																																																																		
Circle desired Seat-to-Floor Height or enter information below:																																																																		
<table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">6" Casters</th> </tr> </thead> <tbody> <tr> <td>12"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td></td> </tr> <tr> <td>16"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>20"</td> </tr> <tr> <td>20"</td> <td>14"</td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td></td> </tr> <tr> <td>22"</td> <td></td> <td>15"</td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>20"</td> </tr> </tbody> </table>	RW	6" Casters						12"	14"	15"	16"	17"	18"		16"			16"	17"	18"	20"	20"	14"	15"	16"	17"	18"		22"		15"	16"	17"	18"	19"	24"			16"	17"	18"	20"	<table border="1"> <thead> <tr> <th>RW</th> <th colspan="6">8" Casters</th> </tr> </thead> <tbody> <tr> <td>16"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> <tr> <td>24"</td> <td></td> <td></td> <td>16"</td> <td>17"</td> <td>18"</td> <td>19"</td> <td>20"</td> </tr> </tbody> </table>	RW	8" Casters						16"			16"	17"	18"	19"	20"	24"			16"	17"	18"	19"	20"
RW	6" Casters																																																																	
12"	14"	15"	16"	17"	18"																																																													
16"			16"	17"	18"	20"																																																												
20"	14"	15"	16"	17"	18"																																																													
22"		15"	16"	17"	18"	19"																																																												
24"			16"	17"	18"	20"																																																												
RW	8" Casters																																																																	
16"			16"	17"	18"	19"	20"																																																											
24"			16"	17"	18"	19"	20"																																																											
Wheel size _____ Caster size _____ STF _____																																																																		
STF Height is adjustable within the range of the selected wheel size.																																																																		

Casters	
	Urethane
	8" Casters are compatible only with 16" and 24" Rear Wheels.
<input type="checkbox"/>	1211061 6" x 1.25" Caster size
<input type="checkbox"/>	1211063 6" x 1.50" Caster size \$ 85
<input type="checkbox"/>	1211081 8" x 1.25" Caster size
<input type="checkbox"/>	1211084 8" x 2.00" Caster size \$ 135
	Pneumatic
	8" Casters are compatible only with 16" and 24" Rear Wheels.
<input type="checkbox"/>	1211161 6" x 1.25" Caster size \$ 110
<input type="checkbox"/>	1211181 8" x 1.25" Caster size \$ 120
<input type="checkbox"/>	1211184 8" x 2.00" Caster size \$ 150
<input type="checkbox"/>	1211260 6" Soft Roll \$ 95
<input type="checkbox"/>	1211400 Ultra Sport Fork \$ 395 Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.

Rear Wheels	
Composite	
<input type="checkbox"/> 1212012	12" Wheel size
<input type="checkbox"/> 1212016	16" Wheel size
<input type="checkbox"/> 1212020	20" Wheel size
<input type="checkbox"/> 1212022	22" Wheel size
<input type="checkbox"/> 1212024	24" Wheel size
Spoke \$ 75	
<input type="checkbox"/> 1212120	20" Wheel size
<input type="checkbox"/> 1212122	22" Wheel size
<input type="checkbox"/> 1212124	24" Wheel size
<input type="checkbox"/> 1212224	24" Heavy Duty Spoke \$ 350 Required for Heavy Duty Package.

Axle Type	
<input type="checkbox"/> 1213000	Quick Release Rear Axles
<input type="checkbox"/> 1213010	Fixed Rear Axles

Rear Tires	
<input type="checkbox"/> 1214000	Urethane
<input type="checkbox"/> 1214010	Pneumatic \$ 60 Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.
<input type="checkbox"/> 1214020	Pneumatic with Flat-Free Inserts \$ 95 Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.

Hand Rims	
<input type="checkbox"/> 1215000	No Hand Rims
<input type="checkbox"/> 1215010	Standard 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/> 1215020	Plastic Coated 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks	
<input type="checkbox"/> 1216010	Push to Lock
<input type="checkbox"/> 1216020	Pull to Lock
<input type="checkbox"/> 1216030	Attendant Wheel Lock (Push to Lock) \$ 85 Not compatible with TTDS (23" & 24").
<input type="checkbox"/> 1216040	User & Attendant Wheel Lock \$ 125 Available only on 20", 22", or 24" Rear Wheels.
<input type="checkbox"/> 1216100	Wheel Lock Extensions \$ 60 Not compatible with Attendant Wheel Lock.
<input type="checkbox"/> 1216200	Hill Holders \$ 350 Requires Pneumatic Rear Tires with Flat-Free Inserts.

Anti-Tippers	
<input checked="" type="checkbox"/> 1217010	Rear Anti-Tippers \$ 60 Required Safety Option.
<input type="checkbox"/> 1217020	Front Anti-Tippers \$ 300
<input type="checkbox"/> 1217030	Heavy Duty Rear Anti-Tipper \$ 600 Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.

Back Post Style	
<input type="checkbox"/> 1218010	Push Bar Back
<input type="checkbox"/> 1218020	Cane Handle Back Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.
Back Post Height Seat Pan to Handle	
<input type="checkbox"/> 1218425	25" Height
<input type="checkbox"/> 1218431	31" Height

Back Options	
<input type="checkbox"/> 1219010	Angle/Height Adjustable Push Bar Handle \$ 150 Not compatible with Headrest Mounting Plate or Cane Handle Back.
<input type="checkbox"/> 1219010	Button on Outside
<input type="checkbox"/> 1219011	Button on Inside
Seat to Back Angle	
<input type="checkbox"/> 1219090	Set at 90°
<input type="checkbox"/> 1219098	Set at 98°
<input type="checkbox"/> 1219106	Set at 106°
<input type="checkbox"/> 1219114	Set at 114°
<input type="checkbox"/> 1219122	Set at 122°
<input type="checkbox"/> 1219210	Backrest Extension 10" high \$ 120 Not compatible with Angle/Height Adjustable Push Bar Handle.
<input type="checkbox"/> 1219220	Fold Down Back Not compatible with Reclining Back.
<input type="checkbox"/> 1219230	Spring Back Suspension \$ 105 Not compatible with Reclining Back.
<input type="checkbox"/> 1219240	Dynamic Rocker Back Assembly \$ 2,995 Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller. Requires longer Frame.

Back Upholstery	
<input type="checkbox"/> 1220000	No Back Upholstery
<input type="checkbox"/> 1220001	Standard \$ 125 Back height measured from Seat Pan to top of Upholstery
<input type="checkbox"/> 1220021	Set at 21"
<input type="checkbox"/> 1220022	Set at 22"
<input type="checkbox"/> 1220100	Short (25" Back only) \$ 125 Back height measured from Seat Pan to top of Upholstery
<input type="checkbox"/> 1220118	Set at 18"
<input type="checkbox"/> 1220119	Set at 19"
<input type="checkbox"/> 1220120	Set at 20"
<input type="checkbox"/> 1220150	Tall (31" Back only) \$ 125 Back height measured from Seat Pan to top of Upholstery
<input type="checkbox"/> 1220224	Set at 24"
<input type="checkbox"/> 1220225	Set at 25"
<input type="checkbox"/> 1220226	Set at 26"
<input type="checkbox"/> 1220300	Adjustable Tension \$ 250 Back height measured from Seat Pan to top of Upholstery
<input type="checkbox"/> 1220321	Set at 21"
<input type="checkbox"/> 1220322	Set at 22"
<input type="checkbox"/> 1220323	Set at 23"

Backrest			
<input type="checkbox"/>	1221000	No Backrest	
nxt Optima HA Active Contour			
<input type="checkbox"/>	1221100	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1221105	22"-24" width, 16"-20" height	\$ 720
nxt Optima HA Plus Contour			
<input type="checkbox"/>	1221101	16"-20" width, 16"-20" height	\$ 690
<input type="checkbox"/>	1221106	22"-24" width, 16"-20" height	\$ 720
Backrest Size			
Available in even number sizes only.			
<input type="checkbox"/>	1221400	Specify Size (Width" x Height") _____	

Tilt/Recline Options			
Attendant Tilt only			
<input type="checkbox"/>	1223011	Installed on Right Back Post	
<input type="checkbox"/>	1223012	Installed on Left Back Post	
User Tilt only \$ 100			
Not available on widths greater than 22". Full Length Arm Pads only.			
<input type="checkbox"/>	1223021	Installed on Right Armrest	
<input type="checkbox"/>	1223022	Installed on Left Armrest	
Attendant & User Tilt \$ 195			
Not available on widths greater than 22". Full Length Arm Pads only.			
<input type="checkbox"/>	1223031	Back Post mount - Right side	
<input type="checkbox"/>	1223032	Back Post mount - Left side	
<input type="checkbox"/>	1223041	Armrest mount - Right side	
<input type="checkbox"/>	1223042	Armrest mount - Left side	
<input type="checkbox"/>	1023110	Trigger Extension	\$ 50
Tilt Limit			
<input type="checkbox"/>	1223220	Set at 20°	
<input type="checkbox"/>	1223225	Set at 25°	
<input type="checkbox"/>	1223230	Set at 30°	
<input type="checkbox"/>	1223235	Set at 35°	
<input type="checkbox"/>	1223240	Set at 40°	
<input type="checkbox"/>	1223245	Set at 45°	
<input type="checkbox"/>	1223410	Reclining Backrest (Manual)	\$ 830
Maximum recline = 30°. Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, 20"-22" seat depth. Requires longer Frame.			

Headrest			
<input type="checkbox"/>	1224000	No Headrest	
<input type="checkbox"/>	1224010	Headrest Mounting Plate	\$ 75
Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			
<input type="checkbox"/>	1224100	nxt MultiFit Head Support Hardware	\$ 155
nxt MultiFit Head Support Pad \$ 195			
<input type="checkbox"/>	1224105	Regular (9" x 5")	
<input type="checkbox"/>	1224106	Large (10" x 5")	
<input type="checkbox"/>	1224107	X-Large (14" x 5")	

Front Rigging			
<input type="checkbox"/>	1225000	No Front Rigging	
<input type="checkbox"/>	1225001	Sector Blocks only	

<input type="checkbox"/>	1225110	Plug In 70°	
<input type="checkbox"/>	1225120	Plug In Adjustable Knee Angle	\$ 975
<input type="checkbox"/>	1225130	Swing Away 60°	
<input type="checkbox"/>	1225140	Swing Away 70°	
<input type="checkbox"/>	1225200	Swing Away Elevating Legrest	\$ 355

Foot Plates			
<input type="checkbox"/>	1226010	Composite	
<input type="checkbox"/>	1226020	Angle Adjustable	\$ 140
<input type="checkbox"/>	1226030	One Piece Fixed	\$ 350
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1226040	One Piece Flip-up	\$ 395
<input type="checkbox"/>	1226050	Interlocking Angle Adjustable	\$ 195
Seat Widths of 14"-21" only.			

Footrest/Legrest Accessories			
<input type="checkbox"/>	1227010	Swing Away 60° & 70° Calf Pad	\$ 135
Not compatible with 60° & 70° Gel Pad.			
<input type="checkbox"/>	1227020	Swing Away 60° & 70° Gel Pad	\$ 100
Not compatible with 60° & 70° Swing Away Calf Pad.			
<input type="checkbox"/>	1227030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1227040	Adjustable Length Heel Loops	\$ 50
<input type="checkbox"/>	1227050	Bodypoint - Padded Calf Strap	\$ 105
Sizes greater than 22" are not padded.			
Miller's Dynamic Footrest Mounts			
<input type="checkbox"/>	1227060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1227065	For Swing Away 70° Front Rigging	\$ 380

Extension Tubes			
Regular Extension Tubes			
<input type="checkbox"/>	1228010	Standard (13"-17")	
<input type="checkbox"/>	1228020	Long (15"-19")	\$ 75
Extra Strength Extension Tubes (Solid) \$ 100			
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1228110	Standard (13"-17")	
<input type="checkbox"/>	1228120	Long (15"-19")	

Positioning Straps			
Bodypoint - Hip Belts			
<input type="checkbox"/>	1229010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	1229015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1229020	4-Point, Padded	\$ 175
<input type="checkbox"/>	1229030	Bodypoint - Chest Strap	\$ 85

Miscellaneous Options			
<input type="checkbox"/>	1230050	Protective Lower Padding	\$ 425
<input type="checkbox"/>	1230110	IV Pole	\$ 295
Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.			
<input type="checkbox"/>	1230120	O2 Holder	\$ 295
Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.			
<input type="checkbox"/>	1230170	Adhesive Back Hook Strip, 2" x 30"	\$ 10
For Cushion mounting.			