

Purchase Order	Tag
Dealer/Organization Name	
Purchaser Name and Contact Information	
Mailing Address	
Shipping Address	

Please send your completed order form by email to info@pdgmobility.com or by fax to 877-323-9093.

Order Notes

Frame Type
<input checked="" type="checkbox"/> 1201000 Stellar 45° Manual Tilt Wheelchair <i>E1161</i> \$ 2,995

Frame Options
<input type="checkbox"/> 1202000 Standard Base and Seat Frame 250 lbs Weight Capacity.
<input type="checkbox"/> 1202050 Weight Capacity Upgrade <i>K0108</i> \$ 200 Increases Weight Capacity to 300 lbs through addition of a secondary Gas Strut. Available for Seat Width of up to 22".
<input type="checkbox"/> 1202100 Heavy Duty Base and Seat Frame <i>K0108</i> \$ 300 450 lbs Weight Capacity. Requires 20" Composite or 24" Heavy Duty Spoked Wheels. Fixed Axles only. Required for 23" & 24" widths.
<input type="checkbox"/> 1202200 Extended Base with Rear Floor Jacks <i>K0108</i> \$ 600 Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Incompatible with Attendant wheel lock, HD Anti-Tipper, IV Pole, O2 Holder.
<input type="checkbox"/> 1202201 Shortened Strut Bracket \$ 125 Increases heel clearance for foot propelling. Available only for 23" & 24" widths.
<input type="checkbox"/> 1202300 Transit Tie Down System (TTDS) TTDS includes 4 factory installed-transport brackets. TTDS has been tested according to ANSI/RENSA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and 23" & 24" widths.

Seat Width
Standard Seat Widths
<input type="checkbox"/> 1203014 14" Seat Width
<input type="checkbox"/> 1203015 15" Seat Width
<input type="checkbox"/> 1203016 16" Seat Width
<input type="checkbox"/> 1203017 17" Seat Width
<input type="checkbox"/> 1203018 18" Seat Width
<input type="checkbox"/> 1203019 19" Seat Width
<input type="checkbox"/> 1203020 20" Seat Width <i>E2201</i> \$ 325
<input type="checkbox"/> 1203021 21" Seat Width <i>E2201</i> \$ 325
<input type="checkbox"/> 1203022 22" Seat Width <i>E2201</i> \$ 325

Custom Seat Widths
Requires 1202100 Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.
<input type="checkbox"/> 1203023 23" Seat Width <i>E2201</i> \$ 325
<input type="checkbox"/> 1203024 24" Seat Width <i>E2202</i> \$ 525
Use Stellar HD Order Form for Seat Widths beyond 24"

Seat Depth
Standard Seat Depths (Adjustable)
<input type="checkbox"/> 1204016 16" Seat Depth
<input type="checkbox"/> 1204017 17" Seat Depth
<input type="checkbox"/> 1204018 18" Seat Depth
<input type="checkbox"/> 1204019 19" Seat Depth
<input type="checkbox"/> 1204020 20" Seat Depth <i>E2203</i> \$ 475
Custom Seat Depths
<input type="checkbox"/> 1204021 21" Seat Depth <i>E2203</i> \$ 475
<input type="checkbox"/> 1204022 22" Seat Depth <i>E2204</i> \$ 750

Frame Color
PDG Primary Colors
<input type="checkbox"/> 1205000 Tracer Red
<input type="checkbox"/> 1205001 Gloss Black
<input type="checkbox"/> 1205002 Silver
<input type="checkbox"/> 1205003 Black Shimmer
<input type="checkbox"/> 1205004 Orange Fusion
<input type="checkbox"/> 1205005 Karl's Imagination
<input type="checkbox"/> 1205006 Blue Shimmer
<input type="checkbox"/> 1205007 Gloss White
<input type="checkbox"/> 1205008 Dark Titanium
<input type="checkbox"/> 1205009 PDG Premium Colors \$ 150
<input type="checkbox"/> 1205010 Geyser Green
<input type="checkbox"/> 1205011 Lollipop Blue
<input type="checkbox"/> 1205011 Lollipop Red
<input type="checkbox"/> 1205200 Custom Color \$ 300 Contact PDG for specifications.

PDG Product Design Group Inc.
318 East Kent Avenue South
Vancouver BC V5X 4N6 Canada
Phone 604-323-9220 | 888-858-4422

Casters		
Urethane		
8" Casters are compatible only with 16" and 24" Rear Wheels.		
<input type="checkbox"/>	1211061	6" x 1.25" Caster size
<input type="checkbox"/>	1211063	6" x 1.50" Caster size \$ 85
<input type="checkbox"/>	1211081	8" x 1.25" Caster size
<input type="checkbox"/>	1211084	8" x 2.00" Caster size \$ 135
Pneumatic		
8" Casters are compatible only with 16" and 24" Rear Wheels.		
<input type="checkbox"/>	1211161	6" x 1.25" Caster size \$ 110
<input type="checkbox"/>	1211181	8" x 1.25" Caster size \$ 120
<input type="checkbox"/>	1211184	8" x 2.00" Caster size \$ 150
<input type="checkbox"/>	1211260	6" Soft Roll \$ 95
<input type="checkbox"/>	1211400	Ultra Sport Fork E1015 \$ 395
Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.		

Rear Wheels		
Composite		
<input type="checkbox"/>	1212012	12" Wheel size
<input type="checkbox"/>	1212016	16" Wheel size
<input type="checkbox"/>	1212020	20" Wheel size
<input type="checkbox"/>	1212022	22" Wheel size
<input type="checkbox"/>	1212024	24" Wheel size
		Spoke \$ 75
<input type="checkbox"/>	1212120	20" Wheel size
<input type="checkbox"/>	1212122	22" Wheel size
<input type="checkbox"/>	1212124	24" Wheel size
<input type="checkbox"/>	1212224	24" Heavy Duty Spoke \$ 350
Required for Heavy Duty Package.		

Axle Type	
<input type="checkbox"/>	1213000 Quick Release Rear Axles
<input type="checkbox"/>	1213010 Fixed Rear Axles

Rear Tires		
<input type="checkbox"/>	1214000	Urethane
<input type="checkbox"/>	1214010	Pneumatic \$ 60
Pneumatic Tires may result in a finished Seat Height 1/2" greater than stated.		
<input type="checkbox"/>	1214020	Pneumatic with Flat-Free Inserts E2213 \$ 95
Pneumatic tires may result in a finished Seat Height 1/2" greater than stated. Not compatible with Spoked Wheels.		

Hand Rims		
<input type="checkbox"/>	1215000	No Hand Rims
<input type="checkbox"/>	1215010	Standard 20", 22", and 24" Rear Wheels only.
<input type="checkbox"/>	1215020	Plastic Coated 20", 22", and 24" Rear Wheels only. \$ 95

Wheel Locks		
<input type="checkbox"/>	1216010	Push to Lock
<input type="checkbox"/>	1216020	Pull to Lock
<input type="checkbox"/>	1216030	Attendant Wheel Lock (Push to Lock) \$ 85
Not compatible with TTDS (23" & 24")		
<input type="checkbox"/>	1216040	User & Attendant Wheel Lock \$ 125
Available only on 20", 22", or 24" Rear Wheels.		
<input type="checkbox"/>	1216100	Wheel Lock Extensions E0961 \$ 60
Not compatible with Attendant Wheel Lock.		
<input type="checkbox"/>	1216200	Hill Holders \$ 350
Requires Pneumatic Rear Tires with Flat-Free Inserts.		

Anti-Tippers		
<input checked="" type="checkbox"/>	1217010	Rear Anti-Tippers E0971 \$ 60
Required Safety Option.		
<input type="checkbox"/>	1217020	Front Anti-Tippers \$ 300
<input type="checkbox"/>	1217030	Heavy Duty Rear Anti-Tipper \$ 600
Provides additional rearward stability for clients with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.		

Back Post Style		
<input type="checkbox"/>	1218010	Push Bar Back
<input type="checkbox"/>	1218020	Cane Handle Back
Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.		
Back Post Height		
Seat Pan to Handle		
<input type="checkbox"/>	1218425	25" Height
<input type="checkbox"/>	1218431	31" Height

Back Options		
		Angle/Height Adjustable Push Bar Handle \$ 150
Not compatible with Headrest Mounting Plate or Cane Handle Back.		
<input type="checkbox"/>	1219010	Button on Outside
<input type="checkbox"/>	1219011	Button on Inside
Seat to Back Angle		
<input type="checkbox"/>	1219090	Set at 90°
<input type="checkbox"/>	1219098	Set at 98°
<input type="checkbox"/>	1219106	Set at 106°
<input type="checkbox"/>	1219114	Set at 114°
<input type="checkbox"/>	1219122	Set at 122°
<input type="checkbox"/>	1219210	Backrest Extension 10" high \$ 120
Not compatible with Angle/Height Adjustable Push Bar Handle.		
<input type="checkbox"/>	1219220	Fold Down Back
Not compatible with Reclining Back.		
<input type="checkbox"/>	1219230	Spring Back Suspension \$ 105
Not compatible with Reclining Back.		
<input type="checkbox"/>	1219240	Dynamic Rocker Back Assembly K0108 \$ 2,995
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Upholstery must be Adjustable Tension. Finished Back Post Height 3" taller. Requires longer Frame.		

Back Upholstery			
<input type="checkbox"/>	1220000	No Back Upholstery	
		Standard	\$ 125
		Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	1220021	Set at 21"	
<input type="checkbox"/>	1220022	Set at 22"	
		Short (25" Back only)	\$ 125
		Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	1220118	Set at 18"	
<input type="checkbox"/>	1220119	Set at 19"	
<input type="checkbox"/>	1220120	Set at 20"	
		Tall (31" Back only)	\$ 125
		Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	1220224	Set at 24"	
<input type="checkbox"/>	1220225	Set at 25"	
<input type="checkbox"/>	1220226	Set at 26"	
		Adjustable Tension	\$ 250
		Back height measured from Seat Pan to top of Upholstery	
<input type="checkbox"/>	1220321	Set at 21"	
<input type="checkbox"/>	1220322	Set at 22"	
<input type="checkbox"/>	1220323	Set at 23"	

Backrest			
<input type="checkbox"/>	1221000	No Backrest	
		nxt Optima HA	
<input type="checkbox"/>	1221100	AC - 16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1221101	PC - 16"-20" width, 16"-20" height	E2620 \$ 690
<input type="checkbox"/>	1221105	AC - 22"-24" width, 16"-20" height	E2621 \$ 720
<input type="checkbox"/>	1221106	PC - 22"-24" width, 16"-20" height	E2621 \$ 720
		Matrx Elite	
<input type="checkbox"/>	1221110	14"-20" width, 10"-20" height	E2620 \$ 825
<input type="checkbox"/>	1221112	HD - 22"-26" width, 16"-20" height	E2621 \$ 1,095
<input type="checkbox"/>	1221114	HD - 26"-30" width, 16"-20" height	E2621 \$ 1,595
<input type="checkbox"/>	1221116	Deep - 14"-20" width, 12"-20" height	E2620 \$ 950
		ROHO Agility	
<input type="checkbox"/>	1221120	Min - 14"-20" width, 10"-18" height	E2620 \$ 825
<input type="checkbox"/>	1221122	Mid - 14"-20" width, 10"-18" height	E2620 \$ 850
<input type="checkbox"/>	1221124	Max - 14"-20" width, 10"-18" height	E2620 \$ 950
		Varilite Icon	
<input type="checkbox"/>	1221130	Tall - 14"-20" width, 14"-24" height	E2615 \$ 850
<input type="checkbox"/>	1221132	Deep - 14"-20" width, 17"-24" height	E2620 \$ 895
		Backrest Size	
	1221400	Specify Size (Width" x Height") _____	

Tilt/Recline Options			
		Attendant Tilt only	
<input type="checkbox"/>	1223011	Installed on Right Back Post	
<input type="checkbox"/>	1223012	Installed on Left Back Post	
		User Tilt only	\$ 100
		Not available on widths greater than 22". Full Length Arm Pads only.	
<input type="checkbox"/>	1223021	Installed on Right Armrest	
<input type="checkbox"/>	1223022	Installed on Left Armrest	
		Attendant & User Tilt	\$ 195
		Not available on widths greater than 22". Full Length Arm Pads only.	
<input type="checkbox"/>	1223031	Back Post mount - Right side	
<input type="checkbox"/>	1223032	Back Post mount - Left side	
<input type="checkbox"/>	1223041	Armrest mount - Right side	
<input type="checkbox"/>	1223042	Armrest mount - Left side	
<input type="checkbox"/>	1023110	Trigger Extension	\$ 50
		Tilt Limit	
<input type="checkbox"/>	1223220	Set at 20°	
<input type="checkbox"/>	1223225	Set at 25°	
<input type="checkbox"/>	1223230	Set at 30°	
<input type="checkbox"/>	1223235	Set at 35°	
<input type="checkbox"/>	1223240	Set at 40°	
<input type="checkbox"/>	1223245	Set at 45°	
<input type="checkbox"/>	1223410	Reclining Backrest (Manual)	E1225 \$ 830
		Maximum recline = 30°. Not compatible with Headrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, 20"-22" seat depth. Requires longer Frame.	

Headrest			
<input type="checkbox"/>	1224000	No Headrest	
<input type="checkbox"/>	1224010	Headrest Mounting Plate	\$ 75
		Not compatible with Head Rest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	
<input type="checkbox"/>	1224100	nxt MultiFit Head Support Hardware	E1028 \$ 155
		nxt MultiFit Head Support Pad	E0955 \$ 195
<input type="checkbox"/>	1224105	Regular (9" x 5")	
<input type="checkbox"/>	1224106	Large (10" x 5")	
<input type="checkbox"/>	1224107	X-Large (14" x 5")	
		Stealth Headrest Hardware	E1028
<input type="checkbox"/>	1224105	Multi-Axis HMO475	\$ 150
<input type="checkbox"/>	1224106	Link TWB480	\$ 325
		Stealth Headrest Pad	E0955 \$ 200
<input type="checkbox"/>	1224105	10" width	
<input type="checkbox"/>	1224106	14" width	

Front Rigging			
<input type="checkbox"/>	1225000	No Front Rigging	
<input type="checkbox"/>	1225001	Sector Blocks only	
<input type="checkbox"/>	1225110	Plug In 70°	
<input type="checkbox"/>	1225120	Plug In Adjustable Knee Angle	\$ 975
<input type="checkbox"/>	1225130	Swing Away 60°	
<input type="checkbox"/>	1225140	Swing Away 70°	
<input type="checkbox"/>	1225200	Swing Away Elevating Legrest	E0990 \$ 355

Foot Plates			
<input type="checkbox"/>	1226010	Composite	
<input type="checkbox"/>	1226020	Angle Adjustable	K0040 \$ 140
<input type="checkbox"/>	1226030	One Piece Fixed	\$ 350
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1226040	One Piece Flip-up	\$ 395
<input type="checkbox"/>	1226050	Interlocking Angle Adjustable	\$ 195
Seat Widths of 14"-21" only.			

Footrest/Legrest Accessories			
<input type="checkbox"/>	1227010	Swing Away 60° & 70° Calf Pad	\$ 135
Not compatible with 60° & 70° Gel Pad.			
<input type="checkbox"/>	1227020	Swing Away 60° & 70° Gel Pad	\$ 100
Not compatible with 60° & 70° Swing Away Calf Pad.			
<input type="checkbox"/>	1227030	Swing Away ELR Gel Knee Pad	\$ 100
<input type="checkbox"/>	1227040	Adjustable Length Heel Loops	E0951 \$ 50
<input type="checkbox"/>	1227050	Bodypoint - Padded Calf Strap	K0038 \$ 105
Sizes greater than 22" are not padded.			
Miller's Dynamic Footrest Mounts			
<input type="checkbox"/>	1227060	For Plug In Front Rigging	\$ 190
<input type="checkbox"/>	1227065	For Swing Away 70° Front Rigging	\$ 380

Extension Tubes			
Regular Extension Tubes			
<input type="checkbox"/>	1228010	Standard (13"-17")	
<input type="checkbox"/>	1228020	Long (15"-19")	\$ 75
Extra Strength Extension Tubes (Solid)			
Available on Plug In Front Rigging only.			
<input type="checkbox"/>	1228110	Standard (13"-17")	\$ 100
<input type="checkbox"/>	1228120	Long (15"-19")	

Positioning Straps			
Bodypoint - Hip Belt E0978			
<input type="checkbox"/>	1229010	Center-Pull Auto Style	\$ 105
<input type="checkbox"/>	1229015	Center-Pull Auto Style, Padded	\$ 155
<input type="checkbox"/>	1229020	4-Point, Padded	\$ 175
<input type="checkbox"/>	1229030	Bodypoint - Chest Strap	E0960 \$ 85

Miscellaneous Options			
<input type="checkbox"/>	1230050	Protective Lower Padding	\$ 425
<input type="checkbox"/>	1230110	IV Pole	E1029 \$ 295
Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.			
<input type="checkbox"/>	1230120	O2 Holder	E2208 \$ 295
Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.			
<input type="checkbox"/>	1230170	Adhesive Back Hook Strip, 2" x 30"	\$ 10
For Cushion mounting.			