

FUZE T20



The *Fuze T20* offers up to 20° of tilt with ease of propelling for occupants and caregivers. The chair has a sleek aluminum frame and is *easy to set up and adjust*. Its forward pivot location provides *low knee-rise* during tilt.



Hand
Propulsion



Foot
Propulsion



Self-Tilt



Desk/Table
Accessible



Transit Tie
Downs

Fuze T20 A sleek 20° knee pivot tilt chair

visit pdgmobility.com/t20
for more information

Available options See order form

- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ Adjustable width frame
- ✓ Can be converted to a Fuze T50



DYNAMIC TILT
Up to 20°



WEIGHT CAPACITY
250 lbs



SEAT WIDTH
14" to 22"



SEAT DEPTH
13" to 22"



SEAT HEIGHT
13" to 20"



BACK SEAT ANGLE
84° to 120°



OVERALL WIDTH
+9.5" to 10.5"



CHAIR WEIGHT
42 lbs

