

FUZE T50n



The *Fuze T50n* has a sleek aluminum frame with an infinitely adjustable tilt range of 50°. The *n* stands for *narrow*, referring to the chair's overall width of *seat width plus only 4"*. This makes the T50n substantially narrower than other tilt chairs.

Available options See order form

- ✓ Transit tie downs WC-19 & ISO-7176-19
- ✓ 30° manual dynamic recline
- ✓ Power tilt up to 300 lbs
- ✓ Adjustable width frame



The overall width of a T50n with a 20" seat is only 24".

By comparison, a standard T50 with the same seat is 28.5" wide (27.5" without handrims).



DYNAMIC TILT
Up to 50°



WEIGHT CAPACITY
Up to 300 lbs



SEAT WIDTH
16" to 24"



SEAT DEPTH
15" to 22"



SEAT HEIGHT
16" to 20"



BACK SEAT ANGLE
90° to 122°



OVERALL WIDTH
+4.0"



CHAIR WEIGHT
44 lbs



“a fully featured 50° tilt chair with the footprint of a transport chair”

Upgrade to a smaller size.
pdgmobility.com/narrow