



Torr Brown, Chief Engineer
PDG Product Design Group Inc.
103 – 318 East Kent Avenue South,
Vancouver, BC Canada
VSX 4N6

4th July 2017

Dear Torr,

Reference: Fatigue testing of the Product Design Group Inc. Stellar HD manual wheelchair to 272 kg (600 lbs) user weight, our Laboratory reference is 492721.

This letter is to advise you that this Laboratory has finished fatigue testing your product to AS/NZS 3695.1:2011 Wheelchairs part 1: Requirements and test methods for manual wheelchairs which includes references to the following standard:
AS/NZS ISO 7176 Section 8 – Requirements and test methods for static, impact and fatigue strength tests,
for 272 kg / 600 lbs user weight.

The Product Design Group Inc. Stellar HD tilting manual wheelchair as represented by the test sample, and detailed in the test report, passed the requirements of 2 drum and drop test fatigue strength testing. Full details of these tests can be found in the official report.

Please note that as this wheelchair has only been tested to the fatigue strength aspects of two drum and drop test, full compliance to the wheelchair standard cannot be granted by this test laboratory or claimed by you as manufacturer.

Please do not hesitate to contact us if you have any enquiries or comments regarding this testing.

Kind regards,

Wayne Wurfel
Senior Test Technician, Authorised signatory

NovitaTech (a division of Novita Children's Services)
171 Days Road, Regency Park, South Australia 5010
PO Box 2438, Regency Park, South Australia 5942
Telephone 1300 855 585 • Fax (08) 8243 8208
enquiries@novita.org.au • www.novitatech.org.au
ABN 53 137 537 636

