



BENTLEY LT-R

Order Form - Canada

Prices in Canadian Dollars. January 1, 2022

Customer and Order Information

Purchase Order *	Tag	v.0524
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

Base Chair

Bentley LT-R 20° Manual Tilt Wheelchair

<input checked="" type="checkbox"/> 2701000	Bentley 20° Manual Tilt Wheelchair	\$ 4,150
<input checked="" type="checkbox"/> 2701009	Lightweight Aluminum Frame Upgrade Reverse Wheel Layout enhances user access for self-propulsion. pdgmobility.com/blt-r	\$ 300
<input type="checkbox"/> ⓘ	Choose the standard Bentley with a steel frame for agitated or heavier users. See separate order form.	

Weight Capacity

<input checked="" type="checkbox"/> 2702000	250 lbs Chair Weight Capacity		
User Weight Range Required for setup of tilt uplift force.			
<input type="checkbox"/> 2702430	Up to 75 lbs	<input type="checkbox"/> 2702431	75 lbs to 100 lbs
<input type="checkbox"/> 2702432	100 lbs to 115 lbs	<input type="checkbox"/> 2702433	115 lbs to 150 lbs
<input type="checkbox"/> 2702434	150 lbs to 200 lbs	<input type="checkbox"/> 2702435	200 lbs to 250 lbs

Frame Options

<input type="checkbox"/> 2702300	Transit Tie Down System (TTDS)
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.	

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths					
<input type="checkbox"/> 2703014	14" Seat Width	<input type="checkbox"/> 2703015	15" Seat Width	<input type="checkbox"/> 2703016	16" Seat Width
<input type="checkbox"/> 2703017	17" Seat Width	<input type="checkbox"/> 2703018	18" Seat Width	<input type="checkbox"/> 2703019	19" Seat Width
<input type="checkbox"/> 2703020	20" Seat Width	<input type="checkbox"/> 2703021	21" Seat Width	<input type="checkbox"/> 2703022	22" Seat Width

Seat Depth

Standard Seat Depths (Adjustable)					
<input type="checkbox"/> 2704016	16" Seat Depth	<input type="checkbox"/> 2704017	17" Seat Depth	<input type="checkbox"/> 2704018	18" Seat Depth
<input type="checkbox"/> 2704019	19" Seat Depth	<input type="checkbox"/> 2704020	20" Seat Depth		
Custom Seat Depths					
<input type="checkbox"/> 2704021	21" Seat Depth	<input type="checkbox"/> 2704022	22" Seat Depth	\$ 400	

Frame Color

PDG Premium Colors					
<input type="checkbox"/> 2705000	Tracer Red	<input type="checkbox"/> 2705001	Gloss Black Shimmer	<input type="checkbox"/> 2705002	Silver
<input type="checkbox"/> 2705003	Fine Texture Black	<input type="checkbox"/> 2705004	Orange Fusion	<input type="checkbox"/> 2705005	Karl's Imagination
<input type="checkbox"/> 2705006	Blue Shimmer	<input type="checkbox"/> 2705007	Gloss White	<input type="checkbox"/> 2705008	Dark Titanium
<input type="checkbox"/> 2705009	Geyser Green	<input type="checkbox"/> 2705010	Lollipop Blue	<input type="checkbox"/> 2705011	Lollipop Red
<input type="checkbox"/> 2705012	Bubble Gum Pink	<input type="checkbox"/> 2705013	Berry Purple	<input type="checkbox"/> 2705014	Daffodil Yellow
<input type="checkbox"/> 2705200	Custom Color Contact PDG for specifications.				\$ 300

Armrests

2706000 **No Armrests**

Height Adjustable Armrests (Single Post) \$ 295

2706111 8"-11" 2706114 10"-14"

2706100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
Seat Pan to top of Armrest: _____ Does not include height of Armpad

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2706401 Rear 2706402 Mid 2706403 Front

Armpads

Regular Armpads

2707100 Desk Length - 2" x 10" 2707101 Full Length - 2" x 14"

Waterfall Armpads \$ 115

2707105 Desk Length - 2" x 10" 2707106 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

2707110 Desk Length - 3" x 10" 2707111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

2707120 Full Length - 2" x 12" 2707121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125

2707130 Right (4.5" x 18") 2707131 Left (4.5" x 18")

Armrest Options

2707210 **Armrest Wideners** Move Armrests outward by 1" on each side. Do not increase seat size. \$ 195

2707310 **Clothing Guards** Available for 10"-14" Height range only. \$ 110

Seat Pan

2708000 **No Seat Pan**

2708001 **Solid Aluminum Seat Pan** \$ 135

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low						High						
	RW	12"	13"	14"	15"	16"	RW	14"	15"	16"	17"	18"	19"
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low						High						
	RW	12"	13"	14"	15"	16"	RW	14"	15"	16"	17"	18"	19"
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Large Wheels.

<input type="checkbox"/>	2711061	6" x 1.25" Caster size	<input type="checkbox"/>	2711081	8" x 1.25" Caster size	
<input type="checkbox"/>	2711063	6" x 1.50" Caster size				\$ 85
<input type="checkbox"/>	2711084	8" x 2.00" Caster size				\$ 135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Large Wheels.

<input type="checkbox"/>	2711161	6" x 1.25" Caster size				\$ 110
<input type="checkbox"/>	2711181	8" x 1.25" Caster size				\$ 120
<input type="checkbox"/>	2711184	8" x 2.00" Caster size				\$ 150

<input type="checkbox"/>	2711260	Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters				\$ 275
--------------------------	---------	--	--	--	--	--------

<input type="checkbox"/>	2711400	Frog Legs Suspension Forks	Requires 6" Soft-Roll Caster.			\$ 595
--------------------------	---------	-----------------------------------	-------------------------------	--	--	--------

Large Wheels

Composite Wheels

<input type="checkbox"/>	2712020	20" Wheel size	<input type="checkbox"/>	2712022	22" Wheel size	<input type="checkbox"/>	2712024	24" Wheel size	
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	--

Spoke Wheels

<input type="checkbox"/>	2712120	20" Wheel size	<input type="checkbox"/>	2712122	22" Wheel size	<input type="checkbox"/>	2712124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------



Large Wheel Position Indicate Large Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/>	2712401	Rear	<input type="checkbox"/>	2712402	Mid	<input type="checkbox"/>	2712403	Front	
--------------------------	---------	------	--------------------------	---------	-----	--------------------------	---------	-------	--

Axle Type

<input type="checkbox"/>	2713000	Quick Release Axles	
--------------------------	---------	----------------------------	--

<input type="checkbox"/>	2713010	Fixed Axles	
--------------------------	---------	--------------------	--

Large Tires

<input type="checkbox"/>	2714000	Urethane Tires	
--------------------------	---------	-----------------------	--

<input type="checkbox"/>	2714010	Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated.	\$ 60
--------------------------	---------	------------------------	--	-------

<input type="checkbox"/>	2714020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$ 95
--------------------------	---------	---	---	-------

Handrims

<input type="checkbox"/>	2715000	No Handrims	
--------------------------	---------	--------------------	--

<input type="checkbox"/>	2715010	Standard Aluminum Handrims	
--------------------------	---------	-----------------------------------	--

<input type="checkbox"/>	2715020	Plastic Coated Handrims	\$ 150
--------------------------	---------	--------------------------------	--------

Wheel Locks & Brakes

User Wheel Locks

<input type="checkbox"/>	2716010	Push to Lock	<input type="checkbox"/>	2716020	Pull to Lock
--------------------------	---------	--------------	--------------------------	---------	--------------

Wheel Lock Options

<input type="checkbox"/>	2716100	Wheel Lock Extensions	\$ 60
--------------------------	---------	------------------------------	-------

<input type="checkbox"/>	2716200	Hill Holders	Requires Pneumatic Large Tires with Flat-Free Inserts.	\$ 450
--------------------------	---------	---------------------	--	--------

Anti-Tippers

<input checked="" type="checkbox"/>	2717010	Front Anti-Tippers	Required Safety Option.	\$ 75
-------------------------------------	---------	---------------------------	-------------------------	-------



BENTLEY LT-R

Order Form - Canada

Prices in Canadian Dollars. January 1, 2022

Back Posts

Back Post Style

2718010 Push Bar Back

Back Post Height

 Measured from the Seat Pan to the Handle.

2718425 Short (22") 2718431 Tall (25")

Seat to Back Angle

 Seat Angles above 98° may cause interference at Seat Depths above 19".

2719090 Set at 90° 2719098 Set at 98° 2719106 Set at 106°
 2719114 Set at 114° 2719122 Set at 122°

Back Options

Angle/Height Adjustable Push Bar Handle

 Not compatible with Headrest Mounting Plate.

\$ 150

2719010 Button on Outside 2719011 Button on Inside

10" Backrest Extension

 Requires the Angle/Height Adjustable Push Bar Handle.

\$ 150

Headrest Mounting Plate

\$ 95

Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar. Not required if mounting to a solid Back.

2719220 **Fold Down Back** Not compatible with Reclining Back.

Back Upholstery

2720000 **No Back Upholstery**

Standard Upholstery

 Not available with Cane Handle Back.

\$ 195

2720418 Short (18" Height) 2720421 Regular (21" Height) 2720424 Tall (24" Height)

2720450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back.

\$ 295

Position of Upholstery

 Select the desired position of the top of the Upholstery relative to the Seat Pan.

Style	22" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Style	25" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation

Attendant Tilt only

 Select trigger mounting location below.

2723011 Right Back Post 2723012 Left Back Post

User Tilt only

 Full Length Arm Pads only. Select trigger mounting location below.

\$ 100

2723021 Right Armrest 2723022 Left Armrest

Attendant & User Tilt

 Full Length Arm Pads only. Select trigger mounting locations below.

\$ 195

2723031 Right Back Post 2723032 Left Back Post

2723041 Right Armrest 2723042 Left Armrest

2723110 **Tilt Trigger Extension**

\$ 50

2723400 **Power Tilt** 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt
 Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both. Not available for 12" STF Height.

\$ 1,795

Toggle Switch 2723421 Installed on Right side 2723422 Installed on Left side
 Button Switch 2723431 Installed on Right side 2723432 Installed on Left side
 Stealth Egg Switch 2723440 Not installed

Optional for Egg Switch: 2723445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625

Specify Headrest to be used: _____

Recline

<input type="checkbox"/>	2723310	Reclining Backrest (Manual)	\$ 895
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.			

Front Rigging

<input type="checkbox"/>	2725000	No Front Rigging	
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out			
<input type="checkbox"/>	2725410	Fixed 70° Not available for 12" STF Height.	\$ 95
<input type="checkbox"/>	2725420	Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) 12" STF Height is only compatible with 0°-60°. Adjust angle by loosening one bolt. Tool-free leg-length adjustment.	\$ 295
<input type="checkbox"/>	2725110	Standard Plug In 70° Front Rigging Legacy Design. Not compatible with Easy-Out.	


Footplates

<input type="checkbox"/>	2726010	Composite Footplates	
<input type="checkbox"/>	2726020	Aluminum Angle Adjustable Footplates	\$ 165
<input type="checkbox"/>	2726030	One Piece Fixed Footplate Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	2726040	One Piece Flip-up Footplate Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	2726050	Interlocking Angle Adjustable Footplates Available for Seat Widths up to 21" only.	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/>	2727012	Calf Pads Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	2727032	Gel Knee Pads	\$ 100
<input type="checkbox"/>	2727040	Length Adjustable Heel Loops	\$ 50
<input type="checkbox"/>	2727050	Bodypoint Padded Calf Strap	\$ 195

Front Rigging Extension Tubes

Regular Extension Tubes			
<input type="checkbox"/>	2728010	Short 14"-18" for Easy-Out Front Rigging 13"-17" for Standard Plug In Front Rigging	
<input type="checkbox"/>	2728020	Long 18"-24" for Easy-Out Front Rigging 15"-19" for Standard Plug In Front Rigging	\$ 75
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used. Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.		

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Large Wheels. Made from durable and lightweight aluminum.			
<input type="checkbox"/>	2730020	Black <input type="checkbox"/> 2730021 Frame Color	\$ 195
<input type="checkbox"/>	2730050	Protective Lower Padding	\$ 795
<input type="checkbox"/>	2730080	Gel Hip Guides Reduces available Seat Width by 2".	\$ 395
<input type="checkbox"/>	2730110	IV Pole	\$ 295
<input type="checkbox"/>	2730120	O2 Holder	\$ 395
<input type="checkbox"/>	2730170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10



BENTLEY LT-R

Order Form - Canada

Prices in Canadian Dollars. January 1, 2022

On-Chair Seating & Positioning



Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.