



BENTLEY LT

Order Form - Canada

Prices in Canadian Dollars. January 1, 2020

Customer and Order Information

Purchase Order *	Tag	v.1206
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in *required

Frame Type

<input checked="" type="checkbox"/> 2401000	Bentley 20° Manual Tilt Wheelchair	\$ 3,995
<input checked="" type="checkbox"/> 2401009	Lightweight Aluminum Frame Upgrade (Bentley LT) Aluminum Frame is not recommended for agitated users.	\$ 300
ⓘ For agitated or heavier users, choose the standard Bentley with steel frame. See separate order form.		

Frame Options

<input checked="" type="checkbox"/> 2402000	Standard Base and Seat Frame 250 lbs Weight Capacity.	
<input type="checkbox"/> 2402300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.	
User Weight Range Required for setup of tilt uplift force.		
<input type="checkbox"/> 2402430 Up to 75 lbs	<input type="checkbox"/> 2402431 75 lbs to 100 lbs	<input type="checkbox"/> 2402432 100 lbs to 115 lbs
<input type="checkbox"/> 2402433 115 lbs to 150 lbs	<input type="checkbox"/> 2402434 150 lbs to 200 lbs	<input type="checkbox"/> 2402435 200 lbs to 250 lbs

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths		
<input type="checkbox"/> 2403014 14" Seat Width	<input type="checkbox"/> 2403015 15" Seat Width	<input type="checkbox"/> 2403016 16" Seat Width
<input type="checkbox"/> 2403017 17" Seat Width	<input type="checkbox"/> 2403018 18" Seat Width	<input type="checkbox"/> 2403019 19" Seat Width
<input type="checkbox"/> 2403020 20" Seat Width	<input type="checkbox"/> 2403021 21" Seat Width	<input type="checkbox"/> 2403022 22" Seat Width

Seat Depth

Standard Seat Depths (Adjustable)		
<input type="checkbox"/> 2404016 16" Seat Depth	<input type="checkbox"/> 2404017 17" Seat Depth	<input type="checkbox"/> 2404018 18" Seat Depth
<input type="checkbox"/> 2404019 19" Seat Depth	<input type="checkbox"/> 2404020 20" Seat Depth	
Custom Seat Depths		
<input type="checkbox"/> 2404021 21" Seat Depth	<input type="checkbox"/> 2404022 22" Seat Depth	\$ 350

Frame Color

PDG Premium Colors		
<input type="checkbox"/> 2405000 Tracer Red	<input type="checkbox"/> 1005001 Gloss Black Shimmer	<input type="checkbox"/> 1005002 Silver
<input type="checkbox"/> 2405003 Fine Texture Black	<input type="checkbox"/> 1005004 Orange Fusion	<input type="checkbox"/> 1005005 Karl's Imagination
<input type="checkbox"/> 2405006 Blue Shimmer	<input type="checkbox"/> 1005007 Gloss White	<input type="checkbox"/> 1005008 Dark Titanium
<input type="checkbox"/> 2405009 Geyser Green	<input type="checkbox"/> 1005010 Lollypop Blue	<input type="checkbox"/> 1005011 Lollypop Red
<input type="checkbox"/> 2405012 Bubble Gum Pink	<input type="checkbox"/> 1005013 Berry Purple	<input type="checkbox"/> 1005014 Daffodil Yellow
<input type="checkbox"/> 2405200 Custom Color Contact PDG for specifications.		\$ 300

Armrest Type

2406000 **No Armrests**

Height Adjustable Armrests (Single Post) \$ 275

2406111 8"-11" 1006114 10"-14"

2406100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2406401 Rear 2406402 Mid 2406403 Front

Armpads & Armrest Options

Regular Armpads

2407100 Desk Length - 2" x 10" 1007101 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

2407110 Desk Length - 3" x 10" 1007111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

2407120 Full Length - 2" x 12" 1007121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125

2407130 Right (4.5" x 18") 1007131 Left (4.5" x 18")

2407310 **Clothing Guards** Available for 10"-14" Height range only. \$ 110

Seat Pan

2408000 **No Seat Pan**

2408001 **Solid Aluminum Seat Pan** \$ 120

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.

6" Casters	RW	Low					High					Extra-High				
	12"	13"	14"	15"			15"	16"	17"			15"	16"	17"	18"	
16"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"
20"	13"	14"	15"			15"	16"	17"			15"	16"	17"	18"		
22"		14"	15"	16"				16"	17"	18"			16"	17"	18"	19"
24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"

8" Casters	RW	Low			High			Extra-High								
	16"			15"	16"	17"			17"	18"	19"			17"	18"	19"
24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"

Wheel size _____
Caster size _____
STF height _____
Range _____

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

2411061 6" x 1.25" Caster size 2411081 8" x 1.25" Caster size

2411063 6" x 1.50" Caster size \$ 85

2411084 8" x 2.00" Caster size \$ 135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

2411161 6" x 1.25" Caster size \$ 110

2411181 8" x 1.25" Caster size \$ 120

2411184 8" x 2.00" Caster size \$ 150

2411260 **Frog Legs 6" x 1.40" Soft Roll Casters** \$ 95

2411400 **Frog Legs Suspension Forks** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. \$ 395

Rear Wheels


Composite Wheels

- | | | |
|---|---|---|
| <input type="checkbox"/> 2412012 12" Wheel size | <input type="checkbox"/> 1012016 16" Wheel size | <input type="checkbox"/> 1012020 20" Wheel size |
| <input type="checkbox"/> 2412022 22" Wheel size | <input type="checkbox"/> 1012024 24" Wheel size | |

Spoke Wheels

- | | | |
|---|---|---|
| <input type="checkbox"/> 2412120 20" Wheel size | <input type="checkbox"/> 1012122 22" Wheel size | <input type="checkbox"/> 1012124 24" Wheel size |
|---|---|---|

\$ 75

 **Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- | | | |
|---------------------------------------|--------------------------------------|--|
| <input type="checkbox"/> 2412401 Rear | <input type="checkbox"/> 2412402 Mid | <input type="checkbox"/> 2412403 Front |
|---------------------------------------|--------------------------------------|--|

Axle Type

- | |
|--|
| <input type="checkbox"/> 2413000 Quick Release Rear Axles |
| <input type="checkbox"/> 2413010 Fixed Rear Axles |

Rear Tires

- | | | |
|--|---|-------|
| <input type="checkbox"/> 2414000 Urethane Tires | | |
| <input type="checkbox"/> 2414010 Pneumatic Tires | May result in a finished Seat Height 1/2" greater than stated. | \$ 60 |
| <input type="checkbox"/> 2414020 Pneumatic Tires with Flat-Free Inserts | May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. | \$ 95 |

Hand Rims

Handrims are only available with 20", 22", and 24" Rear Wheels.

- | | | |
|--|--|--------|
| <input type="checkbox"/> 2415000 No Handrims | | |
| <input type="checkbox"/> 2415010 Standard Aluminum Handrims | | |
| <input type="checkbox"/> 2415020 Plastic Coated Handrims | | \$ 135 |

Wheel Locks

User Wheel Locks

- | | | |
|--|---|--------|
| <input type="checkbox"/> 2416010 Push to Lock | <input type="checkbox"/> 2416020 Pull to Lock | |
| <input type="checkbox"/> 2416030 Attendant Wheel Locks (Push to Lock) | Not compatible with TTDS. | \$ 85 |
| <input type="checkbox"/> 2416040 User & Attendant Wheel Locks | Available with 20", 22", or 24" Rear Wheels only. | \$ 125 |
| <input type="checkbox"/> 2416100 Wheel Lock Extensions | Not compatible with Attendant Wheel Lock. | \$ 60 |
| <input type="checkbox"/> 2416200 Hill Holders | Requires Pneumatic Rear Tires with Flat-Free Inserts. | \$ 450 |
| <input type="checkbox"/> 2416300 Drum Brakes | Available with 20" Rear Wheels only. | \$ 995 |

Anti-Tippers

- | | | |
|--|-------------------------|-------|
| <input checked="" type="checkbox"/> 2417010 Rear Anti-Tippers | Required Safety Option. | \$ 60 |
|--|-------------------------|-------|

Back Post Style

- | | | |
|--|---|---|
| <input checked="" type="checkbox"/> 2418010 Push Bar Back | | |
| Back Post Height Please note: We have updated the descriptions/measurements below. The actual back post heights remain unchanged. | | |
| <input type="checkbox"/> 2418425 Short (22") <i>Previously 25"</i> | <input type="checkbox"/> 2418431 Tall (25") <i>Previously 31"</i> | Measured from the Seat Pan to the Handle. |

Back Options

- | | | |
|---|---|--|
| Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate. | | \$ 150 |
| <input type="checkbox"/> 2419010 Button on Outside | <input type="checkbox"/> 1019011 Button on Inside | |
| Seat to Back Angle | | |
| <input type="checkbox"/> 2419090 Set at 90° | <input type="checkbox"/> 1019098 Set at 98° | <input type="checkbox"/> 1019106 Set at 106° |
| <input type="checkbox"/> 2419114 Set at 114° | <input type="checkbox"/> 1019122 Set at 122° | |

See next page for more options >

<input type="checkbox"/>	2419210	10" Backrest Extension	Requires the Angle/Height Adjustable Push Bar Handle.	\$ 150
<input type="checkbox"/>	2419212	Headrest Mounting Plate	Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar. Not required if mounting to a solid Back.	\$ 95
<input checked="" type="checkbox"/>	2419220	Fold Down Back	Not compatible with Reclining Back.	

Back Upholstery

<input type="checkbox"/>	2420000	No Back Upholstery		\$ 195																																																																																					
Standard Upholstery Not available with Cane Handle Back.																																																																																									
<input type="checkbox"/>	2420418	Short (18" Height)	<input type="checkbox"/> 2420421 Regular (21" Height)	<input type="checkbox"/> 2420424 Tall (24" Height)																																																																																					
<input type="checkbox"/>	2420450	Tension Adjustable Upholstery	Height of Upholstery is 16". Not available with Cane Handle Back.	\$ 295																																																																																					
Position of Upholstery Position of top of Upholstery relative to the Seat Pan. Circle below of select from the drop-down menu on the right.																																																																																									
		<table border="1"> <thead> <tr> <th>Upholstery</th> <th colspan="8">22" Back Post</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td>17"</td><td>18"</td><td>19"</td><td>20"</td><td>21"</td><td>22"</td><td></td><td></td> </tr> <tr> <td>Regular</td> <td></td><td></td><td></td><td>20"</td><td>21"</td><td>22"</td><td></td><td></td> </tr> <tr> <td>Tension Adj</td> <td></td><td></td><td>19"</td><td>20"</td><td>21"</td><td>22"</td><td></td><td></td> </tr> </tbody> </table>	Upholstery	22" Back Post								Short	17"	18"	19"	20"	21"	22"			Regular				20"	21"	22"			Tension Adj			19"	20"	21"	22"			<table border="1"> <thead> <tr> <th>Upholstery</th> <th colspan="8">25" Back Post</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td>17"</td><td>18"</td><td>19"</td><td>20"</td><td>21"</td><td>22"</td><td>23"</td><td>24"</td><td>25"</td> </tr> <tr> <td>Regular</td> <td></td><td></td><td></td><td>20"</td><td>21"</td><td>22"</td><td>23"</td><td>24"</td><td>25"</td> </tr> <tr> <td>Tall</td> <td></td><td></td><td></td><td></td><td></td><td></td><td>23"</td><td>24"</td><td>25"</td> </tr> <tr> <td>Tension Adj</td> <td></td><td></td><td>19"</td><td>20"</td><td>21"</td><td>22"</td><td>23"</td><td>24"</td><td>25"</td> </tr> </tbody> </table>	Upholstery	25" Back Post								Short	17"	18"	19"	20"	21"	22"	23"	24"	25"	Regular				20"	21"	22"	23"	24"	25"	Tall							23"	24"	25"	Tension Adj			19"	20"	21"	22"	23"	24"	25"	Height _____
Upholstery	22" Back Post																																																																																								
Short	17"	18"	19"	20"	21"	22"																																																																																			
Regular				20"	21"	22"																																																																																			
Tension Adj			19"	20"	21"	22"																																																																																			
Upholstery	25" Back Post																																																																																								
Short	17"	18"	19"	20"	21"	22"	23"	24"	25"																																																																																
Regular				20"	21"	22"	23"	24"	25"																																																																																
Tall							23"	24"	25"																																																																																
Tension Adj			19"	20"	21"	22"	23"	24"	25"																																																																																

Tilt/Recline Options

Attendant Tilt only Select trigger mounting location below.				
<input type="checkbox"/>	2423011	Right Back Post	<input type="checkbox"/> 2423012 Left Back Post	
User Tilt only Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below.				
<input type="checkbox"/>	2423021	Right Armrest	<input type="checkbox"/> 2423022 Left Armrest	\$ 100
Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below.				
<input type="checkbox"/>	2423031	Right Back Post	<input type="checkbox"/> 2423032 Left Back Post	\$ 195
<input type="checkbox"/>	2423041	Right Armrest	<input type="checkbox"/> 2423042 Left Armrest	
<input type="checkbox"/>	2423110	Trigger Extension		\$ 50
<input type="checkbox"/>	2423310	Reclining Backrest (Manual)	30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19".	\$ 895
<input type="checkbox"/>	2423400	Power Tilt	250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.	\$ 1,795
		Toggle Switch	<input type="checkbox"/> 2423421 Installed on Right side	<input type="checkbox"/> 2423422 Installed on Left side
		Button Switch	<input type="checkbox"/> 2423431 Installed on Right side	<input type="checkbox"/> 2423432 Installed on Left side
		Stealth Egg Switch	<input type="checkbox"/> 2423440 Not installed	
		Optional for Egg Switch:	<input type="checkbox"/> 2423445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models. \$ 625
<i>Specify Headrest to be used: _____</i>				

Front Rigging

<input type="checkbox"/>	2425000	No Front Rigging	<input type="checkbox"/> 2425001 Sector Blocks only	
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out				
<input type="checkbox"/>	2425410	Fixed 70°		\$ 95
<input type="checkbox"/>	2425420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°) Adjust angle by loosening one bolt. Tool-free leg-length adjustment.	\$ 295
Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.				
<input type="checkbox"/>	2425130	Swing Away 60°	<input type="checkbox"/> 2425140 Swing Away 70°	
<input type="checkbox"/>	2425200	Swing Away Elevating Legrests	Includes Calf Pads.	\$ 375
<input type="checkbox"/>	2425300	1" Offset Hangers for Swing Away		\$ 75
<input type="checkbox"/>	2425110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.	


Footplates

<input type="checkbox"/>	2426010	Composite Footplates		
<input type="checkbox"/>	2426020	Aluminum Angle Adjustable Footplates		\$ 155
<input type="checkbox"/>	2426030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	2426040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	2426050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.	\$ 195

Footrest/Legrest Accessories

Calf Pads Not compatible with Gel Padding.				
<input type="checkbox"/>	2427012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	2427010	For Swing Away Front Rigging		\$ 135
<input type="checkbox"/>	2427020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pads				
<input type="checkbox"/>	2427032	For Easy-Out™ Angle Adjustable Front Rigging		\$ 200
<input type="checkbox"/>	2427030	For Swing Away ELR Front Rigging		\$ 100
<input type="checkbox"/>	2427040	Length Adjustable Heel Loops		\$ 50
<input type="checkbox"/>	2427050	Bodypoint Padded Calf Strap		\$ 195


Front Rigging Extension Tubes

Regular Extension Tubes				
<input type="checkbox"/>	2428010	Short	14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	2428020	Long	18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.			
	Seat Pan to Foot Plates: _____		Measured along the length of the Front Rigging.	

Miscellaneous Options

One-Arm Drive Requires 24" Composite Wheels with Urethane Tires. 18" STF minimum. Limits tilt. Not compatible with User & Attendant Wheel Locks, Quick Release Rear Axles.				\$ 3,500
<input type="checkbox"/>	2430010	Installed on Right side	<input type="checkbox"/> 1030011 Installed on Left side	
<input type="checkbox"/>	2430050	Protective Lower Padding		\$ 675
<input type="checkbox"/>	2430080	Gel Hip Guides	Reduces available Seat Width by 2".	\$ 375
<input type="checkbox"/>	2430110	IV Pole	Not compatible with Attendant Wheel Locks.	\$ 295
<input type="checkbox"/>	2430120	O2 Holder	Not compatible with Attendant Wheel Locks.	\$ 395
<input type="checkbox"/>	2430170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$ 10

On-Chair Seating & Positioning

	Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating			
---	---	--	--	--

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.