



BENTLEY LT

Order Form - Canada

Prices in Canadian Dollars. January 1, 2020

Customer and Order Information

| | | |
|--|------------------|---|
| Purchase Order * | Tag | v.0801 |
| Dealer/Provider Name * | | |
| Purchaser Name and Contact Information * | | |
| Shipping Address * | | |
| User Information (if available) | Weight _____ lbs | Height _____ ft/in *required |

Frame Type

| | | |
|---|--|----------|
| <input checked="" type="checkbox"/> 2401000 | Bentley 20° Manual Tilt Wheelchair | \$ 3,995 |
| <input checked="" type="checkbox"/> 2401009 | Lightweight Aluminum Frame Upgrade (Bentley LT) Aluminum Frame is not recommended for agitated users. | \$ 300 |
| ⓘ For agitated or heavier users, choose the standard Bentley with steel frame. See separate order form. | | |

Frame Options

| | | |
|---|--|---|
| <input checked="" type="checkbox"/> 2402000 | Standard Base and Seat Frame 250 lbs Weight Capacity. | |
| <input type="checkbox"/> 2402300 | Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks. | |
| User Weight Range Required for setup of tilt uplift force. | | |
| <input type="checkbox"/> 2402430 Up to 75 lbs | <input type="checkbox"/> 2402431 75 lbs to 100 lbs | <input type="checkbox"/> 2402432 100 lbs to 115 lbs |
| <input type="checkbox"/> 2402433 115 lbs to 150 lbs | <input type="checkbox"/> 2402434 150 lbs to 200 lbs | <input type="checkbox"/> 2402435 200 lbs to 250 lbs |

Seat Width Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

| | | |
|---|---|---|
| Standard Seat Widths | | |
| <input type="checkbox"/> 2403014 14" Seat Width | <input type="checkbox"/> 2403015 15" Seat Width | <input type="checkbox"/> 2403016 16" Seat Width |
| <input type="checkbox"/> 2403017 17" Seat Width | <input type="checkbox"/> 2403018 18" Seat Width | <input type="checkbox"/> 2403019 19" Seat Width |
| <input type="checkbox"/> 2403020 20" Seat Width | <input type="checkbox"/> 2403021 21" Seat Width | <input type="checkbox"/> 2403022 22" Seat Width |

Seat Depth

| | | |
|---|---|---|
| Standard Seat Depths (Adjustable) | | |
| <input type="checkbox"/> 2404016 16" Seat Depth | <input type="checkbox"/> 2404017 17" Seat Depth | <input type="checkbox"/> 2404018 18" Seat Depth |
| <input type="checkbox"/> 2404019 19" Seat Depth | <input type="checkbox"/> 2404020 20" Seat Depth | |
| Custom Seat Depths | | |
| <input type="checkbox"/> 2404021 21" Seat Depth | <input type="checkbox"/> 2404022 22" Seat Depth | \$ 350 |

Frame Color

| | | |
|--|--|---|
| PDG Premium Colors | | |
| <input type="checkbox"/> 2405000 Tracer Red | <input type="checkbox"/> 1005001 Gloss Black Shimmer | <input type="checkbox"/> 1005002 Silver |
| <input type="checkbox"/> 2405003 Fine Texture Black | <input type="checkbox"/> 1005004 Orange Fusion | <input type="checkbox"/> 1005005 Karl's Imagination |
| <input type="checkbox"/> 2405006 Blue Shimmer | <input type="checkbox"/> 1005007 Gloss White | <input type="checkbox"/> 1005008 Dark Titanium |
| <input type="checkbox"/> 2405009 Geyser Green | <input type="checkbox"/> 1005010 Lollypop Blue | <input type="checkbox"/> 1005011 Lollypop Red |
| <input type="checkbox"/> 2405012 Bubble Gum Pink | <input type="checkbox"/> 1005013 Berry Purple | <input type="checkbox"/> 1005014 Daffodil Yellow |
| <input type="checkbox"/> 2405200 Custom Color Contact PDG for specifications. | | \$ 300 |

Armrest Type

2406000 **No Armrests**

Height Adjustable Armrests (Single Post) \$ 275

2406111 8"-11" 1006114 10"-14"

2406100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2406401 Rear 2406402 Mid 2406403 Front

Armpads & Armrest Options

Regular Armpads

2407100 Desk Length - 2" x 10" 1007101 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

2407110 Desk Length - 3" x 10" 1007111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

2407120 Full Length - 2" x 12" 1007121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125

2407130 Right (4.5" x 18") 1007131 Left (4.5" x 18")

2407310 **Clothing Guards** Available for 10"-14" Height range only. \$ 110

Seat Pan

2408000 **No Seat Pan**

2408001 **Solid Aluminum Seat Pan** \$ 120

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.

| 6" Casters | RW | Low | | | | | High | | | | | Extra-High | | | | | |
|------------|-----|-----|-----|-----|-----|-----|------|-----|-----|-----|-----|------------|-----|-----|-----|-----|-----|
| | 12" | 13" | 14" | 15" | 16" | 17" | 15" | 16" | 17" | 18" | 19" | 15" | 16" | 17" | 18" | 19" | 20" |
| 16" | | | | 15" | 16" | 17" | | | 17" | 18" | 19" | | | 17" | 18" | 19" | 20" |
| 20" | 13" | 14" | 15" | | | | 15" | 16" | 17" | | | 15" | 16" | 17" | 18" | | |
| 22" | | 14" | 15" | 16" | | | | 16" | 17" | 18" | | | 16" | 17" | 18" | 19" | |
| 24" | | | 15" | 16" | 17" | | | 17" | 18" | 19" | | | 17" | 18" | 19" | 20" | |

| 8" Casters | RW | Low | | | High | | | Extra-High | | | |
|------------|-----|-----|-----|-----|------|-----|-----|------------|-----|-----|-----|
| | 16" | 15" | 16" | 17" | 17" | 18" | 19" | 17" | 18" | 19" | 20" |
| 24" | | 15" | 16" | 17" | 17" | 18" | 19" | 17" | 18" | 19" | 20" |

Wheel size _____

Caster size _____

STF height _____

Range _____

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

2411061 6" x 1.25" Caster size 2411081 8" x 1.25" Caster size

2411063 6" x 1.50" Caster size \$ 85

2411084 8" x 2.00" Caster size \$ 135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

2411161 6" x 1.25" Caster size \$ 110

2411181 8" x 1.25" Caster size \$ 120

2411184 8" x 2.00" Caster size \$ 150

2411260 **Frog Legs 6" x 1.40" Soft Roll Casters** \$ 95

2411400 **Frog Legs Suspension Forks** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. \$ 395

Rear Wheels

Composite Wheels

- | | | |
|---|---|---|
| <input type="checkbox"/> 2412012 12" Wheel size | <input type="checkbox"/> 1012016 16" Wheel size | <input type="checkbox"/> 1012020 20" Wheel size |
| <input type="checkbox"/> 2412022 22" Wheel size | <input type="checkbox"/> 1012024 24" Wheel size | |

Spoke Wheels

- | | | |
|---|---|---|
| <input type="checkbox"/> 2412120 20" Wheel size | <input type="checkbox"/> 1012122 22" Wheel size | <input type="checkbox"/> 1012124 24" Wheel size |
|---|---|---|

\$ 75

Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- | | | |
|---------------------------------------|--------------------------------------|--|
| <input type="checkbox"/> 2412401 Rear | <input type="checkbox"/> 2412402 Mid | <input type="checkbox"/> 2412403 Front |
|---------------------------------------|--------------------------------------|--|

Axle Type

- | |
|--|
| <input type="checkbox"/> 2413000 Quick Release Rear Axles |
| <input type="checkbox"/> 2413010 Fixed Rear Axles |

Rear Tires

- | | | |
|--|---|-------|
| <input type="checkbox"/> 2414000 Urethane Tires | | |
| <input type="checkbox"/> 2414010 Pneumatic Tires | May result in a finished Seat Height 1/2" greater than stated. | \$ 60 |
| <input type="checkbox"/> 2414020 Pneumatic Tires with Flat-Free Inserts | May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. | \$ 95 |

Hand Rims

Handrims are only available with 20", 22", and 24" Rear Wheels.

- | | | |
|--|--|--------|
| <input type="checkbox"/> 2415000 No Handrims | | |
| <input type="checkbox"/> 2415010 Standard Aluminum Handrims | | |
| <input type="checkbox"/> 2415020 Plastic Coated Handrims | | \$ 135 |

Wheel Locks

User Wheel Locks

- | | | |
|--|---|--------|
| <input type="checkbox"/> 2416010 Push to Lock | <input type="checkbox"/> 2416020 Pull to Lock | |
| <input type="checkbox"/> 2416030 Attendant Wheel Locks (Push to Lock) | Not compatible with TTDS. | \$ 85 |
| <input type="checkbox"/> 2416040 User & Attendant Wheel Locks | Available with 20", 22", or 24" Rear Wheels only. | \$ 125 |
| <input type="checkbox"/> 2416100 Wheel Lock Extensions | Not compatible with Attendant Wheel Lock. | \$ 60 |
| <input type="checkbox"/> 2416200 Hill Holders | Requires Pneumatic Rear Tires with Flat-Free Inserts. | \$ 450 |
| <input type="checkbox"/> 2416300 Drum Brakes | Available with 20" Rear Wheels only. <i>COMING SOON</i> | \$ TBD |

Anti-Tippers

- | | | |
|--|-------------------------|-------|
| <input checked="" type="checkbox"/> 2417010 Rear Anti-Tippers | Required Safety Option. | \$ 60 |
|--|-------------------------|-------|

Back Post Style

- | | | |
|--|---|---|
| <input checked="" type="checkbox"/> 2418010 Push Bar Back | | |
| Back Post Height Please note: We have updated the descriptions/measurements below. The actual back post heights remain unchanged. | | |
| <input type="checkbox"/> 2418425 Short (22") <i>Previously 25"</i> | <input type="checkbox"/> 2418431 Tall (25") <i>Previously 31"</i> | Measured from the Seat Pan to the Handle. |

Back Options

- | | | |
|---|---|--|
| Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate. | | \$ 150 |
| <input type="checkbox"/> 2419010 Button on Outside | <input type="checkbox"/> 1019011 Button on Inside | |
| Seat to Back Angle | | |
| <input type="checkbox"/> 2419090 Set at 90° | <input type="checkbox"/> 1019098 Set at 98° | <input type="checkbox"/> 1019106 Set at 106° |
| <input type="checkbox"/> 2419114 Set at 114° | <input type="checkbox"/> 1019122 Set at 122° | |

See next page for more options >



| | | | | | |
|---|---------|--------------------------------|---|----|-----|
| <input type="checkbox"/> | 2419210 | 10" Backrest Extension | Requires the Angle/Height Adjustable Push Bar Handle. | \$ | 150 |
| <input type="checkbox"/> | 2419212 | Headrest Mounting Plate | | \$ | 95 |
| Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar. Not required if mounting to a solid Back. | | | | | |
| <input checked="" type="checkbox"/> | 2419220 | Fold Down Back | Not compatible with Reclining Back. | | |

Back Upholstery

| <input type="checkbox"/> | 2420000 | No Back Upholstery | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
|--|---------------|--------------------------------------|---|---------------|----------------------|-----|-----|-----|-----|--|--|-------|-----|-----|-----|-----|-----|-----|--|--|--|---------|--|--|--|-----|-----|-----|--|--|--|-------------|--|--|-----|-----|-----|-----|--|--|--|---|--|--|------------|---------------|--|--|--|--|--|--|--|-------|-----|-----|-----|-----|-----|-----|-----|-----|-----|---------|--|--|--|-----|-----|-----|-----|-----|-----|------|--|--|--|--|--|--|-----|-----|-----|-------------|--|--|-----|-----|-----|-----|-----|-----|-----|--------|-------|
| Standard Upholstery Not available with Cane Handle Back. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | 2420418 | Short (18" Height) | <input type="checkbox"/> | 2420421 | Regular (21" Height) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | 2420424 | Tall (24" Height) | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <input type="checkbox"/> | 2420450 | Tension Adjustable Upholstery | Height of Upholstery is 16". Not available with Cane Handle Back. | \$ | 295 | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Position of Upholstery Position of top of Upholstery relative to the Seat Pan. Circle below of select from the drop-down menu on the right. | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| <table border="1"> <thead> <tr> <th>Upholstery</th> <th colspan="8">22" Back Post</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td>17"</td><td>18"</td><td>19"</td><td>20"</td><td>21"</td><td>22"</td><td></td><td></td><td></td> </tr> <tr> <td>Regular</td> <td></td><td></td><td></td><td>20"</td><td>21"</td><td>22"</td><td></td><td></td><td></td> </tr> <tr> <td>Tension Adj</td> <td></td><td></td><td>19"</td><td>20"</td><td>21"</td><td>22"</td><td></td><td></td><td></td> </tr> </tbody> </table> | | | Upholstery | 22" Back Post | | | | | | | | Short | 17" | 18" | 19" | 20" | 21" | 22" | | | | Regular | | | | 20" | 21" | 22" | | | | Tension Adj | | | 19" | 20" | 21" | 22" | | | | <table border="1"> <thead> <tr> <th>Upholstery</th> <th colspan="8">25" Back Post</th> </tr> </thead> <tbody> <tr> <td>Short</td> <td>17"</td><td>18"</td><td>19"</td><td>20"</td><td>21"</td><td>22"</td><td>23"</td><td>24"</td><td>25"</td> </tr> <tr> <td>Regular</td> <td></td><td></td><td></td><td>20"</td><td>21"</td><td>22"</td><td>23"</td><td>24"</td><td>25"</td> </tr> <tr> <td>Tall</td> <td></td><td></td><td></td><td></td><td></td><td></td><td>23"</td><td>24"</td><td>25"</td> </tr> <tr> <td>Tension Adj</td> <td></td><td></td><td>19"</td><td>20"</td><td>21"</td><td>22"</td><td>23"</td><td>24"</td><td>25"</td> </tr> </tbody> </table> | | | Upholstery | 25" Back Post | | | | | | | | Short | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" | Regular | | | | 20" | 21" | 22" | 23" | 24" | 25" | Tall | | | | | | | 23" | 24" | 25" | Tension Adj | | | 19" | 20" | 21" | 22" | 23" | 24" | 25" | Height | _____ |
| Upholstery | 22" Back Post | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Short | 17" | 18" | 19" | 20" | 21" | 22" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Regular | | | | 20" | 21" | 22" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Tension Adj | | | 19" | 20" | 21" | 22" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Upholstery | 25" Back Post | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Short | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Regular | | | | 20" | 21" | 22" | 23" | 24" | 25" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Tall | | | | | | | 23" | 24" | 25" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Tension Adj | | | 19" | 20" | 21" | 22" | 23" | 24" | 25" | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Tilt/Recline Options

| | | | | | | | |
|--|--------------------------|------------------------------------|--|------------------------------|--|---------|------------------------|
| Attendant Tilt only Select trigger mounting location below. | | | | | | | |
| <input type="checkbox"/> | 2423011 | Right Back Post | <input type="checkbox"/> | 2423012 | Left Back Post | | |
| User Tilt only Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below. | | | | | | | |
| <input type="checkbox"/> | 2423021 | Right Armrest | <input type="checkbox"/> | 2423022 | Left Armrest | | |
| Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below. | | | | | | | |
| <input type="checkbox"/> | 2423031 | Right Back Post | <input type="checkbox"/> | 2423032 | Left Back Post | | |
| <input type="checkbox"/> | 2423041 | Right Armrest | <input type="checkbox"/> | 2423042 | Left Armrest | | |
| <input type="checkbox"/> | 2423110 | Trigger Extension | | | \$ 50 | | |
| <input type="checkbox"/> | 2423310 | Reclining Backrest (Manual) | | | \$ 895 | | |
| 30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19". | | | | | | | |
| <input type="checkbox"/> | 2423400 | Power Tilt | 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt | | \$ 1,795 | | |
| Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both. | | | | | | | |
| | Toggle Switch | <input type="checkbox"/> | 2423421 | Installed on Right side | <input type="checkbox"/> | 2423422 | Installed on Left side |
| | Button Switch | <input type="checkbox"/> | 2423431 | Installed on Right side | <input type="checkbox"/> | 2423432 | Installed on Left side |
| | Stealth Egg Switch | <input type="checkbox"/> | 2423440 | Not installed | | | |
| | Optional for Egg Switch: | <input type="checkbox"/> | 2423445 | Swing Away Mounting Hardware | Not compatible with all Headrest makes and models. | \$ | 625 |
| Specify Headrest to be used: _____ | | | | | | | |

Front Rigging

| | | | | | |
|---|---------|---|--|---------|---------------------------|
| <input type="checkbox"/> | 2425000 | No Front Rigging | <input type="checkbox"/> | 2425001 | Sector Blocks only |
| Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out | | | | | |
| <input type="checkbox"/> | 2425410 | Fixed 70° | | | \$ 95 |
| <input type="checkbox"/> | 2425420 | Universal Angle Adjustable (0°-110°) | Set at angle: _____ (use increments of 10°) | | \$ 295 |
| Adjust angle by loosening one bolt. Tool-free leg-length adjustment. | | | | | |
| Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out. | | | | | |
| <input type="checkbox"/> | 2425130 | Swing Away 60° | <input type="checkbox"/> | 2425140 | Swing Away 70° |
| <input type="checkbox"/> | 2425200 | Swing Away Elevating Legrests | Includes Calf Pads. | | \$ 375 |
| <input type="checkbox"/> | 2425300 | 1" Offset Hangers for Swing Away | | | \$ 75 |
| <input type="checkbox"/> | 2425110 | Standard Plug In 70° Front Rigging | Legacy Design. Not compatible with Easy-Out. | | |


Footplates

| | | | | |
|--------------------------|---------|---|--|--------|
| <input type="checkbox"/> | 2426010 | Composite Footplates | | |
| <input type="checkbox"/> | 2426020 | Aluminum Angle Adjustable Footplates | | \$ 155 |
| <input type="checkbox"/> | 2426030 | One Piece Fixed Footplate | Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. | \$ 350 |
| <input type="checkbox"/> | 2426040 | One Piece Flip-up Footplate | Requires selection of Composite or Aluminum Angle Adjustable. | \$ 395 |
| <input type="checkbox"/> | 2426050 | Interlocking Angle Adjustable Footplates | Available for Seat Widths up to 21" only. | \$ 195 |

Footrest/Legrest Accessories

| | | | | |
|---|---------|--|---|--------|
| Calf Pads Not compatible with Gel Padding. | | | | |
| <input type="checkbox"/> | 2427012 | For Easy-Out™ Front Rigging | Includes reinforced premium Calf Pad mount. | \$ 200 |
| <input type="checkbox"/> | 2427010 | For Swing Away Front Rigging | | \$ 135 |
| <input type="checkbox"/> | 2427020 | Gel Padding | Not compatible with Calf Pad or Easy-Out Front Rigging. | \$ 125 |
| Gel Knee Pads | | | | |
| <input type="checkbox"/> | 2427032 | For Easy-Out™ Angle Adjustable Front Rigging | | \$ 200 |
| <input type="checkbox"/> | 2427030 | For Swing Away ELR Front Rigging | | \$ 100 |
| <input type="checkbox"/> | 2427040 | Length Adjustable Heel Loops | | \$ 50 |
| <input type="checkbox"/> | 2427050 | Bodypoint Padded Calf Strap | | \$ 195 |


Front Rigging Extension Tubes

| | | | | |
|--|--|-------|--|-------|
| Regular Extension Tubes | | | | |
| <input type="checkbox"/> | 2428010 | Short | 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging | |
| <input type="checkbox"/> | 2428020 | Long | 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging | \$ 75 |
|  | Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used. | | | |
| | Seat Pan to Foot Plates: _____ | | Measured along the length of the Front Rigging. | |

Miscellaneous Options

| | | | | |
|--|---------|---|---|----------|
| One-Arm Drive Requires 24" Composite Wheels with Urethane Tires. 18" STF minimum. Limits tilt. Not compatible with User & Attendant Wheel Locks, Quick Release Rear Axles. | | | | \$ 3,500 |
| <input type="checkbox"/> | 2430010 | Installed on Right side | <input type="checkbox"/> 1030011 Installed on Left side | |
| <input type="checkbox"/> | 2430050 | Protective Lower Padding | | \$ 675 |
| <input type="checkbox"/> | 2430080 | Gel Hip Guides | Reduces available Seat Width by 2". | \$ 375 |
| <input type="checkbox"/> | 2430110 | IV Pole | Not compatible with Attendant Wheel Locks. | \$ 295 |
| <input type="checkbox"/> | 2430120 | O2 Holder | Not compatible with Attendant Wheel Locks. | \$ 395 |
| <input type="checkbox"/> | 2430170 | Adhesive Back Hook Strip, 2" x 30" | For Cushion mounting. | \$ 10 |

On-Chair Seating & Positioning

| | | | | |
|---|---|--|--|--|
|  | Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating | | | |
|---|---|--|--|--|

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.