



Customer and Order Information

| | | |
|--|------------------|-------------------------------------|
| Purchase Order * | Tag | v.0801 |
| Dealer/Provider Name * | | |
| Purchaser Name and Contact Information * | | |
| Shipping Address * | | |
| User Information (if available) | Weight _____ lbs | Height _____ ft/in <i>*required</i> |

Frame Type

| | | |
|---|--|----------|
| <input checked="" type="checkbox"/> 1101000 | Eclipse Bariatric Manual Wheelchair Extended delivery lead times may apply. | \$ 3,295 |
|---|--|----------|

Frame Options

| | | |
|--|---|--------|
| Heavy Duty Base and Seat Frame | | |
| <input type="checkbox"/> 1102000 | 600 lbs Weight Capacity | |
| <input type="checkbox"/> 1102100 | 1000 lbs Weight Capacity Not compatible with Transit Tie Down System. | \$ 300 |
| <input type="checkbox"/> 1102300 | Transit Tie Down System (TTDS) | \$ 395 |
| TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 600 lbs. Not compatible with O2 holder or IV pole. | | |

Seat Width Overall chair width is +8"

| | | |
|---|--|--|
| Seat Widths Seat Width is adjustable with the optional Adjustment Kit (sold separately). | | |
| <input type="checkbox"/> 1103520 | 20" Seat Width | <input type="checkbox"/> 1103522 22" Seat Width \$ 650 |
| <input type="checkbox"/> 1103524 | 24" Seat Width | <input type="checkbox"/> 1103526 26" Seat Width \$ 650 |
| <input type="checkbox"/> 1103528 | 28" Seat Width | <input type="checkbox"/> 1103530 30" Seat Width \$ 850 |
| <input type="checkbox"/> 1103532 | 32" Seat Width | <input type="checkbox"/> 1103536 36" Seat Width \$ 1,050 |
| <input type="checkbox"/> 1103540 | 40" Seat Width | <input type="checkbox"/> 1103544 44" Seat Width \$ 1,250 |
| <input type="checkbox"/> 1103599 | Unlisted Seat Width (between sizes) Seat Width: _____ (use increments of 1") | <input type="checkbox"/> 1103548 48" Seat Width \$ 1,250 |
| <i>Add this fee to the price of the next-bigger size</i> | | |
| <input type="checkbox"/> 1103599 Unlisted Seat Width (between sizes) Seat Width: _____ (use increments of 1") <i>Add this fee to the price of the next-bigger size</i> \$ 100 | | |
| i Seat Widths greater than 48" are available for Custom Chairs | | |

Seat Depth

| | | |
|----------------------------------|----------------|--|
| Seat Depths (Adjustable) | | |
| <input type="checkbox"/> 1104016 | 16" Seat Depth | <input type="checkbox"/> 1104017 17" Seat Depth |
| <input type="checkbox"/> 1104018 | 18" Seat Depth | <input type="checkbox"/> 1104019 19" Seat Depth |
| <input type="checkbox"/> 1104020 | 20" Seat Depth | <input type="checkbox"/> 1104021 21" Seat Depth |
| <input type="checkbox"/> 1104022 | 22" Seat Depth | |
| Custom Seat Depths \$ 350 | | |
| <input type="checkbox"/> 1104023 | 23" Seat Depth | <input type="checkbox"/> 1104024 24" Seat Depth Not available with TTDS. |

Frame Color

| | | |
|----------------------------------|---|--|
| PDG Premium Colors | | |
| <input type="checkbox"/> 1105000 | Tracer Red | <input type="checkbox"/> 1105001 Gloss Black Shimmer |
| <input type="checkbox"/> 1105002 | Silver | <input type="checkbox"/> 1105002 Silver |
| <input type="checkbox"/> 1105003 | Fine Texture Black | <input type="checkbox"/> 1105004 Orange Fusion |
| <input type="checkbox"/> 1105005 | Karl's Imagination | <input type="checkbox"/> 1105005 Karl's Imagination |
| <input type="checkbox"/> 1105006 | Blue Shimmer | <input type="checkbox"/> 1105007 Gloss White |
| <input type="checkbox"/> 1105008 | Dark Titanium | <input type="checkbox"/> 1105007 Gloss White |
| <input type="checkbox"/> 1105009 | Geyser Green | <input type="checkbox"/> 1105010 Lollipop Blue |
| <input type="checkbox"/> 1105011 | Lollipop Red | <input type="checkbox"/> 1105010 Lollipop Blue |
| <input type="checkbox"/> 1105012 | Bubble Gum Pink | <input type="checkbox"/> 1105013 Berry Purple |
| <input type="checkbox"/> 1105014 | Daffodil Yellow | <input type="checkbox"/> 1105013 Berry Purple |
| <input type="checkbox"/> 1105014 | Daffodil Yellow | |
| <input type="checkbox"/> 1105200 | Custom Color Contact PDG for specifications. | \$ 300 |

| Armrest Type | |
|---|--|
| <input type="checkbox"/> 1106000 | No Armrests |
| Dual Post Height Adjustable Armrests \$ 295 | |
| <input type="checkbox"/> 1106214 | 10"-14" <input type="checkbox"/> 1106216 12"-16" |
| <input type="checkbox"/> 1106200 | Dual Post Fixed Height Armrests (12") |

| Armpads & Armrest Options | |
|--|--|
| Regular Armpads | |
| <input type="checkbox"/> 1107100 | Desk Length - 2" x 10" <input type="checkbox"/> 1107101 Full Length - 2" x 14" |
| Wide Armpads May interfere with wheels in some tilt configurations. \$ 200 | |
| <input type="checkbox"/> 1107110 | Desk Length - 3" x 10" <input type="checkbox"/> 1107111 Full Length - 3" x 12" |
| Gel Armpads May interfere with wheels in some tilt configurations. \$ 135 | |
| <input type="checkbox"/> 1107120 | Full Length - 2" x 12" <input type="checkbox"/> 1107121 Full Length - 3.5" x 12" |
| <input type="checkbox"/> 1107250 | Forward Guards \$ 200 |
| <input type="checkbox"/> 1107400 | 3" Armrest Extension - Forward \$ 75 |

| Seat Pan | |
|---|--|
| <input checked="" type="checkbox"/> 1108001 | Solid Aluminum Seat Pan \$ 155 |

| Seat Height | |
|---|---|
| Front Seat Height Angle Adjustable Foot Plates required for 14" and 15". Front and Rear Height can be max. 1" apart. | |
| <input type="checkbox"/> 1110014 | 14" STF Front <input type="checkbox"/> 1110015 15" STF Front <input type="checkbox"/> 1110016 16" STF Front |
| <input type="checkbox"/> 1110017 | 17" STF Front <input type="checkbox"/> 1110018 18" STF Front |
| <input type="checkbox"/> 1110019 | 19" STF Front <input type="checkbox"/> 1110020 20" STF Front \$ 200 |
| Rear Seat Height Angle Adjustable Foot Plates required for 14" and 15". Front and Rear Height can be max. 1" apart. | |
| <input type="checkbox"/> 1110114 | 14" STF Rear <input type="checkbox"/> 1110115 15" STF Rear <input type="checkbox"/> 1110116 16" STF Rear |
| <input type="checkbox"/> 1110117 | 17" STF Rear <input type="checkbox"/> 1110118 18" STF Rear |
| <input type="checkbox"/> 1110119 | 19" STF Rear <input type="checkbox"/> 1110120 20" STF Rear \$ 200 |

| Casters | |
|---|--|
| <input checked="" type="checkbox"/> 1111050 | 5" Casters With Heavy Duty Steel Fork |

| Rear Wheels | |
|---|------------------------------------|
| <input checked="" type="checkbox"/> 1112224 | 24" Heavy Duty Spoke Wheels |

| Axle Type | |
|---|--|
| <input checked="" type="checkbox"/> 1113000 | Quick Release Rear Axles |
| Rear Axle Position Positions 3 and 4 not compatible with TTDS. | |
| <input type="checkbox"/> 1113101 | Set at 1 (Rear) <input type="checkbox"/> 1113102 Set at 2 |
| <input type="checkbox"/> 1113103 | Set at 3 <input type="checkbox"/> 1113104 Set at 4 (Front) |

| Rear Tires | |
|---|-----------------------|
| <input checked="" type="checkbox"/> 1114000 | Urethane Tires |

| Hand Rims | |
|----------------------------------|--|
| <input type="checkbox"/> 1115000 | No Handrims |
| <input type="checkbox"/> 1115010 | Standard Aluminum Handrims |
| <input type="checkbox"/> 1115020 | Plastic Coated Handrims \$ 250 |

Wheel Locks

User Wheel Locks

1116010 Push to Lock 1116020 Pull to Lock

1116100 **Wheel Lock Extensions** \$ 60

Anti-Tippers

1117010 **Rear Anti-Tippers** Required Safety Option. \$ 60

Back Post Style

1118030 **Push Bar Back, 8° Bend**

1118040 **Cane Handle Back, 8° Bend** Not compatible with Headrest Mounting Plate.

Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. \$ 395

1118102 1" Offset per side 1118103 1.5" Offset per side 1118104 2" Offset per side

Back Post Height Please note: We have updated the descriptions/measurements below. The actual back post heights remain unchanged.

1118425 Regular (22") *Previously 25"* 1118431 Tall (25") *Previously 30"* Measured from the Seat Pan to the Handle.

Back Options

Seat to Back Angle

1119085 Set at 85° 1119090 Set at 90°

1119095 Set at 95° 1119100 Set at 100°

1119210 **10" Backrest Extension** \$ 150

1119212 **Headrest Mounting Plate** Not required if mounting to a solid Back. \$ 95

Back Upholstery

1120000 **No Back Upholstery**

Tension Adjustable \$ 375

1120421 Regular (18" Height) 1120424 Tall (22" Height)

Position of Upholstery Position of top of Upholstery relative to the Seat Pan.

| Upholstery | 22" Back Post | | | | | | | |
|------------|---------------|-----|-----|-----|-----|-----|--|--|
| Regular | 17" | 18" | 19" | 20" | 21" | 22" | | |
| Tall | | | | 20" | 21" | 22" | | |

| Upholstery | 25" Back Post | | | | | | | | |
|------------|---------------|-----|-----|-----|-----|-----|-----|-----|-----|
| Regular | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" |
| Tall | | | | | 21" | 22" | 23" | 24" | 25" |

Position _____

Recline Options

1123310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Backrest Extension. Reclining Backrest is 1" narrower than Seat Width. Back Frame is not removable.

Front Rigging

1125000 **No Front Rigging** 1125001 **Sector Blocks only**

1125002 **Center Mount Front Rigging Bracket** Compatible with Swing Away 60° only. \$ 300

PDG Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out

1125410 Fixed 60° \$ 95

1125420 Universal Angle Adjustable (0°-60°) Set at angle: _____ (use increments of 10°) \$ 295

Adjust angle by loosening one bolt. Tool-free leg-length adjustment.

Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.

1125140 Swing Away 60°

1125200 Swing Away Elevating Legrests Includes Calf Pads. \$ 375

Footplates

| | | | |
|-------------------------------------|---------|---|--------|
| <input checked="" type="checkbox"/> | 1126020 | Aluminum Angle Adjustable Footplates | \$ 175 |
|-------------------------------------|---------|---|--------|

Footrest/Legrest Accessories

| | | | |
|--------------------------|---------|---|--------|
| <input type="checkbox"/> | 1127015 | Oversized 9" x 10" Calf Pad Not compatible with Seat Width of less than 25". | \$ 695 |
| <input type="checkbox"/> | 1127030 | Swing Away ELR Gel Knee Pad | \$ 100 |
| <input type="checkbox"/> | 1127040 | Length Adjustable Heel Loops | \$ 50 |
| <input type="checkbox"/> | 1127051 | Bodypoint Calf Strap | \$ 325 |

Front Rigging Extension Tubes

Regular Extension Tubes

| | | | |
|--------------------------|---------|---|-------|
| <input type="checkbox"/> | 1128010 | Short <u>14"-18"</u> for Easy-Out Front Rigging <u>13"-17"</u> for Swing Away Front Rigging | |
| <input type="checkbox"/> | 1128020 | Long <u>18"-24"</u> for Easy-Out Front Rigging <u>15"-19"</u> for Swing Away Front Rigging | \$ 75 |



Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging.

Miscellaneous Options

| | | | |
|--------------------------|---------|---|--------|
| <input type="checkbox"/> | 1130110 | IV Pole Not compatible with TTDS. | \$ 295 |
| <input type="checkbox"/> | 1130120 | O2 Holder Not compatible with TTDS. | \$ 395 |
| <input type="checkbox"/> | 1130170 | Adhesive Back Hook Strip, 2" x 30" For Cushion mounting. | \$ 10 |

On-Chair Seating & Positioning

| | | | |
|---|---|--|--|
|  | Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating | | |
|---|---|--|--|

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.