

Customer and Order Information		
Purchase Order *	Tag	v.0910
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in * required

Order Notes/Instructions

Base Chair	
Fuze M50 50° Manual Tilt Wheelchair with Mid-Wheel Drive	
<input checked="" type="checkbox"/> 2801000	Fuze 50° Manual Tilt Wheelchair \$ 4,300
<input checked="" type="checkbox"/> 2801008	Mid-Wheel Drive Frame Upgrade Adds articulating rear casters and places drive wheels in the middle of the chair at the center of gravity. \$ 500

Frame Type & Weight Capacity	
Fixed Width Frame Select a Fixed Seat Width below.	
<input type="checkbox"/> 2802010	300 lbs Weight Capacity
<input type="checkbox"/> 2802011	400 lbs Weight Capacity Not compatible with Transit Tie Down System. Requires Flat-free Tires or Heavy Duty Spoke Drive Wheels.
<input type="checkbox"/> 2802025	Adjustable Width Frame 250 lbs Weight Capacity. Select an Adjustable Seat Width below. \$ 100

Frame Options	
<input type="checkbox"/> 2802300	Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.

Seat Width		Overall chair width is +8.5".
Fixed Seat Widths Order with Fixed Frame.		
<input type="checkbox"/> 2803015	15" Seat Width	<input type="checkbox"/> 2803016
<input type="checkbox"/> 2803018	18" Seat Width	<input type="checkbox"/> 2803017
<input type="checkbox"/> 2803020	20" Seat Width	<input type="checkbox"/> 2803021
<input type="checkbox"/> 2803023	23" Seat Width	<input type="checkbox"/> 2803022
<input type="checkbox"/> 2803024	24" Seat Width	\$ 375
Adjustable Seat Widths Order with Adjustable Frame.		
<input type="checkbox"/> 2803215	15"-19" Range - Set at 15"	<input type="checkbox"/> 2803216
<input type="checkbox"/> 2803218	15"-19" Range - Set at 18"	<input type="checkbox"/> 2803217
		\$ 375



Seat Depth	
Standard Seat Depths (Adjustable)	
<input type="checkbox"/> 2804015	15" Seat Depth
<input type="checkbox"/> 2804018	18" Seat Depth
<input type="checkbox"/> 2804020	20" Seat Depth
<input type="checkbox"/> 2804016	16" Seat Depth
<input type="checkbox"/> 2804019	19" Seat Depth
<input type="checkbox"/> 2804017	17" Seat Depth
Custom Seat Depths Requires Stroller Handle Back - Straight. Not compatible with Reclining or Fold-down Back.	
<input type="checkbox"/> 2804021	21" Seat Depth
<input type="checkbox"/> 2804022	22" Seat Depth
	\$ 400

Frame Colour

PDG Premium Colours View colour swatches on PDG's website. pdgmobility.com/colours

- | | | |
|---|---|---|
| <input type="checkbox"/> 2805001 Black Shimmer | <input type="checkbox"/> 2805015 Sky Blue | <input type="checkbox"/> 2805016 Poppy Red |
| <input type="checkbox"/> 2805003 Fine Texture Black | <input type="checkbox"/> 2805010 Lollypop Blue | <input type="checkbox"/> 2805011 Lollypop Red |
| <input type="checkbox"/> 2805017 Arctic Grey | <input type="checkbox"/> 2805018 Tropical Teal | <input type="checkbox"/> 2805019 Tangerine Orange |
| <input type="checkbox"/> 2805020 Olive Green | <input type="checkbox"/> 2805021 Lush Lilac | <input type="checkbox"/> 2805012 Bubble Gum Pink |
| <input type="checkbox"/> 2805022 Leaf Green | <input type="checkbox"/> 2805023 Boysenberry Purple | <input type="checkbox"/> 2805014 Daffodil Yellow |
| <input type="checkbox"/> 2805200 Custom Colour Contact PDG for specifications. \$ 300 | | |

Armrests

- 2806000 **No Armrests**
- Height Adjustable Armrests (Single Post)** \$ 295
- 2806114 10"-14" Range 2806116 12"-16" Range
- 2806100 **Fixed Height Armrests (12")** Full Length Armpad only.
-  **Armrest Height** If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
Seat Pan to top of Armrest: _____ Does not include height of Armpad
-  **Armrest Position** Indicate preference for Armrest position. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.
- 2806401 Rear Default Position 2806403 Front
- Locking Flip-up Armrests** Not recommended for users with excessive strength or tone. \$ 295
Not compatible with Fold-down Back. Only available with Regular Full Length and Desk Length Armpads.
- 2806310 Set at 10" 2806311 Set at 11" 2806312 Set at 12"
 2806313 Set at 13" 2806314 Set at 14"

Armpads

- Regular Armpads**
- 2807100 Desk Length - 2" x 10" 2807101 Full Length - 2" x 14" \$ 115
- Waterfall Armpads**
- 2807105 Desk Length - 2" x 10" 2807106 Full Length - 2" x 14" \$ 200
- Wide Armpads** May interfere with wheels in some tilt configurations.
- 2807110 Desk Length - 3" x 10" 2807111 Full Length - 3" x 12" \$ 135
- Gel Armpads** May interfere with wheels in some tilt configurations.
- 2807120 Full Length - 2" x 12" 2807121 Full Length - 3.5" x 12"


Armrest Options

- 2807210 **Armrest Wideners** Move Armrests outward by 0.5" on each side for 1" total additional width. Seat size does not change. \$ 95

Seat Pan

- 2808000 **No Seat Pan**
- 2808001 **Solid Aluminum Seat Pan** \$ 135

Seat Height

-  **Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Axle Plate (Low/Standard). RC refers to the size of the Rear Casters.
Chairs with Swing-Away Front Rigging are set up with up to 5° of pre-tilt to avoid interference. This increases the actual STF Height at the front of the seat by up to 1" compared to the selected value. Easy-Out is not affected.

5" Front Casters									
RW	13"	14"	15"	16"	17"	18"	19"	20"	RC
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				5"
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Std. STF		6"
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					6"

6" Front Casters & 6" Forks									
RW	13"	14"	15"	16"	17"	18"	19"	20"	RC
20"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				5"
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Std. STF		6"
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				6"

6" Front Casters & 8" Forks									
RW	13"	14"	15"	16"	17"	18"	19"	20"	RC
20"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Std. STF		5"
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				6"
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				6"

6" Front Casters & Frog Legs Suspension Forks									
RW	13"	14"	15"	16"	17"	18"	19"	20"	RC
20"				<input type="checkbox"/>	<input type="checkbox"/>				5"
22"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		Std. STF		6"
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				6"

Front Casters

Urethane Casters

 Select Front Caster size here. Rear Caster size as per STF tables.

<input type="checkbox"/> 2811051	5" x 1.25" Caster size	<input type="checkbox"/> 2811061	6" x 1.25" Caster size	
<input type="checkbox"/> 2811063	6" x 1.50" Caster size			\$ 85
<input type="checkbox"/> 2811260	Soft-Roll Casters	6" x 1.50" Caster size. Urethane tire & Silver Anodized Aluminum hub. Available as Front Casters only.		\$ 275
<input type="checkbox"/> 2811400	Frog Legs Suspension Forks	Requires 6" x 1.50" Soft-Roll Casters. Suspension Forks used for Front Casters only.		\$ 595

Drive Wheels

Composite Wheels

<input type="checkbox"/> 2812020	20" Wheel size	<input type="checkbox"/> 2812022	22" Wheel size	<input type="checkbox"/> 2812024	24" Wheel size	
<input type="checkbox"/> 2812120	20" Wheel size	<input type="checkbox"/> 2812122	22" Wheel size	<input type="checkbox"/> 2812124	24" Wheel size	\$ 75
<input type="checkbox"/> 2812224	24" Heavy Duty Spoke Wheels					\$ 350

Axle Type

<input type="checkbox"/> 2813000	Quick Release Drive Wheel Axles	
<input type="checkbox"/> 2813010	Fixed Drive Wheel Axles	

Drive Wheel Tires

<input type="checkbox"/> 2814000	Urethane Tires		
<input type="checkbox"/> 2814010	Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated.	\$ 60
<input type="checkbox"/> 2814020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$ 95

Handrims

<input type="checkbox"/> 2815000	No Handrims		
<input type="checkbox"/> 2815010	Standard Aluminum Handrims		
<input type="checkbox"/> 2815020	Plastic Coated Handrims		\$ 150
<input type="checkbox"/> 2815030	Projection Handrims	Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$ 495

Wheel Locks

User Wheel Locks

<input type="checkbox"/> 2816010	Push to Lock	<input type="checkbox"/> 2816020	Pull to Lock	
<input type="checkbox"/> 2816310	Left Activation	<input type="checkbox"/> 2816311	Right Activation	
	User-activated Unilateral Wheel Lock	Pull to Lock. Available up to 22" Seat Width only.		\$ 795

Wheel Lock Options

<input type="checkbox"/> 2816100	Wheel Lock Extensions		\$ 60
----------------------------------	------------------------------	--	-------

Back Posts

Back Post Style

<input type="checkbox"/> 2818100	Stroller Handle Back	<input type="checkbox"/> 2818105	Cane Handle Back	Not compatible with Headrest Mounting Plate.	
<input type="checkbox"/> 2818200	Straight	<input type="checkbox"/> 2818205	8° Bend		
<input type="checkbox"/> 2818420	Short (17")	<input type="checkbox"/> 2818425	Tall (23")		
<input type="checkbox"/> 1519220	Fold-down	<input type="checkbox"/> 1519221	Non-folding		
<input type="checkbox"/> 1519090	Set at 90°	<input type="checkbox"/> 1519096	Set at 98°	<input type="checkbox"/> 1519102	Set at 106°
<input type="checkbox"/> 1519108	Set at 114°	<input type="checkbox"/> 1519114	Set at 122°		

Back Options

<input type="checkbox"/> 2819010 Angle/Height Adjustable Stroller Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 2819010 Button on Outside	<input type="checkbox"/> 2819011 Button on Inside	
<input type="checkbox"/> 2818102 Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12". Not compatible with Flip-up Armrests, Reclining Backrest.		\$ 395
<input type="checkbox"/> 2818102 1" Offset per side	<input type="checkbox"/> 15118103 1.5" Offset per side	<input type="checkbox"/> 2818104 2" Offset per side
<input type="checkbox"/> 2819212 Headrest Mounting Plate Not compatible with Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back.		\$ 95

Back Upholstery

<input type="checkbox"/> 2820000 No Back Upholstery																			
<input type="checkbox"/> 2820418 Standard Upholstery Not available with Cane Handle Back.		\$ 195																	
<input type="checkbox"/> 2820418 Short (18" Height)	<input type="checkbox"/> 2820421 Regular (21" Height)																		
<input type="checkbox"/> 2820450 Tension Adjustable Upholstery Height of Upholstery is 16". Not available with Cane Handle Back.		\$ 295																	
Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.																			
17" Back Post					23" Back Post														
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"	Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>									Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
										Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
										Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation

<input type="checkbox"/> 2823011 Attendant Tilt only Select trigger mounting location below.		
<input type="checkbox"/> 2823011 Right Back Post	<input type="checkbox"/> 2823012 Left Back Post	
<input type="checkbox"/> 2823220 Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.		
<input type="checkbox"/> 2823220 Set at 20°	<input type="checkbox"/> 2823230 Set at 30°	<input type="checkbox"/> 2823240 Set at 40°

Tip Configuration (Front Caster Lift)

<input type="checkbox"/> 2823451 Attendant Tip only Select trigger mounting location below.		
<input type="checkbox"/> 2823451 Right Back Post	<input type="checkbox"/> 2823452 Left Back Post	<input type="checkbox"/> 2823453 Omit Trigger No Front Caster Lift
<input type="checkbox"/> 2823455 Tip Clearance Limit (Front Caster Lift) Actual tip may vary by up to 0.5". No selection required if Trigger is omitted.		
<input type="checkbox"/> 2823455 4" Caster Lift Default	<input type="checkbox"/> 2823456 2" Caster Lift	

Recline

<input type="checkbox"/> 2823310 Reclining Backrest (Manual)	\$ 895
30° maximum recline. Not compatible with Cane Handle Back, Fold-down Back. Not available with 21" & 22" Seat Depth. Finished Back Post Height is 3" taller than selected Back Post Height.	

Front Rigging

<input type="checkbox"/> 2825000 No Front Rigging	<input type="checkbox"/> 2825001 Sector Blocks only	
<input type="checkbox"/> 2825410 Easy-Out Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out		
<input type="checkbox"/> 2825410 Fixed 70°		\$ 95
<input type="checkbox"/> 2825420 Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings. Optional for Universal AA: <input type="checkbox"/> 2825421 AA High-mount Option For shorter Lower Leg Lengths of 8"-13".		
<input type="checkbox"/> 2825130 Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out, Push to Lock Wheel Locks. Requires pre-tilt of up to 5°.		
<input type="checkbox"/> 2825130 Swing Away 60°	<input type="checkbox"/> 2825140 Swing Away 70°	
<input type="checkbox"/> 2825200 Swing Away Elevating Legrests	Includes Calf Pads.	\$ 375


Footplates

<input type="checkbox"/> 2826010 Composite Footplates	
<input type="checkbox"/> 2826020 Aluminum Angle Adjustable Footplates	\$ 165
<input type="checkbox"/> 2826040 One Piece Flip-up Footplate Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/> 2826050 Interlocking Angle Adjustable Footplates Available for Seat Widths up to 21" only.	\$ 195

Footrest/Legrest Accessories

Calf Pad Not compatible with Gel Padding.		
<input type="checkbox"/>	2827012 For Easy-Out Front Rigging Includes reinforced premium Calf Pad mount. Requires pre-tilt of up to 5°.	\$ 200
<input type="checkbox"/>	2827010 For Swing Away Front Rigging	\$ 135
<input type="checkbox"/>	2827020 Gel Padding Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pad		
<input type="checkbox"/>	2827032 For Easy-Out Angle Adjustable Front Rigging	\$ 200
<input type="checkbox"/>	2827030 For Swing Away ELR Front Rigging	\$ 100
<input type="checkbox"/>	2827040 Length Adjustable Heel Loops	\$ 50
<input type="checkbox"/>	2827050 Bodypoint Padded Calf Strap Sizes greater than 22" are not padded.	\$ 195

Front Rigging Extension Tubes

Regular Extension Tubes No selection required for Easy-Out AA High-mount Option.		
<input type="checkbox"/>	2828010 Short <u>14"-18"</u> for Easy-Out Front Rigging <u>13"-17"</u> for Swing Away Front Rigging	
<input type="checkbox"/>	2828020 Long <u>18"-24"</u> for Easy-Out Front Rigging <u>15"-19"</u> for Swing Away Front Rigging	\$ 75
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used. Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.	

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Drive Wheels. Made from durable and lightweight aluminum.		\$ 195
<input type="checkbox"/>	2830020 Black <input type="checkbox"/> 2830021 Frame Color	
Stealth Flip Down Abductor Includes wedge-style pad. Requires Solid Seat Pan.		\$ 650
<input type="checkbox"/>	2830060 Small Pad 3"x4"x3"x2" size <input type="checkbox"/> 2830062 Large Pad 5"x5"x4"x3" size	
<input type="checkbox"/>	2830080 Gel Hip Guides Reduces available Seat Width by 2".	\$ 495
<input type="checkbox"/>	2830090 Hip Belt Mount Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.	\$ 55
<input type="checkbox"/>	2830152 Ventilator Bracket Compatible with LTV1 (laptop) ventilator.	\$ 295
<input type="checkbox"/>	2830160 Clear Lap Tray (Slide On) Available for Seat Widths of 17"-22" only. Will interfere with Drive Wheel when in tilt.	\$ 490
<input type="checkbox"/>	2830200 Q'Straint Lap Belt Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.	\$ 245
<input type="checkbox"/>	2830222 Medical Necessities Backpack Black with Permobil Logo.	\$ 95
<input type="checkbox"/>	2830170 Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10

On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Please send your completed order form by email to sales.pdg@permobil.com.