

### Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

### Base Chair

<input checked="" type="checkbox"/> 1501000	<b>Fuze T50 50° Manual Tilt Wheelchair</b>	\$ 4,150
	Also available: the narrow <b>Fuze T50n</b> with an overall chair width of seat plus only 4". See separate order form. <a href="http://pdgmobility.com/narrow">pdgmobility.com/narrow</a>	
	<p><b>FUZE T50n</b></p> <p>Overall chair width of seat plus only 4". Full 50° of tilt. The new Fuze T50n. Upgrade to a smaller size.</p> <p><a href="http://pdgmobility.com/narrow">pdgmobility.com/narrow</a></p> 	

### Frame Type & Weight Capacity

<b>Fixed Width Frame</b> Select a Fixed Seat Width below.		
<input type="checkbox"/> 1502010	300 lbs Weight Capacity	
<input type="checkbox"/> 1502011	400 lbs Weight Capacity Not compatible with Transit Tie Down System. Requires Flat-free Rear Tires or Heavy Duty Spoke Rear Wheels.	
<input type="checkbox"/> 1502025	<b>Adjustable Width Frame</b> 250 lbs Weight Capacity. Select an Adjustable Seat Width below.	\$ 100

### Frame Options

<input type="checkbox"/> 1502300	<b>Transit Tie Down System (TTDS)</b>
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.	

### Seat Width

Overall chair width is +7.5" for 12/16" RW and +8.5" for 20/22/24" RW

<b>Fixed Seat Widths</b> Order with Fixed Frame.					
<input type="checkbox"/> 1503015	15" Seat Width	<input type="checkbox"/> 1503016	16" Seat Width	<input type="checkbox"/> 1503017	17" Seat Width
<input type="checkbox"/> 1503018	18" Seat Width	<input type="checkbox"/> 1503019	19" Seat Width	<input type="checkbox"/> 1503020	20" Seat Width
<input type="checkbox"/> 1503021	21" Seat Width	<input type="checkbox"/> 1503022	22" Seat Width		
<input type="checkbox"/> 1503023	23" Seat Width	<input type="checkbox"/> 1503024	24" Seat Width		
\$ 375					
<b>Adjustable Seat Widths</b> Order with Adjustable Frame.					
<input type="checkbox"/> 1503113	13"-15" Range - Set at 13"	<input type="checkbox"/> 1503114	13"-15" Range - Set at 14"	<input type="checkbox"/> 1503115	13"-15" Range - Set at 15"
<input type="checkbox"/> 1503215	15"-20" Range - Set at 15"	<input type="checkbox"/> 1503216	15"-20" Range - Set at 16"	<input type="checkbox"/> 1503217	15"-20" Range - Set at 17"
<input type="checkbox"/> 1503218	15"-20" Range - Set at 18"	<input type="checkbox"/> 1503219	15"-20" Range - Set at 19"	<input type="checkbox"/> 1503220	15"-20" Range - Set at 20"

### Seat Depth

#### Standard Seat Depths (Adjustable)

- |   |   |   |
|---|---|---|
| <input type="checkbox"/> 1504015 15" Seat Depth | <input type="checkbox"/> 1504016 16" Seat Depth | <input type="checkbox"/> 1504017 17" Seat Depth |
| <input type="checkbox"/> 1504018 18" Seat Depth | <input type="checkbox"/> 1504019 19" Seat Depth | <input type="checkbox"/> 1504020 20" Seat Depth |

#### Custom Seat Depths Requires Push Bar Back - Straight. Not compatible with Reclining or Fold Down Back.

- |   |   |
|---|---|
| <input type="checkbox"/> 1504021 21" Seat Depth | <input type="checkbox"/> 1504022 22" Seat Depth |
|---|---|

\$ 400

### Frame Color

#### PDG Premium Colors

- |   |  |   |
|---|--|---|
| <input type="checkbox"/> 1505000 Tracer Red         | <input type="checkbox"/> 1505001 Gloss Black Shimmer | <input type="checkbox"/> 1505002 Silver             |
| <input type="checkbox"/> 1505003 Fine Texture Black | <input type="checkbox"/> 1505004 Orange Fusion       | <input type="checkbox"/> 1505005 Karl's Imagination |
| <input type="checkbox"/> 1505006 Blue Shimmer       | <input type="checkbox"/> 1505007 Gloss White         | <input type="checkbox"/> 1505008 Dark Titanium      |
| <input type="checkbox"/> 1505009 Geyser Green       | <input type="checkbox"/> 1505010 Lollypop Blue       | <input type="checkbox"/> 1505011 Lollypop Red       |
| <input type="checkbox"/> 1505012 Bubble Gum Pink    | <input type="checkbox"/> 1505013 Berry Purple        | <input type="checkbox"/> 1505014 Daffodil Yellow    |

- 1505200 **Custom Color** Contact PDG for specifications.

\$ 300

### Armrests

- 1506000 **No Armrests**

- 1506114 **Single Post Height Adjustable Armrests (10"-14")**

\$ 295

- 1506100 **Single Post Fixed Height Armrests (12")** Full Length Armpad only.



#### Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.

Seat Pan to top of Armrest: \_\_\_\_\_ Does not include height of Armpad



#### Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.

- |                                      |  |
|--------------------------------------|--|
| <input type="checkbox"/> 1506402 Mid | <input type="checkbox"/> 1506403 Front |
|--------------------------------------|--|

#### Locking Flip-up Armrests Not recommended for users with excessive strength or tone.

\$ 295

Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.

- |   |   |   |
|---|---|---|
| <input type="checkbox"/> 1506310 Set at 10" | <input type="checkbox"/> 1506311 Set at 11" | <input type="checkbox"/> 1506312 Set at 12" |
| <input type="checkbox"/> 1506313 Set at 13" | <input type="checkbox"/> 1506314 Set at 14" |   |

### Armpads

#### Regular Armpads

- |   |   |
|---|---|
| <input type="checkbox"/> 1507100 Desk Length - 2" x 10" | <input type="checkbox"/> 1507101 Full Length - 2" x 14" |
|---|---|

#### Waterfall Armpads

\$ 115

- |   |   |
|---|---|
| <input type="checkbox"/> 1507105 Desk Length - 2" x 10" | <input type="checkbox"/> 1507106 Full Length - 2" x 14" |
|---|---|

#### Wide Armpads May interfere with wheels in some tilt configurations.

\$ 200

- |   |   |
|---|---|
| <input type="checkbox"/> 1507110 Desk Length - 3" x 10" | <input type="checkbox"/> 1507111 Full Length - 3" x 12" |
|---|---|

#### Gel Armpads May interfere with wheels in some tilt configurations.

\$ 135

- |   |   |
|---|---|
| <input type="checkbox"/> 1507120 Full Length - 2" x 12" | <input type="checkbox"/> 1507121 Full Length - 3.5" x 12" |
|---|---|

#### Contour Arm Trough with Hand Pad

\$ 125

Not compatible with Lap Tray, Locking Flip-up Armrests, Dual Post Armrests. May interfere with wheels in some tilt configurations.

- |   |  |
|---|--|
| <input type="checkbox"/> 1507130 Right (4.5" x 18") | <input type="checkbox"/> 1507131 Left (4.5" x 18") |
|---|--|

### Armrest Options

- 1507210 **Armrest Wideners** Move Armrests outward by 1" on each side. Do not increase seat size.

\$ 195

### Seat Pan

- 1508000 **No Seat Pan**

- 1508001 **Solid Aluminum Seat Pan**

\$ 135

### Seat Height



**Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Axle Plate (Low/Standard/High).

5" Casters								
RW	13"	14"	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	High STF
16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Std. STF
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

6" Casters & 6" Forks								
RW	13"	14"	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	High STF
16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Std. STF
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

6" Casters & 8" Forks								
RW	13"	14"	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	Std. STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	High STF
16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters & 8" Forks								
RW	13"	14"	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	High STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16"	<input type="checkbox"/>	Std. STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Frog Legs 6" Soft Roll & Suspension Forks								
RW	13"	14"	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	High STF
16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Std. STF
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Casters

#### Urethane Casters

- 1511051 5" x 1.25" Caster size  1511061 6" x 1.25" Caster size  1511081 8" x 1.25" Caster size
- 1511063 6" x 1.50" Caster size \$ 85
- 1511084 8" x 2.00" Caster size \$ 135

#### Pneumatic Casters

- 1511161 6" x 1.25" Caster size \$ 110
- 1511181 8" x 1.25" Caster size \$ 120
- 1511184 8" x 2.00" Caster size \$ 150

- 1511260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275

- 1511400 **Frog Legs Suspension Forks** Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster. \$ 595

### Rear Wheels

#### Composite Wheels

- 1512012 12" Wheel size  1512016 16" Wheel size  1512020 20" Wheel size
- 1512022 22" Wheel size  1512024 24" Wheel size

#### Spoke Wheels

- 1512120 20" Wheel size  1512122 22" Wheel size  1512124 24" Wheel size \$ 75

- 1512224 **24" Heavy Duty Spoke Wheels** \$ 350



#### Rear Wheel Position

 Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- 1512401 Rear More stable  1512402 Mid  1512403 Front Less stable

### Axle Type

- 1513000 **Quick Release Rear Axles**

- 1513010 **Fixed Rear Axles**

### Rear Tires

<input type="checkbox"/>	1514000	<b>Urethane Tires</b>		
<input type="checkbox"/>	1514010	<b>Pneumatic Tires</b>	May result in a finished Seat Height 1/2" greater than stated.	\$ 60
<input type="checkbox"/>	1514020	<b>Pneumatic Tires with Flat-Free Inserts</b>	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$ 95

### Handrims

Handrims are available only with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/>	1515000	<b>No Handrims</b>		
<input type="checkbox"/>	1515010	<b>Standard Aluminum Handrims</b>		
<input type="checkbox"/>	1515020	<b>Plastic Coated Handrims</b>		\$ 150

### Wheel Locks & Brakes

#### User Wheel Locks

<input type="checkbox"/>	1516010	Push to Lock	<input type="checkbox"/>	1516020	Pull to Lock	
<input type="checkbox"/>	1516030	<b>Attendant Wheel Locks (Push to Lock)</b>	Not compatible with Drum Brakes, Wheel Lock Extensions.			\$ 85
<input type="checkbox"/>	1516040	<b>User &amp; Attendant Wheel Locks</b>				\$ 125
<input type="checkbox"/>	1516005	<b>Foot Operated Wheel Lock</b>	Rear-mounted unilateral single-activation lock. Not compatible with O2 Holder, IV Pole, Ventilator Tray, One-Arm Drive.			\$ 595
<input type="checkbox"/>	1516300	<b>Drum Brakes</b>	Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.			\$ 795

### Wheel Lock Options

<input type="checkbox"/>	1516100	<b>Wheel Lock Extensions</b>	Not compatible with Attendant Wheel Lock.	\$ 60
--------------------------	---------	------------------------------	---	-------

### Anti-Tippers

<input checked="" type="checkbox"/>	1517010	<b>Rear Anti-Tippers</b>	Required Safety Option.	\$ 75
-------------------------------------	---------	--------------------------	-------------------------	-------

### Back Posts

#### Back Post Style

<input type="checkbox"/>	1518100	Push Bar Back	<input type="checkbox"/>	1518105	Cane Handle Back	Not compatible with Headrest Mounting Plate.		
<b>Back Post Type</b>								
<input type="checkbox"/>	1518200	Straight	<input type="checkbox"/>	1518205	8° Bend			
<b>Back Post Height</b> Measured from the Seat Pan to the Handle.								
<input type="checkbox"/>	1518420	Short (17")	<input type="checkbox"/>	1518425	Tall (23")			
<b>Seat to Back Angle</b>								
<input type="checkbox"/>	1519090	Set at 90°	<input type="checkbox"/>	1519096	Set at 98°	<input type="checkbox"/>	1519102	Set at 106°
<input type="checkbox"/>	1519108	Set at 114°	<input type="checkbox"/>	1519114	Set at 122°			

### Back Options

<b>Angle/Height Adjustable Push Bar Handle</b>				Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/>	1519010	Button on Outside	<input type="checkbox"/>	1519011	Button on Inside	
<input type="checkbox"/>	1519050	<b>Removable Cane Handle Push Bar</b>	Not compatible with Push Bar Back.			\$ 140
<input type="checkbox"/>	1519212	<b>Headrest Mounting Plate</b>	Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.			\$ 95
<input type="checkbox"/>	1519220	<b>Fold Down Back</b>	<input type="checkbox"/>	1519221	<b>Non-Folding Back</b>	
<input type="checkbox"/>	1519240	<b>Seating Dynamics Dynamic Rocker Back Interface</b>				\$ 2,995
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle, Standard Upholstery. Finished Back Post Height is 5" taller than selected Back Post Height.						



# FUZE T50

**Back Upholstery**

1520000 **No Back Upholstery**

**Standard Upholstery** Not available with Cane Handle Back. \$ 195

1520418 Short (18" Height)  1520421 Regular (21" Height)

1520450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. \$ 295

**Position of Upholstery** Select the desired position of the top of the Upholstery relative to the Seat Pan.

17" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								

23" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

**Tilt Activation**

**Attendant Tilt only** Select trigger mounting location below.

1523011 Right Back Post  1523012 Left Back Post

**Power Tilt** 40° tilt range. Includes 2 batteries for off-chair charging. [pdgmobility.com/power-tilt](http://pdgmobility.com/power-tilt)

1523400 250 lbs Weight Capacity \$ 1,795

1523401 350 lbs Weight Capacity Requires Fixed Frame. \$ 2,095

Toggle Switch  1523421 Installed on Right side  1523422 Installed on Left side

Button Switch  1523431 Installed on Right side  1523432 Installed on Left side

Stealth Egg Switch  1523440 Not installed

Optional for Egg Switch:  1523445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625

Specify Headrest to be used: \_\_\_\_\_

**Tilt Limit** Optional. If no tilt limit is selected, the full tilt range will be available.

1523220 Set at 20°  1523230 Set at 30°  1523240 Set at 40°

**Recline**

1523310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back. Not available with 22" Seat Depth. Finished Back Post Height is 3" taller than selected Back Post Height.

**Front Rigging**

1525000 **No Front Rigging**  1525001 **Sector Blocks only**

**Easy-Out™ Front Rigging** Effortless insertion and removal. No Hanger Bracket remains on frame. [pdgmobility.com/easy-out](http://pdgmobility.com/easy-out)

1525410 Fixed 70° \$ 95

1525420 Universal Angle Adjustable (0°-110°) Set at angle: \_\_\_\_\_ (use increments of 10°) \$ 295

Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.

**Swing Away Front Rigging** Legacy Design. Not compatible with Easy-Out. Requires pre-tilt of 5°.

1525140 Swing Away 70°

1525160 Swing Away 90° Includes 11"-16" Extension Tubes. Not compatible with 8" Casters. \$ 75

1525200 Swing Away Elevating Legrests Includes Calf Pads. \$ 375

1525300 1" Offset Hangers for Swing Away \$ 75

1525110 **Standard Plug In 70° Front Rigging** Legacy Design. Not compatible with Easy-Out.

**Footplates**

1526010 **Composite Footplates**

1526020 **Aluminum Angle Adjustable Footplates** \$ 165

1526030 **One Piece Fixed Footplate** Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. \$ 350


1526040 **One Piece Flip-up Footplate** Requires selection of Composite or Aluminum Angle Adjustable. \$ 395

1526050 **Interlocking Angle Adjustable Footplates** Available for Seat Widths up to 21" only. \$ 195

## Footrest/Legrest Accessories

<b>Calf Pad</b> Not compatible with Gel Padding.		
<input type="checkbox"/>	1527012 For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	1527010 For Swing Away Front Rigging	\$ 135
<input type="checkbox"/>	1527020 <b>Gel Padding</b> Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
<b>Gel Knee Pad</b>		
<input type="checkbox"/>	1527032 For Easy-Out™ Angle Adjustable Front Rigging	\$ 200
<input type="checkbox"/>	1527030 For Swing Away ELR Front Rigging	\$ 100
<input type="checkbox"/>	1527040 <b>Length Adjustable Heel Loops</b>	\$ 50
<input type="checkbox"/>	1527050 <b>Bodypoint Padded Calf Strap</b> Sizes greater than 22" are not padded.	\$ 195

## Front Rigging Extension Tubes

<b>Regular Extension Tubes</b> No selection required for 90° Front Rigging		
<input type="checkbox"/>	1528010 Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	1528020 Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used. Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.	

## Miscellaneous Options

<b>One-Arm Drive</b> Requires 24" Composite Rear Wheels and forward wheel position. Available with 18" and 20" STF only. Not compatible with Attendant or Foot Operated Wheel Locks, Drum Brakes.		\$ 3,500	
<input type="checkbox"/>	1530010 Installed on Right side	<input type="checkbox"/>	1530011 Installed on Left side
<b>Spoke Guards</b> Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.		\$ 195	
<input type="checkbox"/>	1530020 Black	<input type="checkbox"/>	1530021 Frame Color
<input type="checkbox"/>	1530060 <b>Stealth Flip Down Abductor</b> Requires Solid Seat Pan.	\$ 595	
<input type="checkbox"/>	1530070 <b>QR Bracket - Accessory Mount</b>	\$ 200	
<input type="checkbox"/>	1530080 <b>Gel Hip Guides</b> Not compatible with Dual Post Arms. Reduces available Seat Width by 2".	\$ 395	
<input type="checkbox"/>	1530090 <b>Hip Belt Mount</b> Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.	\$ 55	
<input type="checkbox"/>	1530110 <b>IV Pole</b>	\$ 295	
<input type="checkbox"/>	1530120 <b>O2 Holder</b> May add to overall width of the Chair.	\$ 395	
<input type="checkbox"/>	1530150 <b>Ventilator Tray</b> Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 16". Not compatible with O2 Holder, Power Tilt, Foot Operated Wheel Locks.	\$ 295	
<input type="checkbox"/>	1530152 <b>Ventilator Bracket</b> Compatible with LTV1 (laptop) ventilator.	\$ 295	
<input type="checkbox"/>	1530154 <b>Ventilator Battery Tray</b> 16" minimum Seat Width. Reduces ground clearance. Fits U1 Battery only. Restricted to 20", 22", and 24" Rear Wheels with STF restrictions.	\$ 195	
<input type="checkbox"/>	1530160 <b>Clear Lap Tray (Slide On)</b> Available for Seat Widths of 16"-22" only. Not compatible with Contour Arm Trough.	\$ 490	
<input type="checkbox"/>	1530170 <b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.	\$ 10	

## On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

## Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.