



FUZE T50n

Order Form - Canada

Prices in Canadian Dollars. January 1, 2022

Customer and Order Information

Purchase Order *	Tag	v.1117
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

Base Chair

Fuze T50n 50° Manual Tilt Wheelchair

<input checked="" type="checkbox"/>	2601000	Fuze T50 50° Manual Tilt Wheelchair	\$ 4,150
<input checked="" type="checkbox"/>	2601008	Narrow Frame Upgrade Reduces overall chair width to seat width plus only 4".	\$ 200
<input type="checkbox"/>	The standard Fuze T50 offers additional options that are not available on the T50n. See separate order form.		

Frame Type & Weight Capacity

Fixed Width Frame Select a Fixed Seat Width below.

<input type="checkbox"/>	2602010	300 lbs Weight Capacity	
<input type="checkbox"/>	2602011	400 lbs Weight Capacity Not compatible with Transit Tie Down System.	
<input type="checkbox"/>	2602025	Adjustable Width Frame 250 lbs Weight Capacity. Select an Adjustable Seat Width below.	\$ 100

Frame Options

<input type="checkbox"/>	2602222	Take-Apart Frame The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. pdgmobility.com/take-apart	\$ 500
This feature reduces the weight of the heaviest piece of the wheelchair to 29 lbs or less. Perfect for transporting it in the trunk of a car or SUV. Not compatible with Power Tilt, Ventilator Tray, Ventilator Battery Tray.			
<input type="checkbox"/>	2602300	Transit Tie Down System (TTDS)	
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.			

Seat Width

Overall chair width is +4.0" for 12" RW

Fixed Seat Widths Order with Fixed Frame.

<input type="checkbox"/>	2603016	16" Seat Width	<input type="checkbox"/>	2603017	17" Seat Width	<input type="checkbox"/>	2603018	18" Seat Width
<input type="checkbox"/>	2603019	19" Seat Width	<input type="checkbox"/>	2603020	20" Seat Width	<input type="checkbox"/>	2603021	21" Seat Width
<input type="checkbox"/>	2603022	22" Seat Width						
<input type="checkbox"/>	2603023	23" Seat Width	<input type="checkbox"/>	2603024	24" Seat Width			\$ 375

Adjustable Seat Widths Order with Adjustable Frame.

<input type="checkbox"/>	2603216	16"-19" Range - Set at 16"	<input type="checkbox"/>	2603217	16"-19" Range - Set at 17"
<input type="checkbox"/>	2603218	16"-19" Range - Set at 18"	<input type="checkbox"/>	2603219	16"-19" Range - Set at 19"
<input type="checkbox"/>	2603319	19"-20" Range - Set at 19"	<input type="checkbox"/>	2603320	19"-20" Range - Set at 20"

Seat Depth

Standard Seat Depths (Adjustable)

<input type="checkbox"/>	2604015	15" Seat Depth	<input type="checkbox"/>	2604016	16" Seat Depth	<input type="checkbox"/>	2604017	17" Seat Depth
<input type="checkbox"/>	2604018	18" Seat Depth	<input type="checkbox"/>	2604019	19" Seat Depth	<input type="checkbox"/>	2604020	20" Seat Depth
<input type="checkbox"/>	2604021	21" Seat Depth	<input type="checkbox"/>	2604022	22" Seat Depth			\$ 400

Custom Seat Depths Requires Push Bar Back - Straight. Not compatible with Reclining or Fold Down Back.

Frame Color

PDG Premium Colors

- | | | |
|---|--|---|
| <input type="checkbox"/> 2605000 Tracer Red | <input type="checkbox"/> 2605001 Gloss Black Shimmer | <input type="checkbox"/> 2605002 Silver |
| <input type="checkbox"/> 2605003 Fine Texture Black | <input type="checkbox"/> 2605004 Orange Fusion | <input type="checkbox"/> 2605005 Karl's Imagination |
| <input type="checkbox"/> 2605006 Blue Shimmer | <input type="checkbox"/> 2605007 Gloss White | <input type="checkbox"/> 2605008 Dark Titanium |
| <input type="checkbox"/> 2605009 Geyser Green | <input type="checkbox"/> 2605010 Lollipop Blue | <input type="checkbox"/> 2605011 Lollipop Red |
| <input type="checkbox"/> 2605012 Bubble Gum Pink | <input type="checkbox"/> 2605013 Berry Purple | <input type="checkbox"/> 2605014 Daffodil Yellow |

- 2605200 **Custom Color** Contact PDG for specifications. \$ 300

Armrests

- 2606000 **No Armrests**
- 2606114 **Single Post Height Adjustable Armrests (10"-14")** \$ 295



Armrest Height

If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.

Seat Pan to top of Armrest: _____ Does not include height of Armpad



Armrest Position

Indicate preference for Armrest position. If left blank, the default setting will be used.

- 2606402 Mid 2606403 Front

Locking Flip-up Armrests

Not recommended for users with excessive strength or tone. \$ 295

Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.

- 2606310 Set at 10" 2606311 Set at 11" 2606312 Set at 12"
- 2606313 Set at 13" 2606314 Set at 14"

Armpads

Regular Armpads

- 2607100 Desk Length - 2" x 10" 2607101 Full Length - 2" x 14"

Waterfall Armpads

- 2607105 Desk Length - 2" x 10" 2607106 Full Length - 2" x 14" \$ 115

Wide Armpads

May interfere with wheels in some tilt configurations. \$ 200

- 2607110 Desk Length - 3" x 10" 2607111 Full Length - 3" x 12"

Gel Armpads

May interfere with wheels in some tilt configurations. \$ 135

- 2607120 Full Length - 2" x 12" 2607121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad

Not compatible with Lap Tray, Locking Flip-up Armrests. May interfere with wheels in some tilt configurations. \$ 125

- 2607130 Right (4.5" x 18") 2607131 Left (4.5" x 18")

Armrest Options

- 2607210 **Armrest Wideners** Move Armrests outward by 1" on each side. Do not increase seat size. May increase overall width of the chair. \$ 195

Seat Pan

- 2608000 **No Seat Pan**
- 2608001 **Solid Aluminum Seat Pan** \$ 135

Seat Height

- Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable.

12" Rear Wheels					
Caster size	16"	17"	18"	19"	20"
5"	<input type="checkbox"/>	<input type="checkbox"/>			
6"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
8"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6" Soft Roll	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Casters

Urethane Casters

- | | | | | | | |
|----------------------------------|------------------------|----------------------------------|------------------------|----------------------------------|------------------------|--------|
| <input type="checkbox"/> 2611051 | 5" x 1.25" Caster size | <input type="checkbox"/> 2611061 | 6" x 1.25" Caster size | <input type="checkbox"/> 2611081 | 8" x 1.25" Caster size | |
| <input type="checkbox"/> 2611063 | 6" x 1.50" Caster size | | | | | \$ 85 |
| <input type="checkbox"/> 2611084 | 8" x 2.00" Caster size | | | | | \$ 135 |

Pneumatic Casters

- | | | | | | | |
|----------------------------------|------------------------|--|--|--|--|--------|
| <input type="checkbox"/> 2611161 | 6" x 1.25" Caster size | | | | | \$ 110 |
| <input type="checkbox"/> 2611181 | 8" x 1.25" Caster size | | | | | \$ 120 |
| <input type="checkbox"/> 2611184 | 8" x 2.00" Caster size | | | | | \$ 150 |

- | | | | | | | |
|----------------------------------|--|--|--|--|--|--------|
| <input type="checkbox"/> 2611260 | Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters | | | | | \$ 275 |
|----------------------------------|--|--|--|--|--|--------|

- | | | | | | | |
|----------------------------------|-----------------------------------|--|--|--|--|--------|
| <input type="checkbox"/> 2611400 | Frog Legs Suspension Forks | Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster. | | | | \$ 595 |
|----------------------------------|-----------------------------------|--|--|--|--|--------|

Rear Wheels

- | | | | | | | |
|---|-----------------------------|--|--|--|--|--|
| <input checked="" type="checkbox"/> 2612012 | 12" Composite Wheels | | | | | |
|---|-----------------------------|--|--|--|--|--|



Rear Wheel Position

Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- | | | | | | | |
|----------------------------------|------------------|----------------------------------|-----|----------------------------------|-------------------|--|
| <input type="checkbox"/> 2612401 | Rear More stable | <input type="checkbox"/> 2612402 | Mid | <input type="checkbox"/> 2612403 | Front Less stable | |
|----------------------------------|------------------|----------------------------------|-----|----------------------------------|-------------------|--|

Axle Type

- | | | | | | | |
|----------------------------------|---------------------------------|--|--|--|--|--|
| <input type="checkbox"/> 2613000 | Quick Release Rear Axles | | | | | |
|----------------------------------|---------------------------------|--|--|--|--|--|

- | | | | | | | |
|----------------------------------|-------------------------|--|--|--|--|--|
| <input type="checkbox"/> 2613010 | Fixed Rear Axles | | | | | |
|----------------------------------|-------------------------|--|--|--|--|--|

Rear Tires

- | | | | | | | |
|----------------------------------|-----------------------|--|--|--|--|--|
| <input type="checkbox"/> 2614000 | Urethane Tires | | | | | |
|----------------------------------|-----------------------|--|--|--|--|--|

- | | | | | | | |
|----------------------------------|------------------------|--|--|--|--|-------|
| <input type="checkbox"/> 2614010 | Pneumatic Tires | May result in a finished Seat Height 1/2" greater than stated. | | | | \$ 60 |
|----------------------------------|------------------------|--|--|--|--|-------|

- | | | | | | | |
|----------------------------------|---|--|--|--|--|-------|
| <input type="checkbox"/> 2614020 | Pneumatic Tires with Flat-Free Inserts | May result in a finished Seat Height 1/2" greater than stated. | | | | \$ 95 |
|----------------------------------|---|--|--|--|--|-------|

Handrims

- | | | | | | | |
|---|--------------------|--|--|--|--|--|
| <input checked="" type="checkbox"/> 2615000 | No Handrims | | | | | |
|---|--------------------|--|--|--|--|--|

Wheel Locks

- | | | | | | | |
|----------------------------------|---|--|--|--|--|--|
| <input type="checkbox"/> 2616030 | Attendant Wheel Locks (Push to Lock) | | | | | |
|----------------------------------|---|--|--|--|--|--|

- | | | | | | | |
|----------------------------------|---------------------------------|---|--|--|--|--------|
| <input type="checkbox"/> 2616005 | Foot Operated Wheel Lock | Rear-mounted unilateral single-activation lock. Compatible only with 19" & 20" STF Height. Not compatible with Ventilator Tray. | | | | \$ 595 |
|----------------------------------|---------------------------------|---|--|--|--|--------|

Wheel Lock Options

- | | | | | | | |
|----------------------------------|------------------------------|--|--|--|--|-------|
| <input type="checkbox"/> 2616100 | Wheel Lock Extensions | | | | | \$ 60 |
|----------------------------------|------------------------------|--|--|--|--|-------|

Anti-Tippers

- | | | | | | | |
|---|--------------------------|-------------------------|--|--|--|-------|
| <input checked="" type="checkbox"/> 2617010 | Rear Anti-Tippers | Required Safety Option. | | | | \$ 75 |
|---|--------------------------|-------------------------|--|--|--|-------|

Back Posts

Back Post Style

- | | | | | | | |
|----------------------------------|---------------|----------------------------------|------------------|--|--|--|
| <input type="checkbox"/> 2618100 | Push Bar Back | <input type="checkbox"/> 2618105 | Cane Handle Back | Not compatible with Headrest Mounting Plate. | | |
|----------------------------------|---------------|----------------------------------|------------------|--|--|--|

Back Post Type

- | | | | | | | |
|----------------------------------|----------|----------------------------------|---------|--|--|--|
| <input type="checkbox"/> 2618200 | Straight | <input type="checkbox"/> 2618205 | 8° Bend | | | |
|----------------------------------|----------|----------------------------------|---------|--|--|--|

Back Post Height

Measured from the Seat Pan to the Handle.

- | | | | | | | |
|----------------------------------|-------------|----------------------------------|------------|--|--|--|
| <input type="checkbox"/> 2618420 | Short (17") | <input type="checkbox"/> 2618425 | Tall (23") | | | |
|----------------------------------|-------------|----------------------------------|------------|--|--|--|

See next page for more options >



FUZE T50n

Order Form - Canada

Prices in Canadian Dollars. January 1, 2022

Seat to Back Angle

- 2619090 Set at 90°
 2619096 Set at 98°
 2619102 Set at 106°
 2619108 Set at 114°
 2619114 Set at 122°

Back Options

- Angle/Height Adjustable Push Bar Handle** Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150
 2619010 Button on Outside 2619011 Button on Inside
 2619050 **Removable Cane Handle Push Bar** Not compatible with Push Bar Back. \$ 140
 2619212 **Headrest Mounting Plate** Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back. \$ 95
 2619220 **Fold Down Back** 2619221 **Non-Folding Back**
 2619240 **Seating Dynamics Dynamic Rocker Back Interface** \$ 2,995
 Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle, Standard Upholstery.
 Finished Back Post Height is 5" taller than selected Back Post Height.

Back Upholstery

- 2620000 **No Back Upholstery**
Standard Upholstery Not available with Cane Handle Back. \$ 195
 2620418 Short (18" Height) 2620421 Regular (21" Height)
 2620450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. \$ 295

Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

		17" Back Post							
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								

		23" Back Post							
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation

- Attendant Tilt only** Select trigger mounting location below.
 2623011 Right Back Post 2623012 Left Back Post
Power Tilt 40° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt
 2623400 250 lbs Weight Capacity \$ 1,795
 2623401 350 lbs Weight Capacity Requires Fixed Frame. \$ 2,095
 Toggle Switch 2623421 Installed on Right side 2623422 Installed on Left side
 Button Switch 2623431 Installed on Right side 2623432 Installed on Left side
 Stealth Egg Switch 2623440 Not installed
 Optional for Egg Switch: 2623445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625
 Specify Headrest to be used: _____

Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.

- 2623220 Set at 20°
 2623230 Set at 30°
 2623240 Set at 40°

Recline

- 2623310 **Reclining Backrest (Manual)** \$ 895
 30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back. Not available with 22" Seat Depth.
 Finished Back Post Height is 3" taller than selected Back Post Height.

Front Rigging

<input type="checkbox"/>	2625000	No Front Rigging	<input type="checkbox"/>	2625001	Sector Blocks only		
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out							
<input type="checkbox"/>	2625410	Fixed 70°					\$ 95
<input type="checkbox"/>	2625420	Universal Angle Adjustable (0°-110°)	Set at angle: _____	(use increments of 10°)			\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.							
Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out. Requires pre-tilt of 5°.							
<input type="checkbox"/>	2625140	Swing Away 70°					
<input type="checkbox"/>	2625160	Swing Away 90°	Includes 11"-16" Extension Tubes. Not compatible with 8" Casters.				\$ 75
<input type="checkbox"/>	2625200	Swing Away Elevating Legrests	Includes Calf Pads.				\$ 375
<input type="checkbox"/>	2625300	1" Offset Hangers for Swing Away					\$ 75
<input type="checkbox"/>	2625110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.				


Footplates

<input type="checkbox"/>	2626010	Composite Footplates					
<input type="checkbox"/>	2626020	Aluminum Angle Adjustable Footplates					\$ 165
<input type="checkbox"/>	2626030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.				\$ 350
<input type="checkbox"/>	2626040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.				\$ 395
<input type="checkbox"/>	2626050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.				\$ 195

Footrest/Legrest Accessories

Calf Pad Not compatible with Gel Padding.							
<input type="checkbox"/>	2627012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.				\$ 200
<input type="checkbox"/>	2627010	For Swing Away Front Rigging					\$ 135
<input type="checkbox"/>	2627020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.				\$ 125
Gel Knee Pad							
<input type="checkbox"/>	2627032	For Easy-Out™ Angle Adjustable Front Rigging					\$ 200
<input type="checkbox"/>	2627030	For Swing Away ELR Front Rigging					\$ 100
<input type="checkbox"/>	2627040	Length Adjustable Heel Loops					\$ 50
<input type="checkbox"/>	2627050	Bodypoint Padded Calf Strap	Sizes greater than 22" are not padded.				\$ 195

Front Rigging Extension Tubes

Regular Extension Tubes No selection required for 90° Front Rigging							
<input type="checkbox"/>	2628010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging			
<input type="checkbox"/>	2628020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging			\$ 75
 Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.							
Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.							



FUZE T50n

Order Form - Canada

Prices in Canadian Dollars. January 1, 2022

Miscellaneous Options

<input type="checkbox"/>	2630060	Stealth Flip Down Abductor	Requires Solid Seat Pan.	\$ 595
<input type="checkbox"/>	2630070	QR Bracket - Accessory Mount		\$ 200
<input type="checkbox"/>	2630080	Gel Hip Guides	Reduces available Seat Width by 2".	\$ 395
<input type="checkbox"/>	2630090	Hip Belt Mount	Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.	\$ 55
<input type="checkbox"/>	2630150	Ventilator Tray	Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 20". Not compatible with Power Tilt, Foot Operated Wheel Locks.	\$ 295
<input type="checkbox"/>	2630152	Ventilator Bracket	Compatible with LTV1 (laptop) ventilator.	\$ 295
<input type="checkbox"/>	2630154	Ventilator Battery Tray	20" min. Seat Width. Reduces ground clearance. Fits U1 Battery only. STF restrictions.	\$ 195
<input type="checkbox"/>	2630160	Clear Lap Tray (Slide On)	Available for Seat Widths of 16"-22" only. Not compatible with Contour Arm Trough.	\$ 490
<input type="checkbox"/>	2630170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$ 10

On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.