

### Cushion

1. Select the type



2. Select the model



General Use	<b>Stealth Simplicity</b> 350 lbs weight capacity.		
	<input type="checkbox"/> ##09141	16"-20" width, 16"-20" depth	E2601 \$ 310
	<input type="checkbox"/> ##09142	22"-24" width, 18"-20" depth	E2602 \$ 370

Positioning	<b>ROHO AirLITE</b> 300 lbs weight capacity.		
	<input type="checkbox"/> ##09120	14"-20" width, 14"-18" depth	E2605 \$ 375

Skin Protection & Positioning	<b>Stealth Spectrum Foam</b> 650 lbs weight capacity.		
	<input type="checkbox"/> ##09143	16"-20" width, 16"-24" depth	E2607 \$ 580
	<input type="checkbox"/> ##09144	22"-32" width, 18"-28" depth	E2608 \$ 630
	<b>Stealth Tru-Comfort 2</b> 450 lbs weight capacity.		
	<input type="checkbox"/> ##09146	16"-20" width, 16"-24" depth	E2607 \$ 515
	<input type="checkbox"/> ##09147	22" width, 18"-24" depth	E2608 \$ 590
	<b>nxt BioFit</b> 300 lbs weight capacity.		
	<input type="checkbox"/> ##09101	14"-20" width, 16"-20" depth	E2607 \$ 595
	<input type="checkbox"/> ##09102	22"-24" width, 16"-20" depth	E2608 \$ 625
	<b>Matrx VI</b> 300 lbs weight capacity.		
	<input type="checkbox"/> ##09111	14"-20" width, 15"-20" depth	E2607 \$ 680
	<input type="checkbox"/> ##09112	22" width, 18"-20" depth	E2608 \$ 745
<b>Matrx VI Heavy Duty</b> 600 lbs weight capacity.			
<input type="checkbox"/> ##09115	18"-21" width, 18"-26" depth	E2607 \$ 1,050	
<input type="checkbox"/> ##09116	22"-30" width, 18"-26" depth	E2608 \$ 1,250	

Skin Protection & Adjustable	<b>ROHO Single Valve High Profile</b> No weight limit if cushion is properly sized.		
	<input type="checkbox"/> ##09125	14"-20" width, 16"-22" depth	E2622 \$ 840
	<input type="checkbox"/> ##09126	22"-24" width, 16"-22" depth	E2623 \$ 960
	<b>STAR Stabil-Air Single Valve</b> No weight limit if cushion is properly sized.		
	<input type="checkbox"/> ##09151	16"-20" width, 16"-20" depth	Standard (even) sizes. E2622 \$ 1,200
<input type="checkbox"/> ##09152	16"-20" width, 16"-20" depth	Custom (odd and larger) sizes. E2622 / E2623 \$ 1,300	

Skin Protection, Positioning & Adjustable	<b>ROHO Quadro Select High Profile</b> No weight limit if cushion is properly sized.		
	<input type="checkbox"/> ##09121	14"-20" width, 16"-22" depth	E2624 \$ 840
	<input type="checkbox"/> ##09122	22"-24" width, 16"-22" depth	E2625 \$ 960

Extra-Wide  No weight limit if cushion is properly sized. Custom sizes available.  Custom-made foam cushion with fluid-proof yet breathable fabric cover.	<b>PDG Solar Basic</b>		<input type="checkbox"/> ##09200	20"-28" width, 16"-24" depth	\$ 395
	High resiliency polyfoam with 2" medium-density base and 1" lower-density top.		<input type="checkbox"/> ##09201	29"-36" width, 16"-24" depth	\$ 495
			<input type="checkbox"/> ##09202	37"-48" width, 16"-24" depth	\$ 595
	<b>PDG Solar Contour</b>		<input type="checkbox"/> ##09210	20"-28" width, 16"-24" depth	\$ 695
	Aids proper pelvic alignment, includes pre-ischial bar, leg shelf, adductor buildups.		<input type="checkbox"/> ##09211	29"-36" width, 16"-24" depth	\$ 795
			<input type="checkbox"/> ##09212	37"-48" width, 16"-24" depth	\$ 895
	<b>PDG Solar Pressure Reduction</b>		<input type="checkbox"/> ##09220	20"-28" width, 16"-24" depth	\$ 1,095
	All features of the Contour, plus added pincore latex rubber and containment foam outline.		<input type="checkbox"/> ##09221	29"-36" width, 16"-24" depth	\$ 1,195
			<input type="checkbox"/> ##09222	37"-48" width, 16"-24" depth	\$ 1,295

3. Specify the size →

**Cushion Size** Sizes based on availability from the cushion manufacturer.

Specify Size (Width" x Depth") \_\_\_\_\_

## Backrest

1. Select the type  
↓

2. Select the contour depth  
↓

3. Select the model  
↓

Positioning, posterior-lateral	2" depth	<input type="checkbox"/> ##21141 <b>Stealth ADI AL Active</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height E2615 \$ 945
	4" depth	<input type="checkbox"/> ##21142 <b>Stealth ADI AL Deep</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height E2615 \$ 1,045
	6" depth	<input type="checkbox"/> ##21143 <b>Stealth ADI AL Extra Deep</b> 250 lbs weight capacity. 14"-20" width, 10"-20" height E2615 \$ 1,145

Positioning, planar back with lateral supports	2" to 3" depth	<input type="checkbox"/> ##21110 <b>Matrx Elite E2</b> 300 lbs weight capacity. 12"-20" width, 10"-20" height E2620 \$ 985
		<input type="checkbox"/> ##21111 <b>Matrx Elite E2 Heavy Duty</b> 500 lbs weight capacity. 21" width, 16"-20" height E2620 \$ 1,225
		<input type="checkbox"/> ##21112 22"-24" width, 16"-20" height E2621 \$ 1,325
		<input type="checkbox"/> ##21114 25"-30" width, 16"-20" height E2621 \$ 1,995
		<input type="checkbox"/> ##21120 <b>ROHO Agility Min</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height E2620 \$ 900
		<input type="checkbox"/> ##21100 <b>nxt Xtend</b> 300 lbs weight capacity. 16"-20" width, 16"-20" height E2620 \$ 690
	3.5" to 4" depth	<input type="checkbox"/> ##21105 22"-24" width, 16"-20" height E2621 \$ 720
		<input type="checkbox"/> ##21122 <b>ROHO Agility Mid</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height E2620 \$ 930
		<input type="checkbox"/> ##21101 <b>nxt Xtend Deep</b> 300 lbs weight capacity. 16"-20" width, 16"-20" height E2620 \$ 690
	6" to 7" depth	<input type="checkbox"/> ##21106 22"-24" width, 16"-20" height E2621 \$ 720
		<input type="checkbox"/> ##21116 <b>Matrx Elite E2 Deep</b> 300 lbs weight capacity. 14"-20" width, 12"-20" height E2620 \$ 1,175
		<input type="checkbox"/> ##21117 22"-24" width, 12"-20" height E2621 \$ 1,425
<input type="checkbox"/> ##21124 <b>ROHO Agility Max</b> 300 lbs weight capacity. 14"-20" width, 10"-18" height E2620 \$ 1,085		

4. Specify the size →

**Backrest Size** Sizes based on availability from the backrest manufacturer.

Specify Size (Width" x Height") \_\_\_\_\_

## Pediatric Seating System

4030400 **Diego Pediatric Seating System**

Select this item to request a customized quotation for delivery of a Diego System with this chair.

Back & Seat Matrix				
in (cm)	10	12	14	16
Width	(25.40)	(30.48)	(35.56)	(40.64)
10 (25.40)	✓	✓	✗	✗
12 (30.48)	✗	✓	✓	✗
14 (35.56)	✗	✗	✓	✓

**A Head System** Head Pad & Hardware

**B Contoured Back System** Base, Cushion & Cover

**C Contoured Seat System** Base, Cushion, Cover & Hardware

**Multi-Axis Hardware** The HMO501, included as part of a package for the Diego Head System, allows growth adjustability and on-demand configurations.

**FlipLock** Configurable as a between carers or full width set up.

**Width Growth** Adapted from our TFB System, the back wings allow for a pediatric growth adjustment.

**Back Hardware Cover** Easy to remove back cover for quick on-demand adjustments and a more aesthetic look.

**Coolcore Covers** Our Diego System uses the Coolcore cover technology which adapt all-year-round for a more comfortable experience.

**Laterals** Modular Laterals either Fixed or Swing Away

**diego** a Stealth Product

**WC-20** \* Size requirements need to be met for safety purposes. Head restraint IT

TBD

### Headrest

Stealth Comfort Plus Headrest	<b>Stealth Hardware</b>			
	<input type="checkbox"/>	##24110 Multi-Axis HMO475-1	E1028	\$ 225
	<input type="checkbox"/>	##24111 UniLink	E1028	\$ 350
	<b>Comfort Plus 2 Pad</b>			
	<input type="checkbox"/>	##24115 Large (10" x 6")	<input type="checkbox"/>	##24116 X-Large (14" x 6")
nxt MultiFit Head Support	<b>nxt Hardware</b>		E1028	\$ 285
	<b>Multi-Adjustable Pad</b>		E0955	\$ 295
	<input type="checkbox"/>	##24105 Contoured (9" x 5")	<input type="checkbox"/>	##24107 X-Large (14" x 5")
<input type="checkbox"/>	##24106 Large (10" x 5")			

### Positioning Belts & Straps

Bodypoint	<b>Hip Belts</b>			
	<input type="checkbox"/>	##29010 Center-Pull Auto Style	E0978	\$ 115
	<input type="checkbox"/>	##29015 Center-Pull Auto Style, Padded	E0978	\$ 215
	<input type="checkbox"/>	##29020 4-Point, Padded	E0978	\$ 235
<input type="checkbox"/>	##29030 <b>Chest Strap</b>	E0960	\$ 165	

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.

PDG is proud to supply on-chair seating and positioning products from these manufacturers:

