



# STELLAR GLT

Order Form - Canada

Prices in Canadian Dollars. January 1, 2021

## Customer and Order Information

Purchase Order *	Tag	v.0901
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

## Base Chair

<b>Stellar GLT 45° Manual Tilt Wheelchair</b>		
<input checked="" type="checkbox"/> 2501000	Stellar GL 45° Manual Tilt Wheelchair	\$ 4,150
<input checked="" type="checkbox"/> 2501009	Lightweight Aluminum Frame Upgrade Aluminum Frame is not recommended for agitated users.	\$ 300
<input type="checkbox"/> ⓘ	Choose the standard <b>Stellar GL</b> with a steel frame for agitated or heavier users. See separate order form.	

## Weight Capacity

<input checked="" type="checkbox"/> 2502000	<b>250 lbs Chair Weight Capacity</b>		
<b>User Weight Range</b> Required for setup of tilt uplift force.			
<input type="checkbox"/> 1902430	Up to 75 lbs	<input type="checkbox"/> 1902431	75 lbs to 100 lbs
<input type="checkbox"/> 1902432	100 lbs to 115 lbs	<input type="checkbox"/> 1902433	115 lbs to 150 lbs
<input type="checkbox"/> 1902434	150 lbs to 200 lbs	<input type="checkbox"/> 1902435	200 lbs to 250 lbs

## Frame Options

<input type="checkbox"/> 2502300	<b>Transit Tie Down System (TTDS)</b>
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.	

## Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

<b>Standard Seat Widths</b>			
<input type="checkbox"/> 2503014	14" Seat Width	<input type="checkbox"/> 1903015	15" Seat Width
<input type="checkbox"/> 1903016	16" Seat Width	<input type="checkbox"/> 2503017	17" Seat Width
<input type="checkbox"/> 1903018	18" Seat Width	<input type="checkbox"/> 2503018	18" Seat Width
<input type="checkbox"/> 1903019	19" Seat Width	<input type="checkbox"/> 2503020	20" Seat Width
<input type="checkbox"/> 1903021	21" Seat Width	<input type="checkbox"/> 1903022	22" Seat Width

## Seat Depth

<b>Standard Seat Depths (Adjustable)</b>			
<input type="checkbox"/> 2504016	16" Seat Depth	<input type="checkbox"/> 1904017	17" Seat Depth
<input type="checkbox"/> 1904018	18" Seat Depth	<input type="checkbox"/> 2504019	19" Seat Depth
<input type="checkbox"/> 1904020	20" Seat Depth	<input type="checkbox"/> 2504021	21" Seat Depth
<input type="checkbox"/> 1904022	22" Seat Depth		
<b>Custom Seat Depths</b>			
<input type="checkbox"/> 2504021	21" Seat Depth	<input type="checkbox"/> 1904022	22" Seat Depth
			\$ 400

## Frame Color

<b>PDG Premium Colors</b>			
<input type="checkbox"/> 2505000	Tracer Red	<input type="checkbox"/> 1905001	Gloss Black Shimmer
<input type="checkbox"/> 1905002	Silver	<input type="checkbox"/> 2505003	Fine Texture Black
<input type="checkbox"/> 1905003	Karl's Imagination	<input type="checkbox"/> 1905004	Orange Fusion
<input type="checkbox"/> 1905004	Dark Titanium	<input type="checkbox"/> 2505006	Blue Shimmer
<input type="checkbox"/> 1905005	Lollypop Red	<input type="checkbox"/> 1905007	Gloss White
<input type="checkbox"/> 1905006	Lollypop Blue	<input type="checkbox"/> 2505009	Geyser Green
<input type="checkbox"/> 1905007	Lollypop Red	<input type="checkbox"/> 1905010	Lollypop Blue
<input type="checkbox"/> 1905008	Daffodil Yellow	<input type="checkbox"/> 2505012	Bubble Gum Pink
<input type="checkbox"/> 1905009	Daffodil Yellow	<input type="checkbox"/> 1905013	Berry Purple
<input type="checkbox"/> 2505200	<b>Custom Color</b> Contact PDG for specifications.		
			\$ 300

### Armrests

2506000 **No Armrests**

**Height Adjustable Armrests (Single Post)** \$ 295

2506111 8"-11" Range  2506114 10"-14" Range

2506100 **Fixed Height Armrests (12")** Full Length Armpad only.

**Armrest Height** If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.  
 Seat Pan to top of Armrest: \_\_\_\_\_ Does not include height of Armpad

**Armrest Position** Indicate preference for Armrest position. If left blank, the default setting will be used.  
 In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2506401 Rear  2506402 Mid  2506403 Front

### Armpads

**Regular Armpads**

2507100 Desk Length - 2" x 10"  1907101 Full Length - 2" x 14"

**Wide Armpads** May interfere with wheels in some tilt configurations. \$ 200

2507110 Desk Length - 3" x 10"  1907111 Full Length - 3" x 12"

**Gel Armpads** May interfere with wheels in some tilt configurations. \$ 135

2507120 Full Length - 2" x 12"  1907121 Full Length - 3.5" x 12"

**Contour Arm Trough with Hand Pad** May interfere with wheels in some tilt configurations. \$ 125

2507130 Right (4.5" x 18")  1907131 Left (4.5" x 18")

### Armrest Options

2507210 **Armrest Wideners** Move Armrests outward by 1" on each side. Do not increase seat size. \$ 195

2507310 **Clothing Guards** Available for 10"-14" Height range only. \$ 110

### Seat Pan

2508000 **No Seat Pan**

2508001 **Solid Aluminum Seat Pan** \$ 135

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Casters

**Urethane Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.

2511061 6" x 1.25" Caster size  2511081 8" x 1.25" Caster size


2511063 6" x 1.50" Caster size \$ 85

2511084 8" x 2.00" Caster size \$ 135

See next page for more options >

<b>Pneumatic Casters</b>		8" Casters are compatible only with 16" and 24" Rear Wheels.	
<input type="checkbox"/>	2511161	6" x 1.25" Caster size	\$ 110
<input type="checkbox"/>	2511181	8" x 1.25" Caster size	\$ 120
<input type="checkbox"/>	2511184	8" x 2.00" Caster size	\$ 150
<input type="checkbox"/>	2511260	<b>Frog Legs 6" x 1.40" Soft Roll Casters</b>	\$ 95
<input type="checkbox"/>	2511400	<b>Frog Legs Suspension Forks</b> Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster.	\$ 395

### Rear Wheels

<b>Composite Wheels</b>			
<input type="checkbox"/>	2512012	12" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	2512022	22" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1912016	16" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1912024	24" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1912020	20" Wheel size	
<b>Spoke Wheels</b>			\$ 75
<input type="checkbox"/>	2512120	20" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1912122	22" Wheel size	<input type="checkbox"/>
<input type="checkbox"/>	1912124	24" Wheel size	
	<b>Rear Wheel Position</b> Indicate Rear Wheel positioning preference. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.		
<input type="checkbox"/>	2512401	Rear	<input type="checkbox"/>
<input type="checkbox"/>	2512402	Mid	<input type="checkbox"/>
<input type="checkbox"/>	2512403	Front	

### Axle Type

<input type="checkbox"/>	2513000	<b>Quick Release Rear Axles</b>
<input type="checkbox"/>	2513010	<b>Fixed Rear Axles</b>

### Rear Tires

<input type="checkbox"/>	2514000	<b>Urethane Tires</b>	
<input type="checkbox"/>	2514010	<b>Pneumatic Tires</b> May result in a finished Seat Height 1/2" greater than stated.	\$ 60
<input type="checkbox"/>	2514020	<b>Pneumatic Tires with Flat-Free Inserts</b> May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$ 95

### Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/>	2515000	<b>No Handrims</b>	
<input type="checkbox"/>	2515010	<b>Standard Aluminum Handrims</b>	
<input type="checkbox"/>	2515020	<b>Plastic Coated Handrims</b>	\$ 135

### Wheel Locks & Brakes

<b>User Wheel Locks</b>			
<input type="checkbox"/>	2516010	Push to Lock	<input type="checkbox"/>
<input type="checkbox"/>	2516020	Pull to Lock	
<input type="checkbox"/>	2516030	<b>Attendant Wheel Locks (Push to Lock)</b>	\$ 85
<input type="checkbox"/>	2516040	<b>User &amp; Attendant Wheel Locks</b> Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes.	\$ 125
<input type="checkbox"/>	2516300	<b>Drum Brakes</b> Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.	\$ 795

### Wheel Lock Options

<input type="checkbox"/>	2516100	<b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	\$ 60
<input type="checkbox"/>	2516200	<b>Hill Holders</b> Requires Pneumatic Rear Tires with Flat-Free Inserts.	\$ 450

### Anti-Tippers

<input checked="" type="checkbox"/>	2517010	<b>Rear Anti-Tippers</b> Required Safety Option.	\$ 75
-------------------------------------	---------	--	-------

### Back Posts

**Back Post Style**

2518010 Push Bar Back

---

**Back Post Height** Measured from the Seat Pan to the Handle.

2518425 Short (22")  2518431 Tall (25")

---

**Seat to Back Angle**

2519090 Set at 90°  2519098 Set at 98°  2519106 Set at 106°

2519114 Set at 114°  2519122 Set at 122°

### Back Options

**Angle/Height Adjustable Push Bar Handle** Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150

2519010 Button on Outside  1919011 Button on Inside

---

2519210 **10" Backrest Extension** Requires the Angle/Height Adjustable Push Bar Handle. \$ 150

---

2519212 **Headrest Mounting Plate** \$ 95  
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.

---

2519220 **Fold Down Back** Not compatible with Reclining Back.

### Back Upholstery

2520000 **No Back Upholstery**

---

**Standard Upholstery** Not available with Cane Handle Back. \$ 195

2520418 Short (18" Height)  2520421 Regular (21" Height)  2520424 Tall (24" Height)

---

2520450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. \$ 295

---

**Position of Upholstery** Select the desired position of the top of the Upholstery relative to the Seat Pan.

		22" Back Post								
Style		17"	18"	19"	20"	21"	22"	23"	24"	25"
Short		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

		25" Back Post								
Style		17"	18"	19"	20"	21"	22"	23"	24"	25"
Short		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall								<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Tilt Activation

**Attendant Tilt only** Select trigger mounting location below.

2523011 Right Back Post  2523012 Left Back Post

---

**User Tilt only** Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below. \$ 100

2523021 Right Armrest  2523022 Left Armrest

---

**Attendant & User Tilt** Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below. \$ 195

2523031 Right Back Post  2523032 Left Back Post

2523041 Right Armrest  2523042 Left Armrest

---

2523110 **Tilt Trigger Extension** \$ 50

---

2523400 **Power Tilt** 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. [pdgmobility.com/power-tilt](http://pdgmobility.com/power-tilt) \$ 1,795  
Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.

Toggle Switch  1923421 Installed on Right side  1923422 Installed on Left side

Button Switch  1923431 Installed on Right side  1923432 Installed on Left side

Stealth Egg Switch  1923440 Not installed

Optional for Egg Switch:  1923445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625  
*Specify Headrest to be used:* \_\_\_\_\_

---

**Tilt Limit** Optional. If no tilt limit is selected, the full tilt range will be available.

2523220 Set at 20°  2523230 Set at 30°  2523240 Set at 40°

### Recline

<input type="checkbox"/>	2523310	<b>Reclining Backrest (Manual)</b>		\$ 895
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.				

### Front Rigging

<input type="checkbox"/>	2525000	<b>No Front Rigging</b>	<input type="checkbox"/>	2525001	<b>Sector Blocks only</b>		
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>							
<input type="checkbox"/>	2525410	Fixed 70°				\$	95
<input type="checkbox"/>	2525420	Universal Angle Adjustable (0°-110°)	Set at angle: _____	(use increments of 10°)		\$	295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.							
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.							
<input type="checkbox"/>	2525130	Swing Away 60°	<input type="checkbox"/>	2525140	Swing Away 70°		
<input type="checkbox"/>	2525200	Swing Away Elevating Legrests	Includes Calf Pads.			\$	375
<input type="checkbox"/>	2525300	1" Offset Hangers for Swing Away				\$	75
<input type="checkbox"/>	2525110	<b>Standard Plug In 70° Front Rigging</b>	Legacy Design. Not compatible with Easy-Out.				


### Footplates

<input type="checkbox"/>	2526010	<b>Composite Footplates</b>					
<input type="checkbox"/>	2526020	<b>Aluminum Angle Adjustable Footplates</b>				\$	165
<input type="checkbox"/>	2526030	<b>One Piece Fixed Footplate</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.			\$	350
<input type="checkbox"/>	2526040	<b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.			\$	395
<input type="checkbox"/>	2526050	<b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.				

### Footrest/Legrest Accessories

<b>Calf Pads</b> Not compatible with Gel Padding.							
<input type="checkbox"/>	2527012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.			\$	200
<input type="checkbox"/>	2527010	For Swing Away Front Rigging				\$	135
<input type="checkbox"/>	2527020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.				
<b>Gel Knee Pads</b>							
<input type="checkbox"/>	2527032	For Easy-Out™ Angle Adjustable Front Rigging				\$	200
<input type="checkbox"/>	2527030	For Swing Away ELR Front Rigging				\$	100
<input type="checkbox"/>	2527040	<b>Length Adjustable Heel Loops</b>				\$	50
<input type="checkbox"/>	2527050	<b>Bodypoint Padded Calf Strap</b>				\$	195

### Front Rigging Extension Tubes

<b>Regular Extension Tubes</b>							
<input type="checkbox"/>	2528010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging			
<input type="checkbox"/>	2528020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging		\$	75
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.						
Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.							



# STELLAR GLT

## Order Form - Canada

Prices in Canadian Dollars. January 1, 2021

### Miscellaneous Options

<input type="checkbox"/>	2530080	<b>Gel Hip Guides</b>	Reduces available Seat Width by 2".	\$	395
<input type="checkbox"/>	2530110	<b>IV Pole</b>	Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later.	\$	295
<input type="checkbox"/>	2530120	<b>O2 Holder</b>	Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later.	\$	395
<input type="checkbox"/>	2530170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.	\$	10

### On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.