



### Customer and Order Information

Purchase Order *	Tag	v.0119
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

### Base Chair

**Stellar LEAP Manual Posterior & Anterior Tilt Wheelchair** Extened lead times may apply.

<input checked="" type="checkbox"/> 2201000	Stellar GL Manual Wheelchair with 20° Posterior Tilt	\$ 4,150
<input checked="" type="checkbox"/> 2201101	25° Anterior Tilt with Linked Back 30° Anterior Tilt available with reduced Posterior Tilt range of 15°.	\$ 100
<input checked="" type="checkbox"/> 2201102	30° Dynamic Back Recline	\$ 895

### Weight Capacity

<input checked="" type="checkbox"/> 2202000	<b>250 lbs Chair Weight Capacity</b>
---	--------------------------------------

### Frame Options

<input type="checkbox"/> 2202300	<b>Transit Tie Down System (TTDS)</b>
<small>TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.</small>	

### Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW



<b>Standard Seat Widths</b>					
<input type="checkbox"/> 2203014	14" Seat Width	<input type="checkbox"/> 2203015	15" Seat Width	<input type="checkbox"/> 2203016	16" Seat Width
<input type="checkbox"/> 2203017	17" Seat Width	<input type="checkbox"/> 2203018	18" Seat Width	<input type="checkbox"/> 2203019	19" Seat Width
<input type="checkbox"/> 2203020	20" Seat Width	<input type="checkbox"/> 2203021	21" Seat Width	<input type="checkbox"/> 2203022	22" Seat Width

### Seat Depth

<b>Standard Seat Depths (Adjustable)</b>					
<input type="checkbox"/> 2204016	16" Seat Depth	<input type="checkbox"/> 2204017	17" Seat Depth	<input type="checkbox"/> 2204018	18" Seat Depth
<input type="checkbox"/> 2204019	19" Seat Depth	<input type="checkbox"/> 2204020	20" Seat Depth		
<b>Custom Seat Depths</b>					
<input type="checkbox"/> 2204021	21" Seat Depth	<input type="checkbox"/> 2204022	22" Seat Depth	\$ 400	

### Frame Color

<b>PDG Premium Colors</b>					
<input type="checkbox"/> 2205000	Tracer Red	<input type="checkbox"/> 2205001	Gloss Black Shimmer	<input type="checkbox"/> 2205002	Silver
<input type="checkbox"/> 2205003	Fine Texture Black	<input type="checkbox"/> 2205004	Orange Fusion	<input type="checkbox"/> 2205005	Karl's Imagination
<input type="checkbox"/> 2205006	Blue Shimmer	<input type="checkbox"/> 2205007	Gloss White	<input type="checkbox"/> 2205008	Dark Titanium
<input type="checkbox"/> 2205009	Geyser Green	<input type="checkbox"/> 2205010	Lollypop Blue	<input type="checkbox"/> 2205011	Lollypop Red
<input type="checkbox"/> 2205012	Bubble Gum Pink	<input type="checkbox"/> 2205013	Berry Purple	<input type="checkbox"/> 2205014	Daffodil Yellow
<input type="checkbox"/> 2205200	<b>Custom Color</b>	<small>Contact PDG for specifications.</small>			\$ 300

Armrests		
<input type="checkbox"/> 2206000	<b>No Armrests</b>	
<b>Single Post Height Adjustable Armrests</b> <span style="float: right;">\$ 295</span>		
<input type="checkbox"/> 2206111 8"-11" Range	<input type="checkbox"/> 2206114 10"-14" Range	
<b>Dual Post Height Adjustable Armrests</b> <span style="float: right;">\$ 295</span>		
<input type="checkbox"/> 2206211 8"-11" Range	<input type="checkbox"/> 2206214 10"-14" Range	<input type="checkbox"/> 2206216 12"-16" Range
<b>Fixed Height Armrests (12")</b> Full Length Armpad only.		
<input type="checkbox"/> 2206100 Single Post	<input type="checkbox"/> 2206200 Dual Post	
 <b>Armrest Height</b> If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used. Seat Pan to top of Armrest: _____ Does not include height of Armpad		
 <b>Armrest Position</b> Indicate preference for Armrest position. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.		
<input type="checkbox"/> 2206401 Rear	<input type="checkbox"/> 2206402 Mid	<input type="checkbox"/> 2206403 Front
<b>Locking Flip-up</b> Not recommended for users with excessive strength or tone. <span style="float: right;">\$ 295</span>		
Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.		
<input type="checkbox"/> 2206310 Set at 10"	<input type="checkbox"/> 2206311 Set at 11"	<input type="checkbox"/> 2206312 Set at 12"
<input type="checkbox"/> 2206313 Set at 13"	<input type="checkbox"/> 2206314 Set at 14"	

Armpads	
<b>Regular Armpads</b>	
<input type="checkbox"/> 2207100 Desk Length - 2" x 10"	<input type="checkbox"/> 2207101 Full Length - 2" x 14"
<b>Waterfall Armpads</b> <span style="float: right;">\$ 115</span>	
<input type="checkbox"/> 2207105 Desk Length - 2" x 10"	<input type="checkbox"/> 2207106 Full Length - 2" x 14"
<b>Wide Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 200</span>	
<input type="checkbox"/> 2207110 Desk Length - 3" x 10"	<input type="checkbox"/> 2207111 Full Length - 3" x 12"
<b>Gel Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 135</span>	
<input type="checkbox"/> 2207120 Full Length - 2" x 12"	<input type="checkbox"/> 2207121 Full Length - 3.5" x 12"
<b>Contour Arm Trough with Hand Pad</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 125</span>	
<input type="checkbox"/> 2207130 Right (4.5" x 18")	<input type="checkbox"/> 2207131 Left (4.5" x 18")

Armrest Options	
<b>Armrest Wideners</b> Move Armrests outward by 1" on each side. Do not increase seat size.	
<input type="checkbox"/> 2207210 For Single Post Armrests	May cause tilt interference at lower Seat Height and/or Seat Depth settings. <span style="float: right;">\$ 195</span>
<input type="checkbox"/> 2207220 For Dual Post Armrests	Limits tilt. Not compatible with 20", 22", or 24" Rear Wheels, Clothing Guards. <span style="float: right;">\$ 375</span>
<b>Armrest Extensions</b> Not compatible with Flip-up or Single Post Armrests. <span style="float: right;">\$ 75</span>	
<input type="checkbox"/> 2207400 3" Extension - Forward	<input type="checkbox"/> 2207401 3" Extension - Rearward

Seat Pan	
<input type="checkbox"/> 2208000	<b>No Seat Pan</b>
<input type="checkbox"/> 2208001	<b>Solid Aluminum Seat Pan</b> <span style="float: right;">\$ 135</span>

## Seat Height



**Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

## Casters

**Urethane Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.

- |   |   |        |
|---|---|--------|
| <input type="checkbox"/> 2211061 6" x 1.25" Caster size | <input type="checkbox"/> 2211081 8" x 1.25" Caster size |        |
| <input type="checkbox"/> 2211063 6" x 1.50" Caster size |   | \$ 85  |
| <input type="checkbox"/> 2211084 8" x 2.00" Caster size |   | \$ 135 |

**Pneumatic Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.

- |   |  |        |
|---|--|--------|
| <input type="checkbox"/> 2211161 6" x 1.25" Caster size |  | \$ 110 |
| <input type="checkbox"/> 2211181 8" x 1.25" Caster size |  | \$ 120 |
| <input type="checkbox"/> 2211184 8" x 2.00" Caster size |  | \$ 150 |

- |   |  |        |
|---|--|--------|
| <input type="checkbox"/> 2211260 <b>Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters</b> |  | \$ 275 |
|---|--|--------|

- |   |  |        |
|---|--|--------|
| <input type="checkbox"/> 2211400 <b>Frog Legs Suspension Forks</b> Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster. |  | \$ 595 |
|---|--|--------|

## Rear Wheels

**Composite Wheels**

- |   |   |   |
|---|---|---|
| <input type="checkbox"/> 2212012 12" Wheel size | <input type="checkbox"/> 2212016 16" Wheel size | <input type="checkbox"/> 2212020 20" Wheel size |
| <input type="checkbox"/> 2212022 22" Wheel size | <input type="checkbox"/> 2212024 24" Wheel size |   |

**Spoke Wheels**

- |   |   |   |       |
|---|---|---|-------|
| <input type="checkbox"/> 2212120 20" Wheel size | <input type="checkbox"/> 2212122 22" Wheel size | <input type="checkbox"/> 2212124 24" Wheel size | \$ 75 |
|---|---|---|-------|

- |   |  |        |
|---|--|--------|
| <input type="checkbox"/> 2212224 <b>24" Heavy Duty Spoke Wheels</b> |  | \$ 350 |
|---|--|--------|



**Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- |                                       |                                      |  |
|---------------------------------------|--------------------------------------|--|
| <input type="checkbox"/> 2212401 Rear | <input type="checkbox"/> 2212402 Mid | <input type="checkbox"/> 2212403 Front |
|---------------------------------------|--------------------------------------|--|

## Axle Type

- |  |
|--|
| <input type="checkbox"/> 2213000 <b>Quick Release Rear Axles</b> |
|--|

- |  |
|--|
| <input type="checkbox"/> 2213010 <b>Fixed Rear Axles</b> |
|--|

## Rear Tires

- |  |
|--|
| <input type="checkbox"/> 2214000 <b>Urethane Tires</b> |
|--|

- |  |       |
|--|-------|
| <input type="checkbox"/> 2214010 <b>Pneumatic Tires</b> May result in a finished Seat Height 1/2" greater than stated. | \$ 60 |
|--|-------|

- |  |       |
|--|-------|
| <input type="checkbox"/> 2214020 <b>Pneumatic Tires with Flat-Free Inserts</b> May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. | \$ 95 |
|--|-------|

## Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

- |   |
|---|
| <input type="checkbox"/> 2215000 <b>No Handrims</b> |
|---|

- |   |
|---|
| <input type="checkbox"/> 2215010 <b>Standard Handrims</b> |
|---|

- |   |        |
|---|--------|
| <input type="checkbox"/> 2215020 <b>Plastic Coated Handrims</b> | \$ 150 |
|---|--------|



## Wheel Locks & Brakes

<b>User Wheel Locks</b>		
<input type="checkbox"/> 2216010	Push to Lock	<input type="checkbox"/> 2216020 Pull to Lock
<input type="checkbox"/> 2216030	<b>Attendant Wheel Locks (Push to Lock)</b>	\$ 85
<input type="checkbox"/> 2216040	<b>User &amp; Attendant Wheel Locks</b> Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes.	\$ 125
<input type="checkbox"/> 2216300	<b>Drum Brakes</b> Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.	\$ 795

## Wheel Lock Options

<input type="checkbox"/> 2216100	<b>Wheel Lock Extensions</b> Not compatible with Attendant Wheel Lock.	\$ 60
<input type="checkbox"/> 2216200	<b>Hill Holders</b> Requires Pneumatic Rear Tires with Flat-Free Inserts.	\$ 450

## Anti-Tippers

<input checked="" type="checkbox"/> 2217010	<b>Rear Anti-Tippers</b> Required Safety Option.	\$ 75
<input type="checkbox"/> 2217020	<b>Front Anti-Tippers</b>	\$ 795
<input type="checkbox"/> 2217030	<b>Heavy Duty Rear Anti-Tipper</b> Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.	\$ 895

## Back Posts

<b>Back Post Style</b>		
<input type="checkbox"/> 2218010	Push Bar Back	<input type="checkbox"/> 2218020 Cane Handle Back Not compatible with Headrest Mounting Plate.
<b>Back Post Height</b> Measured from the Seat Pan to the Handle.		
<input type="checkbox"/> 2218420	Short (17")	<input type="checkbox"/> 1518425 Tall (23")
<b>Seat to Back Angle</b>		
<input type="checkbox"/> 2219090	Set at 90°	<input type="checkbox"/> 2219098 Set at 98° <input type="checkbox"/> 2219106 Set at 106°
<input type="checkbox"/> 2219114	Set at 114°	<input type="checkbox"/> 2219122 Set at 122°

## Back Options

<b>Angle/Height Adjustable Push Bar Handle</b> Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/> 2219010	Button on Outside	<input type="checkbox"/> 2219011 Button on Inside
<input type="checkbox"/> 2219210	<b>10" Backrest Extension</b> Not compatible with Angle/Height Adjustable Push Bar Handle.	\$ 150
<input type="checkbox"/> 2219212	<b>Headrest Mounting Plate</b> Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 95
<input checked="" type="checkbox"/> 2219220	<b>Fold Down Back</b>	

## Back Upholstery

<input type="checkbox"/> 2220000	<b>No Back Upholstery</b>																		
<b>Standard Upholstery</b> Not available with Cane Handle Back.		\$ 195																	
<input type="checkbox"/> 2220418	Short (18" Height)	<input type="checkbox"/> 2220421 Regular (21" Height)																	
<input type="checkbox"/> 2220450	<b>Tension Adjustable Upholstery</b> Height of Upholstery is 16". Not available with Cane Handle Back.	\$ 295																	
<b>Position of Upholstery</b> Select the desired position of the top of the Upholstery relative to the Seat Pan.																			
17" Back Post										23" Back Post									
<b>Style</b>	<b>17"</b>	<b>18"</b>	<b>19"</b>	<b>20"</b>	<b>21"</b>	<b>22"</b>	<b>23"</b>	<b>24"</b>	<b>25"</b>	<b>Style</b>	<b>17"</b>	<b>18"</b>	<b>19"</b>	<b>20"</b>	<b>21"</b>	<b>22"</b>	<b>23"</b>	<b>24"</b>	<b>25"</b>
Short	<input type="checkbox"/>									Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Regular										Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Tension Adj										Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

### Tilt Activation

<b>Attendant Tilt only</b> Select trigger mounting location below.		
<input type="checkbox"/> 2223011 Right Back Post	<input type="checkbox"/> 2223012 Left Back Post	
<b>User Tilt only</b> Full Length Arm Pads only. Select trigger mounting location below.		\$ 100
<input type="checkbox"/> 2223021 Right Armrest	<input type="checkbox"/> 2223022 Left Armrest	
<b>Attendant &amp; User Tilt</b> Full Length Arm Pads only. Select trigger mounting locations below.		\$ 195
<input type="checkbox"/> 2223031 Right Back Post	<input type="checkbox"/> 2223032 Left Back Post	
<input type="checkbox"/> 2223041 Right Armrest	<input type="checkbox"/> 2223042 Left Armrest	
<input type="checkbox"/> 2223400 <b>Power Tilt</b> 250 lbs Weight Capacity. 30° Anterior to 20° Posterior tilt range. Includes 2 batteries for off-chair charging. Minimum Seat Width is 16".		\$ 1,795
Toggle Switch	<input type="checkbox"/> 2223421 Installed on Right side	<input type="checkbox"/> 2223422 Installed on Left side
Button Switch	<input type="checkbox"/> 2223431 Installed on Right side	<input type="checkbox"/> 2223432 Installed on Left side
Stealth Egg Switch	<input type="checkbox"/> 2223440 Not installed	
Optional for Egg Switch:	<input type="checkbox"/> 2223445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models. \$ 625
Specify Headrest to be used: _____		
<b>Anterior Tilt Limit</b> Total tilt range is 45°. Posterior Tilt equals 45° less the value selected below.		
<input type="checkbox"/> 2223220 Set at 20°	<input type="checkbox"/> 2223225 Set at 25°	<input type="checkbox"/> 2223230 Set at 30°

### Front Rigging

<input type="checkbox"/> 2225000 <b>No Front Rigging</b>		
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>		
<input type="checkbox"/> 2225410 Fixed 70°		\$ 95
<input type="checkbox"/> 2225420 Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.		

### Footplates

<input type="checkbox"/> 2226010 <b>Composite Footplates</b>		
<input type="checkbox"/> 2226020 <b>Aluminum Angle Adjustable Footplates</b>		\$ 165
<input type="checkbox"/> 2226030 <b>One Piece Fixed Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/> 2226040 <b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/> 2226050 <b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.	\$ 195

### Footrest/Legrest Accessories

<input type="checkbox"/> 2227040 <b>Length Adjustable Heel Loops</b>		\$ 50
<input type="checkbox"/> 2227050 <b>Bodypoint Padded Calf Strap</b>		\$ 195

### Front Rigging Extension Tubes

<b>Regular Extension Tubes</b>		
<input type="checkbox"/> 2228010 Short 14"-18" for Easy-Out Front Rigging		
<input type="checkbox"/> 2228020 Long 18"-24" for Easy-Out Front Rigging		\$ 75
<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.		
Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.		

### Miscellaneous Options

<b>Spoke Guards</b> Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.		\$ 195
<input type="checkbox"/> 2230020 Black	<input type="checkbox"/> 2230021 Frame Color	
<input type="checkbox"/> 2230110 <b>IV Pole</b> Not compatible with HD Anti Tipper, Attendant Wheel Locks.		\$ 295
<input type="checkbox"/> 2230120 <b>O2 Holder</b> Not compatible with HD Anti Tipper, Attendant Wheel Locks.		\$ 395
<input type="checkbox"/> 2230170 <b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.		\$ 10



# STELLAR LEAP

**Order Form - Canada**

Prices in Canadian Dollars. January 1, 2022

## On-Chair Seating & Positioning



Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

## Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.