



Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

Base Chair

<input checked="" type="checkbox"/> 1001000	Bentley 20° Manual Tilt Wheelchair	\$ 4,200
i	Also available: the lightweight Bentley LT with aluminum frame. See separate order form. pdgmobility.com/blt	

Weight Capacity

Chair Weight Capacity		
<input type="checkbox"/> 1002000	250 lbs Capacity Available for 14" to 22" Seat Width.	
<input type="checkbox"/> 1002050	300 lbs Capacity Available for 14" to 22" Seat Width.	\$ 200
<input type="checkbox"/> 1005006	450 lbs Capacity Available for 20" to 32" Seat Width.	\$ 500
User Weight Range Required for setup of tilt uplift force.		
<input type="checkbox"/> 1002430	Up to 75 lbs	
<input type="checkbox"/> 1002431	75 lbs to 100 lbs	
<input type="checkbox"/> 1002432	100 lbs to 115 lbs	
<input type="checkbox"/> 1002433	115 lbs to 150 lbs	
<input type="checkbox"/> 1002434	150 lbs to 200 lbs	
<input type="checkbox"/> 1002435	200 lbs to 250 lbs	
<input type="checkbox"/> 1002436	250 lbs to 300 lbs	
<input type="checkbox"/> 1002437	300 lbs to 375 lbs	
<input type="checkbox"/> 1002438	375 lbs to 450 lbs	

Frame Options

<input type="checkbox"/> 1002200	Extended Base with Rear Floor Jacks	\$ 900
Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.		
<input type="checkbox"/> 1002300	Transit Tie Down System (TTDS)	
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Seat Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.		

Seat Width Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths		
<input type="checkbox"/> 1003014	14" Seat Width	
<input type="checkbox"/> 1003015	15" Seat Width	
<input type="checkbox"/> 1003016	16" Seat Width	
<input type="checkbox"/> 1003017	17" Seat Width	
<input type="checkbox"/> 1003018	18" Seat Width	
<input type="checkbox"/> 1003019	19" Seat Width	
<input type="checkbox"/> 1003020	20" Seat Width	
<input type="checkbox"/> 1003021	21" Seat Width	
<input type="checkbox"/> 1003022	22" Seat Width	
Custom Seat Widths Requires 450 lbs Weight Capacity Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1003023	23" Seat Width	\$ 375
<input type="checkbox"/> 1003024	24" Seat Width	\$ 375
<input type="checkbox"/> 1003025	25" Seat Width	\$ 700
<input type="checkbox"/> 1003026	26" Seat Width	\$ 700
<input type="checkbox"/> 1003027	27" Seat Width	\$ 700
<input type="checkbox"/> 1003028	28" Seat Width	\$ 700
<input type="checkbox"/> 1003029	29" Seat Width	\$ 700
<input type="checkbox"/> 1003030	30" Seat Width	\$ 700
<input type="checkbox"/> 1003031	31" Seat Width	\$ 700
<input type="checkbox"/> 1003032	32" Seat Width	\$ 700

Seat Depth



Standard Seat Depths (Adjustable)		
<input type="checkbox"/> 1004016	16" Seat Depth	
<input type="checkbox"/> 1004017	17" Seat Depth	
<input type="checkbox"/> 1004018	18" Seat Depth	
<input type="checkbox"/> 1004019	19" Seat Depth	
<input type="checkbox"/> 1004020	20" Seat Depth	
Custom Seat Depths		
<input type="checkbox"/> 1004021	21" Seat Depth	\$ 400
<input type="checkbox"/> 1004022	22" Seat Depth	\$ 400

Frame Color

PDG Premium Colors

<input type="checkbox"/> 1005000 Tracer Red	<input type="checkbox"/> 1005001 Gloss Black Shimmer	<input type="checkbox"/> 1005002 Silver
<input type="checkbox"/> 1005003 Fine Texture Black	<input type="checkbox"/> 1005004 Orange Fusion	<input type="checkbox"/> 1005005 Karl's Imagination
<input type="checkbox"/> 1005006 Blue Shimmer	<input type="checkbox"/> 1005007 Gloss White	<input type="checkbox"/> 1005008 Dark Titanium
<input type="checkbox"/> 1005009 Geyser Green	<input type="checkbox"/> 1005010 Lollipop Blue	<input type="checkbox"/> 1005011 Lollipop Red
<input type="checkbox"/> 1005012 Bubble Gum Pink	<input type="checkbox"/> 1005013 Berry Purple	<input type="checkbox"/> 1005014 Daffodil Yellow
<input type="checkbox"/> 1005200 Custom Color Contact PDG for specifications. \$ 300		

Armrests

<input type="checkbox"/> 1006000 No Armrests			\$ 295
Single Post Height Adjustable Armrests			
<input type="checkbox"/> 1006111 8"-11" Range	<input type="checkbox"/> 1006114 10"-14" Range		
Dual Post Height Adjustable Armrests			
<input type="checkbox"/> 1006211 8"-11" Range	<input type="checkbox"/> 1006214 10"-14" Range	<input type="checkbox"/> 1006216 12"-16" Range	\$ 295
Fixed Height Armrests (12") Full Length Armpad only.			
<input type="checkbox"/> 1006100 Single Post	<input type="checkbox"/> 1006200 Dual Post		
 Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used. Seat Pan to top of Armrest: _____ Does not include height of Armpad			
 Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 1006401 Rear	<input type="checkbox"/> 1006402 Mid	<input type="checkbox"/> 1006403 Front	
Locking Flip-up Armrests Not recommended for users with excessive strength or tone. \$ 295 Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.			
<input type="checkbox"/> 1006310 Set at 10" Height	<input type="checkbox"/> 1006311 Set at 11" Height	<input type="checkbox"/> 1006312 Set at 12" Height	
<input type="checkbox"/> 1006313 Set at 13" Height	<input type="checkbox"/> 1006314 Set at 14" Height		

Armpads

Regular Armpads			
<input type="checkbox"/> 1007100 Desk Length - 2" x 10"	<input type="checkbox"/> 1007101 Full Length - 2" x 14"		
Waterfall Armpads			
<input type="checkbox"/> 1007105 Desk Length - 2" x 10"	<input type="checkbox"/> 1007106 Full Length - 2" x 14"		
Wide Armpads May interfere with wheels in some tilt configurations. \$ 200			
<input type="checkbox"/> 1007110 Desk Length - 3" x 10"	<input type="checkbox"/> 1007111 Full Length - 3" x 12"		
Gel Armpads May interfere with wheels in some tilt configurations. \$ 135			
<input type="checkbox"/> 1007120 Full Length - 2" x 12"	<input type="checkbox"/> 1007121 Full Length - 3.5" x 12"		
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125			
<input type="checkbox"/> 1007130 Right (4.5" x 18")	<input type="checkbox"/> 1007131 Left (4.5" x 18")		

Armrest Options

Armrest Wideners Move Armrests outward by 1" on each side. Do not increase seat size.			
<input type="checkbox"/> 1007210 For Single Post Armrests	May cause tilt interference at lower Seat Height and/or Seat Depth settings.		
<input type="checkbox"/> 1007220 For Dual Post Armrests	Limits tilt. Not compatible with 20", 22", or 24" Rear Wheels, Clothing Guards.		
Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.			
<input type="checkbox"/> 1007310 For Single Post Armrests	Available for 10"-14" Height range only.		
<input type="checkbox"/> 1007320 For Dual Post Armrests			
Armrest Extensions Not compatible with Flip-up or Single Post Armrests. \$ 75			
<input type="checkbox"/> 1007400 3" Extension - Forward	<input type="checkbox"/> 1007401 3" Extension - Rearward		

Seat Pan

- 1008000 **No Seat Pan**
- 1008001 **Solid Aluminum Seat Pan** Required for 22" Seat Depth. \$ 135

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

	Low						High						
	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
6" Casters	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	22"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8" Casters	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

- Urethane Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1011061 6" x 1.25" Caster size 1011081 8" x 1.25" Caster size
 - 1011063 6" x 1.50" Caster size \$ 85
 - 1011084 8" x 2.00" Caster size \$ 135
- Pneumatic Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1011161 6" x 1.25" Caster size \$ 110
 - 1011181 8" x 1.25" Caster size \$ 120
 - 1011184 8" x 2.00" Caster size \$ 150
 - 1011260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275
 - 1011400 **Frog Legs Suspension Forks** Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster. \$ 695

Rear Wheels

- Composite Wheels**
- 1012012 12" Wheel size 1012016 16" Wheel size 1012020 20" Wheel size
 - 1012022 22" Wheel size 1012024 24" Wheel size
- Spoke Wheels** \$ 75
- 1012120 20" Wheel size 1012122 22" Wheel size 1012124 24" Wheel size
 - 1012224 **24" Heavy Duty Spoke Wheels** \$ 350
- Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.
- 1012401 Rear 1012402 Mid 1012403 Front

Axle Type

- 1013000 **Quick Release Rear Axles**
- 1013010 **Fixed Rear Axles**

Rear Tires

- 1014000 **Urethane Tires**
- 1014010 **Pneumatic Tires** May result in a finished Seat Height 1/2" greater than stated. \$ 60
- 1014020 **Pneumatic Tires with Flat-Free Inserts** May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. \$ 95

Handrims		Handrims are only available with 20", 22", and 24" Rear Wheels.
<input type="checkbox"/>	1015000 No Handrims	
<input type="checkbox"/>	1015010 Standard Aluminum Handrims	
<input type="checkbox"/>	1015020 Plastic Coated Handrims	\$ 150

Wheel Locks & Brakes				
User Wheel Locks				
<input type="checkbox"/>	1016010 Push to Lock	<input type="checkbox"/>	1016020 Pull to Lock	
<input type="checkbox"/>	1016030 Attendant Wheel Locks (Push to Lock)			\$ 85
<input type="checkbox"/>	1016040 User & Attendant Wheel Locks	Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes.		\$ 125
<input type="checkbox"/>	1016300 Drum Brakes	Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.		\$ 995

Wheel Lock Options			
<input type="checkbox"/>	1016100 Wheel Lock Extensions	Not compatible with Attendant Wheel Lock.	\$ 60
<input type="checkbox"/>	1016200 Hill Holders	Requires Pneumatic Rear Tires with Flat-Free Inserts.	\$ 450

Anti-Tippers			
<input checked="" type="checkbox"/>	1017010 Rear Anti-Tippers	Required Safety Option.	\$ 75
<input type="checkbox"/>	1017020 Front Anti-Tippers		\$ 795
<input type="checkbox"/>	1017030 Heavy Duty Rear Anti-Tipper	Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder.	\$ 895

Back Posts					
Back Post Style					
<input type="checkbox"/>	1018010 Push Bar Back	<input type="checkbox"/>	1018020 Cane Handle Back	Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.	
Back Post Height		Measured from the Seat Pan to the Handle.			
<input type="checkbox"/>	1018425 Short (22")	<input type="checkbox"/>	1018431 Tall (25")		
Seat to Back Angle					
<input type="checkbox"/>	1019090 Set at 90°	<input type="checkbox"/>	1019098 Set at 98°	<input type="checkbox"/>	1019106 Set at 106°
<input type="checkbox"/>	1019114 Set at 114°	<input type="checkbox"/>	1019122 Set at 122°		

Back Options				
Angle/Height Adjustable Push Bar Handle		Not compatible with Headrest Mounting Plate or Cane Handle Back.	\$ 150	
<input type="checkbox"/>	1019010 Button on Outside	<input type="checkbox"/>	1019011 Button on Inside	
<input type="checkbox"/>	1019210 10" Backrest Extension	Requires the Angle/Height Adjustable Push Bar Handle.	\$ 195	
<input type="checkbox"/>	1019212 Headrest Mounting Plate	Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 95	
<input type="checkbox"/>	1019220 Fold Down Back	Not compatible with Reclining Back.		
<input checked="" type="checkbox"/>	1019230 Spring Back Suspension	Not compatible with Reclining Back.		
<input type="checkbox"/>	1019240 Seating Dynamics Dynamic Rocker Back Interface	Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle, Standard Upholstery. Finished Back Post Height is 4" taller than selected Back Post Height.	\$ 2,995	

Back Upholstery					
<input type="checkbox"/>	1020000 No Back Upholstery				
Standard Upholstery		Not available with Cane Handle Back.	\$ 265		
<input type="checkbox"/>	1020418 Short (18" Height)	<input type="checkbox"/>	1020421 Regular (21" Height)	<input type="checkbox"/>	1020424 Tall (24" Height)

See next page for more options >



Tension Adjustable Upholstery

Height of Upholstery is 16". Not available with Cane Handle Back.

- 1020450 For 14" to 22" Seat Width \$ 365
- 1020451 For 24" to 32" Seath Width \$ 525

Position of Upholstery

Select the desired position of the top of the Upholstery relative to the Seat Pan.

Style	22" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Style	25" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation

Attendant Tilt only

Select trigger mounting location below. Chairs with 450 lbs Weight Capacity require two tilt triggers (L and R) for tilt activation.

- 1023011 Right Back Post 1023012 Left Back Post

User Tilt only

Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below.

- 1023021 Right Armrest 1023022 Left Armrest \$ 100

Attendant & User Tilt

Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below.

- 1023031 Right Back Post 1023032 Left Back Post
- 1023041 Right Armrest 1023042 Left Armrest \$ 195

- 1023110 **Tilt Trigger Extension** \$ 50

- 1023400 **Power Tilt** 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt \$ 2,095
Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.

- Toggle Switch 1023421 Installed on Right side 1023422 Installed on Left side
- Button Switch 1023431 Installed on Right side 1023432 Installed on Left side
- Stealth Egg Switch 1023440 Not installed

- Optional for Egg Switch: 1023445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625

Specify Headrest to be used: _____

Recline

- 1023310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.

Front Rigging

- 1025000 **No Front Rigging** 1025001 **Sector Blocks only**

Easy-Out™ Front Rigging

Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out

- 1025410 Fixed 70° \$ 95
- 1025420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) \$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.

Swing Away Front Rigging

Legacy Design. Not compatible with Easy-Out.

- 1025130 Swing Away 60° 1025140 Swing Away 70°
- 1025200 Swing Away Elevating Legrests Includes Calf Pads. \$ 375
- 1025300 1" Offset Hangers for Swing Away \$ 75

- 1025110 **Standard Plug In 70° Front Rigging** Legacy Design. Not compatible with Easy-Out.


Footplates

- 1026010 **Composite Footplates**
- 1026020 **Aluminum Angle Adjustable Footplates** \$ 165
- 1026030 **One Piece Fixed Footplate** Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable. \$ 350
- 1026040 **One Piece Flip-up Footplate** Requires selection of Composite or Aluminum Angle Adjustable. \$ 395
- 1026050 **Interlocking Angle Adjustable Footplates** Available for Seat Widths up to 21" only. \$ 195

Footrest/Legrest Accessories

Calf Pads Not compatible with Gel Padding.			
<input type="checkbox"/>	1027012	For Easy-Out™ Front Rigging Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	1027010	For Swing Away Front Rigging	\$ 135
<input type="checkbox"/>	1027020	Gel Padding Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pads			
<input type="checkbox"/>	1027032	For Easy-Out™ Angle Adjustable Front Rigging	\$ 200
<input type="checkbox"/>	1027030	For Swing Away ELR Front Rigging	\$ 100
<input type="checkbox"/>	1027040	Length Adjustable Heel Loops	\$ 50
<input type="checkbox"/>	1027050	Bodypoint Padded Calf Strap Sizes greater than 22" are not padded.	\$ 195
Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.			
<input type="checkbox"/>	1027060	For Standard Plug In Front Rigging	\$ 445
<input type="checkbox"/>	1027065	For Swing Away 70° Front Rigging Not available on 13" STF Height.	\$ 575


Front Rigging Extension Tubes

Regular Extension Tubes			
<input type="checkbox"/>	1028010	Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	1028020	Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only. Not compatible with Swing Away or Easy-Out Front Rigging.			\$ 100
<input type="checkbox"/>	1028110	Standard (13"-17")	<input type="checkbox"/>
		1028120	Long (15"-19")
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.		
	Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.		

Miscellaneous Options

One-Arm Drive Maximum Seat Width of 22". 24" Composite Wheels with Urethane Tires required. 18" STF minimum. Limits tilt. Not compatible with User & Attendant Wheel Locks, Drum Brakes, Quick Release Rear Axles.			\$ 3,500
<input type="checkbox"/>	1030010	Installed on Right side	<input type="checkbox"/>
		1030011	Installed on Left side
Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.			\$ 235
<input type="checkbox"/>	1030020	Black	<input type="checkbox"/>
		1030021	Frame Color
<input type="checkbox"/>	1030050	Protective Lower Padding	\$ 995
<input type="checkbox"/>	1030080	Gel Hip Guides Reduces available Seat Width by 2".	\$ 495
<input type="checkbox"/>	1030110	IV Pole Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$ 345
<input type="checkbox"/>	1030120	O2 Holder Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	\$ 495
<input type="checkbox"/>	1030170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$ 10

On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.