

Cushion

1. Select the type
↓

2. Select the model
↓

General Use	Stealth Simplicity 350 lbs weight capacity.	
	<input type="checkbox"/> ##09141 16"-20" width, 16"-20" depth	\$ 390
	<input type="checkbox"/> ##09142 22"-24" width, 18"-20" depth	\$ 505

Positioning	ROHO AirLITE 300 lbs weight capacity.	
	<input type="checkbox"/> ##09120 14"-20" width, 14"-18" depth	\$ 415

Skin Protection & Positioning	Stealth Spectrum Foam 650 lbs weight capacity.	
	<input type="checkbox"/> ##09143 16"-20" width, 16"-24" depth	\$ 1,140
	<input type="checkbox"/> ##09144 22"-32" width, 18"-28" depth	\$ 1,420
	Stealth Tru-Comfort 2 450 lbs weight capacity.	
	<input type="checkbox"/> ##09146 16"-20" width, 16"-24" depth	\$ 1,010
	<input type="checkbox"/> ##09147 22" width, 18"-24" depth	\$ 1,225
	nxt BioFit 300 lbs weight capacity.	
	<input type="checkbox"/> ##09101 14"-20" width, 16"-20" depth	\$ 655
	<input type="checkbox"/> ##09102 22"-24" width, 16"-20" depth	\$ 690
	Matrx VI 300 lbs weight capacity.	
	<input type="checkbox"/> ##09111 14"-20" width, 15"-20" depth	\$ 750
	<input type="checkbox"/> ##09112 22" width, 18"-20" depth	\$ 820
Matrx VI Heavy Duty 600 lbs weight capacity.		
<input type="checkbox"/> ##09115 18"-21" width, 18"-26" depth	\$ 1,155	
<input type="checkbox"/> ##09116 22"-30" width, 18"-26" depth	\$ 1,375	

Skin Protection & Adjustable	ROHO Single Valve High Profile No weight limit if cushion is properly sized.	
	<input type="checkbox"/> ##09125 14"-20" width, 16"-22" depth	\$ 925
	<input type="checkbox"/> ##09126 22"-24" width, 16"-22" depth	\$ 1,055
	STAR Stabil-Air Single Valve No weight limit if cushion is properly sized.	
	<input type="checkbox"/> ##09151 16"-20" width, 16"-20" depth Standard (even) sizes.	\$ 1,320
<input type="checkbox"/> ##09152 16"-20" width, 16"-20" depth Custom (odd and larger) sizes.	\$ 1,430	

Skin Protection, Positioning & Adjustable	ROHO Quadro Select High Profile No weight limit if cushion is properly sized.	
	<input type="checkbox"/> ##09121 14"-20" width, 16"-22" depth	\$ 925
	<input type="checkbox"/> ##09122 22"-24" width, 16"-22" depth	\$ 1,055

Extra-Wide No weight limit if cushion is properly sized. Custom sizes available. Custom-made foam cushion with fluid-proof yet breathable fabric cover.	PDG Solar Basic High resiliency polyfoam with 2" medium-density base and 1" lower-density top.	<input type="checkbox"/> ##09200 20"-28" width, 16"-24" depth	\$ 435
		<input type="checkbox"/> ##09201 29"-36" width, 16"-24" depth	\$ 545
		<input type="checkbox"/> ##09202 37"-48" width, 16"-24" depth	\$ 655
	PDG Solar Contour Aids proper pelvic alignment, includes pre-ischial bar, leg shelf, adductor buildups.	<input type="checkbox"/> ##09210 20"-28" width, 16"-24" depth	\$ 765
		<input type="checkbox"/> ##09211 29"-36" width, 16"-24" depth	\$ 875
		<input type="checkbox"/> ##09212 37"-48" width, 16"-24" depth	\$ 985
	PDG Solar Pressure Reduction All features of the Contour, plus added pincore latex rubber and containment foam outline.	<input type="checkbox"/> ##09220 20"-28" width, 16"-24" depth	\$ 1,205
		<input type="checkbox"/> ##09221 29"-36" width, 16"-24" depth	\$ 1,315
		<input type="checkbox"/> ##09222 37"-48" width, 16"-24" depth	\$ 1,425

3. Specify the size →

Cushion Size Sizes based on availability from the cushion manufacturer.

Specify Size (Width" x Depth") _____

Backrest

1. Select the type



2. Select the contour depth



3. Select the model



Positioning, posterior-lateral	2" depth	<input type="checkbox"/> #21141 Stealth ADI AL Active 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,365
	4" depth	<input type="checkbox"/> #21142 Stealth ADI AL Deep 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,585
	6" depth	<input type="checkbox"/> #21143 Stealth ADI AL Extra Deep 250 lbs weight capacity. 14"-20" width, 10"-20" height \$ 1,755

Positioning, planar back with lateral supports	2" to 3" depth	<input type="checkbox"/> #21110 Matrx Elite E2 300 lbs weight capacity. 12"-20" width, 10"-20" height \$ 1,080
		<input type="checkbox"/> #21111 Matrx Elite E2 Heavy Duty 500 lbs weight capacity. 21" width, 16"-20" height \$ 1,345
		<input type="checkbox"/> #21112 22"-24" width, 16"-20" height \$ 1,455
		<input type="checkbox"/> #21114 25"-30" width, 16"-20" height \$ 2,195
	3.5" to 4" depth	<input type="checkbox"/> #21120 ROHO Agility Min 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 990
		<input type="checkbox"/> #21100 nxt Xtend 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 760
		<input type="checkbox"/> #21105 22"-24" width, 16"-20" height \$ 875
		<input type="checkbox"/> #21122 ROHO Agility Mid 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 1,025
	6" to 7" depth	<input type="checkbox"/> #21101 nxt Xtend Deep 300 lbs weight capacity. 16"-20" width, 16"-20" height \$ 760
		<input type="checkbox"/> #21106 22"-24" width, 16"-20" height \$ 875
		<input type="checkbox"/> #21116 Matrx Elite E2 Deep 300 lbs weight capacity. 14"-20" width, 12"-20" height \$ 1,290
		<input type="checkbox"/> #21117 22"-24" width, 12"-20" height \$ 1,565
		<input type="checkbox"/> #21124 ROHO Agility Max 300 lbs weight capacity. 14"-20" width, 10"-18" height \$ 1,195

4. Specify the size →

Backrest Size Sizes based on availability from the backrest manufacturer.

Specify Size (Width" x Height") _____

Pediatric Seating System

4030400 **Diego Pediatric Seating System**

Select this item to request a customized quotation for delivery of a Diego System with this chair.

Back & Seat Matrix				
in (cm)	10	12	14	16
Width	(25.40)	(30.48)	(35.56)	(40.64)
Height / Depth				
10 (25.40)	✓	✓	✗	✗
12 (30.48)	✗	✓	✓	✗
14 (35.56)	✗	✗	✓	✓



TBD

Headrest

Stealth Comfort Plus Headrest	Stealth Hardware		
	<input type="checkbox"/> ##24110	Multi-Axis HMO475-1	\$ 490
	<input type="checkbox"/> ##24111	UniLink	\$ 760
	Comfort Plus 2 Pad		\$ 500
	<input type="checkbox"/> ##24115	Large (10" x 6")	<input type="checkbox"/> ##24116 X-Large (14" x 6")
nxt MultiFit Head Support	<input type="checkbox"/> ##24100 nxt Hardware		\$ 315
	Multi-Adjustable Pad		\$ 325
	<input type="checkbox"/> ##24105	Contoured (9" x 5")	<input type="checkbox"/> ##24107 X-Large (14" x 5")
	<input type="checkbox"/> ##24106	Large (10" x 5")	

Positioning Belts & Straps

Bodypoint	Hip Belts		
	<input type="checkbox"/> ##29010	Center-Pull Auto Style	\$ 130
	<input type="checkbox"/> ##29015	Center-Pull Auto Style, Padded	\$ 240
	<input type="checkbox"/> ##29020	4-Point, Padded	\$ 260
	<input type="checkbox"/> ##29030	Chest Strap	\$ 185

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.

PDG is proud to supply on-chair seating and positioning products from these manufacturers:

