



### Customer and Order Information

Purchase Order *	Tag	v.0119
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

### Base Chair

<input checked="" type="checkbox"/> 1901000	<b>Stellar GL 45° Manual Tilt Wheelchair</b>	\$ 4,200
<input type="checkbox"/> ⓘ	Also available: the lightweight <b>Stellar GLT</b> with aluminum frame. See separate order form. <a href="http://pdgmobility.com/glt">pdgmobility.com/glt</a>	

### Weight Capacity

<b>Chair Weight Capacity</b>		
<input type="checkbox"/> 1902000	250 lbs Capacity Available for 14" to 22" Seat Width.	
<input type="checkbox"/> 1902050	300 lbs Capacity Available for 14" to 22" Seat Width.	\$ 200
<input type="checkbox"/> 1905006	450 lbs Capacity Available for 20" to 32" Seat Width. The <b>Stellar HD</b> offers superior strength and performance. See separate order form.	\$ 500
<b>User Weight Range</b> Required for setup of tilt uplift force.		
<input type="checkbox"/> 1902430	Up to 75 lbs	
<input type="checkbox"/> 1902431	75 lbs to 100 lbs	
<input type="checkbox"/> 1902432	100 lbs to 115 lbs	
<input type="checkbox"/> 1902433	115 lbs to 150 lbs	
<input type="checkbox"/> 1902434	150 lbs to 200 lbs	
<input type="checkbox"/> 1902435	200 lbs to 250 lbs	
<input type="checkbox"/> 1902436	250 lbs to 300 lbs	
<input type="checkbox"/> 1902437	300 lbs to 375 lbs	
<input type="checkbox"/> 1902438	375 lbs to 450 lbs	

### Frame Options

<input type="checkbox"/> 1902200	<b>Extended Base with Rear Floor Jacks</b>	\$ 900
Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.		
<input type="checkbox"/> 1902300	<b>Transit Tie Down System (TTDS)</b>	
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Seat Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.		

### Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

<b>Standard Seat Widths</b>		
<input type="checkbox"/> 1903014	14" Seat Width	
<input type="checkbox"/> 1903015	15" Seat Width	
<input type="checkbox"/> 1903016	16" Seat Width	
<input type="checkbox"/> 1903017	17" Seat Width	
<input type="checkbox"/> 1903018	18" Seat Width	
<input type="checkbox"/> 1903019	19" Seat Width	
<input type="checkbox"/> 1903020	20" Seat Width	
<input type="checkbox"/> 1903021	21" Seat Width	
<input type="checkbox"/> 1903022	22" Seat Width	
<b>Custom Seat Widths</b> Requires 450 lbs Weight Capacity Frame. Extended delivery lead times may apply.		
<input type="checkbox"/> 1903023	23" Seat Width	\$ 375
<input type="checkbox"/> 1903024	24" Seat Width	\$ 375
<input type="checkbox"/> 1903025	25" Seat Width	\$ 700
<input type="checkbox"/> 1903026	26" Seat Width	\$ 700
<input type="checkbox"/> 1903027	27" Seat Width	
<input type="checkbox"/> 1903028	28" Seat Width	
<input type="checkbox"/> 1903029	29" Seat Width	
<input type="checkbox"/> 1903030	30" Seat Width	
<input type="checkbox"/> 1903031	31" Seat Width	
<input type="checkbox"/> 1903032	32" Seat Width	

### Seat Depth



<b>Standard Seat Depths (Adjustable)</b>		
<input type="checkbox"/> 1904016	16" Seat Depth	
<input type="checkbox"/> 1904017	17" Seat Depth	
<input type="checkbox"/> 1904018	18" Seat Depth	
<input type="checkbox"/> 1904019	19" Seat Depth	
<input type="checkbox"/> 1904020	20" Seat Depth	
<b>Custom Seat Depths</b>		
<input type="checkbox"/> 1904021	21" Seat Depth	
<input type="checkbox"/> 1904022	22" Seat Depth	\$ 400

### Frame Color

#### PDG Premium Colors

<input type="checkbox"/> 1905000 Tracer Red	<input type="checkbox"/> 1905001 Gloss Black Shimmer	<input type="checkbox"/> 1905002 Silver
<input type="checkbox"/> 1905003 Fine Texture Black	<input type="checkbox"/> 1905004 Orange Fusion	<input type="checkbox"/> 1905005 Karl's Imagination
<input type="checkbox"/> 1905006 Blue Shimmer	<input type="checkbox"/> 1905007 Gloss White	<input type="checkbox"/> 1905008 Dark Titanium
<input type="checkbox"/> 1905009 Geyser Green	<input type="checkbox"/> 1905010 Lollipop Blue	<input type="checkbox"/> 1905011 Lollipop Red
<input type="checkbox"/> 1905012 Bubble Gum Pink	<input type="checkbox"/> 1905013 Berry Purple	<input type="checkbox"/> 1905014 Daffodil Yellow
<input type="checkbox"/> 1905200 <b>Custom Color</b> Contact PDG for specifications. <span style="float: right;">\$ 300</span>		

### Armrests

<input type="checkbox"/> 1906000 <b>No Armrests</b>			\$ 295
<b>Single Post Height Adjustable Armrests</b>			
<input type="checkbox"/> 1906111 8"-11" Range	<input type="checkbox"/> 1906114 10"-14" Range		
<b>Dual Post Height Adjustable Armrests</b>			
<input type="checkbox"/> 1906211 8"-11" Range	<input type="checkbox"/> 1906214 10"-14" Range	<input type="checkbox"/> 1906216 12"-16" Range	\$ 295
<b>Fixed Height Armrests (12")</b> Full Length Armpad only.			
<input type="checkbox"/> 1906100 Single Post	<input type="checkbox"/> 1906200 Dual Post		
 <b>Armrest Height</b> If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used. Seat Pan to top of Armrest: _____ Does not include height of Armpad			
 <b>Armrest Position</b> Indicate preference for Armrest position. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 1906401 Rear	<input type="checkbox"/> 1906402 Mid	<input type="checkbox"/> 1906403 Front	
<b>Locking Flip-up</b> Not recommended for users with excessive strength or tone. <span style="float: right;">\$ 295</span> Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.			
<input type="checkbox"/> 1906310 Set at 10"	<input type="checkbox"/> 1906311 Set at 11"	<input type="checkbox"/> 1906312 Set at 12"	
<input type="checkbox"/> 1906313 Set at 13"	<input type="checkbox"/> 1906314 Set at 14"		

### Armpads

<b>Regular Armpads</b>			
<input type="checkbox"/> 1907100 Desk Length - 2" x 10"	<input type="checkbox"/> 1907101 Full Length - 2" x 14"		
<b>Waterfall Armpads</b>			
<input type="checkbox"/> 1907105 Desk Length - 2" x 10"	<input type="checkbox"/> 1907106 Full Length - 2" x 14"		
<b>Wide Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 200</span>			
<input type="checkbox"/> 1907110 Desk Length - 3" x 10"	<input type="checkbox"/> 1907111 Full Length - 3" x 12"		
<b>Gel Armpads</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 135</span>			
<input type="checkbox"/> 1907120 Full Length - 2" x 12"	<input type="checkbox"/> 1907121 Full Length - 3.5" x 12"		
<b>Contour Arm Trough with Hand Pad</b> May interfere with wheels in some tilt configurations. <span style="float: right;">\$ 125</span>			
<input type="checkbox"/> 1907130 Right (4.5" x 18")	<input type="checkbox"/> 1907131 Left (4.5" x 18")		

### Armrest Options

<b>Armrest Wideners</b> Move Armrests outward by 1" on each side. Do not increase seat size.			
<input type="checkbox"/> 1907210 For Single Post Armrests	May cause tilt interference at lower Seat Height and/or Seat Depth settings.	\$ 195	
<input type="checkbox"/> 1907220 For Dual Post Armrests	Limits tilt. Not compatible with 20", 22", or 24" Rear Wheels, Clothing Guards.	\$ 375	
<b>Clothing Guards</b> Not compatible with Flip-up Armrests, Armrest Wideners.			
<input type="checkbox"/> 1907310 For Single Post Armrests	Available for 10"-14" Height range only.	\$ 110	
<input type="checkbox"/> 1907320 For Dual Post Armrests		\$ 165	
<b>Armrest Extensions</b> Not compatible with Flip-up or Single Post Armrests. <span style="float: right;">\$ 75</span>			
<input type="checkbox"/> 1907400 3" Extension - Forward	<input type="checkbox"/> 1907401 3" Extension - Rearward		

### Seat Pan

- 1908000 **No Seat Pan**
- 1908001 **Solid Aluminum Seat Pan** \$ 135

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
6" Casters	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8" Casters	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Casters

- Urethane Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1911061 6" x 1.25" Caster size  1911081 8" x 1.25" Caster size
  - 1911063 6" x 1.50" Caster size \$ 85
  - 1911084 8" x 2.00" Caster size \$ 135
- Pneumatic Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1911161 6" x 1.25" Caster size \$ 110
  - 1911181 8" x 1.25" Caster size \$ 120
  - 1911184 8" x 2.00" Caster size \$ 150
  - 1911260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275
  - 1911400 **Frog Legs Suspension Forks** Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster. \$ 695

### Rear Wheels

- Composite Wheels**
- 1912012 12" Wheel size  1912016 16" Wheel size  1912020 20" Wheel size
  - 1912022 22" Wheel size  1912024 24" Wheel size
- Spoke Wheels** \$ 75
- 1912120 20" Wheel size  1912122 22" Wheel size  1912124 24" Wheel size
  - 1912224 **24" Heavy Duty Spoke Wheels** \$ 350
- Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.  
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.
- 1912401 Rear  1912402 Mid  1912403 Front

### Axle Type

- 1913000 **Quick Release Rear Axles**
- 1913010 **Fixed Rear Axles**

### Rear Tires

- 1914000 **Urethane Tires**
- 1914010 **Pneumatic Tires** May result in a finished Seat Height 1/2" greater than stated. \$ 60
- 1914020 **Pneumatic Tires with Flat-Free Inserts** May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. \$ 95

Handrims		Handrims are only available with 20", 22", and 24" Rear Wheels.
<input type="checkbox"/>	1915000 <b>No Handrims</b>	
<input type="checkbox"/>	1915010 <b>Standard Aluminum Handrims</b>	
<input type="checkbox"/>	1915020 <b>Plastic Coated Handrims</b>	\$ 150

Wheel Locks & Brakes				
<b>User Wheel Locks</b>				
<input type="checkbox"/>	1916010 Push to Lock	<input type="checkbox"/>	1916020 Pull to Lock	
<input type="checkbox"/>	1916030 <b>Attendant Wheel Locks (Push to Lock)</b>			\$ 85
<input type="checkbox"/>	1916040 <b>User &amp; Attendant Wheel Locks</b>	Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes.		\$ 125
<input type="checkbox"/>	1916300 <b>Drum Brakes</b>	Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.		\$ 995

Wheel Lock Options			
<input type="checkbox"/>	1916100 <b>Wheel Lock Extensions</b>	Not compatible with Attendant Wheel Lock.	\$ 60
<input type="checkbox"/>	1916200 <b>Hill Holders</b>	Requires Pneumatic Rear Tires with Flat-Free Inserts.	\$ 450

Anti-Tippers			
<input checked="" type="checkbox"/>	1917010 <b>Rear Anti-Tippers</b>	Required Safety Option.	\$ 75
<input type="checkbox"/>	1917020 <b>Front Anti-Tippers</b>		\$ 795
<input type="checkbox"/>	1917030 <b>Heavy Duty Rear Anti-Tipper</b>		\$ 895
Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.			

Back Posts					
<b>Back Post Style</b>					
<input type="checkbox"/>	1918010 Push Bar Back	<input type="checkbox"/>	1918020 Cane Handle Back	Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.	
<b>Back Post Height</b>		Measured from the Seat Pan to the Handle.			
<input type="checkbox"/>	1918425 Short (22")	<input type="checkbox"/>	1918431 Tall (25")		
<b>Seat to Back Angle</b>					
<input type="checkbox"/>	1919090 Set at 90°	<input type="checkbox"/>	1919098 Set at 98°	<input type="checkbox"/>	1919106 Set at 106°
<input type="checkbox"/>	1919114 Set at 114°	<input type="checkbox"/>	1919122 Set at 122°		

Back Options				
<b>Angle/Height Adjustable Push Bar Handle</b>		Not compatible with Headrest Mounting Plate or Cane Handle Back.	\$ 150	
<input type="checkbox"/>	1919010 Button on Outside	<input type="checkbox"/>	1919011 Button on Inside	
<input type="checkbox"/>	1919210 <b>10" Backrest Extension</b>	Requires the Angle/Height Adjustable Push Bar Handle.	\$ 195	
<input type="checkbox"/>	1919212 <b>Headrest Mounting Plate</b>	Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 95	
<input type="checkbox"/>	1919220 <b>Fold Down Back</b>	Not compatible with Reclining Back.		
<input type="checkbox"/>	1919230 <b>Spring Back Suspension</b>	Not compatible with Reclining Back.	\$ 105	
<input type="checkbox"/>	1919240 <b>Seating Dynamics Dynamic Rocker Back Interface</b>	Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle, Standard Upholstery. Finished Back Post Height is 4" taller than selected Back Post Height.	\$ 2,995	

Back Upholstery					
<input type="checkbox"/>	1920000 <b>No Back Upholstery</b>				
<b>Standard Upholstery</b>		Not available with Cane Handle Back.	\$ 265		
<input type="checkbox"/>	1920418 Short (18" Height)	<input type="checkbox"/>	1920421 Regular (21" Height)	<input type="checkbox"/>	1920424 Tall (24" Height)

See next page for more options >

### Tension Adjustable Upholstery

Height of Upholstery is 16". Not available with Cane Handle Back.

- 1920450 For 14" to 22" Seat Width \$ 365
- 1920451 For 24" to 32" Seath Width \$ 525

### Position of Upholstery

Select the desired position of the top of the Upholstery relative to the Seat Pan.

Style	22" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Style	25" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

## Tilt Activation

### Attendant Tilt only

Select trigger mounting location below. Chairs with 450 lbs Weight Capacity require two tilt triggers (L and R) for tilt activation.

- 1923011 Right Back Post  1923012 Left Back Post

### User Tilt only

Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below.

- 1923021 Right Armrest  1923022 Left Armrest \$ 100

### Attendant & User Tilt

Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below.

- 1923031 Right Back Post  1923032 Left Back Post
- 1923041 Right Armrest  1923042 Left Armrest \$ 195

### Tilt Trigger Extension

- 1923110 \$ 50

### Power Tilt

250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. [pdgmobility.com/power-tilt](http://pdgmobility.com/power-tilt)

Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.

- Toggle Switch  1923421 Installed on Right side  1923422 Installed on Left side
- Button Switch  1923431 Installed on Right side  1923432 Installed on Left side
- Stealth Egg Switch  1923440 Not installed

Optional for Egg Switch:  1923445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625

*Specify Headrest to be used:* \_\_\_\_\_

### Tilt Limit

Optional. If no tilt limit is selected, the full tilt range will be available.

- 1923220 Set at 20°  1923230 Set at 30°  1923240 Set at 40°

## Recline

- 1923310 **Reclining Backrest (Manual)** \$ 895

30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.  
Finished Back Post Height is 3" taller than selected Back Post Height.

## Front Rigging

- 1925000 **No Front Rigging**  1925001 **Sector Blocks only**

### Easy-Out™ Front Rigging

Effortless insertion and removal. No Hanger Bracket remains on frame. [pdgmobility.com/easy-out](http://pdgmobility.com/easy-out)

- 1925410 Fixed 70° \$ 95
- 1925420 Universal Angle Adjustable (0°-110°) Set at angle: \_\_\_\_\_ (use increments of 10°) \$ 295

Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.

### Swing Away Front Rigging

Legacy Design. Not compatible with Easy-Out.

- 1925130 Swing Away 60°  1925140 Swing Away 70°
- 1925200 Swing Away Elevating Legrests Includes Calf Pads. \$ 375
- 1925300 1" Offset Hangers for Swing Away \$ 75

- 1925110 **Standard Plug In 70° Front Rigging** Legacy Design. Not compatible with Easy-Out.

## Footplates

- 1926010 **Composite Footplates**
- 1926020 **Aluminum Angle Adjustable Footplates** \$ 165


See next page for more options >

<input type="checkbox"/>	1926030	<b>One Piece Fixed Footplate</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1926040	<b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1926050	<b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.	\$ 195

### Footrest/Legrest Accessories

<b>Calf Pads</b> Not compatible with Gel Padding.				
<input type="checkbox"/>	1927012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	1927010	For Swing Away Front Rigging		\$ 135
<input type="checkbox"/>	1927020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
<b>Gel Knee Pads</b>				
<input type="checkbox"/>	1927032	For Easy-Out™ Angle Adjustable Front Rigging		\$ 200
<input type="checkbox"/>	1927030	For Swing Away ELR Front Rigging		\$ 100
<input type="checkbox"/>	1927040	<b>Length Adjustable Heel Loops</b>		\$ 50
<input type="checkbox"/>	1927050	<b>Bodypoint Padded Calf Strap</b>	Sizes greater than 22" are not padded.	\$ 195
<b>Miller's Dynamic Footrest Mounts</b> Not compatible with Easy-Out Front Rigging.				
<input type="checkbox"/>	1927060	For Standard Plug In Front Rigging		\$ 445
<input type="checkbox"/>	1927065	For Swing Away 70° Front Rigging		\$ 575

### Front Rigging Extension Tubes

<b>Regular Extension Tubes</b>				
<input type="checkbox"/>	1928010	Short	14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	1928020	Long	18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75
<b>Extra Strength Extension Tubes (Solid)</b> Available on Plug In Front Rigging only. Not compatible with Swing Away or Easy-Out Front Rigging.				
<input type="checkbox"/>	1928110	Standard (13"-17")	<input type="checkbox"/>	1928120 Long (15"-19")
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.			
	Seat Pan to Foot Plates: _____		Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.	

### Miscellaneous Options

<b>Spoke Guards</b> Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.				
<input type="checkbox"/>	1930020	Black	<input type="checkbox"/>	1930021 Frame Color
<input type="checkbox"/>	1930080	<b>Gel Hip Guides</b>	Reduces available Seat Width by 2".	\$ 495
<input type="checkbox"/>	1930110	<b>IV Pole</b>	Not compatible with HD Anti Tipper, Attendant Wheel Locks.	\$ 345
<input type="checkbox"/>	1930120	<b>O2 Holder</b>	Not compatible with HD Anti Tipper, Attendant Wheel Locks.	\$ 495
<input type="checkbox"/>	1930170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.	\$ 10

### On-Chair Seating & Positioning

	Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. <a href="http://pdgmobility.com/seating">pdgmobility.com/seating</a>
---	---

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.