



# STELLAR GLT

## Order Form - Distributor

Prices in Canadian Dollars. January 1, 2022

### Customer and Order Information

Purchase Order *	Tag	v.0420
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

### Base Chair

#### Stellar GLT 45° Manual Tilt Wheelchair

- 2501000 Stellar GL 45° Manual Tilt Wheelchair \$ 4,200
- 2501009 Lightweight Aluminum Frame Upgrade Aluminum Frame is not recommended for agitated users. \$ 300
- Choose the standard **Stellar GL** with a steel frame for agitated or heavier users. See separate order form.

### Weight Capacity

- 2502000 **250 lbs Chair Weight Capacity**
- User Weight Range** Required for setup of tilt uplift force.
- 1902430 Up to 75 lbs
- 1902433 115 lbs to 150 lbs
- 2502431 75 lbs to 100 lbs
- 2502434 150 lbs to 200 lbs
- 2502432 100 lbs to 115 lbs
- 2502435 200 lbs to 250 lbs

### Frame Options

- 2502300 **Transit Tie Down System (TTDS)**  
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

### Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

#### Standard Seat Widths

- 2503014 14" Seat Width
- 2503017 17" Seat Width
- 2503020 20" Seat Width
- 2503015 15" Seat Width
- 2503018 18" Seat Width
- 2503021 21" Seat Width
- 2503016 16" Seat Width
- 2503019 19" Seat Width
- 2503022 22" Seat Width

### Seat Depth

#### Standard Seat Depths (Adjustable)

- 2504016 16" Seat Depth
- 2504019 19" Seat Depth
- 2504017 17" Seat Depth
- 2504020 20" Seat Depth
- 2504018 18" Seat Depth

#### Custom Seat Depths

- 2504021 21" Seat Depth
  - 2504022 22" Seat Depth
- \$ 400

### Frame Color

#### PDG Premium Colors

- 2505000 Tracer Red
- 2505003 Fine Texture Black
- 2505006 Blue Shimmer
- 2505009 Geyser Green
- 2505012 Bubble Gum Pink
- 2505001 Gloss Black Shimmer
- 2505004 Orange Fusion
- 2505007 Gloss White
- 2505010 Lollipop Blue
- 2505013 Berry Purple
- 2505002 Silver
- 2505005 Karl's Imagination
- 2505008 Dark Titanium
- 2505011 Lollipop Red
- 2505014 Daffodil Yellow

- 2505200 **Custom Color** Contact PDG for specifications. \$ 300

### Armrests

2506000 **No Armrests**

**Height Adjustable Armrests (Single Post)** \$ 295

2506111 8"-11" Range  2506114 10"-14" Range

2506100 **Fixed Height Armrests (12")** Full Length Armpad only.

**Armrest Height** If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.  
 Seat Pan to top of Armrest: \_\_\_\_\_ Does not include height of Armpad

**Armrest Position** Indicate preference for Armrest position. If left blank, the default setting will be used.  
 In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2506401 Rear  2506402 Mid  2506403 Front

### Armpads

**Regular Armpads**

2507100 Desk Length - 2" x 10"  2507101 Full Length - 2" x 14"

**Waterfall Armpads** \$ 115

2507105 Desk Length - 2" x 10"  2507106 Full Length - 2" x 14"

**Wide Armpads** May interfere with wheels in some tilt configurations. \$ 200

2507110 Desk Length - 3" x 10"  2507111 Full Length - 3" x 12"

**Gel Armpads** May interfere with wheels in some tilt configurations. \$ 135

2507120 Full Length - 2" x 12"  2507121 Full Length - 3.5" x 12"

**Contour Arm Trough with Hand Pad** May interfere with wheels in some tilt configurations. \$ 125

2507130 Right (4.5" x 18")  2507131 Left (4.5" x 18")

### Armrest Options

2507210 **Armrest Wideners** Move Armrests outward by 1" on each side. Do not increase seat size. \$ 195

2507310 **Clothing Guards** Available for 10"-14" Height range only. \$ 110

### Seat Pan

2508000 **No Seat Pan**

2508001 **Solid Aluminum Seat Pan** \$ 135

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

Casters	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
6" Casters	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8" Casters	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Casters

#### Urethane Casters

8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2511061	6" x 1.25" Caster size	<input type="checkbox"/>	2511081	8" x 1.25" Caster size	
<input type="checkbox"/>	2511063	6" x 1.50" Caster size				\$ 85
<input type="checkbox"/>	2511084	8" x 2.00" Caster size				\$ 135

#### Pneumatic Casters

8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2511161	6" x 1.25" Caster size				\$ 110
<input type="checkbox"/>	2511181	8" x 1.25" Caster size				\$ 120
<input type="checkbox"/>	2511184	8" x 2.00" Caster size				\$ 150

<input type="checkbox"/>	2511260	<b>Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters</b>				\$ 275
--------------------------	---------	--	--	--	--	--------

<input type="checkbox"/>	2511400	<b>Frog Legs Suspension Forks</b>	Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster.			\$ 695
--------------------------	---------	-----------------------------------	--	--	--	--------

### Rear Wheels

#### Composite Wheels

<input type="checkbox"/>	2512012	12" Wheel size	<input type="checkbox"/>	2512016	16" Wheel size	<input type="checkbox"/>	2512020	20" Wheel size
<input type="checkbox"/>	2512022	22" Wheel size	<input type="checkbox"/>	2512024	24" Wheel size			

#### Spoke Wheels

<input type="checkbox"/>	2512120	20" Wheel size	<input type="checkbox"/>	2512122	22" Wheel size	<input type="checkbox"/>	2512124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------



#### Rear Wheel Position

Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/>	2512401	Rear	<input type="checkbox"/>	2512402	Mid	<input type="checkbox"/>	2512403	Front
--------------------------	---------	------	--------------------------	---------	-----	--------------------------	---------	-------

### Axle Type

<input type="checkbox"/>	2513000	<b>Quick Release Rear Axles</b>
<input type="checkbox"/>	2513010	<b>Fixed Rear Axles</b>

### Rear Tires

<input type="checkbox"/>	2514000	<b>Urethane Tires</b>			
<input type="checkbox"/>	2514010	<b>Pneumatic Tires</b>	May result in a finished Seat Height 1/2" greater than stated.		\$ 60
<input type="checkbox"/>	2514020	<b>Pneumatic Tires with Flat-Free Inserts</b>	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.		\$ 95

### Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/>	2515000	<b>No Handrims</b>	
<input type="checkbox"/>	2515010	<b>Standard Aluminum Handrims</b>	
<input type="checkbox"/>	2515020	<b>Plastic Coated Handrims</b>	\$ 150

### Wheel Locks & Brakes

#### User Wheel Locks

<input type="checkbox"/>	2516010	Push to Lock	<input type="checkbox"/>	2516020	Pull to Lock	
<input type="checkbox"/>	2516030	<b>Attendant Wheel Locks (Push to Lock)</b>			\$ 85	
<input type="checkbox"/>	2516040	<b>User &amp; Attendant Wheel Locks</b>	Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes.			\$ 125
<input type="checkbox"/>	2516300	<b>Drum Brakes</b>	Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.			\$ 995

### Wheel Lock Options

<input type="checkbox"/>	2516100	<b>Wheel Lock Extensions</b>	Not compatible with Attendant Wheel Lock.		\$ 60
<input type="checkbox"/>	2516200	<b>Hill Holders</b>	Requires Pneumatic Rear Tires with Flat-Free Inserts.		\$ 450

### Anti-Tippers

<input checked="" type="checkbox"/>	2517010	<b>Rear Anti-Tippers</b>	Required Safety Option.		\$ 75
-------------------------------------	---------	--------------------------	-------------------------	--	-------



### Back Posts

**Back Post Style**

2518010 Push Bar Back

---

**Back Post Height** Measured from the Seat Pan to the Handle.

2518425 Short (22")  2518431 Tall (25")

---

**Seat to Back Angle**

2519090 Set at 90°  2519098 Set at 98°  2519106 Set at 106°

2519114 Set at 114°  2519122 Set at 122°

### Back Options

**Angle/Height Adjustable Push Bar Handle** Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150

2519010 Button on Outside  2519011 Button on Inside

---

2519210 **10" Backrest Extension** Requires the Angle/Height Adjustable Push Bar Handle. \$ 195

---

2519212 **Headrest Mounting Plate** \$ 95  
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.

---

2519220 **Fold Down Back** Not compatible with Reclining Back.

### Back Upholstery

2520000 **No Back Upholstery**

---

**Standard Upholstery** Not available with Cane Handle Back. \$ 265

2520418 Short (18" Height)  2520421 Regular (21" Height)  2520424 Tall (24" Height)

---

2520450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. \$ 365

---

**Position of Upholstery** Select the desired position of the top of the Upholstery relative to the Seat Pan.

22" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

25" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Tilt Activation

**Attendant Tilt only** Select trigger mounting location below.

2523011 Right Back Post  2523012 Left Back Post

---

**User Tilt only** Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below. \$ 100

2523021 Right Armrest  2523022 Left Armrest

---

**Attendant & User Tilt** Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below. \$ 195

2523031 Right Back Post  2523032 Left Back Post

2523041 Right Armrest  2523042 Left Armrest

---

2523110 **Tilt Trigger Extension** \$ 50

---

2523400 **Power Tilt** 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. [pdgmobility.com/power-tilt](http://pdgmobility.com/power-tilt) \$ 2,095  
Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.

Toggle Switch  2523421 Installed on Right side  2523422 Installed on Left side

Button Switch  2523431 Installed on Right side  2523432 Installed on Left side

Stealth Egg Switch  2523440 Not installed

Optional for Egg Switch:  2523445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. \$ 625  
Specify Headrest to be used: \_\_\_\_\_

---

**Tilt Limit** Optional. If no tilt limit is selected, the full tilt range will be available.

2523220 Set at 20°  2523230 Set at 30°  2523240 Set at 40°

### Recline

<input type="checkbox"/>	2523310	<b>Reclining Backrest (Manual)</b>		\$ 895
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post Height.				

### Front Rigging

<input type="checkbox"/>	2525000	<b>No Front Rigging</b>	<input type="checkbox"/>	2525001	<b>Sector Blocks only</b>		
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>							
<input type="checkbox"/>	2525410	Fixed 70°					\$ 95
<input type="checkbox"/>	2525420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)				\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.							
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.							
<input type="checkbox"/>	2525130	Swing Away 60°	<input type="checkbox"/>	2525140	Swing Away 70°		
<input type="checkbox"/>	2525200	Swing Away Elevating Legrests	Includes Calf Pads.				\$ 375
<input type="checkbox"/>	2525300	1" Offset Hangers for Swing Away					\$ 75
<input type="checkbox"/>	2525110	<b>Standard Plug In 70° Front Rigging</b>	Legacy Design. Not compatible with Easy-Out.				


### Footplates

<input type="checkbox"/>	2526010	<b>Composite Footplates</b>					
<input type="checkbox"/>	2526020	<b>Aluminum Angle Adjustable Footplates</b>					\$ 165
<input type="checkbox"/>	2526030	<b>One Piece Fixed Footplate</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.				\$ 350
<input type="checkbox"/>	2526040	<b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.				\$ 395
<input type="checkbox"/>	2526050	<b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.				\$ 195

### Footrest/Legrest Accessories

<b>Calf Pads</b> Not compatible with Gel Padding.							
<input type="checkbox"/>	2527012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.				\$ 200
<input type="checkbox"/>	2527010	For Swing Away Front Rigging					\$ 135
<input type="checkbox"/>	2527020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.				\$ 125
<b>Gel Knee Pads</b>							
<input type="checkbox"/>	2527032	For Easy-Out™ Angle Adjustable Front Rigging					\$ 200
<input type="checkbox"/>	2527030	For Swing Away ELR Front Rigging					\$ 100
<input type="checkbox"/>	2527040	<b>Length Adjustable Heel Loops</b>					\$ 50
<input type="checkbox"/>	2527050	<b>Bodypoint Padded Calf Strap</b>					\$ 195

### Front Rigging Extension Tubes

<b>Regular Extension Tubes</b>							
<input type="checkbox"/>	2528010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging			
<input type="checkbox"/>	2528020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging			\$ 75
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.						
Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.							



# STELLAR GLT

## Order Form - Distributor

Prices in Canadian Dollars. January 1, 2022

### Miscellaneous Options

<b>Spoke Guards</b> Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.		\$ 235
<input type="checkbox"/> 2530020	Black	<input type="checkbox"/> 2530021 Frame Color
<input type="checkbox"/> 2530080	<b>Gel Hip Guides</b> Reduces available Seat Width by 2".	\$ 495
<input type="checkbox"/> 2530110	<b>IV Pole</b> Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later.	\$ 345
<input type="checkbox"/> 2530120	<b>O2 Holder</b> Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later.	\$ 495
<input type="checkbox"/> 2530170	<b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.	\$ 10

### On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

### Order Notes

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.