

Customer and Order Information		
Purchase Order *	Tag	v.0318
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in *required

Base Chair		
Stellar HD 45° Manual Tilt Wheelchair Extended delivery lead times may apply.		
<input checked="" type="checkbox"/> 2101000	Stellar 45° Manual Tilt Wheelchair	\$ 4,200
<input checked="" type="checkbox"/> 2101201	Reinforced Frame Package Heavy-gauge steel reinforcement of frame for bariatric users.	\$ 1,395
Choose the Stellar Impact for highly agitated users up to 250 lbs. See separate order form.		

Weight Capacity		
Chair Weight Capacity		
<input type="checkbox"/> 2102045	450 lbs Capacity	
<input type="checkbox"/> 2102060	600 lbs Capacity Tilt is limited to maximum 25°.	\$ 300
User Weight Range Required for setup of tilt uplift force.		
<input type="checkbox"/> 2102212	Up to 250 lbs	<input type="checkbox"/> 2102308 250 lbs to 350 lbs
<input type="checkbox"/> 2102310	350 lbs to 425 lbs	<input type="checkbox"/> 2102312 425 lbs to 600 lbs

Seat Width			Overall chair width is +11"
Standard Seat Widths			
<input type="checkbox"/> 2103022	22" Seat Width	<input type="checkbox"/> 2103024	24" Seat Width
<input type="checkbox"/> 2103028	28" Seat Width	<input type="checkbox"/> 2103030	30" Seat Width
<input type="checkbox"/> 2103026	26" Seat Width	<input type="checkbox"/> 2103032	32" Seat Width
			\$ 375
			\$ 800

Seat Depth		
Standard Seat Depths (Adjustable)		
<input type="checkbox"/> 2104018	18" Seat Depth	<input type="checkbox"/> 2104019 19" Seat Depth
<input type="checkbox"/> 2104021	21" Seat Depth	<input type="checkbox"/> 2104020 20" Seat Depth
<input type="checkbox"/> 2104022	22" Seat Depth	
		\$ 300

Frame Color		
PDG Premium Colors		
<input type="checkbox"/> 2105000	Tracer Red	<input type="checkbox"/> 2105001 Gloss Black Shimmer
<input type="checkbox"/> 2105003	Fine Texture Black	<input type="checkbox"/> 2105002 Silver
<input type="checkbox"/> 2105006	Blue Shimmer	<input type="checkbox"/> 2105004 Orange Fusion
<input type="checkbox"/> 2105009	Geyser Green	<input type="checkbox"/> 2105005 Karl's Imagination
<input type="checkbox"/> 2105012	Bubble Gum Pink	<input type="checkbox"/> 2105007 Gloss White
		<input type="checkbox"/> 2105008 Dark Titanium
		<input type="checkbox"/> 2105010 Lollypop Blue
		<input type="checkbox"/> 2105011 Lollypop Red
		<input type="checkbox"/> 2105012 Berry Purple
		<input type="checkbox"/> 2105014 Daffodil Yellow
<input type="checkbox"/> 2105200	Custom Color Contact PDG for specifications.	\$ 300

Armrests		
<input type="checkbox"/> 2106000	No Armrests	
<input type="checkbox"/> 2106214	Dual Post Height Adjustable Armrests (10"-14")	\$ 295
Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.		
	Seat Pan to top of Armrest: _____	Does not include height of Armpad


Armpads		
Regular Armpads		
<input type="checkbox"/> 2107100	Desk Length - 2" x 10"	<input type="checkbox"/> 2107101 Full Length - 2" x 14"
Waterfall Armpads		
<input type="checkbox"/> 2107105	Desk Length - 2" x 10"	<input type="checkbox"/> 2107106 Full Length - 2" x 14"
Wide Armpads May interfere with wheels in some tilt configurations.		
<input type="checkbox"/> 2107110	Desk Length - 3" x 10"	<input type="checkbox"/> 2107111 Full Length - 3" x 12"
Gel Armpads May interfere with wheels in some tilt configurations.		
<input type="checkbox"/> 2107120	Full Length - 2" x 12"	<input type="checkbox"/> 2107121 Full Length - 3.5" x 12"
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.		
<input type="checkbox"/> 2107130	Right (4.5" x 18")	<input type="checkbox"/> 2107131 Left (4.5" x 18")

Armrest Options		
<input type="checkbox"/> 2107320	Clothing Guards	\$ 165

Seat Pan		
<input checked="" type="checkbox"/> 2108001	Solid Aluminum Seat Pan	\$ 175

Seat Height		
Seat to Floor Height		
<input type="checkbox"/> 2110016	Set at 16" STF	<input type="checkbox"/> 2110017 Set at 17" STF
		<input type="checkbox"/> 2110018 Set at 18" STF

Casters		
<input checked="" type="checkbox"/> 2111050	5" Casters with Heavy Duty Steel Fork	

Rear Wheels		
<input checked="" type="checkbox"/> 2112224	24" Heavy Duty Spoke Wheels	
	Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.	
<input type="checkbox"/> 2112401	Rear	<input type="checkbox"/> 2112402 Mid
		<input type="checkbox"/> 2112403 Front

Axle Type		
<input checked="" type="checkbox"/> 2113010	Fixed Rear Axles	

Rear Tires		
<input checked="" type="checkbox"/> 2114000	Urethane Tires	

Handrims		
<input type="checkbox"/> 2115000	No Handrims	
<input type="checkbox"/> 2115010	Standard Aluminum Handrims	
<input type="checkbox"/> 2115020	Plastic Coated Handrims	\$ 250

Wheel Locks		
User Wheel Locks		
<input type="checkbox"/> 2116010	Push to Lock	<input type="checkbox"/> 2116020 Pull to Lock

Wheel Lock Options		
<input type="checkbox"/> 2116100	Wheel Lock Extensions	\$ 60

Anti-Tippers

- 2117010 **Rear Anti-Tippers** Required Safety Option. \$ 75

Back Post

- 2118010 **Push Bar Back**
- Back Post Height** Measured from the Seat Pan to the Handle.
- 2118425 Short (22") 2118431 Tall (25")
- 2119090 **Seat to Back Angle set at 90°**

Back Options

- 2119212 **Headrest Mounting Plate** Not required if mounting to a solid Back. \$ 95
- 2119240 **Seating Dynamics Dynamic Rocker Back Interface** \$ 2,995
- Not compatible with 30° Recline, Standard Upholstery. Finished Back Post Height is 4" taller than selected Back Post Height.

Back Upholstery

- 2120000 **No Back Upholstery**
- Tension Adjustable Upholstery** \$ 525
- 2120421 Regular (18" Height) 2120424 Tall (22" Height)
- Position of Upholstery** Select the desired position of the top of the Upholstery relative to the Seat Pan.
- | | | 22" Back Post | | | | | | | | 25" Back Post | | | | | | | | | | | |
|--|---------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|-----|-----|---------------|--|--|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|--------------------------|
| | | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" | | | 17" | 18" | 19" | 20" | 21" | 22" | 23" | 24" | 25" |
| | Regular | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| | Tall | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | | | | | | | | | | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Tilt Activation

- Attendant Tilt only** Select trigger mounting location below.
- 2123011 Right Back Post 2123012 Left Back Post
- Tilt Limit** Optional. If no tilt limit is selected, the full tilt range will be available. Maximum Tilt Limit for 600 lbs Weight Capacity is 25°.
- 2123220 Set at 20° 2123230 Set at 30° 2123240 Set at 40°

Recline

- 2123310 **Reclining Backrest (Manual)** \$ 895
- 30° maximum recline. Not compatible with Tension Adjustable Upholstery. Finished Back Post Height is 3" taller than selected Back Post Height.

Front Rigging

- 2125000 **No Front Rigging** 2125001 **Sector Blocks only**
- 2125130 **Swing Away 60° Front Rigging**
- 2125200 **Swing Away Elevating Legrests** Includes Calf Pads. \$ 375

Footplates

- 2126020 **Aluminum Angle Adjustable Footplates** \$ 175

Footrest/Legrest Accessories

- 2127030 **Swing Away ELR Gel Knee Pads** \$ 100
- 2127040 **Length Adjustable Heel Loops** \$ 50
- 2127051 **Calf Strap** \$ 195



STELLAR HD

Order Form - Distributor

Prices in Canadian Dollars. January 1, 2022

Front Rigging Extension Tubes

2128010 **Regular Extension Tubes (13"-17")**



Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

2130120 **O2 Holder** \$ 495

2130170 **Adhesive Back Hook Strip, 2" x 30"** For Cushion mounting. \$ 10

On-Chair Seating & Positioning

Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.