



Customer and Order Information

Purchase Order *	Tag	v.0119
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in *required

Base Chair

Stellar LEAP Manual Posterior & Anterior Tilt Wheelchair Extened lead times may apply.

<input checked="" type="checkbox"/> 2201000	Stellar GL Manual Wheelchair with 20° Posterior Tilt	\$ 4,200
<input checked="" type="checkbox"/> 2201101	25° Anterior Tilt with Linked Back 30° Anterior Tilt available with reduced Posterior Tilt range of 15°.	\$ 100
<input checked="" type="checkbox"/> 2201102	30° Dynamic Back Recline	\$ 895

Weight Capacity

<input checked="" type="checkbox"/> 2202000	250 lbs Chair Weight Capacity
---	--------------------------------------

Frame Options

<input type="checkbox"/> 2202300	Transit Tie Down System (TTDS) <small>TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.</small>
----------------------------------	---

Seat Width Overall chair width is +9.5" for 12/16" RW and +10.5" for 20/22/24" RW

Standard Seat Widths

<input type="checkbox"/> 2203014	14" Seat Width	<input type="checkbox"/> 2203015	15" Seat Width	<input type="checkbox"/> 2203016	16" Seat Width
<input type="checkbox"/> 2203017	17" Seat Width	<input type="checkbox"/> 2203018	18" Seat Width	<input type="checkbox"/> 2203019	19" Seat Width
<input type="checkbox"/> 2203020	20" Seat Width	<input type="checkbox"/> 2203021	21" Seat Width	<input type="checkbox"/> 2203022	22" Seat Width

Seat Depth

Standard Seat Depths (Adjustable)

<input type="checkbox"/> 2204016	16" Seat Depth	<input type="checkbox"/> 2204017	17" Seat Depth	<input type="checkbox"/> 2204018	18" Seat Depth
<input type="checkbox"/> 2204019	19" Seat Depth	<input type="checkbox"/> 2204020	20" Seat Depth		

Custom Seat Depths



<input type="checkbox"/> 2204021	21" Seat Depth	<input type="checkbox"/> 2204022	22" Seat Depth	\$ 400
----------------------------------	----------------	----------------------------------	----------------	--------

Frame Color

PDG Premium Colors

<input type="checkbox"/> 2205000	Tracer Red	<input type="checkbox"/> 2205001	Gloss Black Shimmer	<input type="checkbox"/> 2205002	Silver
<input type="checkbox"/> 2205003	Fine Texture Black	<input type="checkbox"/> 2205004	Orange Fusion	<input type="checkbox"/> 2205005	Karl's Imagination
<input type="checkbox"/> 2205006	Blue Shimmer	<input type="checkbox"/> 2205007	Gloss White	<input type="checkbox"/> 2205008	Dark Titanium
<input type="checkbox"/> 2205009	Geyser Green	<input type="checkbox"/> 2205010	Lollypop Blue	<input type="checkbox"/> 2205011	Lollypop Red
<input type="checkbox"/> 2205012	Bubble Gum Pink	<input type="checkbox"/> 2205013	Berry Purple	<input type="checkbox"/> 2205014	Daffodil Yellow

2205200 **Custom Color** Contact PDG for specifications. \$ 300

Armrests	
<input type="checkbox"/> 2206000	No Armrests
Single Post Height Adjustable Armrests \$ 295	
<input type="checkbox"/> 2206111	8"-11" Range <input type="checkbox"/> 2206114 10"-14" Range
Dual Post Height Adjustable Armrests \$ 295	
<input type="checkbox"/> 2206211	8"-11" Range <input type="checkbox"/> 2206214 10"-14" Range <input type="checkbox"/> 2206216 12"-16" Range
Fixed Height Armrests (12") Full Length Armpad only.	
<input type="checkbox"/> 2206100	Single Post <input type="checkbox"/> 2206200 Dual Post
	Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used. Seat Pan to top of Armrest: _____ Does not include height of Armpad
	Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.
<input type="checkbox"/> 2206401	Rear <input type="checkbox"/> 2206402 Mid <input type="checkbox"/> 2206403 Front
Locking Flip-up Not recommended for users with excessive strength or tone. \$ 295 Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.	
<input type="checkbox"/> 2206310	Set at 10" <input type="checkbox"/> 2206311 Set at 11" <input type="checkbox"/> 2206312 Set at 12"
<input type="checkbox"/> 2206313	Set at 13" <input type="checkbox"/> 2206314 Set at 14"

Armpads	
Regular Armpads	
<input type="checkbox"/> 2207100	Desk Length - 2" x 10" <input type="checkbox"/> 2207101 Full Length - 2" x 14"
Waterfall Armpads \$ 115	
<input type="checkbox"/> 2207105	Desk Length - 2" x 10" <input type="checkbox"/> 2207106 Full Length - 2" x 14"
Wide Armpads May interfere with wheels in some tilt configurations. \$ 200	
<input type="checkbox"/> 2207110	Desk Length - 3" x 10" <input type="checkbox"/> 2207111 Full Length - 3" x 12"
Gel Armpads May interfere with wheels in some tilt configurations. \$ 135	
<input type="checkbox"/> 2207120	Full Length - 2" x 12" <input type="checkbox"/> 2207121 Full Length - 3.5" x 12"
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. \$ 125	
<input type="checkbox"/> 2207130	Right (4.5" x 18") <input type="checkbox"/> 2207131 Left (4.5" x 18")

Armrest Options	
Armrest Wideners Move Armrests outward by 1" on each side. Do not increase seat size.	
<input type="checkbox"/> 2207210	For Single Post Armrests May cause tilt interference at lower Seat Height and/or Seat Depth settings. \$ 195
<input type="checkbox"/> 2207220	For Dual Post Armrests Limits tilt. Not compatible with 20", 22", or 24" Rear Wheels, Clothing Guards. \$ 375
Armrest Extensions Not compatible with Flip-up or Single Post Armrests. \$ 75	
<input type="checkbox"/> 2207400	3" Extension - Forward <input type="checkbox"/> 2207401 3" Extension - Rearward

Seat Pan	
<input type="checkbox"/> 2208000	No Seat Pan
<input type="checkbox"/> 2208001	Solid Aluminum Seat Pan \$ 135

Seat Height



Select desired Seat-to-Floor Height STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low						High							
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low						High							
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2211061 6" x 1.25" Caster size 2211081 8" x 1.25" Caster size
- 2211063 6" x 1.50" Caster size \$ 85
- 2211084 8" x 2.00" Caster size \$ 135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Rear Wheels.

- 2211161 6" x 1.25" Caster size \$ 110
- 2211181 8" x 1.25" Caster size \$ 120
- 2211184 8" x 2.00" Caster size \$ 150

- 2211260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275

- 2211400 **Frog Legs Suspension Forks** Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster. \$ 695

Rear Wheels

Composite Wheels

- 2212012 12" Wheel size 2212016 16" Wheel size 2212020 20" Wheel size
- 2212022 22" Wheel size 2212024 24" Wheel size

Spoke Wheels \$ 75

- 2212120 20" Wheel size 2212122 22" Wheel size 2212124 24" Wheel size
- 2212224 **24" Heavy Duty Spoke Wheels** \$ 350



Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- 2212401 Rear 2212402 Mid 2212403 Front

Axle Type

- 2213000 **Quick Release Rear Axles**
- 2213010 **Fixed Rear Axles**

Rear Tires

- 2214000 **Urethane Tires**
- 2214010 **Pneumatic Tires** May result in a finished Seat Height 1/2" greater than stated. \$ 60
- 2214020 **Pneumatic Tires with Flat-Free Inserts** May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. \$ 95

Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

- 2215000 **No Handrims**
- 2215010 **Standard Handrims**
- 2215020 **Plastic Coated Handrims** \$ 150



Wheel Locks & Brakes

User Wheel Locks

- 2216010 Push to Lock 2216020 Pull to Lock

- 2216030 **Attendant Wheel Locks (Push to Lock)** \$ 85

- 2216040 **User & Attendant Wheel Locks** Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes. \$ 125

- 2216300 **Drum Brakes** Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only. \$ 995

Wheel Lock Options

- 2216100 **Wheel Lock Extensions** Not compatible with Attendant Wheel Lock. \$ 60

- 2216200 **Hill Holders** Requires Pneumatic Rear Tires with Flat-Free Inserts. \$ 450

Anti-Tippers

- 2217010 **Rear Anti-Tippers** Required Safety Option. \$ 75

- 2217020 **Front Anti-Tippers** \$ 795

- 2217030 **Heavy Duty Rear Anti-Tipper** \$ 895

Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.

Back Posts

Back Post Style

- 2218010 Push Bar Back 2218020 Cane Handle Back Not compatible with Headrest Mounting Plate.

Back Post Height

 Measured from the Seat Pan to the Handle.

- 2218420 Short (17") 1518425 Tall (23")

Seat to Back Angle

- 2219090 Set at 90° 2219098 Set at 98° 2219106 Set at 106°

- 2219114 Set at 114° 2219122 Set at 122°

Back Options

- Angle/Height Adjustable Push Bar Handle** Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150

- 2219010 Button on Outside 2219011 Button on Inside

- 2219210 **10" Backrest Extension** Not compatible with Angle/Height Adjustable Push Bar Handle. \$ 195

- 2219212 **Headrest Mounting Plate** \$ 95

Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.

- 2219220 **Fold Down Back**

Back Upholstery

- 2220000 **No Back Upholstery**

- Standard Upholstery** Not available with Cane Handle Back. \$ 265

- 2220418 Short (18" Height) 2220421 Regular (21" Height)

- 2220450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. \$ 365

Position of Upholstery

 Select the desired position of the top of the Upholstery relative to the Seat Pan.

17" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								

23" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

Tilt Activation

Attendant Tilt only Select trigger mounting location below.		
<input type="checkbox"/> 2223011 Right Back Post	<input type="checkbox"/> 2223012 Left Back Post	
User Tilt only Full Length Arm Pads only. Select trigger mounting location below.		\$ 100
<input type="checkbox"/> 2223021 Right Armrest	<input type="checkbox"/> 2223022 Left Armrest	
Attendant & User Tilt Full Length Arm Pads only. Select trigger mounting locations below.		\$ 195
<input type="checkbox"/> 2223031 Right Back Post	<input type="checkbox"/> 2223032 Left Back Post	
<input type="checkbox"/> 2223041 Right Armrest	<input type="checkbox"/> 2223042 Left Armrest	
<input type="checkbox"/> 2223400 Power Tilt 250 lbs Weight Capacity. 30° Anterior to 20° Posterior tilt range. Includes 2 batteries for off-chair charging. Minimum Seat Width is 16".		\$ 2,095
Toggle Switch	<input type="checkbox"/> 2223421 Installed on Right side	<input type="checkbox"/> 2223422 Installed on Left side
Button Switch	<input type="checkbox"/> 2223431 Installed on Right side	<input type="checkbox"/> 2223432 Installed on Left side
Stealth Egg Switch	<input type="checkbox"/> 2223440 Not installed	
Optional for Egg Switch:	<input type="checkbox"/> 2223445 Swing Away Mounting Hardware	Not compatible with all Headrest makes and models. \$ 625
Specify Headrest to be used: _____		
Anterior Tilt Limit Total tilt range is 45°. Posterior Tilt equals 45° less the value selected below.		
<input type="checkbox"/> 2223220 Set at 20°	<input type="checkbox"/> 2223225 Set at 25°	<input type="checkbox"/> 2223230 Set at 30°

Front Rigging

<input type="checkbox"/> 2225000 No Front Rigging	
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out	
<input type="checkbox"/> 2225410 Fixed 70°	\$ 95
<input type="checkbox"/> 2225420 Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°) \$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.	


Footplates

<input type="checkbox"/> 2226010 Composite Footplates	
<input type="checkbox"/> 2226020 Aluminum Angle Adjustable Footplates	\$ 165
<input type="checkbox"/> 2226030 One Piece Fixed Footplate Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/> 2226040 One Piece Flip-up Footplate Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/> 2226050 Interlocking Angle Adjustable Footplates Available for Seat Widths up to 21" only.	\$ 195

Footrest/Legrest Accessories

<input type="checkbox"/> 2227040 Length Adjustable Heel Loops	\$ 50
<input type="checkbox"/> 2227050 Bodypoint Padded Calf Strap	\$ 195

Front Rigging Extension Tubes

Regular Extension Tubes	
<input type="checkbox"/> 2228010 Short 14"-18" for Easy-Out Front Rigging	
<input type="checkbox"/> 2228020 Long 18"-24" for Easy-Out Front Rigging	\$ 75
Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.	
 Seat Pan to Foot Plates: _____	Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.		\$ 235
<input type="checkbox"/> 2230020 Black	<input type="checkbox"/> 2230021 Frame Color	
<input type="checkbox"/> 2230110 IV Pole Not compatible with HD Anti Tipper, Attendant Wheel Locks.		\$ 345
<input type="checkbox"/> 2230120 O2 Holder Not compatible with HD Anti Tipper, Attendant Wheel Locks.		\$ 495
<input type="checkbox"/> 2230170 Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.		\$ 10



STELLAR LEAP

Order Form - Distributor

Prices in Canadian Dollars. January 1, 2022

On-Chair Seating & Positioning



Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.