



BENTLEY LT-R

Order Form - USA

Prices in U.S. Dollars. January 1, 2022

Customer and Order Information

Purchase Order *	Tag	v.0524
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

Base Chair

Bentley LT-R 20° Manual Tilt Wheelchair

- 2701000 Bentley 20° Manual Tilt Wheelchair E1161 \$ 3,150
- 2701009 Lightweight Aluminum Frame Upgrade Reverse Wheel Layout enhances user access for self-propulsion. pdgmobility.com/blt-r K0108 \$ 300
- Choose the standard **Bentley** with a steel frame for agitated or heavier users. See separate order form.

Weight Capacity

- 2702000 **250 lbs Chair Weight Capacity**
- User Weight Range** Required for setup of tilt uplift force.
- 2702430 Up to 75 lbs
- 2702431 75 lbs to 100 lbs
- 2702432 100 lbs to 115 lbs
- 2702433 115 lbs to 150 lbs
- 2702434 150 lbs to 200 lbs
- 2702435 200 lbs to 250 lbs

Frame Options

- 2702300 **Transit Tie Down System (TTDS)**
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths

- 2703014 14" Seat Width
- 2703015 15" Seat Width
- 2703016 16" Seat Width
- 2703017 17" Seat Width
- 2703018 18" Seat Width
- 2703019 19" Seat Width
- 2703020 20" Seat Width
- 2703021 21" Seat Width
- 2703022 22" Seat Width E2201 \$ 375

Seat Depth

Standard Seat Depths (Adjustable)

- 2704016 16" Seat Depth
- 2704017 17" Seat Depth
- 2704018 18" Seat Depth
- 2704019 19" Seat Depth
- 2704020 20" Seat Depth E2203 \$ 500

Custom Seat Depths

- 2704021 21" Seat Depth E2203 \$ 525
- 2704022 22" Seat Depth E2204 \$ 800

Frame Color

PDG Premium Colors

- 2705000 Tracer Red
- 2705001 Gloss Black Shimmer
- 2705002 Silver
- 2705003 Fine Texture Black
- 2705004 Orange Fusion
- 2705005 Karl's Imagination
- 2705006 Blue Shimmer
- 2705007 Gloss White
- 2705008 Dark Titanium
- 2705009 Geyser Green
- 2705010 Lollipop Blue
- 2705011 Lollipop Red
- 2705012 Bubble Gum Pink
- 2705013 Berry Purple
- 2705014 Daffodil Yellow
- 2705200 **Custom Color** Contact PDG for specifications. \$ 300

Armrests

2706000 **No Armrests**

Height Adjustable Armrests (Single Post) E0973 \$ 295

2706111 8"-11" 2706114 10"-14"

2706100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.
 Seat Pan to top of Armrest: _____ Does not include height of Armpad

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
 In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2706401 Rear 2706402 Mid 2706403 Front

Armpads

Regular Armpads

2707100 Desk Length - 2" x 10" 2707101 Full Length - 2" x 14"

Waterfall Armpads \$ 115

2707105 Desk Length - 2" x 10" 2707106 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations. \$ 200

2707110 Desk Length - 3" x 10" 2707111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations. \$ 135

2707120 Full Length - 2" x 12" 2707121 Full Length - 3.5" x 12"

Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations. E2209 \$ 125

2707130 Right (4.5" x 18") 2707131 Left (4.5" x 18")

Armrest Options

2707210 **Armrest Wideners** Move Armrests outward by 1" on each side. Do not increase seat size. \$ 195

2707310 **Clothing Guards** Available for 10"-14" Height range only. \$ 110

Seat Pan

2708000 **No Seat Pan**

2708001 **Solid Aluminum Seat Pan** E2231 \$ 135

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low						High						
	RW	12"	13"	14"	15"	16"	RW	14"	15"	16"	17"	18"	19"
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low						High						
	RW	12"	13"	14"	15"	16"	RW	14"	15"	16"	17"	18"	19"
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

Urethane Casters 8" Casters are compatible only with 16" and 24" Large Wheels.

<input type="checkbox"/>	2711061	6" x 1.25" Caster size	<input type="checkbox"/>	2711081	8" x 1.25" Caster size		
<input type="checkbox"/>	2711063	6" x 1.50" Caster size				\$	85
<input type="checkbox"/>	2711084	8" x 2.00" Caster size				\$	135

Pneumatic Casters 8" Casters are compatible only with 16" and 24" Large Wheels.

<input type="checkbox"/>	2711161	6" x 1.25" Caster size				E2214 & E2215	\$	110
<input type="checkbox"/>	2711181	8" x 1.25" Caster size				E2214 & E2215	\$	120
<input type="checkbox"/>	2711184	8" x 2.00" Caster size				E2214 & E2215	\$	150

<input type="checkbox"/>	2711260	Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters					\$	275
--------------------------	---------	--	--	--	--	--	----	-----

<input type="checkbox"/>	2711400	Frog Legs Suspension Forks	Requires 6" Soft-Roll Caster.			E1015	\$	595
--------------------------	---------	-----------------------------------	-------------------------------	--	--	-------	----	-----

Large Wheels

Composite Wheels

<input type="checkbox"/>	2712020	20" Wheel size	<input type="checkbox"/>	2712022	22" Wheel size	<input type="checkbox"/>	2712024	24" Wheel size
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------

Spoke Wheels

<input type="checkbox"/>	2712120	20" Wheel size	<input type="checkbox"/>	2712122	22" Wheel size	<input type="checkbox"/>	2712124	24" Wheel size	\$	75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	----	----



Large Wheel Position Indicate Large Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/>	2712401	Rear	<input type="checkbox"/>	2712402	Mid	<input type="checkbox"/>	2712403	Front
--------------------------	---------	------	--------------------------	---------	-----	--------------------------	---------	-------

Axle Type

<input type="checkbox"/>	2713000	Quick Release Axles
--------------------------	---------	----------------------------

<input type="checkbox"/>	2713010	Fixed Axles
--------------------------	---------	--------------------

Large Tires

<input type="checkbox"/>	2714000	Urethane Tires					
--------------------------	---------	-----------------------	--	--	--	--	--

<input type="checkbox"/>	2714010	Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated.			E2211 & E2212	\$	60
--------------------------	---------	------------------------	--	--	--	---------------	----	----

<input type="checkbox"/>	2714020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.			E2211 & E2213	\$	95
--------------------------	---------	---	---	--	--	---------------	----	----

Handrims

<input type="checkbox"/>	2715000	No Handrims
--------------------------	---------	--------------------

<input type="checkbox"/>	2715010	Standard Aluminum Handrims
--------------------------	---------	-----------------------------------

<input type="checkbox"/>	2715020	Plastic Coated Handrims	\$	150
--------------------------	---------	--------------------------------	----	-----

Wheel Locks & Brakes

User Wheel Locks

<input type="checkbox"/>	2716010	Push to Lock	<input type="checkbox"/>	2716020	Pull to Lock
--------------------------	---------	--------------	--------------------------	---------	--------------

Wheel Lock Options

<input type="checkbox"/>	2716100	Wheel Lock Extensions		E0961	\$	60
--------------------------	---------	------------------------------	--	-------	----	----

<input type="checkbox"/>	2716200	Hill Holders	Requires Pneumatic Large Tires with Flat-Free Inserts.		E0974	\$	450
--------------------------	---------	---------------------	--	--	-------	----	-----

Anti-Tippers

<input checked="" type="checkbox"/>	2717010	Front Anti-Tippers	Required Safety Option.		E0971	\$	75
-------------------------------------	---------	---------------------------	-------------------------	--	-------	----	----



BENTLEY LT-R

Order Form - USA

Prices in U.S. Dollars. January 1, 2022

Back Posts

Back Post Style

2718010 Push Bar Back

Back Post Height

Measured from the Seat Pan to the Handle.

2718425 Short (22") 2718431 Tall (25")

Seat to Back Angle

Seat Angles above 98° may cause interference at Seat Depths above 19".

2719090 Set at 90° 2719098 Set at 98° 2719106 Set at 106°
 2719114 Set at 114° 2719122 Set at 122°

Back Options

Angle/Height Adjustable Push Bar Handle

Not compatible with Headrest Mounting Plate.

2719010 Button on Outside 2719011 Button on Inside \$ 150

2719210 **10" Backrest Extension** Requires the Angle/Height Adjustable Push Bar Handle. E0966 \$ 150

2719212 **Headrest Mounting Plate** \$ 95
 Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar. Not required if mounting to a solid Back.

2719220 **Fold Down Back** Not compatible with Reclining Back.

Back Upholstery

2720000 **No Back Upholstery**

Standard Upholstery

Not available with Cane Handle Back.

2720418 Short (18" Height) 2720421 Regular (21" Height) 2720424 Tall (24" Height) E2611 / E2612 \$ 195

2720450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. E2611 / E2612 \$ 295

Position of Upholstery

Select the desired position of the top of the Upholstery relative to the Seat Pan.

Style	22" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Style	25" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation

Attendant Tilt only

Select trigger mounting location below.

2723011 Right Back Post 2723012 Left Back Post

User Tilt only

Full Length Arm Pads only. Select trigger mounting location below.

2723021 Right Armrest 2723022 Left Armrest \$ 100

Attendant & User Tilt

Full Length Arm Pads only. Select trigger mounting locations below.

2723031 Right Back Post 2723032 Left Back Post \$ 195
 2723041 Right Armrest 2723042 Left Armrest

2723110 **Tilt Trigger Extension** \$ 50

2723400 **Power Tilt** 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt E1002 \$ 2,295
 Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both. Not available for 12" STF Height.

Toggle Switch 2723421 Installed on Right side 2723422 Installed on Left side
 Button Switch 2723431 Installed on Right side 2723432 Installed on Left side
 Stealth Egg Switch 2723440 Not installed

Optional for Egg Switch: 2723445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. E1028 \$ 625

Specify Headrest to be used: _____

Recline

<input type="checkbox"/>	2723310	Reclining Backrest (Manual)		E1225	\$	895
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.						

Front Rigging

<input type="checkbox"/>	2725000	No Front Rigging				
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out						
<input type="checkbox"/>	2725410	Fixed 70°	Not available for 12" STF Height.			\$ 95
<input type="checkbox"/>	2725420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)	12" STF Height is only compatible with 0°-60°.	K0108	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment.						
<input type="checkbox"/>	2725110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.			


Footplates

<input type="checkbox"/>	2726010	Composite Footplates				
<input type="checkbox"/>	2726020	Aluminum Angle Adjustable Footplates		K0040	\$	165
<input type="checkbox"/>	2726030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.			\$ 350
<input type="checkbox"/>	2726040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.			\$ 395
<input type="checkbox"/>	2726050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.	K0040	\$	195

Footrest/Legrest Accessories

<input type="checkbox"/>	2727012	Calf Pads	Includes reinforced premium Calf Pad mount.			\$ 200
<input type="checkbox"/>	2727032	Gel Knee Pads				\$ 100
<input type="checkbox"/>	2727040	Length Adjustable Heel Loops		E0951	\$	50
<input type="checkbox"/>	2727050	Bodypoint Padded Calf Strap		K0038	\$	195

Front Rigging Extension Tubes

Regular Extension Tubes						
<input type="checkbox"/>	2728010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Standard Plug In Front Rigging		
<input type="checkbox"/>	2728020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Standard Plug In Front Rigging		\$ 75
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.					
Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.						

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Large Wheels. Made from durable and lightweight aluminum. K0065 \$ 195						
<input type="checkbox"/>	2730020	Black	<input type="checkbox"/>	2730021	Frame Color	
<input type="checkbox"/>	2730050	Protective Lower Padding				\$ 795
<input type="checkbox"/>	2730080	Gel Hip Guides	Reduces available Seat Width by 2".	E0956 & E1028	\$	395
<input type="checkbox"/>	2730110	IV Pole		K0105	\$	295
<input type="checkbox"/>	2730120	O2 Holder		E2208	\$	395
<input type="checkbox"/>	2730170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.			\$ 10



BENTLEY LT-R

Order Form - USA

Prices in U.S. Dollars. January 1, 2022

On-Chair Seating & Positioning



Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.