



BENTLEY LT

Order Form - USA

Prices in U.S. Dollars. January 1, 2020

Customer and Order Information

Purchase Order *	Tag	v.1206
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in <i>*required</i>

Frame Type

<input checked="" type="checkbox"/> 2401000	Bentley 20° Manual Tilt Wheelchair	<i>E1161</i>	\$ 2,995
<input checked="" type="checkbox"/> 2401009	Lightweight Aluminum Frame Upgrade (Bentley LT) Aluminum Frame is not recommended for agitated users.	<i>K0108</i>	\$ 300
ⓘ For agitated or heavier users, choose the standard Bentley with steel frame. See separate order form.			

Frame Options

<input checked="" type="checkbox"/> 2402000	Standard Base and Seat Frame	250 lbs Weight Capacity.	
<input type="checkbox"/> 2402300	Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks.	
User Weight Range Required for setup of tilt uplift force.			
<input type="checkbox"/> 2402430	Up to 75 lbs	<input type="checkbox"/> 2402431	75 lbs to 100 lbs
<input type="checkbox"/> 2402432	100 lbs to 115 lbs	<input type="checkbox"/> 2402433	115 lbs to 150 lbs
<input type="checkbox"/> 2402434	150 lbs to 200 lbs	<input type="checkbox"/> 2402435	200 lbs to 250 lbs

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths			
<input type="checkbox"/> 2403014	14" Seat Width	<input type="checkbox"/> 2403015	15" Seat Width
<input type="checkbox"/> 2403016	16" Seat Width	<input type="checkbox"/> 2403017	17" Seat Width
<input type="checkbox"/> 2403018	18" Seat Width	<input type="checkbox"/> 2403019	19" Seat Width
<input type="checkbox"/> 2403020	20" Seat Width	<input type="checkbox"/> 2403021	21" Seat Width
<input type="checkbox"/> 2403022	22" Seat Width	<i>E2201</i>	\$ 350

Seat Depth

Standard Seat Depths (Adjustable)			
<input type="checkbox"/> 2404016	16" Seat Depth	<input type="checkbox"/> 2404017	17" Seat Depth
<input type="checkbox"/> 2404018	18" Seat Depth	<input type="checkbox"/> 2404019	19" Seat Depth
<input type="checkbox"/> 2404020	20" Seat Depth	<i>E2203</i>	\$ 500
Custom Seat Depths			
<input type="checkbox"/> 2404021	21" Seat Depth	<i>E2203</i>	\$ 500
<input type="checkbox"/> 2404022	22" Seat Depth	<i>E2204</i>	\$ 750

Frame Color

PDG Premium Colors			
<input type="checkbox"/> 2405000	Tracer Red	<input type="checkbox"/> 1005001	Gloss Black Shimmer
<input type="checkbox"/> 1005002	Silver	<input type="checkbox"/> 2405003	Fine Texture Black
<input type="checkbox"/> 1005004	Orange Fusion	<input type="checkbox"/> 2405006	Blue Shimmer
<input type="checkbox"/> 1005007	Gloss White	<input type="checkbox"/> 2405009	Geyser Green
<input type="checkbox"/> 1005010	Lollypop Blue	<input type="checkbox"/> 2405012	Bubble Gum Pink
<input type="checkbox"/> 1005011	Lollypop Red	<input type="checkbox"/> 2405200	Custom Color Contact PDG for specifications.
<input type="checkbox"/> 1005014	Daffodil Yellow		\$ 300

Armrest Type

2406000 **No Armrests**

2406111 **Height Adjustable Armrests (Single Post)** E0973 \$ 275

2406111 8"-11" 1006114 10"-14"

2406100 **Fixed Height Armrests (12")** Full Length Armpad only.

Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2406401 Rear 2406402 Mid 2406403 Front

Armpads & Armrest Options

2407100 **Regular Armpads**

2407100 Desk Length - 2" x 10" 1007101 Full Length - 2" x 14"

2407110 **Wide Armpads** May interfere with wheels in some tilt configurations. \$ 200

2407110 Desk Length - 3" x 10" 1007111 Full Length - 3" x 12"

2407120 **Gel Armpads** May interfere with wheels in some tilt configurations. \$ 135

2407120 Full Length - 2" x 12" 1007121 Full Length - 3.5" x 12"

2407130 **Contour Arm Trough with Hand Pad** May interfere with wheels in some tilt configurations. E2209 \$ 125

2407130 Right (4.5" x 18") 1007131 Left (4.5" x 18")

2407310 **Clothing Guards** Available for 10"-14" Height range only. \$ 110

Seat Pan

2408000 **No Seat Pan**

2408001 **Solid Aluminum Seat Pan** E2231 \$ 120

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.

6" Casters	RW	Low					High					Extra-High					
	12"	13"	14"	15"	16"	17"	15"	16"	17"	18"	19"	15"	16"	17"	18"	19"	20"
16"				15"	16"	17"			17"	18"	19"			17"	18"	19"	20"
20"	13"	14"	15"				15"	16"	17"			15"	16"	17"	18"		
22"		14"	15"	16"				16"	17"	18"			16"	17"	18"	19"	
24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"	

8" Casters	RW	Low			High			Extra-High			
	16"	15"	16"	17"	17"	18"	19"	17"	18"	19"	20"
24"		15"	16"	17"	17"	18"	19"	17"	18"	19"	20"

Wheel size _____
Caster size _____
STF height _____
Range _____

Casters

2411061 **Urethane Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.

2411061 6" x 1.25" Caster size 2411081 8" x 1.25" Caster size

2411063 6" x 1.50" Caster size \$ 85

2411084 8" x 2.00" Caster size \$ 135

2411161 **Pneumatic Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels. E2214 & E2215 \$ 110


2411181 6" x 1.25" Caster size E2214 & E2215 \$ 120

2411184 8" x 1.25" Caster size E2214 & E2215 \$ 150

2411260 **Frog Legs 6" x 1.40" Soft Roll Casters** \$ 95

2411400 **Frog Legs Suspension Forks** Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. E1015 \$ 395

Rear Wheels

Composite Wheels			
<input type="checkbox"/> 2412012 12" Wheel size	<input type="checkbox"/> 1012016 16" Wheel size	<input type="checkbox"/> 1012020 20" Wheel size	
<input type="checkbox"/> 2412022 22" Wheel size	<input type="checkbox"/> 1012024 24" Wheel size		
Spoke Wheels			\$ 75
<input type="checkbox"/> 2412120 20" Wheel size	<input type="checkbox"/> 1012122 22" Wheel size	<input type="checkbox"/> 1012124 24" Wheel size	
 Rear Wheel Position Indicate Rear Wheel positioning preference. If left blank, the default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 2412401 Rear	<input type="checkbox"/> 2412402 Mid	<input type="checkbox"/> 2412403 Front	

Axle Type

<input type="checkbox"/> 2413000 Quick Release Rear Axles
<input type="checkbox"/> 2413010 Fixed Rear Axles

Rear Tires

<input type="checkbox"/> 2414000 Urethane Tires			
<input type="checkbox"/> 2414010 Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212	\$ 60
<input type="checkbox"/> 2414020 Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	E2211 & E2213	\$ 95

Hand Rims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/> 2415000 No Handrims	
<input type="checkbox"/> 2415010 Standard Aluminum Handrims	
<input type="checkbox"/> 2415020 Plastic Coated Handrims	\$ 135

Wheel Locks

User Wheel Locks			
<input type="checkbox"/> 2416010 Push to Lock	<input type="checkbox"/> 2416020 Pull to Lock		
<input type="checkbox"/> 2416030 Attendant Wheel Locks (Push to Lock)	Not compatible with TTDS.		\$ 85
<input type="checkbox"/> 2416040 User & Attendant Wheel Locks	Available with 20", 22", or 24" Rear Wheels only.		\$ 125
<input type="checkbox"/> 2416100 Wheel Lock Extensions	Not compatible with Attendant Wheel Lock.	E0961	\$ 60
<input type="checkbox"/> 2416200 Hill Holders	Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974	\$ 450
<input type="checkbox"/> 2416300 Drum Brakes	Available with 20" Rear Wheels only.	E2228	\$ 995

Anti-Tippers

<input checked="" type="checkbox"/> 2417010 Rear Anti-Tippers	Required Safety Option.	E0971	\$ 60
--	-------------------------	-------	-------

Back Post Style

<input checked="" type="checkbox"/> 2418010 Push Bar Back	
Back Post Height Please note: We have updated the descriptions/measurements below. The actual back post heights remain unchanged.	
<input type="checkbox"/> 2418425 Short (22") <i>Previously 25"</i>	<input type="checkbox"/> 2418431 Tall (25") <i>Previously 31"</i> Measured from the Seat Pan to the Handle.

Back Options

Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate.			\$ 150
<input type="checkbox"/> 2419010 Button on Outside	<input type="checkbox"/> 1019011 Button on Inside		
Seat to Back Angle			
<input type="checkbox"/> 2419090 Set at 90°	<input type="checkbox"/> 1019098 Set at 98°	<input type="checkbox"/> 1019106 Set at 106°	
<input type="checkbox"/> 2419114 Set at 114°	<input type="checkbox"/> 1019122 Set at 122°		

See next page for more options >



BENTLEY LT

Order Form - USA

Prices in U.S. Dollars. January 1, 2020

<input type="checkbox"/>	2419210	10" Backrest Extension	Requires the Angle/Height Adjustable Push Bar Handle.	E0966	\$	150
<input type="checkbox"/>	2419212	Headrest Mounting Plate	Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar. Not required if mounting to a solid Back.		\$	95
<input checked="" type="checkbox"/>	2419220	Fold Down Back	Not compatible with Reclining Back.			

Back Upholstery

<input type="checkbox"/>	2420000	No Back Upholstery									
		Standard Upholstery	Not available with Cane Handle Back.	E2611 / E2612	\$	195					
<input type="checkbox"/>	2420418	Short (18" Height)	<input type="checkbox"/>	2420421	Regular (21" Height)	<input type="checkbox"/>	2420424	Tall (24" Height)			
<input type="checkbox"/>	2420450	Tension Adjustable Upholstery	Height of Upholstery is 16". Not available with Cane Handle Back.	E2611 / E2612	\$	295					
		Position of Upholstery	Position of top of Upholstery relative to the Seat Pan. Circle below of select from the drop-down menu on the right.								
		Upholstery	22" Back Post								
		Short	17"	18"	19"	20"	21"	22"			
		Regular				20"	21"	22"			
		Tension Adj			19"	20"	21"	22"			
		Upholstery	25" Back Post								
		Short	17"	18"	19"	20"	21"	22"	23"	24"	25"
		Regular				20"	21"	22"	23"	24"	25"
		Tall							23"	24"	25"
		Tension Adj			19"	20"	21"	22"	23"	24"	25"
									Height	_____	

Tilt/Recline Options

		Attendant Tilt only	Select trigger mounting location below.					
<input type="checkbox"/>	2423011	Right Back Post	<input type="checkbox"/>	2423012	Left Back Post			
		User Tilt only	Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below.			\$ 100		
<input type="checkbox"/>	2423021	Right Armrest	<input type="checkbox"/>	2423022	Left Armrest			
		Attendant & User Tilt	Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below.			\$ 195		
<input type="checkbox"/>	2423031	Right Back Post	<input type="checkbox"/>	2423032	Left Back Post			
<input type="checkbox"/>	2423041	Right Armrest	<input type="checkbox"/>	2423042	Left Armrest			
<input type="checkbox"/>	2423110	Trigger Extension				\$ 50		
<input type="checkbox"/>	2423310	Reclining Backrest (Manual)				E1225 \$ 895		
			30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19".					
<input type="checkbox"/>	2423400	Power Tilt	250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt		E1002	\$ 2,295		
			Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.					
		Toggle Switch	<input type="checkbox"/>	2423421	Installed on Right side	<input type="checkbox"/>	2423422	Installed on Left side
		Button Switch	<input type="checkbox"/>	2423431	Installed on Right side	<input type="checkbox"/>	2423432	Installed on Left side
		Stealth Egg Switch	<input type="checkbox"/>	2423440	Not installed			
		Optional for Egg Switch:	<input type="checkbox"/>	2423445	Swing Away Mounting Hardware	Not compatible with all Headrest makes and models.		E1028 \$ 625
					Specify Headrest to be used: _____			

Front Rigging

<input type="checkbox"/>	2425000	No Front Rigging	<input type="checkbox"/>	2425001	Sector Blocks only	
		Easy-Out™ Front Rigging	Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out			
<input type="checkbox"/>	2425410	Fixed 70°				\$ 95
<input type="checkbox"/>	2425420	Universal Angle Adjustable (0°-110°)	Set at angle: _____	(use increments of 10°)	K0108	\$ 295
			Adjust angle by loosening one bolt. Tool-free leg-length adjustment.			
		Swing Away Front Rigging	Legacy Design. Not compatible with Easy-Out.			
<input type="checkbox"/>	2425130	Swing Away 60°	<input type="checkbox"/>	2425140	Swing Away 70°	
<input type="checkbox"/>	2425200	Swing Away Elevating Legrests	Includes Calf Pads.		E0990	\$ 375
<input type="checkbox"/>	2425300	1" Offset Hangers for Swing Away				\$ 75
<input type="checkbox"/>	2425110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.			


Footplates

<input type="checkbox"/>	2426010	Composite Footplates		
<input type="checkbox"/>	2426020	Aluminum Angle Adjustable Footplates	K0040	\$ 155
<input type="checkbox"/>	2426030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	2426040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	2426050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.	K0040 \$ 195

Footrest/Legrest Accessories

Calf Pads Not compatible with Gel Padding.				
<input type="checkbox"/>	2427012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	2427010	For Swing Away Front Rigging		\$ 135
<input type="checkbox"/>	2427020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pads				
<input type="checkbox"/>	2427032	For Easy-Out™ Angle Adjustable Front Rigging		\$ 200
<input type="checkbox"/>	2427030	For Swing Away ELR Front Rigging		\$ 100
<input type="checkbox"/>	2427040	Length Adjustable Heel Loops		E0951 \$ 50
<input type="checkbox"/>	2427050	Bodypoint Padded Calf Strap		K0038 \$ 195


Front Rigging Extension Tubes

Regular Extension Tubes				
<input type="checkbox"/>	2428010	Short	14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging	
<input type="checkbox"/>	2428020	Long	18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.			
	Seat Pan to Foot Plates: _____		Measured along the length of the Front Rigging.	

Miscellaneous Options

One-Arm Drive Requires 24" Composite Wheels with Urethane Tires. 18" STF minimum. Limits tilt.				
Not compatible with User & Attendant Wheel Locks, Quick Release Rear Axles.				
<input type="checkbox"/>	2430010	Installed on Right side	<input type="checkbox"/> 1030011 Installed on Left side	E0958 \$ 3,500
<input type="checkbox"/>	2430050	Protective Lower Padding		\$ 675
<input type="checkbox"/>	2430080	Gel Hip Guides	Reduces available Seat Width by 2".	E0956 & E1028 \$ 375
<input type="checkbox"/>	2430110	IV Pole	Not compatible with Attendant Wheel Locks.	E1029 \$ 295
<input type="checkbox"/>	2430120	O2 Holder	Not compatible with Attendant Wheel Locks.	E2208 \$ 395
<input type="checkbox"/>	2430170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.	\$ 10

On-Chair Seating & Positioning

	Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating			
---	---	--	--	--

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.