



# BENTLEY LT

Order Form - USA

Prices in U.S. Dollars. January 1, 2022

## Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in
		*required

## Base Chair

### Bentley LT 20° Manual Tilt Wheelchair

<input checked="" type="checkbox"/>	2401000	Bentley 20° Manual Tilt Wheelchair	E1161	\$ 3,150
<input checked="" type="checkbox"/>	2401009	Lightweight Aluminum Frame Upgrade Aluminum Frame is not recommended for agitated users.	K0108	\$ 300
<input type="checkbox"/>	Choose the standard <b>Bentley</b> with a steel frame for agitated or heavier users. See separate order form.			

## Weight Capacity

<input checked="" type="checkbox"/>	2402000	<b>250 lbs Chair Weight Capacity</b>		
<b>User Weight Range</b> Required for setup of tilt uplift force.				
<input type="checkbox"/>	2402430	Up to 75 lbs	<input type="checkbox"/>	2402431 75 lbs to 100 lbs
<input type="checkbox"/>	2402432	100 lbs to 115 lbs	<input type="checkbox"/>	2402433 115 lbs to 150 lbs
<input type="checkbox"/>	2402434	150 lbs to 200 lbs	<input type="checkbox"/>	2402435 200 lbs to 250 lbs

## Frame Options

<input type="checkbox"/>	2402300	<b>Transit Tie Down System (TTDS)</b>
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details.		

## Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

### Standard Seat Widths

<input type="checkbox"/>	2403014	14" Seat Width	<input type="checkbox"/>	2403015	15" Seat Width	<input type="checkbox"/>	2403016	16" Seat Width	
<input type="checkbox"/>	2403017	17" Seat Width	<input type="checkbox"/>	2403018	18" Seat Width	<input type="checkbox"/>	2403019	19" Seat Width	
<input type="checkbox"/>	2403020	20" Seat Width	<input type="checkbox"/>	2403021	21" Seat Width	<input type="checkbox"/>	2403022	22" Seat Width	
								E2201	\$ 375

## Seat Depth

### Standard Seat Depths (Adjustable)

<input type="checkbox"/>	2404016	16" Seat Depth	<input type="checkbox"/>	2404017	17" Seat Depth		
<input type="checkbox"/>	2404018	18" Seat Depth	<input type="checkbox"/>	2404019	19" Seat Depth		
<input type="checkbox"/>	2404020	20" Seat Depth					
						E2203	\$ 500

### Custom Seat Depths

<input type="checkbox"/>	2404021	21" Seat Depth	E2203	\$ 525
<input type="checkbox"/>	2404022	22" Seat Depth	E2204	\$ 800

## Frame Color

### PDG Premium Colors

<input type="checkbox"/>	2405000	Tracer Red	<input type="checkbox"/>	2405001	Gloss Black Shimmer	<input type="checkbox"/>	2405002	Silver
<input type="checkbox"/>	2405003	Fine Texture Black	<input type="checkbox"/>	2405004	Orange Fusion	<input type="checkbox"/>	2405005	Karl's Imagination
<input type="checkbox"/>	2405006	Blue Shimmer	<input type="checkbox"/>	2405007	Gloss White	<input type="checkbox"/>	2405008	Dark Titanium
<input type="checkbox"/>	2405009	Geyser Green	<input type="checkbox"/>	2405010	Lollypop Blue	<input type="checkbox"/>	2405011	Lollypop Red
<input type="checkbox"/>	2405012	Bubble Gum Pink	<input type="checkbox"/>	2405013	Berry Purple	<input type="checkbox"/>	2405014	Daffodil Yellow
<input type="checkbox"/>	2405200	<b>Custom Color</b> Contact PDG for specifications.						\$ 300

## Armrests

2406000 **No Armrests**

**Height Adjustable Armrests (Single Post)** E0973 \$ 295

2406111 8"-11"  2406114 10"-14"

2406100 **Fixed Height Armrests (12")** Full Length Armpad only.

**Armrest Height** If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.  
Seat Pan to top of Armrest: \_\_\_\_\_ Does not include height of Armpad

**Armrest Position** Indicate preference for Armrest position. If left blank, the default setting will be used.  
In the interest of maintaining chair stably and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

2406401 Rear  2406402 Mid  2406403 Front

## Armpads

**Regular Armpads**

2407100 Desk Length - 2" x 10"  2407101 Full Length - 2" x 14" \$ 115

**Waterfall Armpads**

2407105 Desk Length - 2" x 10"  2407106 Full Length - 2" x 14" \$ 200

**Wide Armpads** May interfere with wheels in some tilt configurations.

2407110 Desk Length - 3" x 10"  2407111 Full Length - 3" x 12" \$ 135

**Gel Armpads** May interfere with wheels in some tilt configurations.

2407120 Full Length - 2" x 12"  2407121 Full Length - 3.5" x 12" \$ 125

**Contour Arm Trough with Hand Pad** May interfere with wheels in some tilt configurations.

2407130 Right (4.5" x 18")  2407131 Left (4.5" x 18") E2209 \$ 125

## Armrest Options

2407210 **Armrest Wideners** Move Armrests outward by 1" on each side. Do not increase seat size. \$ 195

2407310 **Clothing Guards** Available for 10"-14" Height range only. \$ 110

## Seat Pan

2408000 **No Seat Pan**

2408001 **Solid Aluminum Seat Pan** E2231 \$ 135

## Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

6" Casters	Low						High						
	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

8" Casters	Low						High						
	RW	13"	14"	15"	16"	17"	RW	15"	16"	17"	18"	19"	20"
16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Casters

#### Urethane Casters

8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2411061	6" x 1.25" Caster size	<input type="checkbox"/>	2411081	8" x 1.25" Caster size		
<input type="checkbox"/>	2411063	6" x 1.50" Caster size					\$ 85
<input type="checkbox"/>	2411084	8" x 2.00" Caster size					\$ 135

#### Pneumatic Casters

8" Casters are compatible only with 16" and 24" Rear Wheels.

<input type="checkbox"/>	2411161	6" x 1.25" Caster size				E2214 & E2215	\$ 110
<input type="checkbox"/>	2411181	8" x 1.25" Caster size				E2214 & E2215	\$ 120
<input type="checkbox"/>	2411184	8" x 2.00" Caster size				E2214 & E2215	\$ 150

<input type="checkbox"/>	2411260	<b>Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters</b>					\$ 275
--------------------------	---------	--	--	--	--	--	--------

<input type="checkbox"/>	2411400	<b>Frog Legs Suspension Forks</b>	Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster.			E1015	\$ 595
--------------------------	---------	-----------------------------------	--	--	--	-------	--------

### Rear Wheels

#### Composite Wheels

<input type="checkbox"/>	2412012	12" Wheel size	<input type="checkbox"/>	2412016	16" Wheel size	<input type="checkbox"/>	2412020	20" Wheel size
<input type="checkbox"/>	2412022	22" Wheel size	<input type="checkbox"/>	2412024	24" Wheel size			

#### Spoke Wheels

<input type="checkbox"/>	2412120	20" Wheel size	<input type="checkbox"/>	2412122	22" Wheel size	<input type="checkbox"/>	2412124	24" Wheel size	\$ 75
--------------------------	---------	----------------	--------------------------	---------	----------------	--------------------------	---------	----------------	-------



#### Rear Wheel Position

Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

<input type="checkbox"/>	2412401	Rear	<input type="checkbox"/>	2412402	Mid	<input type="checkbox"/>	2412403	Front
--------------------------	---------	------	--------------------------	---------	-----	--------------------------	---------	-------

### Axle Type

<input type="checkbox"/>	2413000	<b>Quick Release Rear Axles</b>
<input type="checkbox"/>	2413010	<b>Fixed Rear Axles</b>

### Rear Tires

<input type="checkbox"/>	2414000	<b>Urethane Tires</b>			
<input type="checkbox"/>	2414010	<b>Pneumatic Tires</b>	May result in a finished Seat Height 1/2" greater than stated.	E2211 & E2212	\$ 60
<input type="checkbox"/>	2414020	<b>Pneumatic Tires with Flat-Free Inserts</b>	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	E2211 & E2213	\$ 95

### Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

<input type="checkbox"/>	2415000	<b>No Handrims</b>	
<input type="checkbox"/>	2415010	<b>Standard Aluminum Handrims</b>	
<input type="checkbox"/>	2415020	<b>Plastic Coated Handrims</b>	\$ 150

### Wheel Locks & Brakes

#### User Wheel Locks

<input type="checkbox"/>	2416010	Push to Lock	<input type="checkbox"/>	2416020	Pull to Lock
<input type="checkbox"/>	2416030	<b>Attendant Wheel Locks (Push to Lock)</b>			\$ 85
<input type="checkbox"/>	2416040	<b>User &amp; Attendant Wheel Locks</b>	Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes.		\$ 125
<input type="checkbox"/>	2416300	<b>Drum Brakes</b>	Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.	E2228	\$ 795

### Wheel Lock Options

<input type="checkbox"/>	2416100	<b>Wheel Lock Extensions</b>	Not compatible with Attendant Wheel Lock.	E0961	\$ 60
<input type="checkbox"/>	2416200	<b>Hill Holders</b>	Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974	\$ 450

### Anti-Tippers

<input checked="" type="checkbox"/>	2417010	<b>Rear Anti-Tippers</b>	Required Safety Option.	E0971	\$ 75
-------------------------------------	---------	--------------------------	-------------------------	-------	-------



### Back Posts

#### Back Post Style

2418010 Push Bar Back

#### Back Post Height Measured from the Seat Pan to the Handle.

2418425 Short (22")  2418431 Tall (25")

#### Seat to Back Angle

2419090 Set at 90°  2419098 Set at 98°  2419106 Set at 106°  
 2419114 Set at 114°  2419122 Set at 122°

### Back Options

#### Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate. \$ 150

2419010 Button on Outside  2419011 Button on Inside

2419210 **10" Backrest Extension** Requires the Angle/Height Adjustable Push Bar Handle. *E0966* \$ 150

2419212 **Headrest Mounting Plate** \$ 95  
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar. Not required if mounting to a solid Back.

2419220 **Fold Down Back** Not compatible with Reclining Back.

### Back Upholstery

2420000 **No Back Upholstery**

#### Standard Upholstery Not available with Cane Handle Back. *E2611 / E2612* \$ 195

2420418 Short (18" Height)  2420421 Regular (21" Height)  2420424 Tall (24" Height)

2420450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. *E2611 / E2612* \$ 295

#### Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

Style	22" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

Style	25" Back Post								
	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

### Tilt Activation

#### Attendant Tilt only Select trigger mounting location below.

2423011 Right Back Post  2423012 Left Back Post

#### User Tilt only Full Length Arm Pads only. Select trigger mounting location below. \$ 100

2423021 Right Armrest  2423022 Left Armrest

#### Attendant & User Tilt Full Length Arm Pads only. Select trigger mounting locations below. \$ 195

2423031 Right Back Post  2423032 Left Back Post

2423041 Right Armrest  2423042 Left Armrest

2423110 **Tilt Trigger Extension** \$ 50

2423400 **Power Tilt** 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. [pdgmobility.com/power-tilt](http://pdgmobility.com/power-tilt) *E1002* \$ 2,295  
Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.

Toggle Switch  2423421 Installed on Right side  2423422 Installed on Left side  
 Button Switch  2423431 Installed on Right side  2423432 Installed on Left side  
 Stealth Egg Switch  2423440 Not installed

Optional for Egg Switch:  2423445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. *E1028* \$ 625

*Specify Headrest to be used:* \_\_\_\_\_

### Recline

<input type="checkbox"/>	2423310	<b>Reclining Backrest (Manual)</b>		E1225	\$	895
30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19". Finished Back Post Height is 3" taller than selected Back Post Height.						

### Front Rigging

<input type="checkbox"/>	2425000	<b>No Front Rigging</b>	<input type="checkbox"/>	2425001	<b>Sector Blocks only</b>		
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>							
<input type="checkbox"/>	2425410	Fixed 70°				\$	95
<input type="checkbox"/>	2425420	Universal Angle Adjustable (0°-110°)	Set at angle: _____ (use increments of 10°)			K0108	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.							
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out.							
<input type="checkbox"/>	2425130	Swing Away 60°	<input type="checkbox"/>	2425140	Swing Away 70°		
<input type="checkbox"/>	2425200	Swing Away Elevating Legrests	Includes Calf Pads.			E0990	\$ 375
<input type="checkbox"/>	2425300	1" Offset Hangers for Swing Away					\$ 75
<input type="checkbox"/>	2425110	<b>Standard Plug In 70° Front Rigging</b>	Legacy Design. Not compatible with Easy-Out.				


### Footplates

<input type="checkbox"/>	2426010	<b>Composite Footplates</b>					
<input type="checkbox"/>	2426020	<b>Aluminum Angle Adjustable Footplates</b>				K0040	\$ 165
<input type="checkbox"/>	2426030	<b>One Piece Fixed Footplate</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.				\$ 350
<input type="checkbox"/>	2426040	<b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.				\$ 395
<input type="checkbox"/>	2426050	<b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.			K0040	\$ 195

### Footrest/Legrest Accessories

<b>Calf Pads</b> Not compatible with Gel Padding.							
<input type="checkbox"/>	2427012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.				\$ 200
<input type="checkbox"/>	2427010	For Swing Away Front Rigging					\$ 135
<input type="checkbox"/>	2427020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.				\$ 125
<b>Gel Knee Pads</b>							
<input type="checkbox"/>	2427032	For Easy-Out™ Angle Adjustable Front Rigging					\$ 200
<input type="checkbox"/>	2427030	For Swing Away ELR Front Rigging					\$ 100
<input type="checkbox"/>	2427040	<b>Length Adjustable Heel Loops</b>				E0951	\$ 50
<input type="checkbox"/>	2427050	<b>Bodypoint Padded Calf Strap</b>				K0038	\$ 195

### Front Rigging Extension Tubes

<b>Regular Extension Tubes</b>							
<input type="checkbox"/>	2428010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging			
<input type="checkbox"/>	2428020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging			\$ 75
	<b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.						
Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.							



# BENTLEY LT

## Order Form - USA

Prices in U.S. Dollars. January 1, 2022

### Miscellaneous Options

<b>One-Arm Drive</b> Requires 24" Composite Wheels with Urethane Tires. 18" STF minimum. Limits tilt. Not compatible with User & Attendant Wheel Locks, Drum Brakes, Quick Release Rear Axles.		<i>E0958</i>	\$ 3,500
<input type="checkbox"/> 2430010	Installed on Right side	<input type="checkbox"/> 2430011	Installed on Left side
<b>Spoke Guards</b> Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.		<i>K0065</i>	\$ 195
<input type="checkbox"/> 2430020	Black	<input type="checkbox"/> 2430021	Frame Color
<input type="checkbox"/> 2430050	<b>Protective Lower Padding</b>		\$ 795
<input type="checkbox"/> 2430080	<b>Gel Hip Guides</b> Reduces available Seat Width by 2".	<i>E0956 &amp; E1028</i>	\$ 395
<input type="checkbox"/> 2430110	<b>IV Pole</b> Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later.	<i>K0105</i>	\$ 295
<input type="checkbox"/> 2430120	<b>O2 Holder</b> Not compatible with Attendant Wheel Locks. Requires frame modification and cannot be added to chair later.	<i>E2208</i>	\$ 395
<input type="checkbox"/> 2430170	<b>Adhesive Back Hook Strip, 2" x 30"</b> For Cushion mounting.		\$ 10

### On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

### Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.  
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.