



Customer and Order Information

Purchase Order *	Tag	v.0801
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in <i>*required</i>

Frame Type

<input checked="" type="checkbox"/> 1001000	Bentley 20° Manual Tilt Wheelchair	E1161	\$ 2,995
<input type="checkbox"/> ⓘ	Consider the lightweight Bentley LT with aluminum frame. See separate order form. pdgmobility.com/light		

Frame Options

<input type="checkbox"/> 1002000	Standard Base and Seat Frame	250 lbs Weight Capacity.	
<input type="checkbox"/> 1002050	Upgraded Base and Seat Frame	300 lbs Weight Capacity. Available for Chairs with 14" to 22" Seat Width.	\$ 200
<input type="checkbox"/> 1002100	Heavy Duty Base and Seat Frame	450 lbs Weight Capacity. Available for Chairs with 16" to 22" Seat Width. Required for 23" or wider. Requires 20" Composite or 24" Heavy Duty Spoke Wheels. Fixed Axles only. Not compatible with Transit Tie Down System.	\$ 500
<input type="checkbox"/> 1002200	Extended Base with Rear Floor Jacks	Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.	K0108 \$ 900
<input type="checkbox"/> 1002300	Transit Tie Down System (TTDS)	TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Attendant Wheel Locks, Extended Base, and Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.	
User Weight Range Required for setup of tilt uplift force.			
<input type="checkbox"/> 1002430	Up to 75 lbs	<input type="checkbox"/> 1002431	75 lbs to 100 lbs
<input type="checkbox"/> 1002433	115 lbs to 150 lbs	<input type="checkbox"/> 1002434	150 lbs to 200 lbs
<input type="checkbox"/> 1002436	250 lbs to 300 lbs	<input type="checkbox"/> 1002437	300 lbs to 375 lbs
<input type="checkbox"/> 1002432	100 lbs to 115 lbs	<input type="checkbox"/> 1002435	200 lbs to 250 lbs
<input type="checkbox"/> 1002438	375 lbs to 450 lbs		

Seat Width

Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths			
<input type="checkbox"/> 1003014	14" Seat Width	<input type="checkbox"/> 1003015	15" Seat Width
<input type="checkbox"/> 1003017	17" Seat Width	<input type="checkbox"/> 1003018	18" Seat Width
<input type="checkbox"/> 1003020	20" Seat Width	<input type="checkbox"/> 1003021	21" Seat Width
<input type="checkbox"/> 1003016	16" Seat Width	<input type="checkbox"/> 1003019	19" Seat Width
<input type="checkbox"/> 1003022	22" Seat Width		
			E2201 \$ 350
Custom Seat Widths Requires Heavy Duty Base and Seat Frame. Extended delivery lead times may apply.			
<input type="checkbox"/> 1003023	23" Seat Width		E2201 \$ 350
<input type="checkbox"/> 1003024	24" Seat Width	<input type="checkbox"/> 1003025	25" Seat Width
<input type="checkbox"/> 1003026	26" Seat Width	<input type="checkbox"/> 1003027	27" Seat Width
<input type="checkbox"/> 1003028	28" Seat Width	<input type="checkbox"/> 1003029	29" Seat Width
<input type="checkbox"/> 1003031	31" Seat Width	<input type="checkbox"/> 1003030	30" Seat Width
			E2202 \$ 550
			E2202 \$ 550
			K0108 \$ 1,050
			K0108 \$ 1,050

Seat Depth

Standard Seat Depths (Adjustable)			
<input type="checkbox"/> 1004016	16" Seat Depth	<input type="checkbox"/> 1004017	17" Seat Depth
<input type="checkbox"/> 1004018	18" Seat Depth	<input type="checkbox"/> 1004019	19" Seat Depth
<input type="checkbox"/> 1004020	20" Seat Depth		E2203 \$ 500
Custom Seat Depths			
<input type="checkbox"/> 1004021	21" Seat Depth		E2203 \$ 500
<input type="checkbox"/> 1004022	22" Seat Depth		E2204 \$ 750

Frame Color

PDG Premium Colors

- | | | |
|---|--|---|
| <input type="checkbox"/> 1005000 Tracer Red | <input type="checkbox"/> 1005001 Gloss Black Shimmer | <input type="checkbox"/> 1005002 Silver |
| <input type="checkbox"/> 1005003 Fine Texture Black | <input type="checkbox"/> 1005004 Orange Fusion | <input type="checkbox"/> 1005005 Karl's Imagination |
| <input type="checkbox"/> 1005006 Blue Shimmer | <input type="checkbox"/> 1005007 Gloss White | <input type="checkbox"/> 1005008 Dark Titanium |
| <input type="checkbox"/> 1005009 Geyser Green | <input type="checkbox"/> 1005010 Lollypop Blue | <input type="checkbox"/> 1005011 Lollypop Red |
| <input type="checkbox"/> 1005012 Bubble Gum Pink | <input type="checkbox"/> 1005013 Berry Purple | <input type="checkbox"/> 1005014 Daffodil Yellow |

- 1005200 **Custom Color** Contact PDG for specifications. \$ 300

Armrest Type

- 1006000 **No Armrests**
- 1006111 **Single Post Height Adjustable Armrests** 8"-11" Range 1006114 10"-14" Range E0973 \$ 275
Not compatible with the Heavy Duty Base and Seat Frame.
- 1006211 **Dual Post Height Adjustable Armrests** 8"-11" Range 1006214 10"-14" Range 1006216 12"-16" Range E0973 \$ 295
Required with the Heavy Duty Base and Seat Frame.
- 1006100 **Fixed Height Armrests (12")** Single Post 1006200 Dual Post
Full Length Armpad only.
- Armrest Position** Indicate preference for Armrest position. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.
- 1006401 Rear 1006402 Mid 1006403 Front
- 1006310 **Locking Flip-up Armrests** Set at 10" Height 1006311 Set at 11" Height 1006312 Set at 12" Height E0973 \$ 295
Not recommended for users with excessive strength or tone. Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Requires Regular Full Length Armpads.
- 1006313 Set at 13" Height 1006314 Set at 14" Height

Armpads & Armrest Options

- 1007100 **Regular Armpads** Desk Length - 2" x 10" 1007101 Full Length - 2" x 14"
- 1007110 **Wide Armpads** Desk Length - 3" x 10" 1007111 Full Length - 3" x 12" \$ 200
May interfere with wheels in some tilt configurations.
- 1007120 **Gel Armpads** Full Length - 2" x 12" 1007121 Full Length - 3.5" x 12" \$ 135
May interfere with wheels in some tilt configurations.
- 1007130 **Contour Arm Trough with Hand Pad** Right (4.5" x 18") 1007131 Left (4.5" x 18") E2209 \$ 125
May interfere with wheels in some tilt configurations.
- 1007220 **Armrest Wideners for Dual Post Armrests** \$ 375
Limits tilt. Not compatible with 20", 22" or 24" Rear Wheels, Clothing Guards.
- 1007310 **Clothing Guards** For Single Post Armrests Available for 10"-14" Height range only. \$ 110
Not compatible with Flip-up Armrests, Armrest Wideners.
- 1007320 For Dual Post Armrests \$ 165
- 1007400 **Armrest Extensions** 3" Extension - Forward 1007401 3" Extension - Rearward \$ 75
Not compatible with Flip-up or Single Post Armrests.

Seat Pan

- 1008000 **No Seat Pan**
- 1008001 **Solid Aluminum Seat Pan** Required for 22" Seat Depth. E2231 \$ 120

Seat Height



Select desired Seat-to-Floor Height

STF Height is adjustable within the Range (Low/High/Extra-High) and Rear Wheel size.

6" Casters	RW	Low					High					Extra-High				
	12"	13"	14"	15"			15"	16"	17"			15"	16"	17"	18"	
16"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"
20"	13"	14"	15"			15"	16"	17"			15"	16"	17"	18"		
22"		14"	15"	16"			16"	17"	18"			16"	17"	18"	19"	
24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"

8" Casters	RW	Low					High					Extra-High				
	16"			15"	16"	17"			17"	18"	19"			17"	18"	19"
24"			15"	16"	17"			17"	18"	19"			17"	18"	19"	20"

Wheel size _____

Caster size _____

STF height _____

Range _____

Casters

Urethane Casters

8" Casters are compatible only with 16" and 24" Rear Wheels.

- | | | | | | | | |
|--------------------------|---------|------------------------|--------------------------|---------|------------------------|--|--------|
| <input type="checkbox"/> | 1011061 | 6" x 1.25" Caster size | <input type="checkbox"/> | 1011081 | 8" x 1.25" Caster size | | |
| <input type="checkbox"/> | 1011063 | 6" x 1.50" Caster size | | | | | \$ 85 |
| <input type="checkbox"/> | 1011084 | 8" x 2.00" Caster size | | | | | \$ 135 |

Pneumatic Casters

8" Casters are compatible only with 16" and 24" Rear Wheels.

- | | | | | | | | |
|--------------------------|---------|------------------------|--|--|--|--------------------------|--------|
| <input type="checkbox"/> | 1011161 | 6" x 1.25" Caster size | | | | <i>E2214 & E2215</i> | \$ 110 |
| <input type="checkbox"/> | 1011181 | 8" x 1.25" Caster size | | | | <i>E2214 & E2215</i> | \$ 120 |
| <input type="checkbox"/> | 1011184 | 8" x 2.00" Caster size | | | | <i>E2214 & E2215</i> | \$ 150 |

- | | | | | | | | |
|--------------------------|---------|---|--|--|--|--|-------|
| <input type="checkbox"/> | 1011260 | Frog Legs 6" x 1.40" Soft Roll Casters | | | | | \$ 95 |
|--------------------------|---------|---|--|--|--|--|-------|

- | | | | | | | | |
|--------------------------|---------|-----------------------------------|---|--|--|--------------|--------|
| <input type="checkbox"/> | 1011400 | Frog Legs Suspension Forks | Not compatible with 12", 20", or 22" Rear Wheels. Requires 6" Soft Roll Caster. | | | <i>E1015</i> | \$ 395 |
|--------------------------|---------|-----------------------------------|---|--|--|--------------|--------|

Rear Wheels

Composite Wheels

- | | | | | | | | | |
|--------------------------|---------|----------------|--------------------------|---------|----------------|--------------------------|---------|----------------|
| <input type="checkbox"/> | 1012012 | 12" Wheel size | <input type="checkbox"/> | 1012016 | 16" Wheel size | <input type="checkbox"/> | 1012020 | 20" Wheel size |
| <input type="checkbox"/> | 1012022 | 22" Wheel size | <input type="checkbox"/> | 1012024 | 24" Wheel size | | | |

Spoke Wheels

- | | | | | | | | | | |
|--------------------------|---------|----------------|--------------------------|---------|----------------|--------------------------|---------|----------------|-------|
| <input type="checkbox"/> | 1012120 | 20" Wheel size | <input type="checkbox"/> | 1012122 | 22" Wheel size | <input type="checkbox"/> | 1012124 | 24" Wheel size | \$ 75 |
|--------------------------|---------|----------------|--------------------------|---------|----------------|--------------------------|---------|----------------|-------|

- | | | | | | | | | |
|--------------------------|---------|------------------------------------|--|--|--|--|--|--------|
| <input type="checkbox"/> | 1012224 | 24" Heavy Duty Spoke Wheels | | | | | | \$ 350 |
|--------------------------|---------|------------------------------------|--|--|--|--|--|--------|



Rear Wheel Position

Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- | | | | | | | | | |
|--------------------------|---------|------|--------------------------|---------|-----|--------------------------|---------|-------|
| <input type="checkbox"/> | 1012401 | Rear | <input type="checkbox"/> | 1012402 | Mid | <input type="checkbox"/> | 1012403 | Front |
|--------------------------|---------|------|--------------------------|---------|-----|--------------------------|---------|-------|

Axle Type

- | | | |
|--------------------------|---------|---------------------------------|
| <input type="checkbox"/> | 1013000 | Quick Release Rear Axles |
|--------------------------|---------|---------------------------------|

- | | | |
|--------------------------|---------|-------------------------|
| <input type="checkbox"/> | 1013010 | Fixed Rear Axles |
|--------------------------|---------|-------------------------|

Rear Tires

- | | | | | | |
|--------------------------|---------|-----------------------|--|--|--|
| <input type="checkbox"/> | 1014000 | Urethane Tires | | | |
|--------------------------|---------|-----------------------|--|--|--|

- | | | | | | |
|--------------------------|---------|------------------------|--|--------------------------|-------|
| <input type="checkbox"/> | 1014010 | Pneumatic Tires | May result in a finished Seat Height 1/2" greater than stated. | <i>E2211 & E2212</i> | \$ 60 |
|--------------------------|---------|------------------------|--|--------------------------|-------|

- | | | | | | |
|--------------------------|---------|---|---|--------------------------|-------|
| <input type="checkbox"/> | 1014020 | Pneumatic Tires with Flat-Free Inserts | May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. | <i>E2211 & E2213</i> | \$ 95 |
|--------------------------|---------|---|---|--------------------------|-------|

Handrims

Handrims are only available with 20", 22", and 24" Rear Wheels.

- | | | |
|--------------------------|---------|--------------------|
| <input type="checkbox"/> | 1015000 | No Handrims |
|--------------------------|---------|--------------------|

- | | | |
|--------------------------|---------|-----------------------------------|
| <input type="checkbox"/> | 1015010 | Standard Aluminum Handrims |
|--------------------------|---------|-----------------------------------|

- | | | | | | |
|--------------------------|---------|--------------------------------|--|--|--------|
| <input type="checkbox"/> | 1015020 | Plastic Coated Handrims | | | \$ 135 |
|--------------------------|---------|--------------------------------|--|--|--------|

Wheel Locks

User Wheel Locks		
<input type="checkbox"/>	1016010 Push to Lock	<input type="checkbox"/> 1016020 Pull to Lock
<input type="checkbox"/>	1016030 Attendant Wheel Locks (Push to Lock)	Not compatible with TTDS or with Upgraded or Heavy Duty Base and Seat Frame.
<input type="checkbox"/>	1016040 User & Attendant Wheel Locks	Available with 20", 22", or 24" Rear Wheels only.
<input type="checkbox"/>	1016100 Wheel Lock Extensions	Not compatible with Attendant Wheel Lock. E0961
<input type="checkbox"/>	1016200 Hill Holders	Requires Pneumatic Rear Tires with Flat-Free Inserts. E0974
<input type="checkbox"/>	1016300 Drum Brakes	Available with 20" Rear Wheels only. COMING SOON E2228

Anti-Tippers

<input checked="" type="checkbox"/>	1017010 Rear Anti-Tippers	Required Safety Option. E0971	\$ 60
<input type="checkbox"/>	1017020 Front Anti-Tippers		E0971 \$ 795
<input type="checkbox"/>	1017030 Heavy Duty Rear Anti-Tipper	Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, Rear Floor Jacks, IV Pole, O2 Holder. E0971	\$ 895

Back Post Style

<input type="checkbox"/>	1018010 Push Bar Back	
<input type="checkbox"/>	1018020 Cane Handle Back	Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.
Back Post Height Please note: We have updated the descriptions/measurements below. The actual back post heights remain unchanged.		
<input type="checkbox"/>	1018425 Short (22") <i>Previously 25"</i>	<input type="checkbox"/> 1018431 Tall (25") <i>Previously 31"</i> Measured from the Seat Pan to the Handle.

Back Options

Angle/Height Adjustable Push Bar Handle		Not compatible with Headrest Mounting Plate or Cane Handle Back.	\$ 150
<input type="checkbox"/>	1019010 Button on Outside	<input type="checkbox"/> 1019011 Button on Inside	
Seat to Back Angle			
<input type="checkbox"/>	1019090 Set at 90°	<input type="checkbox"/> 1019098 Set at 98°	<input type="checkbox"/> 1019106 Set at 106°
<input type="checkbox"/>	1019114 Set at 114°	<input type="checkbox"/> 1019122 Set at 122°	
<input type="checkbox"/>	1019210 10" Backrest Extension	Requires the Angle/Height Adjustable Push Bar Handle. E0966	\$ 150
<input type="checkbox"/>	1019212 Headrest Mounting Plate	Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.	\$ 95
<input type="checkbox"/>	1019220 Fold Down Back	Not compatible with Reclining Back.	
<input checked="" type="checkbox"/>	1019230 Spring Back Suspension	Not compatible with Reclining Back.	
<input type="checkbox"/>	1019240 Seating Dynamics Dynamic Rocker Back Interface	Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle. Requires Adjustable Tension Upholstery or rigid Backrest. Finished Back Post Height 3" taller than selected Back Post Height. K0108	\$ 2,995

Back Upholstery

<input type="checkbox"/>	1020000 No Back Upholstery		
Standard Upholstery		Not available with Cane Handle Back. E2611 / E2612	
<input type="checkbox"/>	1020418 Short (18" Height)	<input type="checkbox"/> 1020421 Regular (21" Height)	<input type="checkbox"/> 1020424 Tall (24" Height)
<input type="checkbox"/>	1020450 Tension Adjustable Upholstery	Height of Upholstery is 16". Not available with Cane Handle Back. E2611 / E2612	\$ 295

Position of Upholstery Position of top of Upholstery relative to the Seat Pan. Circle below of select from the drop-down menu on the right.

Upholstery	22" Back Post							
Short	17"	18"	19"	20"	21"	22"		
Regular				20"	21"	22"		
Tension Adj			19"	20"	21"	22"		

Upholstery	25" Back Post								
Short	17"	18"	19"	20"	21"	22"	23"	24"	25"
Regular				20"	21"	22"	23"	24"	25"
Tall							23"	24"	25"
Tension Adj			19"	20"	21"	22"	23"	24"	25"

Height _____

Tilt/Recline Options

Attendant Tilt only Select trigger mounting location below.			
<input type="checkbox"/> 1023011	Right Back Post	<input type="checkbox"/> 1023012	Left Back Post
User Tilt only Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below.			\$ 100
<input type="checkbox"/> 1023021	Right Armrest	<input type="checkbox"/> 1023022	Left Armrest
Attendant & User Tilt Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below.			\$ 195
<input type="checkbox"/> 1023031	Right Back Post	<input type="checkbox"/> 1023032	Left Back Post
<input type="checkbox"/> 1023041	Right Armrest	<input type="checkbox"/> 1023042	Left Armrest
<input type="checkbox"/> 1023110	Trigger Extension		\$ 50
<input type="checkbox"/> 1023310	Reclining Backrest (Manual) 30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery, Dynamic Rocker Back. Requires longer Frame. Maximum Seat Depth is 19".		E1225 \$ 895
<input type="checkbox"/> 1023400	Power Tilt 250 lbs Weight Capacity. 20° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt Minimum Seat Width is 15". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.		E1002 \$ 2,295
	Toggle Switch	<input type="checkbox"/> 1023421	Installed on Right side
	Button Switch	<input type="checkbox"/> 1023431	Installed on Right side
	Stealth Egg Switch	<input type="checkbox"/> 1023440	Not installed
	Optional for Egg Switch:	<input type="checkbox"/> 1023445	Swing Away Mounting Hardware Not compatible with all Headrest makes and models.
			E1028 \$ 625
<i>Specify Headrest to be used: _____</i>			

Front Rigging

<input type="checkbox"/> 1025000	No Front Rigging	<input type="checkbox"/> 1025001	Sector Blocks only	
Easy-Out™ Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out				
<input type="checkbox"/> 1025410	Fixed 70°			\$ 95
<input type="checkbox"/> 1025420	Universal Angle Adjustable (0°-110°)	Set at angle: _____	(use increments of 10°)	\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment.				
Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out.				
<input type="checkbox"/> 1025130	Swing Away 60°	<input type="checkbox"/> 1025140	Swing Away 70°	
<input type="checkbox"/> 1025200	Swing Away Elevating Legrests	Includes Calf Pads.		E0990 \$ 375
<input type="checkbox"/> 1025300	1" Offset Hangers for Swing Away			\$ 75
<input type="checkbox"/> 1025110	Standard Plug In 70° Front Rigging	Legacy Design. Not compatible with Easy-Out.		

Footplates

<input type="checkbox"/> 1026010	Composite Footplates		
<input type="checkbox"/> 1026020	Aluminum Angle Adjustable Footplates		K0040 \$ 155
<input type="checkbox"/> 1026030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/> 1026040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/> 1026050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.	K0040 \$ 195

Footrest/Legrest Accessories

Calf Pads Not compatible with Gel Padding.		
<input type="checkbox"/> 1027012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount. \$ 200
<input type="checkbox"/> 1027010	For Swing Away Front Rigging	\$ 135
<input type="checkbox"/> 1027020	Gel Padding Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pads		
<input type="checkbox"/> 1027032	For Easy-Out™ Angle Adjustable Front Rigging	\$ 200
<input type="checkbox"/> 1027030	For Swing Away ELR Front Rigging	\$ 100

See next page for more options >

<input type="checkbox"/>	1027040	Length Adjustable Heel Loops	<i>E0951</i>	\$ 50
<input type="checkbox"/>	1027050	Bodypoint Padded Calf Strap Sizes greater than 22" are not padded.	<i>K0038</i>	\$ 195
Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.				
<input type="checkbox"/>	1027060	For Standard Plug In Front Rigging		\$ 445
<input type="checkbox"/>	1027065	For Swing Away 70° Front Rigging Not available on 13" STF Height.		\$ 575

Front Rigging Extension Tubes

Regular Extension Tubes

<input type="checkbox"/>	1028010	Short 14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging		
<input type="checkbox"/>	1028020	Long 18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging		\$ 75

Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only. Not compatible with Easy-Out Front Rigging.

<input type="checkbox"/>	1028110	Standard (13"-17")	<input type="checkbox"/>	1028120	Long (15"-19")	\$ 100
--------------------------	---------	--------------------	--------------------------	---------	----------------	--------



Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging.

Miscellaneous Options

One-Arm Drive Maximum Seat Width of 22". 24" Composite Wheels with Urethane Tires required. 18" STF minimum. Limits tilt. *E0958* \$ 3,500

Not compatible with User & Attendant Wheel Locks, Quick Release Rear Axles.

<input type="checkbox"/>	1030010	Installed on Right side	<input type="checkbox"/>	1030011	Installed on Left side	
--------------------------	---------	-------------------------	--------------------------	---------	------------------------	--

<input type="checkbox"/>	1030050	Protective Lower Padding		\$ 675
--------------------------	---------	---------------------------------	--	--------

<input type="checkbox"/>	1030080	Gel Hip Guides Reduces available Seat Width by 2".	<i>E0956 & E1028</i>	\$ 375
--------------------------	---------	---	--------------------------	--------

<input type="checkbox"/>	1030110	IV Pole Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	<i>E1029</i>	\$ 295
--------------------------	---------	---	--------------	--------

<input type="checkbox"/>	1030120	O2 Holder Not compatible with HD Anti Tipper, Rear Floor Jacks, Attendant Wheel Locks.	<i>E2208</i>	\$ 395
--------------------------	---------	---	--------------	--------

<input type="checkbox"/>	1030170	Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.		\$ 10
--------------------------	---------	---	--	-------

On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.