



# FUZE T50n

**Order Form - USA**

Prices in U.S. Dollars. January 1, 2022

## Customer and Order Information

Purchase Order *	Tag	v.0101
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in <span style="float: right;">*required</span>

## Base Chair

**Fuze T50n 50° Manual Tilt Wheelchair**

<input checked="" type="checkbox"/> 2601000	Fuze T50 50° Manual Tilt Wheelchair	E1161	\$ 3,150
<input checked="" type="checkbox"/> 2601008	Narrow Frame Upgrade Reduces overall chair width to seat width plus only 4".	K0108	\$ 200
<b>i</b> The standard <b>Fuze T50</b> offers additional options that are not available on the T50n. See separate order form.			

## Frame Type & Weight Capacity

**Fixed Width Frame** Select a Fixed Seat Width below.

<input type="checkbox"/> 2602010	300 lbs Weight Capacity		
<input type="checkbox"/> 2602011	400 lbs Weight Capacity Not compatible with Transit Tie Down System.		
<input type="checkbox"/> 2602025	<b>Adjustable Width Frame</b> 250 lbs Weight Capacity. Select an Adjustable Seat Width below.		\$ 100

## Frame Options

<input type="checkbox"/> 2602300	<b>Transit Tie Down System (TTDS)</b>
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 300 lbs.	

## Seat Width Overall chair width is +4.0" for 12" RW

**Fixed Seat Widths** Order with Fixed Frame.

<input type="checkbox"/> 2603016	16" Seat Width	<input type="checkbox"/> 2603017	17" Seat Width		
<input type="checkbox"/> 2603018	18" Seat Width	<input type="checkbox"/> 2603019	19" Seat Width		
<input type="checkbox"/> 2603020	20" Seat Width	<input type="checkbox"/> 2603021	21" Seat Width	E2201	\$ 375
<input type="checkbox"/> 2603022	22" Seat Width	<input type="checkbox"/> 2603023	23" Seat Width	E2201	\$ 375
<input type="checkbox"/> 2603024	24" Seat Width			E2202	\$ 575

**Adjustable Seat Widths** Order with Adjustable Frame.

<input type="checkbox"/> 2603216	16"-19" Range - Set at 16"	<input type="checkbox"/> 2603217	16"-19" Range - Set at 17"		
<input type="checkbox"/> 2603218	16"-19" Range - Set at 18"	<input type="checkbox"/> 2603219	16"-19" Range - Set at 19"		
<input type="checkbox"/> 2603319	19"-20" Range - Set at 19"				
<input type="checkbox"/> 2603320	19"-20" Range - Set at 20"			E2201	\$ 375

## Seat Depth

**Standard Seat Depths (Adjustable)**

<input type="checkbox"/> 2604015	15" Seat Depth	<input type="checkbox"/> 2604016	16" Seat Depth	<input type="checkbox"/> 2604017	17" Seat Depth		
<input type="checkbox"/> 2604018	18" Seat Depth	<input type="checkbox"/> 2604019	19" Seat Depth				
<input type="checkbox"/> 2604020	20" Seat Depth					E2203	\$ 500

**Custom Seat Depths** Requires Push Bar Back - Straight. Not compatible with Reclining or Fold Down Back.

<input type="checkbox"/> 2604021	21" Seat Depth					E2203	\$ 525
<input type="checkbox"/> 2604022	22" Seat Depth					E2204	\$ 800

### Frame Color

#### PDG Premium Colors

- |   |  |   |
|---|--|---|
| <input type="checkbox"/> 2605000 Tracer Red         | <input type="checkbox"/> 2605001 Gloss Black Shimmer | <input type="checkbox"/> 2605002 Silver             |
| <input type="checkbox"/> 2605003 Fine Texture Black | <input type="checkbox"/> 2605004 Orange Fusion       | <input type="checkbox"/> 2605005 Karl's Imagination |
| <input type="checkbox"/> 2605006 Blue Shimmer       | <input type="checkbox"/> 2605007 Gloss White         | <input type="checkbox"/> 2605008 Dark Titanium      |
| <input type="checkbox"/> 2605009 Geyser Green       | <input type="checkbox"/> 2605010 Lollipop Blue       | <input type="checkbox"/> 2605011 Lollipop Red       |
| <input type="checkbox"/> 2605012 Bubble Gum Pink    | <input type="checkbox"/> 2605013 Berry Purple        | <input type="checkbox"/> 2605014 Daffodil Yellow    |

2605200 **Custom Color** Contact PDG for specifications.

\$ 300

### Armrests

2606000 **No Armrests**

2606114 **Single Post Height Adjustable Armrests (10"-14")**

E0973 \$ 295



**Armrest Height** If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.

Seat Pan to top of Armrest: \_\_\_\_\_ Does not include height of Armpad



**Armrest Position** Indicate preference for Armrest position. If left blank, the default setting will be used.

- 2606402 Mid  2606403 Front

**Locking Flip-up Armrests** Not recommended for users with excessive strength or tone.

E0973 \$ 295

Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Arm pads.

- |   |   |   |
|---|---|---|
| <input type="checkbox"/> 2606310 Set at 10" | <input type="checkbox"/> 2606311 Set at 11" | <input type="checkbox"/> 2606312 Set at 12" |
| <input type="checkbox"/> 2606313 Set at 13" | <input type="checkbox"/> 2606314 Set at 14" |   |

### Arm pads

#### Regular Arm pads

- 2607100 Desk Length - 2" x 10"  2607101 Full Length - 2" x 14"

#### Waterfall Arm pads

- 2607105 Desk Length - 2" x 10"  2607106 Full Length - 2" x 14"

\$ 115

**Wide Arm pads** May interfere with wheels in some tilt configurations.

\$ 200

- 2607110 Desk Length - 3" x 10"  2607111 Full Length - 3" x 12"

**Gel Arm pads** May interfere with wheels in some tilt configurations.

\$ 135

- 2607120 Full Length - 2" x 12"  2607121 Full Length - 3.5" x 12"

#### Contour Arm Trough with Hand Pad

E2209 \$ 125

Not compatible with Lap Tray, Locking Flip-up Armrests. May interfere with wheels in some tilt configurations.

- 2607130 Right (4.5" x 18")  2607131 Left (4.5" x 18")

### Armrest Options

2607210 **Armrest Wideners** Move Armrests outward by 1" on each side. Do not increase seat size. May increase overall width of the chair.

\$ 195

### Seat Pan

2608000 **No Seat Pan**

2608001 **Solid Aluminum Seat Pan**

E2231 \$ 135

### Seat Height

**Don't forget** **Select desired Seat-to-Floor Height** STF Height is adjustable.

12" Rear Wheels					
Caster size	16"	17"	18"	19"	20"
5"	<input type="checkbox"/>	<input type="checkbox"/>			
6"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
8"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
6" Soft Roll	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

### Casters

#### Urethane Casters

- |                                  |                        |                                  |                        |                                  |                        |        |
|----------------------------------|------------------------|----------------------------------|------------------------|----------------------------------|------------------------|--------|
| <input type="checkbox"/> 2611051 | 5" x 1.25" Caster size | <input type="checkbox"/> 2611061 | 6" x 1.25" Caster size | <input type="checkbox"/> 2611081 | 8" x 1.25" Caster size |        |
| <input type="checkbox"/> 2611063 | 6" x 1.50" Caster size |                                  |                        |                                  |                        | \$ 85  |
| <input type="checkbox"/> 2611084 | 8" x 2.00" Caster size |                                  |                        |                                  |                        | \$ 135 |

#### Pneumatic Casters

- |                                  |                        |  |  |               |        |
|----------------------------------|------------------------|--|--|---------------|--------|
| <input type="checkbox"/> 2611161 | 6" x 1.25" Caster size |  |  | E2214 & E2215 | \$ 110 |
| <input type="checkbox"/> 2611181 | 8" x 1.25" Caster size |  |  | E2214 & E2215 | \$ 120 |
| <input type="checkbox"/> 2611184 | 8" x 2.00" Caster size |  |  | E2214 & E2215 | \$ 150 |

- |                                  |  |  |  |  |        |
|----------------------------------|--|--|--|--|--------|
| <input type="checkbox"/> 2611260 | <b>Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters</b> |  |  |  | \$ 275 |
|----------------------------------|--|--|--|--|--------|

- |                                  |                                   |  |  |       |        |
|----------------------------------|-----------------------------------|--|--|-------|--------|
| <input type="checkbox"/> 2611400 | <b>Frog Legs Suspension Forks</b> | Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster. |  | E1015 | \$ 595 |
|----------------------------------|-----------------------------------|--|--|-------|--------|

### Rear Wheels

- |   |                             |  |  |  |  |
|---|-----------------------------|--|--|--|--|
| <input checked="" type="checkbox"/> 2612012 | <b>12" Composite Wheels</b> |  |  |  |  |
|---|-----------------------------|--|--|--|--|



#### Rear Wheel Position

Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- |                                  |                  |                                  |     |                                  |                   |
|----------------------------------|------------------|----------------------------------|-----|----------------------------------|-------------------|
| <input type="checkbox"/> 2612401 | Rear More stable | <input type="checkbox"/> 2612402 | Mid | <input type="checkbox"/> 2612403 | Front Less stable |
|----------------------------------|------------------|----------------------------------|-----|----------------------------------|-------------------|

### Axle Type

- |                                  |                                 |  |
|----------------------------------|---------------------------------|--|
| <input type="checkbox"/> 2613000 | <b>Quick Release Rear Axles</b> |  |
|----------------------------------|---------------------------------|--|

- |                                  |                         |  |
|----------------------------------|-------------------------|--|
| <input type="checkbox"/> 2613010 | <b>Fixed Rear Axles</b> |  |
|----------------------------------|-------------------------|--|

### Rear Tires

- |                                  |                       |  |  |
|----------------------------------|-----------------------|--|--|
| <input type="checkbox"/> 2614000 | <b>Urethane Tires</b> |  |  |
|----------------------------------|-----------------------|--|--|

- |                                  |                        |  |               |       |
|----------------------------------|------------------------|--|---------------|-------|
| <input type="checkbox"/> 2614010 | <b>Pneumatic Tires</b> | May result in a finished Seat Height 1/2" greater than stated. | E2211 & E2212 | \$ 60 |
|----------------------------------|------------------------|--|---------------|-------|

- |                                  |   |  |               |       |
|----------------------------------|---|--|---------------|-------|
| <input type="checkbox"/> 2614020 | <b>Pneumatic Tires with Flat-Free Inserts</b> | May result in a finished Seat Height 1/2" greater than stated. | E2211 & E2213 | \$ 95 |
|----------------------------------|---|--|---------------|-------|

### Handrims

- |   |                    |  |
|---|--------------------|--|
| <input checked="" type="checkbox"/> 2615000 | <b>No Handrims</b> |  |
|---|--------------------|--|

### Wheel Locks

- |                                  |   |  |
|----------------------------------|---|--|
| <input type="checkbox"/> 2616030 | <b>Attendant Wheel Locks (Push to Lock)</b> |  |
|----------------------------------|---|--|

- |                                  |                                 |   |        |
|----------------------------------|---------------------------------|---|--------|
| <input type="checkbox"/> 2616005 | <b>Foot Operated Wheel Lock</b> | Rear-mounted unilateral single-activation lock. Compatible only with 19" & 20" STF Height. Not compatible with Ventilator Tray. | \$ 595 |
|----------------------------------|---------------------------------|---|--------|

### Wheel Lock Options

- |                                  |                              |       |       |
|----------------------------------|------------------------------|-------|-------|
| <input type="checkbox"/> 2616100 | <b>Wheel Lock Extensions</b> | E0961 | \$ 60 |
|----------------------------------|------------------------------|-------|-------|

### Anti-Tippers

- |   |                          |                         |       |       |
|---|--------------------------|-------------------------|-------|-------|
| <input checked="" type="checkbox"/> 2617010 | <b>Rear Anti-Tippers</b> | Required Safety Option. | E0971 | \$ 75 |
|---|--------------------------|-------------------------|-------|-------|

### Back Posts

#### Back Post Style

- |                                  |               |                                  |                  |  |
|----------------------------------|---------------|----------------------------------|------------------|--|
| <input type="checkbox"/> 2618100 | Push Bar Back | <input type="checkbox"/> 2618105 | Cane Handle Back | Not compatible with Headrest Mounting Plate. |
|----------------------------------|---------------|----------------------------------|------------------|--|

#### Back Post Type

- |                                  |          |                                  |         |
|----------------------------------|----------|----------------------------------|---------|
| <input type="checkbox"/> 2618200 | Straight | <input type="checkbox"/> 2618205 | 8° Bend |
|----------------------------------|----------|----------------------------------|---------|

#### Back Post Height

Measured from the Seat Pan to the Handle.

- |                                  |             |                                  |            |
|----------------------------------|-------------|----------------------------------|------------|
| <input type="checkbox"/> 2618420 | Short (17") | <input type="checkbox"/> 2618425 | Tall (23") |
|----------------------------------|-------------|----------------------------------|------------|

See next page for more options >



# FUZE T50n

**Order Form - USA**

Prices in U.S. Dollars. January 1, 2022

### Seat to Back Angle

- 2619090 Set at 90°     
  2619096 Set at 98°     
  2619102 Set at 106°  
 2619108 Set at 114°     
  2619114 Set at 122°

### Back Options

- Angle/Height Adjustable Push Bar Handle** Not compatible with Headrest Mounting Plate or Cane Handle Back. \$ 150  
 2619010 Button on Outside       2619011 Button on Inside  
 2619050 **Removable Cane Handle Push Bar** Not compatible with Push Bar Back. \$ 140  
 2619212 **Headrest Mounting Plate** Not compatible with Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back. \$ 95  
 2619220 **Fold Down Back**       2619221 **Non-Folding Back**  
 2619240 **Seating Dynamics Dynamic Rocker Back Interface** K0108 \$ 2,995  
 Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle, Standard Upholstery.  
 Finished Back Post Height is 5" taller than selected Back Post Height.

### Back Upholstery

- 2620000 **No Back Upholstery**  
 **Standard Upholstery** Not available with Cane Handle Back. E2611 / E2612 \$ 195  
 2620418 Short (18" Height)       2620421 Regular (21" Height)  
 2620450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. E2611 / E2612 \$ 295

**Position of Upholstery** Select the desired position of the top of the Upholstery relative to the Seat Pan.

17" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								

23" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		

### Tilt Activation

- Attendant Tilt only** Select trigger mounting location below.  
 2623011 Right Back Post       2623012 Left Back Post  
 **Power Tilt** 40° tilt range. Includes 2 batteries for off-chair charging. [pdgmobility.com/power-tilt](http://pdgmobility.com/power-tilt)  
 2623400 250 lbs Weight Capacity E1002 \$ 2,295  
 2623401 350 lbs Weight Capacity Requires Fixed Frame. E1002 \$ 2,595  
 Toggle Switch       2623421 Installed on Right side       2623422 Installed on Left side  
 Button Switch       2623431 Installed on Right side       2623432 Installed on Left side  
 Stealth Egg Switch       2623440 Not installed  
 Optional for Egg Switch:  2623445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. E1028 \$ 625  
 Specify Headrest to be used: \_\_\_\_\_

**Tilt Limit** Optional. If no tilt limit is selected, the full tilt range will be available.

- 2623220 Set at 20°     
  2623230 Set at 30°     
  2623240 Set at 40°

### Recline

- 2623310 **Reclining Backrest (Manual)** E1225 \$ 895  
 30° maximum recline. Not compatible with Cane Handle Back, Fold Down Back. Not available with 22" Seat Depth.  
 Finished Back Post Height is 3" taller than selected Back Post Height.

### Front Rigging

<input type="checkbox"/>	2625000	<b>No Front Rigging</b>	<input type="checkbox"/>	2625001	<b>Sector Blocks only</b>		
<b>Easy-Out™ Front Rigging</b> Effortless insertion and removal. No Hanger Bracket remains on frame. <a href="http://pdgmobility.com/easy-out">pdgmobility.com/easy-out</a>							
<input type="checkbox"/>	2625410	Fixed 70°					\$ 95
<input type="checkbox"/>	2625420	Universal Angle Adjustable (0°-110°)	Set at angle: _____	(use increments of 10°)			\$ 295
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.							
<b>Swing Away Front Rigging</b> Legacy Design. Not compatible with Easy-Out. Requires pre-tilt of 5°.							
<input type="checkbox"/>	2625140	Swing Away 70°					
<input type="checkbox"/>	2625160	Swing Away 90°	Includes 11"-16" Extension Tubes. Not compatible with 8" Casters.				\$ 75
<input type="checkbox"/>	2625200	Swing Away Elevating Legrests	Includes Calf Pads.			E0990	\$ 375
<input type="checkbox"/>	2625300	1" Offset Hangers for Swing Away					\$ 75
<input type="checkbox"/>	2625110	<b>Standard Plug In 70° Front Rigging</b>	Legacy Design. Not compatible with Easy-Out.				


### Footplates

<input type="checkbox"/>	2626010	<b>Composite Footplates</b>					
<input type="checkbox"/>	2626020	<b>Aluminum Angle Adjustable Footplates</b>				K0040	\$ 165
<input type="checkbox"/>	2626030	<b>One Piece Fixed Footplate</b>	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.				\$ 350
<input type="checkbox"/>	2626040	<b>One Piece Flip-up Footplate</b>	Requires selection of Composite or Aluminum Angle Adjustable.				\$ 395
<input type="checkbox"/>	2626050	<b>Interlocking Angle Adjustable Footplates</b>	Available for Seat Widths up to 21" only.			K0040	\$ 195

### Footrest/Legrest Accessories

<b>Calf Pad</b> Not compatible with Gel Padding.							
<input type="checkbox"/>	2627012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.				\$ 200
<input type="checkbox"/>	2627010	For Swing Away Front Rigging					\$ 135
<input type="checkbox"/>	2627020	<b>Gel Padding</b>	Not compatible with Calf Pad or Easy-Out Front Rigging.				\$ 125
<b>Gel Knee Pad</b>							
<input type="checkbox"/>	2627032	For Easy-Out™ Angle Adjustable Front Rigging					\$ 200
<input type="checkbox"/>	2627030	For Swing Away ELR Front Rigging					\$ 100
<input type="checkbox"/>	2627040	<b>Length Adjustable Heel Loops</b>				E0951	\$ 50
<input type="checkbox"/>	2627050	<b>Bodypoint Padded Calf Strap</b>	Sizes greater than 22" are not padded.			K0038	\$ 195

### Front Rigging Extension Tubes

<b>Regular Extension Tubes</b> No selection required for 90° Front Rigging							
<input type="checkbox"/>	2628010	Short	14"-18" for Easy-Out Front Rigging	13"-17" for Swing Away and Standard Plug In Front Rigging			
<input type="checkbox"/>	2628020	Long	18"-24" for Easy-Out Front Rigging	15"-19" for Swing Away and Standard Plug In Front Rigging			\$ 75
 <b>Lower Leg Length</b> If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.							
Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.							



# FUZE T50n

## Order Form - USA

Prices in U.S. Dollars. January 1, 2022

### Miscellaneous Options

<input type="checkbox"/>	2630060	<b>Stealth Flip Down Abductor</b>	Requires Solid Seat Pan.	E0957 & E1028	\$ 595
<input type="checkbox"/>	2630070	<b>QR Bracket - Accessory Mount</b>		E1028	\$ 200
<input type="checkbox"/>	2630080	<b>Gel Hip Guides</b>	Reduces available Seat Width by 2".	E0956 & E1028	\$ 395
<input type="checkbox"/>	2630090	<b>Hip Belt Mount</b>	Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.	E1028	\$ 55
<input type="checkbox"/>	2630150	<b>Ventilator Tray</b>	Reduces ground clearance. Chair may be unstable when unoccupied. Minimum Seat Width is 20". Not compatible with Power Tilt, Foot Operated Wheel Locks.	E1029	\$ 295
<input type="checkbox"/>	2630152	<b>Ventilator Bracket</b>	Compatible with LTV1 (laptop) ventilator.	E1029	\$ 295
<input type="checkbox"/>	2630154	<b>Ventilator Battery Tray</b>	20" min. Seat Width. Reduces ground clearance. Fits U1 Battery only. STF restrictions.	E1029	\$ 195
<input type="checkbox"/>	2630160	<b>Clear Lap Tray (Slide On)</b>	Available for Seat Widths of 16"-22" only. Not compatible with Contour Arm Trough.	E0950	\$ 490
<input type="checkbox"/>	2630170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.		\$ 10

### On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

### Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.  
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.