



Customer and Order Information

Purchase Order *	Tag	v.0119
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in

Base Chair

<input checked="" type="checkbox"/> 1901000	Stellar GL 45° Manual Tilt Wheelchair	E1161	\$ 3,150
<input type="checkbox"/> ⓘ	Also available: the lightweight Stellar GLT with aluminum frame. See separate order form. pdgmobility.com/glt		

Weight Capacity

Chair Weight Capacity			
<input type="checkbox"/> 1902000	250 lbs Capacity	Available for 14" to 22" Seat Width.	
<input type="checkbox"/> 1902050	300 lbs Capacity	Available for 14" to 22" Seat Width.	\$ 200
<input type="checkbox"/> 1905006	450 lbs Capacity	Available for 20" to 32" Seat Width. The Stellar HD offers superior strength and performance. See separate order form.	\$ 500
User Weight Range Required for setup of tilt uplift force.			
<input type="checkbox"/> 1902430	Up to 75 lbs	<input type="checkbox"/> 1902431	75 lbs to 100 lbs
<input type="checkbox"/> 1902433	115 lbs to 150 lbs	<input type="checkbox"/> 1902434	150 lbs to 200 lbs
<input type="checkbox"/> 1902436	250 lbs to 300 lbs	<input type="checkbox"/> 1902437	300 lbs to 375 lbs
		<input type="checkbox"/> 1902432	100 lbs to 115 lbs
		<input type="checkbox"/> 1902435	200 lbs to 250 lbs
		<input type="checkbox"/> 1902438	375 lbs to 450 lbs

Frame Options

<input type="checkbox"/> 1902200	Extended Base with Rear Floor Jacks	K0108	\$ 900
Lifts Rear Wheels off ground for increased stability. Limited ground clearance. Not compatible with Attendant Wheel Lock, HD Anti-Tipper, IV Pole, O2 Holder.			
<input type="checkbox"/> 1902300	Transit Tie Down System (TTDS)		
TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Not compatible with Seat Widths greater than 22". Tested to max. Weight Capacity of 300 lbs.			

Seat Width



Overall chair width is +9.5" with 12/16" Rear Wheels and +10.5" with 20/22/24" Rear Wheels.

Standard Seat Widths			
<input type="checkbox"/> 1903014	14" Seat Width	<input type="checkbox"/> 1903015	15" Seat Width
<input type="checkbox"/> 1903017	17" Seat Width	<input type="checkbox"/> 1903018	18" Seat Width
<input type="checkbox"/> 1903020	20" Seat Width	<input type="checkbox"/> 1903021	21" Seat Width
		<input type="checkbox"/> 1903016	16" Seat Width
		<input type="checkbox"/> 1903019	19" Seat Width
		<input type="checkbox"/> 1903022	22" Seat Width
			E2201 \$ 375
Custom Seat Widths Requires 450 lbs Weight Capacity Frame. Extended delivery lead times may apply.			
<input type="checkbox"/> 1903023	23" Seat Width		E2201 \$ 375
<input type="checkbox"/> 1903024	24" Seat Width	<input type="checkbox"/> 1903025	25" Seat Width
<input type="checkbox"/> 1903026	26" Seat Width	<input type="checkbox"/> 1903027	27" Seat Width
<input type="checkbox"/> 1903028	28" Seat Width	<input type="checkbox"/> 1903029	29" Seat Width
<input type="checkbox"/> 1903031	31" Seat Width	<input type="checkbox"/> 1903030	30" Seat Width
			K0108 \$ 1,100
			K0108 \$ 1,100

Seat Depth

Standard Seat Depths (Adjustable)			
<input type="checkbox"/> 1904016	16" Seat Depth	<input type="checkbox"/> 1904017	17" Seat Depth
<input type="checkbox"/> 1904018	18" Seat Depth	<input type="checkbox"/> 1904019	19" Seat Depth
<input type="checkbox"/> 1904020	20" Seat Depth		E2203 \$ 500
Custom Seat Depths			
<input type="checkbox"/> 1904021	21" Seat Depth		E2203 \$ 525
<input type="checkbox"/> 1904022	22" Seat Depth		E2204 \$ 800

Frame Color			
PDG Premium Colors			
<input type="checkbox"/> 1905000	Tracer Red	<input type="checkbox"/> 1905001	Gloss Black Shimmer
<input type="checkbox"/> 1905003	Fine Texture Black	<input type="checkbox"/> 1905004	Orange Fusion
<input type="checkbox"/> 1905006	Blue Shimmer	<input type="checkbox"/> 1905007	Gloss White
<input type="checkbox"/> 1905009	Geyser Green	<input type="checkbox"/> 1905010	Lollypop Blue
<input type="checkbox"/> 1905012	Bubble Gum Pink	<input type="checkbox"/> 1905013	Berry Purple
<input type="checkbox"/> 1905002	Silver	<input type="checkbox"/> 1905005	Karl's Imagination
		<input type="checkbox"/> 1905008	Dark Titanium
		<input type="checkbox"/> 1905011	Lollypop Red
		<input type="checkbox"/> 1905014	Daffodil Yellow
<input type="checkbox"/> 1905200	Custom Color	Contact PDG for specifications.	
			\$ 300

Armrests			
<input type="checkbox"/> 1906000	No Armrests		
Single Post Height Adjustable Armrests			E0973 \$ 295
<input type="checkbox"/> 1906111	8"-11" Range	<input type="checkbox"/> 1906114	10"-14" Range
Dual Post Height Adjustable Armrests			E0973 \$ 295
<input type="checkbox"/> 1906211	8"-11" Range	<input type="checkbox"/> 1906214	10"-14" Range
		<input type="checkbox"/> 1906216	12"-16" Range
Fixed Height Armrests (12") Full Length Armpad only.			
<input type="checkbox"/> 1906100	Single Post	<input type="checkbox"/> 1906200	Dual Post
	Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.		
	Seat Pan to top of Armrest: _____ Does not include height of Armpad		
	Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.		
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.			
<input type="checkbox"/> 1906401	Rear	<input type="checkbox"/> 1906402	Mid
		<input type="checkbox"/> 1906403	Front
Locking Flip-up Not recommended for users with excessive strength or tone.			E0973 \$ 295
Not compatible with Fold Down Back, Clothing Guards, Back Upholstery. Available with Full Length and Desk Length Armpads.			
<input type="checkbox"/> 1906310	Set at 10"	<input type="checkbox"/> 1906311	Set at 11"
		<input type="checkbox"/> 1906312	Set at 12"
<input type="checkbox"/> 1906313	Set at 13"	<input type="checkbox"/> 1906314	Set at 14"

Armpads			
Regular Armpads			
<input type="checkbox"/> 1907100	Desk Length - 2" x 10"	<input type="checkbox"/> 1907101	Full Length - 2" x 14"
Waterfall Armpads			\$ 115
<input type="checkbox"/> 1907105	Desk Length - 2" x 10"	<input type="checkbox"/> 1907106	Full Length - 2" x 14"
Wide Armpads May interfere with wheels in some tilt configurations.			\$ 200
<input type="checkbox"/> 1907110	Desk Length - 3" x 10"	<input type="checkbox"/> 1907111	Full Length - 3" x 12"
Gel Armpads May interfere with wheels in some tilt configurations.			\$ 135
<input type="checkbox"/> 1907120	Full Length - 2" x 12"	<input type="checkbox"/> 1907121	Full Length - 3.5" x 12"
Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.			E2209 \$ 125
<input type="checkbox"/> 1907130	Right (4.5" x 18")	<input type="checkbox"/> 1907131	Left (4.5" x 18")

Armrest Options			
Armrest Wideners Move Armrests outward by 1" on each side. Do not increase seat size.			
<input type="checkbox"/> 1907210	For Single Post Armrests	May cause tilt interference at lower Seat Height and/or Seat Depth settings.	
<input type="checkbox"/> 1907220	For Dual Post Armrests	Limits tilt. Not compatible with 20", 22", or 24" Rear Wheels, Clothing Guards.	
Clothing Guards Not compatible with Flip-up Armrests, Armrest Wideners.			
<input type="checkbox"/> 1907310	For Single Post Armrests	Available for 10"-14" Height range only.	
<input type="checkbox"/> 1907320	For Dual Post Armrests		
Armrest Extensions Not compatible with Flip-up or Single Post Armrests.			\$ 75
<input type="checkbox"/> 1907400	3" Extension - Forward	<input type="checkbox"/> 1907401	3" Extension - Rearward

Seat Pan

- 1908000 **No Seat Pan**
- 1908001 **Solid Aluminum Seat Pan** E2231 \$ 135

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Subframe Plate (Low/High).

	Low							High						
	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
6" Casters	12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			12"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
8" Casters	RW	13"	14"	15"	16"	17"	18"	RW	15"	16"	17"	18"	19"	20"
	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	16"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	24"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Casters

- Urethane Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1911061 6" x 1.25" Caster size 1911081 8" x 1.25" Caster size
 - 1911063 6" x 1.50" Caster size \$ 85
 - 1911084 8" x 2.00" Caster size \$ 135
- Pneumatic Casters** 8" Casters are compatible only with 16" and 24" Rear Wheels.
- 1911161 6" x 1.25" Caster size E2214 & E2215 \$ 110
 - 1911181 8" x 1.25" Caster size E2214 & E2215 \$ 120
 - 1911184 8" x 2.00" Caster size E2214 & E2215 \$ 150
 - 1911260 **Frog Legs 6" x 1.50" Aluminum Soft-Roll Casters** \$ 275
 - 1911400 **Frog Legs Suspension Forks** Not compatible with 12" Rear Wheels. Requires 6" Soft-Roll Caster. E1015 \$ 595

Rear Wheels

- Composite Wheels**
- 1912012 12" Wheel size 1912016 16" Wheel size 1912020 20" Wheel size
 - 1912022 22" Wheel size 1912024 24" Wheel size
- Spoke Wheels** \$ 75
- 1912120 20" Wheel size 1912122 22" Wheel size 1912124 24" Wheel size
 - 1912224 **24" Heavy Duty Spoke Wheels** \$ 350
- Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.
- 1912401 Rear 1912402 Mid 1912403 Front

Axle Type

- 1913000 **Quick Release Rear Axles**
- 1913010 **Fixed Rear Axles**

Rear Tires

- 1914000 **Urethane Tires**
- 1914010 **Pneumatic Tires** May result in a finished Seat Height 1/2" greater than stated. E2211 & E2212 \$ 60
- 1914020 **Pneumatic Tires with Flat-Free Inserts** May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels. E2211 & E2213 \$ 95



Handrims		Handrims are only available with 20", 22", and 24" Rear Wheels.
<input type="checkbox"/>	1915000 No Handrims	
<input type="checkbox"/>	1915010 Standard Aluminum Handrims	
<input type="checkbox"/>	1915020 Plastic Coated Handrims	\$ 150

Wheel Locks & Brakes		
User Wheel Locks		
<input type="checkbox"/>	1916010 Push to Lock <input type="checkbox"/> 1916020 Pull to Lock	
<input type="checkbox"/>	1916030 Attendant Wheel Locks (Push to Lock)	\$ 85
<input type="checkbox"/>	1916040 User & Attendant Wheel Locks Available with 20", 22", or 24" Rear Wheels only. Not compatible with Drum Brakes.	\$ 125
<input type="checkbox"/>	1916300 Drum Brakes Dual trigger activation (right/left side activated independently). Available with 20" Rear Wheels only.	E2228 \$ 795

Wheel Lock Options		
<input type="checkbox"/>	1916100 Wheel Lock Extensions Not compatible with Attendant Wheel Lock.	E0961 \$ 60
<input type="checkbox"/>	1916200 Hill Holders Requires Pneumatic Rear Tires with Flat-Free Inserts.	E0974 \$ 450

Anti-Tippers		
<input checked="" type="checkbox"/>	1917010 Rear Anti-Tippers Required Safety Option.	E0971 \$ 75
<input type="checkbox"/>	1917020 Front Anti-Tippers	E0971 \$ 795
<input type="checkbox"/>	1917030 Heavy Duty Rear Anti-Tipper	E0971 \$ 895
Provides additional rearward stability for users with excessive movement. Not compatible with Attendant Wheel Locks, IV Pole, O2 Holder.		

Back Posts	
Back Post Style	
<input type="checkbox"/>	1918010 Push Bar Back <input type="checkbox"/> 1918020 Cane Handle Back Seat Widths of 14"-22" only. Not compatible with Headrest Mounting Plate.
Back Post Height Measured from the Seat Pan to the Handle.	
<input type="checkbox"/>	1918425 Short (22") <input type="checkbox"/> 1918431 Tall (25")
Seat to Back Angle	
<input type="checkbox"/>	1919090 Set at 90° <input type="checkbox"/> 1919098 Set at 98° <input type="checkbox"/> 1919106 Set at 106°
<input type="checkbox"/>	1919114 Set at 114° <input type="checkbox"/> 1919122 Set at 122°

Back Options		
Angle/Height Adjustable Push Bar Handle Not compatible with Headrest Mounting Plate or Cane Handle Back.		\$ 150
<input type="checkbox"/>	1919010 Button on Outside <input type="checkbox"/> 1919011 Button on Inside	
<input type="checkbox"/>	1919210 10" Backrest Extension Requires the Angle/Height Adjustable Push Bar Handle.	E0966 \$ 150
<input type="checkbox"/>	1919212 Headrest Mounting Plate	\$ 95
Not compatible with Backrest Extension, Angle/Height Adjustable Push Bar, Cane Handle Back. Not required if mounting to a solid Back.		
<input type="checkbox"/>	1919220 Fold Down Back Not compatible with Reclining Back.	
<input type="checkbox"/>	1919230 Spring Back Suspension Not compatible with Reclining Back.	\$ 105
<input type="checkbox"/>	1919240 Seating Dynamics Dynamic Rocker Back Interface	K0108 \$ 2,995
Not compatible with Fold Down Back, Spring Back Suspension, 30° Recline, Cane Handle, Standard Upholstery. Finished Back Post Height is 4" taller than selected Back Post Height.		

Back Upholstery		
<input type="checkbox"/>	1920000 No Back Upholstery	
Standard Upholstery Not available with Cane Handle Back.		E2611 / E2612 \$ 195
<input type="checkbox"/>	1920418 Short (18" Height) <input type="checkbox"/> 1920421 Regular (21" Height) <input type="checkbox"/> 1920424 Tall (24" Height)	

See next page for more options >



Tension Adjustable Upholstery

Height of Upholstery is 16". Not available with Cane Handle Back.

- 1920450 For 14" to 22" Seat Width E2611 / E2612 \$ 295
- 1920451 For 24" to 32" Seath Width E2611 / E2612 \$ 375

Position of Upholstery

Select the desired position of the top of the Upholstery relative to the Seat Pan.

		22" Back Post							
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

		25" Back Post							
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tall							<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation

Attendant Tilt only

Select trigger mounting location below. Chairs with 450 lbs Weight Capacity require two tilt triggers (L and R) for tilt activation.

- 1923011 Right Back Post 1923012 Left Back Post

User Tilt only

Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting location below.

- 1923021 Right Armrest 1923022 Left Armrest \$ 100

Attendant & User Tilt

Not available on widths greater than 22". Full Length Arm Pads only. Select trigger mounting locations below.

- 1923031 Right Back Post 1923032 Left Back Post
- 1923041 Right Armrest 1923042 Left Armrest \$ 195

- 1923110 Tilt Trigger Extension \$ 50

- 1923400 **Power Tilt** 250 lbs Weight Capacity. 45° tilt range. Includes 2 batteries for off-chair charging. pdgmobility.com/power-tilt E1002 \$ 2,295

Minimum Seat Width is 16". Power Tilt is compatible with either O2 Holder or IV Pole on the same chair, not both.

- Toggle Switch 1923421 Installed on Right side 1923422 Installed on Left side
- Button Switch 1923431 Installed on Right side 1923432 Installed on Left side
- Stealth Egg Switch 1923440 Not installed

- Optional for Egg Switch: 1923445 Swing Away Mounting Hardware Not compatible with all Headrest makes and models. E1028 \$ 625

Specify Headrest to be used: _____

Tilt Limit

Optional. If no tilt limit is selected, the full tilt range will be available.

- 1923220 Set at 20° 1923230 Set at 30° 1923240 Set at 40°

Recline

- 1923310 **Reclining Backrest (Manual)** \$ 895
- 30° maximum recline. Not compatible with Backrest Extension, Cane Handle Back, Spring Back Suspension, Fold Down Back, Tension Adjustable Upholstery.
Finished Back Post Height is 3" taller than selected Back Post Height.

Front Rigging

- 1925000 **No Front Rigging** 1925001 **Sector Blocks only**

Easy-Out™ Front Rigging

Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out

- 1925410 Fixed 70° \$ 95
 - 1925420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) \$ 295
- Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.

Swing Away Front Rigging

Legacy Design. Not compatible with Easy-Out.

- 1925130 Swing Away 60° 1925140 Swing Away 70°
- 1925200 Swing Away Elevating Legrests Includes Calf Pads. E0990 \$ 375
- 1925300 1" Offset Hangers for Swing Away \$ 75

- 1925110 **Standard Plug In 70° Front Rigging** Legacy Design. Not compatible with Easy-Out.

Footplates

- 1926010 **Composite Footplates**
- 1926020 **Aluminum Angle Adjustable Footplates** K0040 \$ 165


See next page for more options >

<input type="checkbox"/>	1926030	One Piece Fixed Footplate	Available on Plug In Front Rigging only. Requires selection of Composite or Aluminum Angle Adjustable.	\$ 350
<input type="checkbox"/>	1926040	One Piece Flip-up Footplate	Requires selection of Composite or Aluminum Angle Adjustable.	\$ 395
<input type="checkbox"/>	1926050	Interlocking Angle Adjustable Footplates	Available for Seat Widths up to 21" only.	K0040 \$ 195

Footrest/Legrest Accessories

Calf Pads Not compatible with Gel Padding.				
<input type="checkbox"/>	1927012	For Easy-Out™ Front Rigging	Includes reinforced premium Calf Pad mount.	\$ 200
<input type="checkbox"/>	1927010	For Swing Away Front Rigging		\$ 135
<input type="checkbox"/>	1927020	Gel Padding	Not compatible with Calf Pad or Easy-Out Front Rigging.	\$ 125
Gel Knee Pads				
<input type="checkbox"/>	1927032	For Easy-Out™ Angle Adjustable Front Rigging		\$ 200
<input type="checkbox"/>	1927030	For Swing Away ELR Front Rigging		\$ 100
<input type="checkbox"/>	1927040	Length Adjustable Heel Loops		E0951 \$ 50
<input type="checkbox"/>	1927050	Bodypoint Padded Calf Strap	Sizes greater than 22" are not padded.	K0038 \$ 195
Miller's Dynamic Footrest Mounts Not compatible with Easy-Out Front Rigging.				
<input type="checkbox"/>	1927060	For Standard Plug In Front Rigging		\$ 445
<input type="checkbox"/>	1927065	For Swing Away 70° Front Rigging		\$ 575


Front Rigging Extension Tubes

Regular Extension Tubes					
<input type="checkbox"/>	1928010	Short	14"-18" for Easy-Out Front Rigging 13"-17" for Swing Away and Standard Plug In Front Rigging		
<input type="checkbox"/>	1928020	Long	18"-24" for Easy-Out Front Rigging 15"-19" for Swing Away and Standard Plug In Front Rigging	\$ 75	
Extra Strength Extension Tubes (Solid) Available on Plug In Front Rigging only. Not compatible with Swing Away or Easy-Out Front Rigging.					
<input type="checkbox"/>	1928110	Standard (13"-17")	<input type="checkbox"/>	1928120 Long (15"-19")	\$ 100
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.				
	Seat Pan to Foot Plates: _____		Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.		

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Rear Wheels. Made from durable and lightweight aluminum.					
<input type="checkbox"/>	1930020	Black	<input type="checkbox"/>	1930021 Frame Color	K0065 \$ 195
<input type="checkbox"/>	1930080	Gel Hip Guides	Reduces available Seat Width by 2".	E0956 & E1028	\$ 395
<input type="checkbox"/>	1930110	IV Pole	Not compatible with HD Anti Tipper, Attendant Wheel Locks.	K0105	\$ 295
<input type="checkbox"/>	1930120	O2 Holder	Not compatible with HD Anti Tipper, Attendant Wheel Locks.	E2208	\$ 395
<input type="checkbox"/>	1930170	Adhesive Back Hook Strip, 2" x 30"	For Cushion mounting.		\$ 10

On-Chair Seating & Positioning

	Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. pdgmobility.com/seating
---	---

Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions. Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to info@pdgmobility.com or by fax to 877 323 9093.