



# STELLAR HD

## Order Form - USA

Prices in U.S. Dollars. January 1, 2020

### Customer and Order Information

Purchase Order *	Tag	v.0801
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ lbs	Height _____ ft/in <i>*required</i>

### Frame Type

<input checked="" type="checkbox"/> 2101000	<b>Stellar 45° Manual Tilt Wheelchair</b>	E1161	\$ 2,995
<input checked="" type="checkbox"/> 2101201	<b>Reinforced Frame Package (Stellar HD)</b> Heavy-gauge steel reinforcement of frame. Extended delivery lead times may apply.	K0108	\$ 1,295
	Consider the <b>Stellar Impact</b> for highly agitated users up to 250 lbs. See separate order form.		

### Frame Options

<b>Heavy Duty Base and Seat Frame</b>			
<input type="checkbox"/> 2102045	450 lbs Weight Capacity (45° of Tilt)		
<input type="checkbox"/> 2102060	600 lbs Weight Capacity (25° of Tilt)		\$ 300
<b>User Weight Range</b> Required for setup of tilt uplift force.			
<input type="checkbox"/> 2102212	Up to 250 lbs	<input type="checkbox"/> 2102308	250 lbs to 350 lbs
<input type="checkbox"/> 2102310	350 lbs to 425 lbs	<input type="checkbox"/> 2102312	425 lbs to 600 lbs

### Seat Width Overall chair width is +11"

<b>Standard Seat Widths</b>			
<input type="checkbox"/> 2103022	22" Seat Width		E2201 \$ 350
<input type="checkbox"/> 2103024	24" Seat Width	<input type="checkbox"/> 2103026	26" Seat Width E2202 \$ 550
<input type="checkbox"/> 2103028	28" Seat Width	<input type="checkbox"/> 2103030	30" Seat Width <input type="checkbox"/> 2103032
			32" Seat Width K0108 \$ 1,050

### Seat Depth

<b>Standard Seat Depths</b>			
<input type="checkbox"/> 2104018	18" Seat Depth	<input type="checkbox"/> 2104019	19" Seat Depth
<input type="checkbox"/> 2104020	20" Seat Depth	<input type="checkbox"/> 2104021	21" Seat Depth E2203 \$ 250
<input type="checkbox"/> 2104022	22" Seat Depth		E2204 \$ 450

### Frame Color

<b>PDG Premium Colors</b>			
<input type="checkbox"/> 2105000	Tracer Red	<input type="checkbox"/> 2105001	Gloss Black Shimmer
<input type="checkbox"/> 2105003	Fine Texture Black	<input type="checkbox"/> 2105004	Orange Fusion
<input type="checkbox"/> 2105006	Blue Shimmer	<input type="checkbox"/> 2105007	Gloss White
<input type="checkbox"/> 2105009	Geyser Green	<input type="checkbox"/> 2105010	Lollypop Blue
<input type="checkbox"/> 2105012	Bubble Gum Pink	<input type="checkbox"/> 2105013	Berry Purple
<input type="checkbox"/> 2105002	Silver	<input type="checkbox"/> 2105005	Karl's Imagination
		<input type="checkbox"/> 2105008	Dark Titanium
		<input type="checkbox"/> 2105011	Lollypop Red
		<input type="checkbox"/> 2105014	Daffodil Yellow
<input type="checkbox"/> 2105200	<b>Custom Color</b> Contact PDG for specifications.		\$ 300

### Armrest Type

<input type="checkbox"/> 2106000	<b>No Armrests</b>		
<input type="checkbox"/> 2106214	<b>Dual Post Height Adjustable Armrests (10"-14")</b>	E0973	\$ 295

### Armpads & Armrest Options

#### Regular Armpads

- 2107100 Desk Length - 2" x 10"  2107101 Full Length - 2" x 14"

#### Wide Armpads May interfere with wheels in some tilt configurations.

\$ 200

- 2107110 Desk Length - 3" x 10"  2107111 Full Length - 3" x 12"

#### Gel Armpads May interfere with wheels in some tilt configurations.

\$ 135

- 2107120 Full Length - 2" x 12"  2107121 Full Length - 3.5" x 12"

#### Contour Arm Trough with Hand Pad May interfere with wheels in some tilt configurations.

*E2209* \$ 125

- 2107130 Right (4.5" x 18")  2107131 Left (4.5" x 18")

- 2107320 **Clothing Guards** \$ 165

### Seat Pan

- 2108001 **Solid Aluminum Seat Pan** *E2231* \$ 155

### Seat Height

#### Seat to Floor Height


- 2110016 Set at 16" STF  2110017 Set at 17" STF  2110018 Set at 18" STF

### Casters

- 2111050 **5" Casters with Heavy Duty Steel Fork**

### Rear Wheels

- 2112224 **24" Heavy Duty Spoke Wheels**

 **Rear Wheel Position** Indicate Rear Wheel positioning preference. If left blank, the default setting will be used.  
In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

- 2112401 Rear  2112402 Mid  2112403 Front

### Axle Type

- 2113010 **Fixed Rear Axles**

### Rear Tires

- 2114000 **Urethane Tires**

### Hand Rims

- 2115000 **No Handrims**
- 2115010 **Standard Aluminum Handrims**
- 2115020 **Plastic Coated Handrims** \$ 250

### Wheel Locks

#### User Wheel Locks

- 2116010 Push to Lock  2116020 Pull to Lock

- 2116100 **Wheel Lock Extensions** *E0961* \$ 60

### Anti-Tippers

- 2117010 **Rear Anti-Tippers** Required Safety Option. *E0971* \$ 60

### Back Post Style

2118010 **Push Bar Back**

**Back Post Height** Please note: We have updated the descriptions/measurements below. The actual back post heights remain unchanged.

2118425 Short (22") *Previously 25"*       2118431 Tall (25") *Previously 31"*      Measured from the Seat Pan to the Handle.

### Back Options

2119090 **Seat to Back Angle set at 90°**

2119212 **Headrest Mounting Plate** Not required if mounting to a solid Back. \$ 95

2119240 **Seating Dynamics Dynamic Rocker Back Interface** K0108 \$ 2,995  
 Not compatible with 30° Recline. Requires Adjustable Tension Upholstery or rigid Backrest. Finished Back Post Height 3" taller than selected Back Post Height.

### Back Upholstery

2120000 **No Back Upholstery**

**Tension Adjustable Upholstery** E2611 / E2612 \$ 375

2120421 Regular (18" Height)       2120424 Tall (22" Height)

**Position of Upholstery** Position of top of Upholstery relative to the Seat Pan. Circle below of select from the drop-down menu on the right.

Upholstery	22" Back Post							
Regular	17"	18"	19"	20"	21"	22"		
Tall				20"	21"	22"		

Upholstery	25" Back Post									
Regular	17"	18"	19"	20"	21"	22"	23"	24"	25"	
Tall					21"	22"	23"	24"	25"	

Position \_\_\_\_\_

### Tilt/Recline Options

**Attendant Tilt only** Select trigger mounting location below.

2123011 Right Back Post       2123012 Left Back Post

**Tilt Limit** Optional. If no tilt limit is selected, the full tilt range will be available. Maximum Tilt Limit for 600 lbs Weight Capacity is 25°.

2123220 Set at 20°       2123230 Set at 30°       2123240 Set at 40°

2123310 **Reclining Backrest (Manual)** 30° maximum recline. Not compatible with Tension Adjustable Upholstery. E1225 \$ 895

### Front Rigging

2125000 **No Front Rigging**       2125001 **Sector Blocks only**

2125130 **Swing Away 60° Front Rigging**

2125200 **Swing Away Elevating Legrests** Includes Calf Pads. E0990 \$ 375

### Footplates

2126020 **Aluminum Angle Adjustable Footplates** K0040 \$ 175

### Footrest/Legrest Accessories

2127030 **Swing Away ELR Gel Knee Pads** \$ 100

2127040 **Length Adjustable Heel Loops** E0951 \$ 50

### Front Rigging Extension Tubes

2128010 **Regular Extension Tubes (13"-17")**

**Lower Leg Length** If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used.

Seat Pan to Foot Plates: \_\_\_\_\_ Measured along the length of the Front Rigging.



# STELLAR HD


**Order Form - USA**

Prices in U.S. Dollars. January 1, 2020

## Miscellaneous Options

<input type="checkbox"/>	2130120	<b>O2 Holder</b>		<i>E2208</i>	\$	395
<input type="checkbox"/>	2130170	<b>Adhesive Back Hook Strip, 2" x 30"</b>	For Cushion mounting.		\$	10

## On-Chair Seating & Positioning

 Refer to the separate Seating & Positioning order form for PDG's on-chair options. Standard PDG discounts apply. [pdgmobility.com/seating](http://pdgmobility.com/seating)

## Order Notes

HCPCS codes provided are based on code verification by PDAC or PDG's interpretation of existing code definitions.  
Use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Please send your completed order form by email to [info@pdgmobility.com](mailto:info@pdgmobility.com) or by fax to 877 323 9093.